



designed for scientists



SY clamp 150

/// Data Sheet

Vessel clamp for the EasySyn reactor system.

Flange size: DN 150.

Suitable for reactor vessels with volumes of 2000 ml (SY 2000 and SY 2000 D) and 5000 ml (SY 5000 and SY 5000 D).



designed for scientists

Technical Data

Height (room dimension) [mm]	25
Width (room dimension) [mm]	243
Depth (room dimension) [mm]	257
Flange dimension	DN150
Flange collar diameter [mm]	159