

Lovibond® Water Testing Tintometer® Group



Nessleriser 2250 with stand

More than 400 different test discs available



- More than 400 different test discs available
- Compensation for colored and turbid samples
- Guaranteed constancy of the colored glass standards
- Integrated prism

Part Number: 172250

Usable in different applications

The Nessleriser 2250 is versatile: it is used in swimming pools as well as in research centers, universities or general drinking water treatment.

Daylight lamp

So that you can always measure safely and stably under unfavorable visibility conditions, this package comes with a mains daylight unit, so the lighting conditions always remain the same no matter where you are measuring.

Integrated prism

The integrated prism visually brings the glass standards of the color disc and the coloured sample into view.

Industry

Chemical Industry | Food and Beverage Industry | Industries Others | Municipalities | NGO | Oil and Gas | Pharmaceutical Industry | Power and Energy

Nessleriser 2250 with stand

The Nessleriser 2250 is a versatile colorimetric system for color measurement and water analysis. Comfortable handling, no compromise on accuracy and reproducibility: The Nessleriser 2250 can be used for many applications. Swimming pools, research or drinking water treatment are just a few examples.

Technical Data

Dimensions "

Accessories

Title	Part Number
Daylight cabinet, mains operated	171020
Stand for Upgrade Kit	172180
Pair Nessler tubes, 113 mm path length, 50 ml, AF 306	353060
Plunger for Nessler tube AF 306 and AF 306/P	353070
Pair Nessler tubes, 113 mm path length with plungers, 50 ml, AF 306/P	353080
AF 315 1 Special Nessler tube (determination of Oxygen with color disc NOE)	353150
Pair Nessler tubes, 250 mm path length with plungers, DB 420	354200
Plunger for Nessler tube DB 420	354229

Tintometer GmbH
Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
Germany

The Tintometer Limited
Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
sales@lovibond.uk
www.lovibond.com
UK

Tintometer China
Room 1001, China Life Tower
16 Chaoyangmenwai Avenue,
Beijing, 100020
Tel.: +86 10 85251111 App. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
China

Tintometer South East Asia
Unit B-3-12, BBT One Boulevard,
Lebuh Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Malaysia

Tintometer Brazil
Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@tintometer.com.br
www.lovibond.com.br
Brazil

Tintometer Inc.
6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.com
USA

Tintometer India Pvt. Ltd.
Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
India

Tintometer Spain
Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
Spain

Technical changes without notice
Printed in Germany

Lovibond® and Tintometer® are Trademarks of the Tintometer Group of Companies