

3500326 Clino® One Way Ultra PW/I Single-Use Mop for Cleanroom

Individualization
on request!

Product description

Gamma irradiated disposable mop for sterile GMP areas.

- o High area performance of up to 25 m² thanks to absorbent material
- o Excellent cleaning performance thanks to special structure
- o Particularly smooth-running and effective thanks to soft bristles
- o Very high-quality, cleanroom-compatible processing



ISO	3	4	5	6	7	8
GMP			A/B		C	D

Product specifications

Article name Article number	Clino® One Way Ultra PW/I 3500326
Dimensions Weight	40 x 12 cm 60 g
Working material (colour)	100 % polyester
Absorbent material	100 % polyester
Backing material (colour)	100 % polyester
Edge material (colour)	100 % polyester
Holding system	Turned Over Pocket
Pocket System material (colour)	70% polyester (microfibre) 30% polyester
Surface coverage	25 m ²
Chemical compatibility	Resistant to a wide range of commercially available cleaning chemicals, with the exception of acids, chlorinated compounds, alkalis and bleaching agents
Shelf life	24 months post-production
Re-use	for single use in cleanrooms. Can be washed and reused in grey areas (non-cleanroom).
Storage	The general storage conditions apply: Protect from harmful influences (liquids, sun or UV light, soiling, pressure) and store at room temperature.
Pieces per PackingUnit Selling unit	5 Pieces 1 PU
Packaging note	The packaging size and the quantity per carton may vary depending on the order quantity and the scope of additional products ordered.
Customs tariff number	63071010
Country of origin	China
EAN-Number	4065237021132
Compatible accessories	Our sales department will be happy to advise you on suitable accessories

Compatible methods Cleaning area

- o MopScoop®
- o MopFloat®
- o EasyMop®
- o Flat press Ringo®



Industry Sectors

Preparation Sterility

This classification is a non-binding recommendation. The user must check the individual suitability and compatibility for his individual area of application. The condition of the material must be checked regularly; damaged products must be replaced immediately. All information without guarantee. Subject to changes and deviations.