

Werksbescheinigung / Declaration of compliance with the order

gemäß / in accordance with DIN EN 10204 - 2.1

Analysemühle / Analytical mill:

Modell / Model:

A 11 B

Produktberührende Materialien / Parts in contact with medium:

Mahlkammer / Mill chamber:	Edelstahl / Stainless steel: AISI 316 L
O-Ring-Dichtung / O-ring seal:	Kunststoff / Plastic: PTFE
Motorwelle / Motor Shaft:	Edelstahl / Stainless steel: AISI 431
Mutter / Nut:	Edelstahl / Stainless steel: AISI 304
Schläger / Hammer: (A11.1)	Edelstahl / Stainless steel: AISI 402
Schläger / Hammer: (A11.2, A11.3)	Edelstahl / Stainless steel: AISI 440B
Schläger / Hammer: (A11.6)	Titan-Legierung / Titanium alloy: UNS R56400 (lasergeschweißt / laser-welded)

Staufen, 07.04.2015

Severin Studer
Corporate Director
Development

Hans-Urs Eckerle
Normenstelle /
Standard Department

Dieses Dokument wurde elektronisch erstellt und ist ohne Unterschrift gültig. (DIN EN 10204, Abschnitt 5) /
This document was created electronically and is valid without signature. (DIN EN 10204, Chapter 5)



designed
to work perfectly

Brazil	India	Malaysia	China	USA	Germany
IKA® do Brasil Rua Matão, 122 Campinas-SP CEP 13030-300, Brazil South and Central America Tel. +55 19 3772-9600 Fax +55 19 3772-9601 info@ika.net.br	IKA® India Private Limited 814/475, Survey No. 129/1 Mysore Road, Kengeri, 560060 Bangalore, Karnataka India Tel. +91 80 26253900 Fax +91 80 26253901 info@ika.in	IKA® Works (Asia) Sdn Bhd Lot PT6445, Jalan Industri 3/4 Rawang Integrated Industrial Park 48000, Rawang, Selangor, Malaysia Tel. +60 3 60995679 Fax +60 3 6092 0240 sales@ika.my	IKA® Works Guangzhou 173 – 175 Friendship Road Guangzhou Economic and Technological Development District 510730 Guangzhou, China Tel. +86 20 8222-6771 Fax +86 20 8222-6776 info@ika.cn	IKA® Works, Inc. 2635 Northchase Pkwy SE Wilmington, NC 28405-7419 USA Tel. +1 910 452-7059 Fax +1 910 452-7693 sales@ika.net	IKA®-Werke GmbH & Co.KG Janke & Kunkel-Str. 10 79219 Staufen Germany Tel. +49 7633 831-0 Fax +49 7633 831-98 sales@ika.de