

LABORHITS

Cell Culture Flask DURAN® TILT

The smart, reusable bottle for cell culture media preparation



Features

- The innovative design enables both, good aseptic and ergonomic practice.
- The GL 56 mouth has a 100 % larger area than the standard GL 45 mouth.
- It easily accommodates even the largest serological pipettes (100 ml), and helps to avoid contamination
- Ergonomically shaped GL 56 screw cap included, made from polypropylene
- Vacuum resistant
- Space saving shape for fridge storage and warming media in a water bath
- Tilted volume indicators in both positions and both sides
- Manufactured from DURAN® borosilicate 3.3

DURAN® TILT is designed for making work in biosafety cabinets and clean hoods safer and easier. The handling of traditional cell culture media bottles in biosafety cabinets and hoods is inconvenient and time-consuming. The innovative features of the DURAN® TILT bottle change everything. The TILT bottle has two positions: upright for filter sterilization or storage, and tilted at 45° for pipetting.

Type	Thread	Height with cap mm	PK	Cat. No.
Cell Culture Flasks DURAN® TILT, incl. screw cap	GL 56	151	1	6.270 514

Accessory available on request