

LABORHITS



Ergonomic GL 45 cap

Retrace Code in cap

4 bottles per pack

Colourful Bottle Tags

Ergonomic bottle shape

Retrace Code on bottle

Cap label

Nominal volume

Bottle label

Proven DURAN® glass

NEW!

DURAN® YOUTILITY Laboratory Bottle system

Complete DURAN® YOUTILITY bottle sets GL 45 and DURAN® YOUTILITY Screw Caps, Pouring Rings, Labels and Bottle Tag selection:

Capacity ml	Ext. diam. mm	Height	Neck thread GL	Colour	PK	Cat. No.
125	55	124	45	clear	4	9.071 679
250	66	158	45	clear	4	9.071 678
500	78	193	45	clear	4	9.071 677
1000	93	253	45	clear	4	9.071 676
125	55	124	45	amber	4	9.071 672
250	66	158	45	amber	4	9.071 673
500	78	193	45	amber	4	9.071 674
1000	93	253	45	amber	4	9.071 675
200 bottle and 200 cap labels in a dispenser box, self-adhesive, white, PES					1	9.209 140