










SELECTION GUIDE

You store chemical products		We propose you		Range	You are in compliance with	Pages
flammables 	Safety cabinets with fire resistance	Cabinets type 105 min	7.90	EN 14470-1	10	
		Cabinets type 60 min	7.60		16	
		Cabinets type 30 min	7030		18	
		Cabinets type 15 min	7015		18	
	Safety cabinets	FM cabinets	2	FM - NFPA - OSHA	28	
		Safety cabinets with extinguisher	1	-	27	
		Safety storage cabinets	3/5	NF EN ISO 1182	30	
		Secure cabinets	14.E	-	43	
		Safety shelvings	Safety shelvings	R	-	90
		Retention tanks	Retention tanks	BR	-	87, 89
Safety cans	Stainless steel cans	S	EPA - UFC	88		
	Steel cans	J	FM - NFPA - OSHA	95		
corrosives 	Safety cabinets built in steel	FM cabinets	4	FM - NFPA - OSHA	35	
		Ventilation cabinets	11.A	-	40	
	High-corrosive products safety cabinets	Corrosion resistant cabinets	8.A & 8.B	-	36, 38	
		Polyethylene cabinets	10	-	37	
	Safety shelvings	Cabinets type 30 min	7030.M	EN 14470-1	39	
	Safety shelvings	Safety shelvings	R	-	90	
	Retention tanks	Retention tanks	BR	-	87, 89	
Safety cans	HPDE cans	J	EPA - UFC	88		
	Steel cans	J	FM	94		
toxics / noxious 	Safety cabinets	Ventilation cabinets	13	-	51	
		Safety cabinets for drums	14.F	-	96	
		Filtering ventilation safety cabinets	12.X	NF X 15-211*	78	
		Filtering ventilation secure cabinets	14.X	-	80	
	Safety shelvings	Safety shelvings	R	-	90	
	Retention tanks	Retention tanks	BR	-	87, 89	
different chemical families (flammables, corrosives, toxics, noxious...) 	Safety cabinets	Under-bench cabinets	6	-	54	
		Multirisk and moduable secure cabinets	14	-	44	
		Compartments cabinets	11.A	-	50	
		Secure cabinets	14.L	-	48	
		Secure cabinets with drawers	14.T	-	49	
		Secure cabinets	14.E	-	43	
		Secure cabinets with glazed doors	18	-	47	
		2 compartments secure cabinets	1.M	-	46	
		Multirisk cabinets type 105 min	7.90M	-	14	
		Multirisk cabinets type 90, 60, 30 and 15 min	11.B	EN 14470-1	52	
		Multirisk cabinets type 20 min	7030M	-	20	
		Filtering ventilation safety cabinets	12.X	NF X 15-211*	78	
		Filtering ventilation secure cabinets	14.X	-	80	
pesticides 	Safety cabinets	Safety cabinets	16.E	-	58	
		Secure cabinets	16.L	-	57	
		Ready to build secure cabinets	16.K	-	58	
	Safety shelvings	HPDE safety shelvings	R	-	90	
		Steel safety shelvings	R	EPA - UFC	90	
Retention tanks	Retention tanks	BR	-	87, 89		
Gas cylinders	Safety cabinets	Antifire safety cabinets	7.90BG	EN 14470-2	22	
filing/ document 	Safe for documents	Fire resistant filing/document cabinets	AD76 & AD79	NT Fire 017	102, 103	
		Safe for documents	AD3 & AD4	NF EN ISO 1182	105, 106	
You need to collect your wastes						
flammables	Safety cabinets	Cabinets type 105 min for recovery of waste	7.90ARD	EN 11470-1	24	
You use chemical products						
flammables, corrosives, toxics, noxious... 	Ductless fume hoods	Labopur ductless fume hoods	H	NF X 15-211 - Classe 2	60, 77	
		PMSF - Ductless fume hood and storage cabinets with recirculating air	H	-	68	
	Laminar flow clean benches	Laminar flow clean benches	H	EN 12469 - ISO 14644-1 classe 3	73	
	Biological safety cabinets	Biological safety cabinets	H	-	74	
You would like to vent your cabinet						
ventilation	Without active charcoal filter	Ventilation boxes	V	-	83	
		ATEX ventilation boxes	V	ATEX	84	
	With active charcoal filter	Recycling air box Labopur	H	NF X 15-211*	82	
		Filtration boxes	V	-	83	
You need equipments						
fire protection 	Safety waste bins	Safety waste bins	P & J	TNO - IBBC - NVBB - FM	99, 100	
	Fire blankets	Fire blankets	A	EN 1869-1997	100	
	Safety ashtrays	Safety ashtrays	C	-	101	
first-aid 	Safety combination units, showers and eye-washers	Combination units	DL	EN 15154-1, EN 15154-2	110, 111	
		Showers connected to water network	DL	EN 15154-1	114, 115	
		Eye-washers connected to water network	DL	EN 15154-2	112, 113	
		Portable shower stations	DL	-	117	
		Autonomous eye-washers	DL	-	117	
		Medicine cabinets	Medicine cabinets	19	-	118

*Active coal filters absorption's capacity and filtration system are tested and certified NF X 15-211 - may 2009.

The manufacturer keeps the right to modify technical characteristics without notice. 01/2014

Flammable cabinets

Corrosive cabinets

Toxic cabinets and multirisk

Pesticide range

Fume hoods, filtering cabinets - Ventilation

Retention and Cans

Anti-fire equipments and Files cabinets

Showers and first-aid equipments