

Key features

- Conical shaped heating element directs radiant heat to the top cavity
- Top cowl deflects heat away from your hand
- Air circulation from the vented housing keeps the base cool enough to hold during operation
- Burner consumes only 400W of power
- Recommended controller, MC5
- One year warranty parts and labour



Electric Bunsen

Electrothermal Electric Bunsens combine the advantages of a regular gas burner with the clean, easy operation of our electric heaters. Corrosion-resistant and cool-to-the-touch base. Radiation from the heater is directed upwards to a focal point. The Electric Bunsen is ideal for heating test tubes, crucibles, small flasks and beakers, regardless of their shape.

MC5 Controller

The MC5 Controller has been designed to provide a comprehensive answer to control the heating of resistive loads such as Electric Bunsen, Heating Mantles, Heating Tapes and Cords for bench top operation delivering up to a maximum of 800 watts. A rod support clamp is provided at the rear of the controller to take a standard 12.5mm diameter rod. An accessory extension mains lead is available where remote operation is required e.g. in a fume extraction unit.

Ordering information

Bunsen Without Controller	Electrical supply	Controller	Electrical supply
BA6101	230V 50/ 60Hz, 480W	MC5	230V
BA6101X1	115v 50/ 60Hz, 430W	MC5X1	115V, US Plug
Bunsen With Controller	Electrical supply		
BA6101/C	230v 50/ 60Hz, 480W		
BA6101/CX1	115v 50/ 60Hz, 430W		



For more information, please visit: www.electrothermal.com

Follow Electrothermal on Twitter! www.twitter.com/electrothermal

