

25% off until September 30, 2023



It's Time to...

Re-stock your lab with your favorite consumables now, and you can benefit from more for your budget. In stock and available for quick delivery!

Until September 30, 2023 discover the offers*;

- > Dive deeper with our high-performance 96- and 384-well deepwell plates
- > Increase purity with unbeatable Eppendorf Combitips advanced®
- > Get done in a snap, with the innovative SnapTec® cap

Your Distributor:



Conor Sweeney
 Technical Sales Consultant
E: csweeney@masontec.ie
M: +353 87 825 3615
masontechnology.ie




*No cash redemption. Errors and mistakes reserved. Promotion is valid in the Eppendorf eShop from April 1, 2023 until June 30, 2023. You can use the promotion code as many times as you like. Cannot be combined with any other offer(s). Promotion is valid for customers in Europe only. Shipping costs may apply. Eppendorf®, the Eppendorf Brand Design, Biopur®, ep Dualfilter T.I.P.S.®, epMotion®, Eppli®, g-Safe®, Mastercycler®, Safe-Lock Tubes®, twin.tec® and Xplorer® are registered trademarks of Eppendorf SE, Hamburg, Germany. Patents are listed on www.eppendorf.com/ip. All rights reserved, including graphics and images. Copyright © 2023 by Eppendorf SE

It's time to dive deeper with our high-performance 96- and 384-well deepwell plates

Eppendorf Deepwell Plates offer high-performance in all applications both manual and automated.

- > RecoverMax® well geometry maximizes sample recovery
- > g-Safe® provides exceptional robustness and centrifugation stability up to $6,000 \times g$
- > OptiTrack® matrix enables you to find your sample up to 30% faster



Until September 30, 2023 get 25% off your order.

Deepwell Plate 96/2000 µL, wells colorless	Order No.
Sterile, White, 20 Plates (5 Bags × 4 Plates)	0030 502 302
Sterile, Yellow, 20 Plates (5 Bags × 4 Plates)	0030 502 310
Sterile, Green, 20 Plates (5 Bags × 4 Plates)	0030 502 337
Sterile, Blue, 20 Plates (5 Bags × 4 Plates)	0030 502 345
Sterile, White, 80 Plates (10 Bags × 8 Plates)	0030 506 308
PCR clean, White, 20 Plates (5 Bags × 4 Plates)	0030 501 306
PCR clean, Yellow, 20 Plates (5 Bags × 4 Plates)	0030 501 314
PCR clean, Green, 20 Plates (5 Bags × 4 Plates)	0030 501 330
PCR clean, Blue, 20 Plates (5 Bags × 4 Plates)	0030 501 349
PCR clean, White, 80 Plates (10 Bags × 8 Plates)	0030 505 301

Deepwell Plate 96/500 µL, wells colorless	Order No.
Sterile, White, 40 Plates (5 Bags × 8 Plates)	0030 502 108
Sterile, Yellow, 40 Plates (5 Bags × 8 Plates)	0030 502 116
Sterile, Green, 40 Plates (5 Bags × 8 Plates)	0030 502 132
Sterile, Blue, 40 Plates (5 Bags × 8 Plates)	0030 502 140
Sterile, White, 120 Plates (10 Bags × 12 Plates)	0030 506 103
PCR clean, White, 40 Plates (5 Bags × 8 Plates)	0030 501 101
PCR clean, Yellow, 40 Plates (5 Bags × 8 Plates)	0030 501 110
PCR clean, Green, 40 Plates (5 Bags × 8 Plates)	0030 501 136
PCR clean, Blue, 40 Plates (5 Bags × 8 Plates)	0030 501 144
PCR clean, White, 120 Plates (10 Bags × 12 Plates)	0030 505 107

Deepwell Plate 96/1000 µL, wells colorless	Order No.
Sterile, White, 20 Plates (5 Bags × 4 Plates)	0030 502 205
Sterile, Yellow, 20 Plates (5 Bags × 4 Plates)	0030 502 213
Sterile, Green, 20 Plates (5 Bags × 4 Plates)	0030 502 230
Sterile, Blue, 20 Plates (5 Bags × 4 Plates)	0030 502 248
Sterile, White, 80 Plates (10 Bags × 8 Plates)	0030 506 200
PCR clean, White, 20 Plates (5 Bags × 4 Plates)	0030 501 209
PCR clean, Yellow, 20 Plates (5 Bags × 4 Plates)	0030 501 217
PCR clean, Green, 20 Plates (5 Bags × 4 Plates)	0030 501 233
PCR clean, Blue, 20 Plates (5 Bags × 4 Plates)	0030 501 241
PCR clean, White, 80 Plates (10 Bags × 8 Plates)	0030 505 204

Deepwell Plate 384/200 µL, wells colorless	Order No.
Sterile, white, 40 plates (5 bags × 8 plates)	0030 522 109
Sterile, white, 120 plates (10 bags × 12 plates)	0030 526 104
PCR clean, white, 40 plates (5 bags × 8 plates)	0030 521 102


Deepwell Plate, DNA LoBind®, wells colorless	Order No.
96/1,000 µL, PCR clean, white, 20 plates (5 bags × 4 plates)	0030 503 201
96/1,000 µL, PCR clean, blue, 20 plates (5 bags × 4 plates)	0030 503 244
96/1,000 µL, PCR clean, white, 80 plates (10 bags × 8 plates)	0030 507 207
96/500 µL, PCR clean, white, 40 plates (5 bags × 8 plates)	0030 503 104
96/500 µL, PCR clean, blue, 40 plates (5 bags × 8 plates)	0030 503 147
96/500 µL, PCR clean, white, 120 plates (10 bags × 12 plates)	0030 507 100
384/200 µL, PCR clean, white, 40 plates (5 bags × 8 plates)	0030 523 105

Deepwell Plate, Protein LoBind®, wells colorless	Order No.
96/2000 µL, PCR clean, white, 20 plates (5 bags × 4 plates)	0030 504 305
96/1000 µL, PCR clean, white, 20 plates (5 bags × 4 plates)	0030 504 208
96/1000 µL, PCR clean, yellow, 20 plates (5 bags × 4 plates)	0030 504 216
96/1000 µL, PCR clean, white, 80 plates (10 bags × 8 plates)	0030 508 203
96/500 µL, PCR clean, white, 40 plates (5 bags × 8 plates)	0030 504 100
96/500 µL, PCR clean, yellow, 40 plates (5 bags × 8 plates)	0030 504 119
96/500 µL, PCR clean, white, 120 plates (10 bags × 12 plates)	0030 508 106
384/200 µL, PCR clean, white, 40 plates (5 bags × 8 plates)	0030 524 101

It's time to increase purity with unbeatable Eppendorf Combitips® advanced

Optimize your workflows involving plates or large series of other vessels with Combitips advanced dispenser tips. They work as a perfect team with the Multipipette®; thanks to positive displacement even viscous or volatile liquids are no problem.

Until September 30, 2023 get 25% off your order.



Combitips® advanced PCR clean	Order No.
0.1 mL, white, colorless tips, 100 pcs. (4 reclosable bags × 25 pcs.)	0030 089 766
0.2 mL, light blue, colorless tips, 100 pcs. (4 reclosable bags × 25 pcs.)	0030 089 774
0.5 mL, violet, colorless tips, 100 pcs. (4 reclosable bags × 25 pcs.)	0030 089 782
1.0 mL, yellow, colorless tips, 100 pcs. (4 reclosable bags × 25 pcs.)	0030 089 790
2.5 mL, green, colorless tips, 100 pcs. (4 reclosable bags × 25 pcs.)	0030 089 804
5.0 mL, blue, colorless tips, 100 pcs. (4 reclosable bags × 25 pcs.)	0030 089 812

In addition to this offer, these tips are available in volumes up to 50 mL and in a range of purity grades. Visit our website for more information.

It's time to snap the cap with our innovative SnapTec® cap

SnapTec caps are firmly fixed to your tube, allowing single-handed operation for quick sample handling and avoiding the contamination risk of mixing lids.

- > Straightforward to integrate into your lab
- > High centrifugation stability (up to 17,000 × g)
- > Tight seal from -86 °C to 100 °C

Until September 30, 2023 get 25% off your order.



Eppendorf Conical Tubes, 25 mL, with SnapTec® cap	Order No.
Starter Pack, 4 Single Tube Stands, 6 adapters for rotors with bore for 50 mL conical tubes, PCR clean, colorless, 200 pcs., 5 bags of 40 tubes each	0030 118 502
Eppendorf Quality®, 200 pcs., 5 bags of 40 tubes each	0030 118 405
PCR Clean, 200 pcs., 5 bags of 40 tubes each	0030 118 413
Sterile, pyrogen-, DNase-, RNase-, and DNA-free, 150 pcs., 6 bags of 25 tubes each	0030 118 421
Sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free, amber (light protection), 150 pcs., 6 bags of 25 tubes each	0030 118 430
DNA LoBind® Tubes 25 mL with SnapTec® cap, PCR clean, colorless, 200 pcs., 5 bags of 40 tubes each	0030 108 523

Eppendorf Conical Tubes SnapTec® 50	Order No.
Eppendorf Quality®, 50 mL, 240 pcs., 8 bags of 30 tubes each	0030 118 650
PCR clean, 50mL 240 pcs., 8 bags of 30 tubes each	0030 118 669
Sterile, pyrogen-, DNase-, RNase-, and DNA-free, 50 mL, 200 pcs., 10 bags of 20 tubes each	0030 118 677

Visit our eShop and filter by cap type to find your SnapTec tube.