

Combination of epT.I.P.S.® Pipette Tips with Eppendorf Research® 3 neo Pipettes

Eppendorf Research® 3 neo	0.1 – 10 µL, dark gray	0.1 – 10 µL M, medium gray	0.1 – 20 µL, medium gray	0.5 – 20 µL L, light gray	2 – 20 µL, yellow	2 – 100 µL, yellow	2 – 200 µL, yellow	20 – 300 µL, orange	50 – 1,000 µL, blue	50 – 1,250 µL, green	50 – 1,250 µL L, dark green	0.1 – 2.5 mL, red	0.1 – 5 mL, violet	0.2 – 5 mL L, violet	0.5 – 10 mL, turquoise	0.5 – 10 mL L, turquoise
epT.I.P.S.® BioBased	■		■	■			■	■	■	■	■	■	■	■		
epT.I.P.S.® Standard	■		■	■			■	■	■	■	■	■	■	■	■	■
epT.I.P.S.® BioBased Reloads			■		■	■	■	■	■	■	■					
ep Dualfilter T.I.P.S.®	■	■		■	■	■	■	■	■		■	■	■	■		■
ep Dualfilter T.I.P.S.® BioBased Reloads		■		■	■	■	■	■	■		■					
ep Dualfilter T.I.P.S.® SealMax®	■			■		■	■	■	■							
ep Dualfilter T.I.P.S.® SealMax® BioBased Reloads				■		■	■	■	■							
Variable volume, 1-channel pipettes																
0.1 – 2.5 µL	■	■	●	●												
0.5 – 10 µL	■	●	■	■	●											
1 – 20 µL	■	◆	◆	●	■											
1 – 20 µL	■				■	■	■	●								
5 – 100 µL	■				◆	■	■	●								
10 – 200 µL	■				◆	◆	■	●								
15 – 300 µL	■				◆	◆	◆	■								
50 – 1,000 µL	■								■	●	●					
0.1 – 2 mL	■											■				
0.25 – 5 mL	■												■	●		
0.5 – 10 mL	■														■	●
■ Test tip (other tips may require adjustment) ● Full volume range ◆ Limited volume range																

Your local distributor: www.eppendorf.com/contact
Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com

Eppendorf®, the Eppendorf Brand Design, epT.I.P.S.®, ep Dualfilter T.I.P.S.®, SealMax® and Eppendorf Research® are registered trademarks of Eppendorf SE, Germany.
All rights reserved, including graphics and images. Copyright © 2025 by Eppendorf SE. Order No: AA01 050 120 – EN1/0825/MCP



For more information go to
www.eppendorf.link/eptips

Combination of epT.I.P.S.® Pipette Tips with Eppendorf Research® plus Pipettes

Eppendorf Research® plus/ Eppendorf Research® plus Move It®	0.1 – 10 µL, dark gray	0.1 – 10 µL M, medium gray	0.1 – 20 µL, medium gray	0.5 – 20 µL L, light gray	0.1 – 20 µL, light pink	2 – 20 µL, yellow	2 – 100 µL, yellow	5 – 100 µL, light yellow	2 – 200 µL, yellow	20 – 300 µL, orange	50 – 1,000 µL, blue	50 – 1,250 µL, green	50 – 1,250 µL L, dark green	0.1 – 2.5 mL, red	0.1 – 5 mL, violet	0.2 – 5 mL L, violet	0.5 – 10 mL, turquoise	0.5 – 10 mL L, turquoise
epT.I.P.S.® BioBased	■		■	■					■	■	■	■	■	■	■	■		
epT.I.P.S.® Standard	■		■	■	■			■	■	■	■	■	■	■	■	■	■	■
epT.I.P.S.® BioBased Reloads			■			■	■		■	■	■	■	■					
ep Dualfilter T.I.P.S.®	■	■		■	■	■	■	■	■	■	■		■	■	■	■		■
ep Dualfilter T.I.P.S.® BioBased Reloads	■	■		■		■	■		■	■	■		■					
ep Dualfilter T.I.P.S.® SealMax®	■			■			■		■	■	■							
ep Dualfilter T.I.P.S.® SealMax® BioBased Reloads				■			■		■	■	■							
Fixed volume, 1-channel pipettes																		
10 µL	■	●	■	■														
20 µL	■			●														
10 µL	■					■	■		■	●								
20 µL	■					■	■		■	●								
25 µL/50 µL/100 µL	■						■		■	●								
200 µL	■								■	●								
200 µL/250 µL/500 µL/1,000 µL	■										■	●	●					
Variable volume, 1-, 8- and 12-channel pipettes																		
0.1 – 2.5 µL	■	■	●	●														
0.5 – 10 µL	■	●	■	■	●													
2 – 20 µL	■	◆	◆	●	■													
2 – 20 µL	■					■	■		■	●								
10 – 100 µL	■					◆	■		■	●								
20 – 200 µL	■					◆	◆		■	●								
30 – 300 µL	■					◆	◆		◆	■								
100 – 1,000 µL	■										■	●	●					
120 – 1,200 µL	■											●	■					
0.25 – 2.5 mL	■													■				
0.5 – 5 mL	■														■	●		
1 – 10 mL	■																■	●
Variable volume, 16- and 24-channel pipettes																		
1 – 20 µL	■				■													
5 – 100 µL	■							■										
Move It®, variable volume, 4-, 6-, 8- and 12-channel pipettes with adjustable tip spacing																		
1 – 20 µL	■				■													
5 – 100 µL	■							■										
30 – 300 µL	■					◆	◆		◆	■								
120 – 1,200 µL	■											●	■					
■ Test tip (other tips may require adjustment) ● Full volume range ◆ Limited volume range																		

Your local distributor: www.eppendorf.com/contact
Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com

Eppendorf®, the Eppendorf Brand Design, epT.I.P.S.®, ep Dualfilter T.I.P.S.®, SealMax®, Eppendorf Research®, and Move It® are registered trademarks of Eppendorf SE, Germany.
All rights reserved, including graphics and images. Copyright © 2025 by Eppendorf SE. Order No: AA01 050 120 – EN1/0825/MCP.



For more information go to
www.eppendorf.link/eptips

Combination of epT.I.P.S.® Pipette Tips with Eppendorf Reference® 2 Pipettes

Eppendorf Reference® 2		0.1 – 10 µL, dark gray	0.1 – 10 µL M, medium gray	0.1 – 20 µL, medium gray	0.5 – 20 µL L, light gray	2 – 20 µL, yellow	2 – 100 µL, yellow	2 – 200 µL, yellow	20 – 300 µL, orange	50 – 1,000 µL, blue	50 – 1,250 µL, green	50 – 1,250 µL L, dark green	0.1 – 2.5 mL, red	0.1 – 5 mL, violet	0.2 – 5 mL L, violet	0.5 – 10 mL, turquoise	0.5 – 10 mL L, turquoise
epT.I.P.S.® BioBased		■		■	■			■	■	■	■	■	■	■	■		
epT.I.P.S.® Standard		■		■	■			■	■	■	■	■	■	■	■	■	■
epT.I.P.S.® BioBased Reloads				■		■	■	■	■	■	■	■					
ep Dualfilter T.I.P.S.®		■	■		■	■	■	■	■	■		■	■		■		
ep Dualfilter T.I.P.S.® BioBased Reloads			■		■		■	■	■	■		■					
ep Dualfilter T.I.P.S.® SealMax®		■			■		■	■	■	■							
ep Dualfilter T.I.P.S.® SealMax® BioBased Reloads					■		■	■	■	■							
Fixed volume, 1-channel pipettes																	
1 µL/2 µL	■	■	●	●													
5 µL/10 µL	■	●	■	■	●												
20 µL	■			●	■												
10 µL	■					■	■	■	●								
20 µL	■					■	■	■	●								
25 µL/50 µL/100 µL	■						■	■	●								
200 µL	■							■	●								
200 µL/250 µL/500 µL/1,000 µL	■									■	●	●					
2 mL/2.5 mL	■										●		■				
Variable volume, 1-, 8- and 12-channel pipettes																	
0.1 – 2.5 µL	■	■	●	●													
0.5 – 10 µL	■	●	■	■	●												
2 – 20 µL	■	◆	◆	●	■												
2 – 20 µL	■					◆	◆	■	●								
10 – 100 µL	■					◆	■	■	●								
20 – 200 µL	■					■	■	■	●								
30 – 300 µL	■					◆	◆	◆	■								
100 – 1,000 µL	■									■	●	●					
0.25 – 2.5 mL	■										●		■				
0.5 – 5 mL	■													■	●		
1 – 10 mL	■															■	●

■ Test tip (other tips may require adjustment) ● Full volume range ◆ Limited volume range

Your local distributor: www.eppendorf.com/contact
Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com

Eppendorf®, the Eppendorf Brand Design, epT.I.P.S.®, ep Dualfilter T.I.P.S.®, SealMax® and Eppendorf Reference® are registered trademarks of Eppendorf SE, Germany.
All rights reserved, including graphics and images. Copyright © 2025 by Eppendorf SE. Order No: AA01 050 120 – EN1/0825/BE



For more information go to
www.eppendorf.link/eptips

Combination of epT.I.P.S.® Pipette Tips with Eppendorf Xplorer®/Xplorer® plus Pipettes

Eppendorf Xplorer®/ Eppendorf Xplorer® plus/ Eppendorf Xplorer® plus Move It®	0.1 – 10 µL, dark gray	0.1 – 10 µL M, medium gray	0.1 – 20 µL, medium gray	0.5 – 20 µL L, light gray	0.1 – 20 µL, light pink	2 – 20 µL, yellow	2 – 100 µL, yellow	5 – 100 µL, light yellow	2 – 200 µL, yellow	20 – 300 µL, orange	50 – 1,000 µL, blue	50 – 1,250 µL, green	50 – 1,250 µL L, dark green	0.1 – 2.5 mL, red	0.1 – 5 mL, violet	0.2 – 5 mL L, violet	0.5 – 10 mL, turquoise	0.5 – 10 mL L, turquoise
epT.I.P.S.® BioBased	■		■	■					■	■	■	■	■	■	■	■		
epT.I.P.S.® Standard	■		■	■	■			■	■	■	■	■	■	■	■	■	■	■
epT.I.P.S.® BioBased Reloads			■			■	■		■	■	■	■	■					
ep Dualfilter T.I.P.S.®	■	■		■	■	■	■	■	■	■	■		■	■	■	■		■
ep Dualfilter T.I.P.S.® BioBased Reloads		■		■		■	■		■	■	■		■					
ep Dualfilter T.I.P.S.® SealMax®	■			■			■		■	■	■							
ep Dualfilter T.I.P.S.® SealMax® BioBased Reloads				■			■		■	■	■							
Variable volume, 1-, 8- and 12-channel pipettes																		
0.5 – 10 µL	■	●	■	■	●													
1 – 20 µL	■	◆	◆	●	■													
5 – 100 µL	■					◆	■		■	●								
10 – 200 µL	■					◆	◆		■	●								
15 – 300 µL	■					◆	◆		◆	■								
50 – 1,000 µL	■										■	●	●					
50 – 1,200 µL	■										◆	■	●					
0.1 – 2.5 mL	■													■				
0.2 – 5 mL	■														■	●		
0.5 – 10 mL	■																■	●
Variable volume, 16- and 24-channel pipettes																		
1 – 20 µL	■				■													
5 – 100 µL	■							■										
Move It®, variable volume, 4-, 6-, 8- and 12-channel pipettes with adjustable tip spacing																		
1 – 20 µL	■				■													
5 – 100 µL	■							■										
15 – 300 µL	■					◆	◆		◆	■								
50 – 1,200 µL	■											■	●					
■ Test tip (other tips may require adjustment) ● Full volume range ◆ Limited volume range																		

Your local distributor: www.eppendorf.com/contact
Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com

Eppendorf®, the Eppendorf Brand Design, epT.I.P.S.®, ep Dualfilter T.I.P.S.®, Eppendorf Xplorer® and Move It® are registered trademarks of Eppendorf SE, Germany.
All rights reserved, including graphics and images. Copyright © 2025 by Eppendorf SE. Order No: AA01 050 120 – EN1/0825



For more information go to
www.eppendorf.link/eptips