


E 106

Insert

for the optimum loading of wide neck flasks, measuring cylinders, etc.



- Equipped with 10 spring hooks x 175 mm height and 16 x 105 mm height
- The distance between the spring hooks is approx. 60 mm


EAN: 4002511115675 / Material number: 03808310 / Old Material Number: 69510601D

Product category	
Laboratory washers	•
Large-chamber laboratory washers	•
SlimLine laboratory washers	•
Cleaning machines for industrial parts cleaning	•
Large-chamber laboratory washers for industrial parts cleaning	•
SlimLine laboratory washers for industrial parts cleaning	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 7111 [IPC]	•
PLW 8505	•
Product category	
Insert/module for beakers	•
Insert/module for wide-necked glassware	•
Area of application	
Reprocessing of laboratory glassware	•
Reprocessing of measuring cylinders	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	180
External dimensions, net width in mm	195

E 106

Insert

for the optimum loading of wide neck flasks, measuring cylinders, etc.



EAN: 4002511115675 / Material number: 03808310 / Old Material Number: 69510601D

External dimensions, net depth in mm	186
External dimensions, gross height in mm	205
External dimensions, gross width in mm	245
External dimensions, gross depth in mm	490
Net weight in kg	0.94
Gross weight in kg	1.42