


E 110

Insert

for the optimum loading of beakers from 250 to 600 ml.



- Suitable for 10 glass beakers between 250–600 ml
- equipped with 10 x 3 stay bars
- For use in lower basket

EAN: 4002511115736 / Material number: 03808390 / Old Material Number: 69511001D

Product category	
Laboratory washers	•
Large-chamber laboratory washers	•
SlimLine laboratory washers	•
Cleaning machines for industrial parts cleaning	•
Large-chamber laboratory washers for industrial parts cleaning	•
SlimLine laboratory washers for industrial parts cleaning	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 7111 [IPC]	•
PLW 8505	•
Product category	
Insert/module for beakers	•
Area of application	
Reprocessing of laboratory glassware	•
Product properties	
Material	Stainless steel, Wire, rilsan-coated
Colour	Stainless Steel, White
Capacity	
Beakers [quantity]	10
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	175

E 110

Insert

for the optimum loading of beakers from 250 to 600 ml.



EAN: 4002511115736 / Material number: 03808390 / Old Material Number: 69511001D

External dimensions, net width in mm	230
External dimensions, net depth in mm	460
External dimensions, gross height in mm	250
External dimensions, gross width in mm	230
External dimensions, gross depth in mm	505
Net weight in kg	0.83
Gross weight in kg	1.35