

E 111

Insert

For the optimum loading of beakers from 600 to 1,000 ml.



- Suitable for 8 glass beakers between 600–1000 ml
- equipped with 8 x 3 stay bars
- For use in lower basket

EAN: 4002511115743 / Material number: 03808420 / Old Material Number: 69511101D


Product category	
Laboratory washers	•
Large-chamber laboratory washers	•
SlimLine laboratory washers	•
Cleaning machines for industrial parts cleaning	•
Large-chamber laboratory washers for industrial parts cleaning	•
SlimLine laboratory washers for industrial parts cleaning	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 7111 [IPC]	•
PLW 8505	•
Product category	
Insert/module for beakers	•
Area of application	
Reprocessing of laboratory glassware	•
Product properties	
Material	Stainless steel, Wire, rilsan-coated
Colour	Stainless Steel, White
Capacity	
Beakers [quantity]	8
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	205

E 111

Insert

For the optimum loading of beakers from 600 to 1,000 ml.

EAN: 4002511115743 / Material number: 03808420 / Old Material Number: 69511101D



External dimensions, net width in mm	230
External dimensions, net depth in mm	460
External dimensions, gross height in mm	230
External dimensions, gross width in mm	250
External dimensions, gross depth in mm	480
Net weight in kg	0.96
Gross weight in kg	1.10