

A 300/3

Injector module

for optimum loading of laboratory glassware – equipped with 8 nozzles.



- Suitable for lab. glassware incl. Erlenmeyer and round flasks.
- Equipped with 4 A 840 nozzles and 4 A 841 nozzles
- Part of the EasyLoad system
- Nozzle spacing 118 mm across the width, 118 mm across the depth

EAN: 4002516188964 / Material number: 11116900 / Old Material Number: 69630003D

Product category	
Laboratory washers	•
SlimLine laboratory washers	•
Cleaning machines for industrial parts cleaning	•
SlimLine laboratory washers for industrial parts cleaning	•
Machine affiliation	
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
PLW 7111 [IPC]	•
Product category	
Insert/module for narrow-necked glassware	•
Insert for various flasks	•
Area of application	
Reprocessing of laboratory glassware	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Number of compartments	8
Part of the EasyLoad system	•
Capacity	
Narrow-necked glassware [number]	8
Dimensions and weight	
External dimensions, net height in mm	228
External dimensions, net width in mm	208
External dimensions, net depth in mm	479
External dimensions, gross height in mm	200
External dimensions, gross width in mm	240
External dimensions, gross depth in mm	515
Net weight in kg	1,47
Gross weight in kg	1,92
Companion pack	
A 840, injector nozzle Ø 6 mm, length 130 mm [number]	4
A 841, injector nozzle Ø 6 mm, length 210 mm [number]	4