


A 302/3

Injector module

for optimum loading of laboratory glassware – equipped with 32 nozzles.



- Suitable for lab. glassware such as Erlenmeyer and round flasks.
- Equipped with 16 A 844 nozzles and 16 A 845 nozzles
- Part of the SimpleLoad system
- Nozzle spacing 2 3/8" across the width, 2 1/4" across the depth

EAN: 4002516188988 / Material number: 11116960 / Old Material Number: 69630203D

Product category

Washer-disinfectors, laboratory	•
SlimLine laboratory glasswashers	•
Machines for cleaning industrial parts	•
SlimLine laboratory glasswashers for cleaning industrial parts	•

Machine affiliation

PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
PLW 7111 [IPC]	•

Product category

Insert/module for narrow-neck glassware	•
Insert for various bottles	•

Application area

Reprocessing of laboratory glassware	•
--------------------------------------	---

Product characteristics

Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Part of the SimpleLoad system	•

Capacity

Narrow-neck glassware [number]	32
--------------------------------	----

Dimensions and weight

External dimensions, net height in inches	5 11/16 (143)
External dimensions, net width in inches	9 3/8 (238)
External dimensions, net depth in inches	18 7/8 (479)
External dimensions, gross height in inches	5 3/8 (135)
External dimensions, gross width in inches	11 7/16 (290)
External dimensions, gross depth in inches	23 1/4 (590)
Net weight in lbs	3 (1.50)
Gross weight in lbs	6 (2.67)

Accessory pack

A 844, injector nozzle Ø 1/8" (2.5 mm), length 3 1/8" (80 mm) [quantity]	16
A 845, injector nozzle Ø 1/8" (2.5 mm), length 5" (125 mm) [quantity]	16