


A 101

Upper basket

Height adjustable, for the optimum loading of inserts.



- Vertical clearance 160 +/- 30 mm, width 475 mm, depth 443 mm
- With integrated spray arm


EAN: 4002516828143 / Material number: 12645590 / Old Material Number: 69610101D

| Product category                                          |                 |
|-----------------------------------------------------------|-----------------|
| Fresh-water dishwashers                                   | •               |
| Washer-disinfectors, laboratory                           | •               |
| Washer-disinfectors, medical                              | •               |
| Cleaning machines for industrial parts cleaning           | •               |
| Washer-disinfectors, fire brigades and emergency services | •               |
| Machine affiliation                                       |                 |
| PG 8063                                                   | •               |
| PLW 8636 [IPC]                                            | •               |
| PLW 8636 [LAB / LAB MON]                                  | •               |
| PG 8562                                                   | •               |
| PG 8582                                                   | •               |
| PG 8582 CD                                                | •               |
| PG 8504                                                   | •               |
| PG 8583                                                   | •               |
| PG 8583 CD                                                | •               |
| PG 8592                                                   | •               |
| PG 8593                                                   | •               |
| PLW 8604                                                  | •               |
| PLW 8683                                                  | •               |
| PLW 8683 CD                                               | •               |
| PLW 8693                                                  | •               |
| PLW 8683 [Safety]                                         | •               |
| PLW 8693 [Safety]                                         | •               |
| PWD 8682                                                  | •               |
| PWD 8682 CD                                               | •               |
| PWD 8692                                                  | •               |
| Product category                                          |                 |
| Upper basket                                              | •               |
| Carrier/basket carrier                                    | •               |
| Upper basket                                              | •               |
| Area of application                                       |                 |
| Reprocessing of laboratory glassware                      | •               |
| Reprocessing of GYN instruments                           | •               |
| Reprocessing of ENT instruments                           | •               |
| Reprocessing of basins                                    | •               |
| Reprocessing of kidney dishes                             | •               |
| Reprocessing of surgical instruments                      | •               |
| Reprocessing of different utensils                        | •               |
| Reprocessing of breathing air regulators                  | •               |
| Product properties                                        |                 |
| Material                                                  | Stainless steel |
| Colour                                                    | Stainless Steel |
| Dimensions and weight                                     |                 |
| External dimensions, net height in mm                     | 206             |
| External dimensions, net width in mm                      | 528             |
| External dimensions, net depth in mm                      | 527             |
| External dimensions, gross height in mm                   | 220             |

A 101

Upper basket

Height adjustable, for the optimum loading of inserts.



EAN: 4002516828143 / Material number: 12645590 / Old Material Number: 69610101D

|                                        |      |
|----------------------------------------|------|
| External dimensions, gross width in mm | 540  |
| External dimensions, gross depth in mm | 550  |
| Net weight in kg                       | 3.03 |
| Gross weight in kg                     | 3.68 |