

A 860
Holder rack
for A 300/2 and A 300/3 injector modules.



- Arrangement: 2 x 4 fields
- For use on A 300/2 and A 300/3 injector modules
- For the optimum positioning of laboratory glassware
- Part of the EasyLoad system

EAN: 4002516855118 / Material number: 12665410 / Old Material Number: 69686001D

Product category

Laboratory washers	•
SlimLine laboratory washers	•
Cleaning machines for industrial parts cleaning	•
SlimLine laboratory washers for industrial parts cleaning	•

Machine affiliation

PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
PLW 7111 [IPC]	•

Product category

Other inserts/other modules	•
-----------------------------	---

Area of application

Reprocessing of laboratory glassware	•
--------------------------------------	---

Product properties

Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Number of compartments	8
Compartment width in mm	114
Compartment length in mm	114
Part of the EasyLoad system	•

Dimensions and weight

External dimensions, net height in mm	30
External dimensions, net width in mm	250
External dimensions, net depth in mm	350
External dimensions, gross height in mm	170
External dimensions, gross width in mm	255
External dimensions, gross depth in mm	490
Net weight in kg	0.48
Gross weight in kg	0.89