

Article number: 14082825

Glove DermaShield 73-711 Size 7.5


[Ansell](#)

DermaShield™ 73-711 is the ideal glove for workers in aseptic environments who are concerned about allergy risks. Its proprietary material formulation is free of latex proteins and accelerators that can cause allergic reactions. A beaded cuff provides a secure fit over the upper arm.



Characteristics

Brand: ANSELL**Sterile:** Yes**Powdered:** No**Antistatic:** No**Silicone-free:** No**Length:** 300 cm**Thickness:** 0.15 mm**Finger Thickness:** 0.16**Colour:** Green**Norm:** EN 556, ISO 10282, EN 421, ISO 11193, ISO 13485, EN ISO 21420, EN ISO 374, EN 420, EN 374**Approvals:** AS/NZS 4179**AQL (Acceptable Quality Level):** 0.65**Cuff Thickness:** 0.12**Cuff Style:** Beaded with SUREFIT™ Technology**Material:** Neoprene**Internal Glove Surface:** Polymer Coated with DERMASHIELD™ Technology**External Glove Surface:** Textured Fingertip**Possible Sensitizer Ingredients:** None**EN 420:2003 + A1:2009:** Yes**EN 421:2010:** Yes**EN ISO 374-1:2016:** KLMNPT**EN ISO 374-5:2016:** Yes

Application

- Recommended in: Chemical

Disclaimer: The content of this document has been composed with the utmost care. However, it is possible that certain information changes over time, becomes inaccurate or incomplete. ERIKS does not guarantee that the information provided on this document is up to date, accurate and complete; the information provided is not intended to be advice. ERIKS shall never be liable for damage resulting from the use of the information provided.