

multiSUB mini



The multiSUB Mini is the smallest unit in the range, designed for low to medium numbers of samples. The small gel size maximises run economy but does not compromise versatility as two tray options are available – 7 x 7cm and 7 x 10cm – and combs ranging from preparative up to 16 samples. Simply by altering the gel tray or comb, this compact unit is capable of resolving up to 64 different samples, prepping 1ml of sample or separating sample bands over a distance of 9cm.

Buffer saver blocks physically reduce the volume of a gel chamber and so reduce buffer requirements, saving cost.

Features

- Available with 7 x 7cm, 7 x 10cm or with both gel trays
- Economic low gel and buffer volumes
- Small lab bench footprint

Ordering Information

6.258 373	multiSUB Mini, 7 x 7cm UV Tray, 2 x 8 sample combs, loading guides and dams
9.584 650	multiSUB Mini, 7 x 10cm UV Tray, 2 x 8 sample combs, loading guides and dams
7.920 043	multiSUB Mini, 7 x 7cm & 7 x 10cm UV Tray, 2 x 8 sample combs, loading guides and dams

Technical Specifications

Gel dimensions (w x L)	7 x 7cm 7 x 10cm
Unit dimensions (w x L x H)	9 x 21 x 9cm
Max. sample	7 x 7cm tray – 32 samples
Buffer volume	225ml
Combs available: No. of samples Thicknesses	1, 2, 4, 8MC, 8, 10, 12MC, 16



Phone: 1300 848 065
Email: info@labfriend.com.au
Web: www.labfriend.com.au

