



#### Package Includes:

- High-Speed Mini-Centrifuge
- 12 place 1.5/2.0 mL tube rotor
- 32 place 0.2 mL tube rotor (singles, or 4 x 8 Strip)
- 12 tube adapters 0.2 mL
- 12 tube adapters 0.5 mL

Maximum speed and flexibility, the Heathrow Scientific gusto® high speed mini centrifuge has a compact design that fits any workstation. Comfortable and easy to operate with tool-free quick rotor exchange. Simply inset tubes, close lid, adjust speed/time and press start. For duplicate spins, simply reload rotor and press start. Quick-spin option included for commonly used 20-second spins.

Safety features help ensure it will not operate when the lid is open, a rotor is jammed, or the unit is out of balance.

Ideal for most protocols requiring fast spins (12,500 rpm / 9,800 x g), and appropriate for PCR, microfilter cell separation, and HPLC protocols.

- Select and set the digital speed display in RCF or RPM
- Back lit display showing time, speed and operational information on the same screen
- UV-Resistant plastic, increased useful life under harsh light conditions
- double insulated, low voltage power supply- push button lid release
- Manual override allows user to open lid in case of power failure

Part No.	Description
<b>HS10050</b>	gusto®, 110–240 VAC, 50/60 Hz

#### Specifications

<b>Speed:</b>	500–12,500 RPM
<b>RCF:</b>	5–9, 800 Xg
<b>Power:</b>	Suitable for worldwide voltages 110 to 240 VAC, 50/60 Hz with 4 interchangeable power cords
<b>Dimensions:</b>	203 x 171 x 114 mm (7.92 x 6.67 x 4.45 in.)
<b>Two-year warranty</b>	



Heathrow Scientific, LLC | 620 Lakeview Parkway, Vernon Hills IL 60061 USA  
Tel.: +1-800-741-4597 | Fax: +1-847-816-5072 | E-mail: info@heatsci.com