

Heathrow Scientific® Biohazard Disposal Bags and Holders



Description Biohazard collection bags printed with biohazard warning symbol and precautions in English, Spanish, and French. For personal or general lab use, sizes to suit all needs. Made of 1.57 mil. thick polypropylene. Includes an indicator patch that darkens when autoclaved. Autoclavable to 250°F (121°C). Disposable.

Epoxy coated steel wire holder (Catalog No. 140-99221) available for 203 x 305 mm bags (Catalog No. 970-91996).

Image(s) are representative of the product group and not necessarily the individual product.

Order		PRICES ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE.						
14-Feb-2026 08:54: AM PST								
Catalog	Item	Availability	Manufacturer	Mfg. Part #	NDC#	List Price	Your Price	Unit
140-99221-EA	Biohazard Bag Holder - 1 ea Size in. (mm): 3.6 x 5 x 0.8 (91 x 128 x 21)	Sign In	Heathrow Scientific	HS10044H		\$25.40		EA
970-91996-PK	Biohazard Bags, 203 x 305 mm (8 x 12 in.), 1.57 mil. Thick - 500 ea Size in. (mm): 8 x 12 (203 x 305)		Heathrow Scientific	HS10320		\$95.35		PK
970-91999-PK	Biohazard Bags, 356 x 483 mm (14 x 19 in.), 1.57 mil. Thick - 200 ea Size in. (mm): 14 x 19 (356 x 483)		Heathrow Scientific	HS10321		\$86.15		PK
970-92002-PK	Biohazard Bags, 483 x 584 mm (19 x 23 in.), 1.57 mil. Thick - 200 ea Size in. (mm): 19 x 23 (483 x 584)		Heathrow Scientific	HS10322		\$133.10		PK
970-92005-PK	Biohazard Bags, 635 x 838 mm (25 x 33 in.), 1.57 mil. Thick - 200 ea Size in. (mm): 25 x 33 (635 x 838)		Heathrow Scientific	HS10323		\$318.45		PK