



Heat Sealer with Cutter

Quickly and efficiently seal polyethylene or polypropylene bags, sheets or flat tubing from 0.01 to 0.8 mm thick

Sliding blade removes excess material 19 mm from the seal.

- Timing knob ensures safe sealing. Simply turn knob to set time (1-8 seconds), press handle down and release when light indicates cycle is complete
- Can seal a bag up to 12 in (30.5 cm) long
- Heating element measures 31 cm long and is only hot when sealing
- Includes one replacement part kit
- Additional replacement part kits can be ordered separately

Replacement Part Kit:

- 2 extra blades
- 2 extra heating elements
- 2 extra PTFE cloths

Item No.	Power Supply		D x W x H in	D x W x H cm	UOM
HS28270	220/230 VAC, U.K	☐			
HS28271	220/230 VAC, Cont. Europe/Israel	◀••▶	3.4 x 17.6 x 6.5	8.6 x 44.6 x 16.6	1 ea, 8 ea/cs
HS28272	110/120 VAC, U.S./Taiwan/Japan	☐☐			
HS28280	Replacement Part Kit				1 ea/pk