

Ausstellungs-Schaugläser, DURAN®, DURAN Group / IDL DURAN® Specimenjar

Schaugläser, zylindrisch, auf Stopfen stehend.

Vol. ml	Ø mm	H mm	Schliff-Ø mm	Best.-Nr.	St.	VE St./Pack
100	52	135	34,5	12 21 58024	1-9	10
300	69	163	45	12 21 58039	1-9	10
600	81	214	50	12 21 58048	1-9	10
750	90	240	60	12 21 58051	1-9	10
1200	100	253	60	12 21 58058	1-9	10

! Preis pro Stück bei Abnahme kompletter Verpackungseinheiten (VE).

