

Heating Mantles, SAF Wärmetechnik

LabHEAT® Metal-cased heating mantles KM-ME for flasks, round bottom

Flexible heater in a plastic-coated, chemical-resistant metal housing with integrated mains and heating zones switch, and additionally with integrated controller for continuous power adjustment. Max. heating element temperature 450°C. 1.5 m supply cable with earth contact type plug and safety switch, nominal voltage 230 V AC.

Vol. ml	Heating cap. W	Quantity heating zones	Order No.	Vol. ml	Heating cap. W	Quantity heating zones	Order No.
50	55	1	45 91 92017	3000	600	2	45 91 92068
100	100	1	45 91 92024	4000	750	2	45 91 92071
250	150	2	45 91 92036	5000	860	2	45 91 92073
500	200	2	45 91 92044	6000	1000	2	45 91 92076
1000	300	2	45 91 92054	10000	1400	2	45 91 92086
2000	500	2	45 91 92063	20000	2000	2	45 91 92091

Metal-cased heating mantles for glasses or other vessel forms available on request.



WWE38018

SAF
WÄRMETECHNIK

