



Kartell LABWARE



Draining Rack



www.kartell.com.au


Kartell LABWARE

NEW IN
2010!


The KARTELL Draining Rack (ART. 213) now comes supplied with a pack of 11 additional small diameter pegs. These small pegs are ideal for items such as test tubes.

Draining Rack

- 72 place support rack
- 4mm rigid High Impact Polystyrene construction
- Peg sockets have closed ends eliminating leakage and consequent bio-hazards
- Seamless drip channel for waste collection
- Drainage tube (provided) attaches to drip channel
- Includes 72 pegs (art. 21308) for items with a neck bore greater than 15mm
- Includes 11 pegs (art. 1213) for items with a neck bore greater than 6mm
- Flexible peg positioning allows for larger objects
- Mounting kit, draining hose and fixing template supplied
- Additional racks can be joined together to make a bank of racks, using the joining plugs provided
- Both packs of spare pegs are available separately



Included:



DRAINING RACK	art.	no. of places	dimensions (mm)
	213	72	450 x 630 x 110

Available From:



SIEPER & Co. Pty Ltd
Wholesale Distributors
Ph: (02) 9352 0600 | www.sieper.com.au