

LABORHITS

LLG-Buffer solution pH

Accuracy $\pm 0,01$ pH

With certificate of analysis

Colour-coded

Advantages:

- Accuracy ± 0.01 pH, tested and certified at 25 °C
- Temperature dependence data table printed on the label
- Traceable to NIST
- Shelf-life: 24 months (pH 7 and 10) or 36 months (pH 4) from date of production



Ready-to-use, for easy and exact calibration of your pH-meters.

The ready-to-use colour coded pH-buffer solutions offer excellent quality, precision and are designed to suit all end user requirements. Supplied in 250 ml or 1 L bottles and conveniently colour coded to eliminate any risk of confusion between the 3 pH-values. Certificates of analysis under www.llg-labware.com.

Twin-neck dispensing bottles with integrated calibration vessel, allows the particularly economical dispensing of solution and prevents contamination.

ph-value at 25 °C	Capacity ml	Colour	Cat. No.
4.00	250	red	9.040 868
7.00	250	yellow	9.040 869
10.00	250	blue	9.040 871

Bottle, 1 liter

ph-value at 25 °C	Capacity ml	Colour	Cat. No.
4.00	1000	red	9.040 865
7.00	1000	yellow	9.040 866
10.00	1000	blue	9.040 867