

LABORHITS



**ergonomic
soft grip**

**colour-coded
volume**

LLG - Fixed Volume Microlitre Pipettes

Ergonomic shape for a very comfortable handling

The rounded shape and the unique soft grip of the LLG pipette ensures a safe and comfortable handling. A minimal effort is required for pipetting so that the risk of repetitive strain injuries (RSI) is reduced.

Maximum Precision and Accuracy

Thermal insulation of the internal components, using the innovative soft grip, improves accuracy.

Universal Blow-out System

The two-step plunger operation allows reverse pipetting. This also facilitates last drop dispensing by soft plunger movements. High accuracy and precision are guaranteed.

Colour-coded push-button for quick operation

The colour code enables a quick choice of the right pipette depending on the volume to be pipetted (see corresponding column in the table below).

Easy Maintenance and Adjustment

The volume setting is easy and precise, using the push-button, even when wearing gloves. Maintenance and adjustment take only a few minutes with just one tool. (Included in delivery)

Contamination-free Pipetting

The pipette is fully autoclavable, ideal for applications where sterility is important.

Satisfaction Guaranteed

Each pipette is individually tested and calibrated twice, conforming to EN ISO 8655 and DIN 12600 Standards. A certificate is included.

Made in Germany, CE certified.

Available with fixed or adjustable volume.



| Capacity µl | Accuracy ≤ ± A% (max. vol.) | Reproducibility ≤ ± CV% (max. vol.) | Colour code | PK | Cat. No. |
|----------------|-----------------------------------|---|-------------|----|-----------|
| 5 | 1.0 | 0.5 | grey | 1 | 9.277 515 |
| 10 | 1.0 | 0.5 | grey | 1 | 9.277 516 |
| 20 | 4.5 | 1.0 | grey | 1 | 9.277 517 |
| 25 | 2.5 | 0.6 | yellow | 1 | 9.277 518 |
| 50 | 2.5 | 0.6 | yellow | 1 | 9.277 519 |
| 100 | 2.5 | 0.6 | yellow | 1 | 9.277 520 |
| 200 | 2.5 | 0.6 | blue | 1 | 9.277 521 |
| 250 | 2.9 | 0.6 | blue | 1 | 9.277 522 |
| 500 | 2.9 | 0.6 | blue | 1 | 9.277 523 |
| 1000 | 2.4 | 0.6 | violet | 1 | 9.277 524 |
| 2000 | 2.4 | 0.6 | violet | 1 | 9.277 525 |
| 2500 | 2.4 | 0.6 | violet | 1 | 9.277 526 |
| 5000 | 2.4 | 0.6 | turquoise | 1 | 9.277 527 |
| 10.000* | 3.0 | 0.6 | turquoise | 1 | 9.277 528 |

* 10.000µl pipette does not have a tip ejector