




## Product data sheet

### LLG-Disposable Visitor gown

Art.-Nr.: 6.282 796

Material : non-woven fabric  
Weight : 30gsm  
Color : white  
Size : universal

	Allowance	
Length (A)	106	± 2 cm
Chest (B)	67	± 2 cm
Upper sleeve L (C)	57	± 2 cm
Armhole height (D)	26	± 2 cm
Cuff (E)	16	± 2 cm
Collar width (F)	26	± 2 cm
Collar height (G)	8	± 2 cm
Neck-cut-out (H)	19	± 2 cm
Buttons space (I)	17	± 2 cm



#### Safety instructions:

1. Identification of Hazards  
No hazards.

2. First-aid measures

In case of rashes or any other discomfort occur, discontinue use immediately and wash hand with plenty of water.  
Consult physician.

3. Fire-fighting measures

Extinguishing media : Water, foam or dry powder

In case of fire, wear respiratory protection equipment.

4. Handling and storage

Store in cool and dry area. Protect from direct sun light and UV light. Avoid prolonged exposure to atmospheric air.

5. Toxicology information

Skin irritation: No irritation (rabbit)

Dermal contact sensitization No sensitization (guinea pig)

#### QUALITY CERTIFY

ISO 13485: 2003

CE MARK

Lab Logistics Group GmbH  
Product Management  
Am Hambuch 1  
D-53340 Meckenheim

2017-04-12