


LABORHITS

NMR tubes for your Lab



NMR tubes, 3 and 5 mm, borosilicate glass 3.3


The NMR tubes have precise tolerances for consistent high accuracy. Most suitable for autosampling NMR systems and high-throughput NMR. Manufactured from borosilicate glass 3.3, conforming to USP Type I and ASTM E438, Type I, Class A

requirements.

- Borosilicate glass 3.3
- Economy quality
- One end closed, rim fire polished
- Applicable until 600 MHz


Please order closing caps separately.

Outer diam. mm	Int. diam. mm	Wall thickness mm	Length mm	PK	Cat. No.
2,95 ± 0,03	2,36 ± 0,03	0,29	178	50	6.281 792
2,95 ± 0,03	2,36 ± 0,03	0,29	203	50	6.281 793
4,95 ± 0,05	4,19 ± 0,05	0,38	178	100	9.400 310
4,95 ± 0,05	4,19 ± 0,05	0,38	203	100	9.400 311



Closing caps for NMR tubes, 3 and 5 mm

Description	Colour	PK	Cat. No.
for 3 mm tubes	red	50	6.282 148
for 3 mm tubes	green	50	6.282 149
for 3 mm tubes	yellow	50	6.282 150
for 3 mm tubes	blue	50	6.282 151
for 5 mm tubes	red	100	9.400 312
for 5 mm tubes	green	100	9.400 313
for 5 mm tubes	yellow	100	9.400 314
for 5 mm tubes	blue	100	9.400 316



NMR Tube Rack, polypropylene

Polypropylene rack with carrying handles features an upper deck with holes and a lower deck support for tube ends. Stackable.

For 72 tubes with 3 or 5mm diameter.

Dimensions: 213 x 115 x 222 mm (LxWxH)

Type	PK	Cat. No.
NMR Tube Rack	1	9.301 065