


LABORHITS

Eluent SafetyCap Starter Kits


Kit I



Kit II



Kit III



Various Kits for Selection

Starter kits for HPLC with up to 4 storage bottles. The kits also contain air valves, fittings for connecting capillaries, and blind fittings for closing connections that are not required.

Description	Thread	PK	Cat. No.
Starter Kit I	GL 45	1	SC000164.1*LE
Starter Kit II	GL 45	1	SC000134.1*LE
Starter Kit III	GL 45	1	SC000135.1*LE

Scope of supply:

Starter Kit I:

4 x Eluent SafetyCaps with 1 PP fitting,
4 x air valves

Starter Kit II:

3 x Eluent SafetyCaps with 1 PP fitting,
1 x Eluent SafetyCap with 2 PP fittings,
4 x air valves,
1 x blind fitting

Starter Kit III:

4 x Eluent SafetyCaps with 2 PP fittings,
4 x air valves,
4 x blind fittings