

	Life Sciences	Affected Site(s): US : NC : Durham	Document Number: LSR00026	
Description: Aspirating Pipets - 5mL			Revision: 7	Page: 1 of 2

APPLICATIONS:

For tissue culture usage or any non-toxic solution that requires aliquoted or aspirating/dispensing in an aseptic environment.

FEATURES:

- Sterility assurance on individually wrapped pipet
 - Contamination control
- Consistent performance
 - Reliability, contamination control
- Aseptic Presentation
 - Contamination control
- Non-pyrogenic
 - Low endotoxins
- Crystal grade polystyrene
 - Can clearly see fluid sample
- Non-wettable
 - Minimum meniscus
- Individually wrapped packaging: Cellulose-based, direct-seal paper on one side and the bottom web is polyethylene plastic material.

REGULATORY COMPLIANCE:

Falcon® products are manufactured under the current ISO 13485 Standards, and the current FDA Quality System Regulation 21 CFR Section 820.

STERILITY:

Product has been sterilized and dosimetrically released per the requirements of ANSI/AAMI/ISO 11137, "Sterilization of Health Care Products – Radiation". Product meets minimum Sterility Assurance Level (SAL) of 10⁻⁶.

Warranty/Disclaimer: Unless otherwise specified, all products are For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.

For a listing of trademarks, visit www.corning.com/clstrademarks. All other trademarks are the property of their respective owners.
© 2022 Corning Incorporated. All rights reserved.

CORNING	Life Sciences	Affected Site(s): US : NC : Durham	Document Number: LSR00026	
Description: Aspirating Pipets - 5mL			Revision: 7	Page: 2 of 2

Falcon® 5 mL Aspirating Pipet

Material: Polystyrene

Related Catalog Numbers:

Sterile


Cat. No.

357501

Packout Description

Individually wrapped 50/Pk,
200/case

DIMENSIONS = XXX mm (XXX inches)			
TIP:		PIPET SIZE: 5 mL	
Length:	= 37.36 (1.471)	INCREMENTS:	N/A
Interior Diameter:	= 1.57 (0.62)		
Outside Diameter:	= 3.35 (.132)		
MOUTHPIECE:		GRADUATIONS: N/A	
Length:	= N/A	OVERFILL CAPACITY:	N/A
Interior Diameter:	= 8.03 (.316)		
Outside Diameter:	= 6.01 (.236)	COLOR CODE:	N/A
Plug Length:	= N/A		
EXTRUSION:		GENERAL DIMENSIONS:	
Interior Diameter:	= 6.01 (.236)	A = 348.51 (13.721)	
Outside Diameter:	= 8.03 (.316)	B = 3.35 (.132) Diameter	
Length:	= 311.15 (12.250)	C = 8.03 (.316) Diameter	
		D = 8.03 (.316) Diameter	



Warranty/Disclaimer: Unless otherwise specified, all products are For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.

For a listing of trademarks, visit www.corning.com/clstrademarks. All other trademarks are the property of their respective owners.
© 2022 Corning Incorporated. All rights reserved.