

## Lab-Jacks MAXI 18/10 stainless steel



Lab-Jack MAXI 18/10 stainless steel

### Properties:

- The MAXI lab jacks offer 30% more height because of an amplified construction
- Scissor-construction made of 18/10-stainless steel
- Adjustment of height with the help of an adjusting wheel.
- Ergonomic knob, rounded and user friendly
- With the operating load we differ between the maximal dynamic operating load (Kg max. dyn.) and the maximal static operating load (Kg max. stat.) of the lab jack. These norms are also according to DIN 12897.

Material	A (Length)	B (Width)	C min	C max	Kg max dyn	Kg max stat	Article no.
18/10 Stainless	160	130	75	400	7	30	11122
18/10 Stainless	200	200	75	400	7	30	11132