

LABNET INTERNATIONAL, INC.

# SPECTRAFUGE™ MINI CENTRIFUGE

JOIN THE MINI REVOLUTION

**Labnet**  
Labnet International, Inc.



- Complete with interchangeable rotors for microtubes and PCR\* strips
- Ideal for microfiltration and quick spins
- Durable lid hinge
- Compact design, small footprint
- On/off switch
- Five translucent lid colors available



## SPECIFICATIONS

Capacity	6 x 1.5/2.0 ml tubes 2 x 0.2 ml strip tubes
Max speed/g-force	6000 rpm/2000 x g <sup>1</sup>
Dimensions (WxDxH)	5.9 x 5.9 x 4.6 in./15 x 15 x 11.7 cm
Weight	1 lb/0.45 kg
Electrical	120V~, 60 Hz or 230V~, 50 Hz

<sup>1</sup>Average max rcf in the strip tube rotor is 1,177 x g

## ORDERING INFORMATION

Catalog No.	Description
C1301*	Spectrafuge Mini with 1.5/2.0 ml rotor, strip tube rotor, 0.5 and 0.4 ml adapters, smoke gray lid
C1301-230V-EU*	Same as above, 230V Europe Power Cord
C1301-230V-UK*	Same as above, 230V UK Power Cord
C1205	Individual adapters for 0.5 ml tubes, pk of 6
C1206	Individual adapters for 0.4 ml tubes, pk of 6
C1222	Individual adapters for 0.2 ml tubes, pk of 6
C1300-RT	Replacement 6 x 1.5 ml rotor
C1300-RTS	Replacement strip tube rotor

\* Add -B for blue lid, -T for teal lid, -P for purple lid, -R for red lid

**Labnet**  
Labnet International, Inc.

P.O. Box 841 • Woodbridge, NJ 07095 • USA • website: <http://www.labnetlink.com>  
US Office • toll free: 888-522-6381 • fax: 732 417-1750 • e-mail: [labnet@labnetlink.com](mailto:labnet@labnetlink.com)  
International Office • phone: +1-732-417-0700 • fax: +1-732-417-1750  
e-mail: [international@labnetlink.com](mailto:international@labnetlink.com)