


LABNET INTERNATIONAL, INC.

SLIDE SPINNER

FOR MICROARRAY SLIDE DRYING

Labnet
Labnet International, Inc.



- Quickly dries 2 slides
- Ideal for removing wash buffer from microarrays prior to scanning
- Removable cassettes conveniently catch all liquid

SPECIFICATIONS

Capacity	2 slides
Maximum Speed	4,800 rpm
Dimensions (WxDxH)	5 x 5 x 4.4 inches/12.7 x 12.7 x 11.2 cm
Weight	1 lbs/0.45 kg
Electrical	120V~ or 230V~, 50/60 Hz

ORDERING INFORMATION

Catalog No.	Description
C1303-T	Slide Spinner for 2 slides, 120V
C1303-T-230V-EU	Slide Spinner for 2 slides, 230V, with European cord
C1303-T-230V-UK	Slide Spinner for 2 slides, 230V, with UK cord
C1303-SC	Additional slide cassettes, pack of 2

Labnet
Labnet International, Inc.

P.O. Box 841 • Woodbridge, NJ 07095 • USA • website: <http://www.labnetlink.com>
US Office • toll free: 888-522-6381 • fax: 732 417-1750 • e-mail: labnet@labnetlink.com
International Office • phone: +1-732-417-0700 • fax: +1-732-417-1750
e-mail: international@labnetlink.com