

LKexv 5400

MediLine



Total gross capacity	554 litre
Ext. dimensions in mm (w/d/h)	750/730/1640
Int. dimensions in mm (w/d/h)	600/560/1452
Energy consumption in 24 h	0,983 kWh
Climate rating ¹	SN-T
Cooling system	forced-air
Defrosting	automatic
Temperature range	+1°C to +15°C
Gastronorm	–
Housing material/colour	steel/white
Door/lid material	Steel
Liner Material	moulded polystyrol white
Type of control	analogue controls
Temperature display	external digital
Interior light	–
Fault: warning signal	–
Interface	–
volt-free alarm contact	no
Shelves	6
Adjustable Shelving	5
Number of drawers	–
Number of baskets	–
Maximum number of baskets	–
Shelf material	glass
Shelf loading	40 kg
Height of adjustable legs	-
Castors/type of castors	–
Door handle	Ergonomic door handle
Foot pedal opener	–
Self-closing door	–
Insulation	–
Defrost water drain	–
Refrigerant	R 600a
Noise level	48 dB(A)
Gross/net weight	84/77 kg
EAN-Nr.	9005382164853
Connection Rating	1.5 A / 150 W

¹ Climate Ratings: **SN** +10°C to +32°C **N** +16°C to +32°C **ST** +16°C to +38°C **T** +16°C to +43°C