

Data Sheet

Page 1/1



CHROMABOND columns Florisil 3 mL, sorbent: 500 mg, PP, with PE filter elements pack of 50

Part Number: MN-730081
Manufacturer: Macherey-Nagel
Quantity: 1 pc

Description

CHROMABOND columns Florisil volume: 3 mL, content of sorbent: 500 mg material: PP, with PE filter elements pack of 50

Technical Information

Brand:	CHROMABOND	Mode of Separation:	Normal Phase Chromatography (NP)
Product Type:	SPE Cartridges		
Particle Size:	150-250 µm		
Chemistry:	Florisil		
Chemistry Producers Term:	Florisil		
Country of Origin:	Germany		
Volume:	3 ml		
Bed Mass [mg]:	500 mg		