



NEW!

LABORHITS



- Optional Plan Apochromatic lenses
- AUTO ON-OFF function
- Light memory function, coded 6-fold nosepiece
- Rackless stage concept
- Ceramic stage insert
- 50W and 100W Halogen version

Clinical & Lab Microscope BA410 ELITE Series

The Elite version of Motic's established BA410 flagship model displays the highest level of optical performance.

A complete set of Plan Apochromatic lenses is now available for maximum demands in colour reproduction and resolution. The BA410E allows a profound diagnosis in pathology, haematology and cytology with their requirements for best colour fidelity. Increased Numerical Apertures provide an expanded resolution power for a more detailed specimen evaluation.

The AUTO ON-OFF function with its energy-saving aspect improves lab safety: a final instrument check before closing time is no longer necessary. A light memory function, based on a coded 6-fold nosepiece, memorizes the light intensity of each nosepiece position. When changing to a specific objective, the last illumination setup of this position is recalled. Once coordinated, there is no need to adjust the illumination when changing magnification. For increased illumination power, a 100W Halogen version is available.

The rackless stage concept, well known from BA210/BA310 Elite models, is now adopted in this flagship model. A ceramic stage insert for highest abrasion resistance guarantees a perfectly flat placing of the glass slide for a lifetime.

Туре	Designation	РК	Cat. No.
BA410E	Binocular, 50W N-WF, 10X/22mm	1	9.727 050
BA410E	Binocular, 100W N-WF, 10X/22mm	1	9.727 051
BA410E	Trinocular, 50W N-WF, 10X/22mm	1	9.727 052
BA410E	Trinocular, 100W N-WF, 10X/22mm	1	9.727 053