

# 2014 EUROPE PRICE LIST

## Digital Microscopes | DM-1802



Reference	Article	Price (€)*
DD99421201	<b>DM-1802</b> Monocular head, 30° inclined with Built-in digital camera CMOS 1/2" chip 3MP - 2048x1536 pixels with USB2.0 output Widefield eyepiece WF10X/18mm Reversed quadruple revolving nosepiece Achromatic objectives EA 4X, 10X, 40X S Coarse and fine focusing system with tension adjustment Built in low position coaxial mechanical stage Focusable 1.25 N.A. Abbe condenser with iris diaphragm and filter holder Halogen illumination 12V/20W with intensity control Power supply 100V-240V (CE) Blue filter and dust cover Motic Images Plus 2.0 for PC & Mac Calibration slide	682,00

## Digital Microscopes | DM-1802



Reference	Article	Price (€)*
<b>DM-1802 accessories</b>		
<b>Eyepieces</b>		
1101001400551	Huygens eyepiece H5X/14.5mm	21,00
12180142	Widefield eyepiece WF10X/18mm	21,00
53100143	Widefield eyepiece WF15X/12mm	21,00
1101001400572	Widefield eyepiece WF20X/11mm	41,00
<b>Micrometer eyepieces</b>		
12180143	Micrometer eyepiece WF10X/18mm, 100 divisions in 10mm and crosshair	49,00
<b>Reticles</b>		
1101001401025	Reticle with 100 divisions in 10mm and crosshair (Ø19mm)	19,00
1101001400075	Reticle holder for 18mm eyepiece	9,00
<b>Achromatic objectives - EA</b>		
12180161	Achromatic objective EA 4X/0.10	16,00
12180162	Achromatic objective EA 10X/0.25	23,00
12180163	Achromatic objective EA 20X/0.40	34,00
12180164	Achromatic objective EA 40X/0.65/S	34,00
12180165	Achromatic objective EA 60X/0.85/S	42,00
12180168	Achromatic objective EA 100X/1.25/S - Oil	42,00
<b>Spare parts</b>		
1101000500172	Blue filter (Ø 32mm)	5,00
1101001201012	Dust cover (size C)	2,00
1101001300011	Immersion oil 5ml	5,00
1101001203491	Aluminium box	68,00
<b>Spare bulbs</b>		
12180258	12V/20W Halogen bulb	7,00