



🇬🇧 Bottle Top Dispenser
FORTUNA® OPTIFIX® BASIC,
 with PTFE encased piston , FIX
 adaptor made of white PP, comes
 with 2 PP adaptors GL 40+45,
 1 PTFE intake tube, single pack

🇩🇪 Flaschenaufsatzdispenser
FORTUNA® OPTIFIX® BASIC,
 mit PTFE ummanteltem Kolben,
 FIX Adapter aus weißem PP,
 inkl. 2 PP-Adapter GL 40+45,
 1 PTFE-Ansaugschlauch, VE 1

ISO 8655

| Volume Volumen | grad. grad. | Accuracy Richtigkeit | Tolerance Toleranz | Art.-No. Art.-Nr. |
|-------------------|----------------|-------------------------|-----------------------|----------------------|
| 0.4 - 2 ml | 0.1 ml | $\leq \pm 0,6 \%$ | $\leq \pm 0.2 \%$ | 101.080-27 |
| 1 - 5 ml | 0.1 ml | $\leq \pm 0,6 \%$ | $\leq \pm 0.2 \%$ | 101.080-33 |
| 2 - 10 ml | 0.2 ml | $\leq \pm 0,6 \%$ | $\leq \pm 0.2 \%$ | 101.080-38 |
| 6 - 30 ml | 0.5 ml | $\leq \pm 0,6 \%$ | $\leq \pm 0.2 \%$ | 101.080-44 |
| 10 - 50 ml | 1.0 ml | $\leq \pm 0,6 \%$ | $\leq \pm 0.2 \%$ | 101.080-48 |
| 20 - 100 ml | 2.0 ml | $\leq \pm 0,6 \%$ | $\leq \pm 0.2 \%$ | 101.080-51 |
| 40 - 200 ml | 5.0 ml | $\leq \pm 0,6 \%$ | $\leq \pm 0.2 \%$ | 101.080-53 |
| 60 - 300 ml* | 5.0 ml | $\leq \pm 0,6 \%$ | $\leq \pm 0.2 \%$ | 101.080-56 |
| 100 - 500 ml* | 5.0 ml | $\leq \pm 0,6 \%$ | $\leq \pm 0.2 \%$ | 101.080-58 |

* additional to ISO 8655
 * in Ergänzung zur ISO 8655

🇬🇧 For dispensing aqueous solutions, acids, crystallizing liquids
🇩🇪 Zur Dosierung von wässrigen Lösungen, Säuren und
auskristallisierenden Flüssigkeiten