



Balances & Scales Pricelist 2024



Ingenious Portfolio

Our rich history and extensive industry knowledge has resulted in a comprehensive line of high-performance balances and scales that continually meet our customers' needs at economical and cost-effective prices.

Look inside for NEW 2024 product and pricing information!

Ingeniously Practical English € MEC

Stay Connected

www.ohaus.com

Don't hesitate to try out our www.ohaus.com.
We added some extra features that will help you get
the information you need with one click.



Dealer Zone

If you are looking for marketing materials visit the
new Dealer Zone to get pictures, videos, brochures,
presentations, latest news and more!

Let us keep you in the loop!



Follow **OHAUS** on **YouTube**, **LinkedIn** and **Facebook** to learn about our latest products,
promotions and success stories with customers around the globe.





Discover our ingenious portfolio

1 Explorer™ Semi-Micro

OHAUS highest performance balances, now available with software update 2.25. Data integrity, enhanced user management system, log file, calibration file and label printing function, all you need in one balance!

► See page 11

2 MB120

MB120 with revolutionary SmartGuide™ feature. Designed to simplify the method creation and development. With MB120, moisture analysis process is faster, easier, and more efficient.

► See page 31

3 Defender™ 6000

The OHAUS IP68/IP69K washdown table scale Defender™ 6000, constructed from 316 stainless steel.

► See page 45

4 Courier™ 7000

New range of scales ideal for general shipping and animal weighing applications.

► See page 84

Quick Finder	6 Industrial scales
Semi-Micro Balances	9 Semi-Micro, Analytical, Precision and Portable Precision balances
	11 Explorer™ Semi-Micro
	12 Adventurer™ Semi-Micro
	13 Pioneer™ Semi-Micro
Analytical Balances	15 Explorer™ Analytical
	16 Adventurer™ Analytical
	17 Pioneer™ Analytical
	18 PR Series Analytical
Precision Balances	20 Explorer™ Precision / High Capacity
	21 Adventurer™ Precision / High Capacity
	22 Pioneer™ Precision / High Capacity
	23 PR Series Precision
Precision Portable Balances	25 Scout™ STX
	26 Scout™ SKX
	27 Navigator™ / Navigator™ NV / Navigator™ NVT
	28 Compass™ CX
	29 Compass™ CR
	29 YA Gold Series
Moisture Analysers	31 MB120 / MB90
	32 MB27 / MB25
	33 MB23
Mechanical Scales & Balances	35 Cent-O-Gram™ 311-00 / Dial-O-Gram™ 310-00 / Triple Beam Pro™ / Triple Beam Junior / Triple Beam 700
	36 Harvard Junior / Harvard Trip
Floor Scales	37 DF 5000 Series Stainless Steel Floor Scales
	38 DF 5000 Series Floor Scales
	39 DF 3000 Series Stainless Steel Floor Scales (i-DF33)
	40 DF 3000 Series Floor Scales (i-DF33)
	41 DF 3000 Series Floor Scales (DF32)
Pallet Scales	43 DFP 5000 Series Pallet Scales
	44 DFP 3000 Series Pallet Scales
Bench Scales	45 Defender™ 6000 Stainless Steel
	46 Defender™ 6000 Hybrid
	47 Defender™ 6000 Front Mount
	49 Defender™ 5000 Stainless Steel (Single Range)
	50 Defender™ 5000 Stainless Steel (Dual Range)
	51 Defender™ 5000 Hybrid
	52 Defender™ 5000 Standard
	53 Defender™ 5000 Front Mount
	54 Defender™ 5000 Front Mount
	56 Defender™ 3000 Stainless Steel
	57 Defender™ 3000 Hybrid
	58 Defender™ 3000 Standard
	108 Defender™ 2500 (Country Specific Items)
	60 Defender™ 2000
Indicators	61 T82XW Indicator
	62 T72XW Indicator
	63 i-DT61XWE Indicator
	64 i-DT61PW Indicator
	65 TD52XW Indicator
	66 TD52P Indicator
	67 i-DT33XW Indicator / i-DT33P Indicator
	68 T32M Indicator / T24PE Indicator
Bench Scale Bases	69 Defender™ 6000 Base
	70 Defender™ 5000 Stainless Steel Bases / Defender™ 3000 Stainless Steel Bases
	71 Defender™ 5000 Bases / Defender™ 3000 Bases
Compact Bench Scales	72 Ranger™ 7000
	73 Ranger™ 4000
	74 Ranger™ 4000 Count
	75 Ranger™ 3000
	76 Ranger™ 3000 Count
	78 Valor™ 7000
	79 Valor™ 4000
	80 Valor™ 3000
	81 Valor™ 2000
	82 Valor™ 1000
Shipping Scales	84 Courier™ 7000
	85 Courier™ 5000
	86 Courier™ 3000
	87 Courier™ 1000
Floor Scale Platforms	89 Defender™ 5000 Lift Platforms / Defender™ 5000 Low Profile Platforms
	90 DF 5000 Series Stainless Steel Floor Scale Platforms / DF 5000 Series Floor Scale Platforms
	91 DF 3000 Series Floor Scale Platforms (i-DF) / DF 3000 Series Floor Scale Platforms (DF)
	93 DFP Series Stainless Steel and Steel Pallet Platforms / DF 3000 Series Stainless Steel and Steel Pallet Platforms
OIML Weights	94 OIML E1 Weights / OIML E1 Weights Sets
	96 OIML E2 Weights / OIML E2 Weights Sets
	98 OIML F1 Weights / OIML F1 Weights Sets
	100 OIML F2 Weights / OIML F2 Weights Sets
	102 OIML M1 Weights / OIML M3 Weights / OIML Weight sets 1 mg / OIML Weight sets 1 g
Printers	103 SF-110B
	104 SF40A
	105 STP103
	106 Printer Cables
Software	106 Serial Port Data Collection Software
Accessories	107 PC Cables
	107 Pipette Adjustment Kit
	107 ION-100A Ionizer

Features

	INTERNAL ADJUSTMENT
	AutoCal™
	USB INTERFACE
	RS232 INTERFACE
	ANALOGUE INPUT
	EC TYPE APPROVED
	IP DUST AND WATER PROTECTED
	BATTERY (NOT SUPPLIED)
	RECHARGEABLE BATTERY
	KNOWLEDGE AND/OR TOOLS NEEDED
	BLUETOOTH
	WI-FI
	ETHERNET
	SD ALIBI MEMORY
	AUTO DOOR
	IONIZER
	NEW

Laboratory Balances
– applications & features

		Semi-Micro, Analytical & Precision										Precision Portable					
		EX 5(d)	EX	EX HC	AX 5(d)	AX	AX HC	PX 5(d)	PX	PX HC	PR	STX	SKX	NV / NVT	CX	CR	YA
Applications	Weighing	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Parts Counting	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Percentage Weighing	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Below Balance Weighing	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Custom unit - Factor entry	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Checkweighing	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Dynamic/Animal Weighing	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Gross/Net/Tare	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Totalization	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Statistics	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Differential Weighing	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Filling	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Density Determination	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Pipette Adjustment	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Display Hold	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Formulation	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Ingredients Costing	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	SQC	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Fill Weight Variation	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	High Point	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Mole weighing	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	User management	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	System log	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Operation	Mains Operation	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Battery Operation	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Rechargeable battery	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Display	LCD	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Backlit LCD	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Backlit 2-lines LCD	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Full colour VGA touch display	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Multilingual	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Output	RS232	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	USB	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Ethernet	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Bluetooth	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	GLP Time/Date	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Other	Auto Door	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Ionizer	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Internal Calibration	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Advanced AutoCal™	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	EC Type Approval	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Print-to-USB drive	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Page No.		11	15	20	12	16	21	13	17	22	18	25	26	27	28	29	29

● standard
 ●* on specific models
 ● option
 ●* accessory

To follow regulations,
 some applications may
 not be available on EC
 Type approved Models

Industry Scales
– applications & features

	Shipping Scales				Food Compact Scales					Compact Scales Ranger					Bench scales & indicators						
	Courier™ 1000	Courier™ 3000	Courier™ 5000	Courier™ 7000	Valor 1000	Valor 2000	Valor 3000	Valor 4000	Valor 7000	3000 Count	3000	4000 Count	4000	7000	D2000	D2500	D3000	D5000 VF	D6000	D7000	D8000
Max Resolution	2000	2000	3000	7500	7500	7500	40000	7500	30000	30000	30000	30000	30000	350000	20000	30000	30000	75000	75000	50000	10000K
Parts counting	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Check Weighing	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Percentage Weighing	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Totalling	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Dynamic weighing	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Weighing below	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Auto Tare	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Numeric Keyboard	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Alphanumeric keyboard	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Touchless sensor	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
LCD Display	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Backlit LCD Display	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
LED display	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Battery Operation	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Rechargeable battery	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
RS232 interface included	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
USB	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Ethernet	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
2nd RS interface	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Time / Date	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
GMP Data Output	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
IP 65 / 68 Protection	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Alibi Memory	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
EC Type Approval	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Bluetooth / Wi-Fi	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Page No.	87	86	85	84	83	82	80	79	78	76	75	74	73	72	68	108	67	65	63	62	61

● standard; ●* on specific models; ●**Approvable by Notified Body; ● Optional 1 - application depends on selected software

To follow regulations, some applications may not be available on EC Type approved models.

Industrial scales

Readability [d]	Capacity max (kg)	Model	Page	M* EC type approval	Internal Adjustment	Counting	RS232 Interface	USB Interface	IP 65 / 68	Rechargeable battery	Battery operation	Bluetooth / Wi-Fi	Ethernet
0.01	0.2	V31XH202	80						•				
0.01	0.4	V31XH402	80								•		
0.01	3	R71MHD3-EU	72		•	•	•	•					•
0.02	6	R71MHD6-EU	72										
0.05	1.5	V71P1502T	78										
0.05	1.5	R31P1502	75			•	•	•					•
0.05	1.5	RC31P1502	76			•	•	•					•
0.05	3	R71MD3-EU	72			•	•	•					•
0.1	0.3	V31XW301	80						•				
0.1	0.5	V31XS01	80								•		
0.1	2	V31XH2	80						•				
0.1	3	V71P3T	78					•				•	
0.1	3	R31P3	75			•	•	•					•
0.1	3	RC31P3	76				•	•					•
0.1	3	R41ME3	73			•	•	•					•
0.1	3	RC41M3	74			•	•	•					•
0.1	3	R71MHD3-EU-M	72	•	•	•	•	•					
0.1	4	V31XH4	80								•		
0.1	6	R71MD6-EU	72			•	•	•					•
0.1	15	R71MHD15-EU	72		•	•	•	•					•
0.1	35	R71MHD35-EU	72										•
0.2	1.5	V22PWE1501T	81						•				
0.2	1.5	V22XWE1501T	81						•				
0.2	1.5	V41PWE1501T	79						•				
0.2	1.5	V41XWE1501T	79						•				
0.2	6	V71P6T	78				•	•					•
0.2	6	R31P6	75			•	•	•					•
0.2	6	RC31P6	76			•	•	•					•
0.2	6	R41ME6	73			•	•	•					•
0.2	6	RC41M6	74			•	•	•					•
0.2	6	R71MHD6-EU-M	72	•	•	•	•	•					•
0.2	15	R71MD15-EU	72			•	•	•					•
0.5	3	I-D61PW3K1S5	47			•	•	•		•			
0.5	3	I-D61XWE3K1S6-EU	45			•	•	•			•		
0.5	3	I-D61PW3K1S6	46			•	•	•		•			
0.5	6	D52P6RQDR5	54			•	•	•					•
0.5	6	D52P6RQDR1	52			•	•	•					•
0.5	6	D52XW6RQDR1-EU	51			•	•	•					•
0.5	6	D52XW6WQDR5-EU	53			•	•	•					•
0.5	6	D52XW6WQDR6-EU	50			•	•	•					•
0.5	6	I-D61PW6K1S5	47			•	•	•		•			
0.5	6	I-D61PW6K1S6	46			•	•	•		•			
0.5	6	I-D61XWE6K1S6-EU	45			•	•	•			•		
0.5	1.5	V41PWE1501T-M	79	•					•				
0.5	1.5	V41XWE1501T-M	79	•					•				
0.5	1.5	V71P1502T-M	78			•	•	•					•
0.5	1.5	R31P1502-M	75			•	•	•					•
0.5	1.5	RC31P1502-M	76			•	•	•					•
0.5	3	V12P3	82										
0.5	3	V12P3T-M	82	•									
0.5	3	V22PWE3T	81						•				
0.5	3	V22XWE3T	81						•				
0.5	3	V31XW3	80								•		
0.5	3	V41PWE3T	79						•				
0.5	3	V41XWE3T	79						•				
0.5	3	R71MD3-EU-M	72	•		•	•	•					•
0.5	15	V71P15T	78										•
0.5	15	R31P15	75			•	•	•					•
0.5	15	RC31P15	76			•	•	•					•
0.5	15	R41ME15	73			•	•	•					•
0.5	15	RC41M15	74			•	•	•					•
0.5	35	R71MD35-EU	72										•
1	3	V31X3	80						•		•		
1	3	V41PWE3T-M	79	•					•				
1	3	V41XWE3T-M	79	•					•				
1	3	V71P3T-M	78			•	•	•					•
1	3	R31P3-M	75			•	•	•					•
1	3	RC31P3-M	76			•	•	•					•
1	3	R41ME3-M	73			•	•	•					•
1	3	RC41M3-M	74			•	•	•					•
1	3	I-D61PW3K1S5-M	47	•		•	•	•		•			
1	3	I-D61PW3K1S6-M	46	•		•	•	•		•			
1	3	I-D61XWE3K1S6-M	45	•		•	•	•			•		
1	6	V12P6	82										
1	6	V12P6T-M	82	•									
1	6	V22PWE6T	81						•				
1	6	V22XWE6T	81						•				
1	6	V31X6	80								•		
1	6	V31XW6	80								•		
1	6	V41PWE6T	79						•				
1	6	V41XWE6T	79						•				
1	6	R71MD6-EU-M	72	•		•	•	•					•
1	15	R71MHD15-EU-M	72	•	•	•	•	•					•
1	15	D52P15RQDR5	54			•	•	•					•
1	15	D52P15RQDR1	52			•	•	•					•
1	15	D52P15RQDL5	54			•	•	•					•
1	15	D52P15RQDL2	52			•	•	•					•
1	15	D52XW15RQDR1-EU	51			•	•	•					•
1	15	D52XW15RQDL2-EU	51			•	•	•					•
1	15	D52XW15WQDR5-EU	53			•	•	•					•
1	15	D52XW15WQDL5-EU	53			•	•	•					•
1	15	D52XW15WQDR6-EU	50			•	•	•					•
1	15	D52XW15WQDL7-EU	50			•	•	•					•
1	15	D52XW15WQR5-EU	53			•	•	•					•
1	15	D52XW15WQR6-EU	49			•	•	•					•
1	15	D52XW15WQL5-EU	53			•	•	•					•
1	15	D52XW15WQL7-EU	49			•	•	•					•
1	15	I-D61PW15K1R5	47			•	•	•		•			
1	15	I-D61PW15K1R6	46			•	•	•		•			
1	15	I-D61XWE15K1R6-EU	45			•	•	•			•		
1	30	V71P30T	78										•
1	30	R31P30	75			•	•	•					•
1	30	RC31P30	76			•	•	•					•
1	30	R41ME30	73			•	•	•					•
1	30	RC41M30	74			•	•	•					•
1	35	R71MHD35-EU-M	72	•	•	•	•	•					•

• Standard; • Optional; *not all features/applications are available on approved models

Readability [d]	Capacity max (kg)	Model	Page	M* EC type approval	Internal Adjustment	Counting	RS232 Interface	USB Interface	IP 65 / 68	Rechargeable battery	Battery operation	Bluetooth / Wi-Fi	Ethernet
1	60	R71MD60-EU	72			•	•	•		•			•
2	6	i-C52M6R EU	85					•			•		
2	6	V41PWE6T-M	79	•									
2	6	V41XWE6T-M	79	•									
2	6	V71P6T-M	78	•									
2	6	R31P6-M	75	•		•	•	•					•
2	6	RC31P6-M	76	•		•	•	•					•
2	6	R41ME6-M	73	•		•	•	•					•
2	6	RC41M6-M	74	•		•	•	•					•
2	6	D52P6RQDR5-M	54	•		•	•	•					•
2	6	D52P6RQDR1-M	52	•		•	•	•					•
2	6	D52XW6RQDR1-M	51	•		•	•	•					•
2	6	D52XW6WQDR5-M	53	•		•	•	•					•
2	6	D52XW6WQDR6-M	50	•		•	•	•					•
2	6	i-D61PW6K1S5-M	47	•		•	•	•			•		•
2	6	i-D61PW6K1S6-M	46	•		•	•	•			•		•
2	6	i-D61XWE6K1S6-M	45	•		•	•	•			•		•
2	15	i-C71M15R EU	84				•	•			•		
2	15	V12P15	82								•		
2	15	V12P15T-M	82	•							•		
2	15	V22PWE15T	81				•				•		
2	15	V22XWE15T	81								•		
2	15	V41PWE15T	79								•		
2	15	V41XWE15T	79								•		
2	15	R71MD15-EU-M	72	•		•		•			•		
2	15	i-D33P15B1R5	58			•	•	•			•		
2	15	i-D33P15B1R1	58			•	•	•			•		
2	15	i-D33XW15B1R5-EU	57			•	•	•			•		
2	15	i-D33XW15B1R1-EU	57			•	•	•			•		
2	15	i-D33XW15C1R5-EU	56			•	•	•			•		
2	15	i-D33XW15C1R6-EU	56			•	•	•			•		
2	30	D52P30RQDL5	54			•	•	•				•	
2	30	D52P30RQDL2	52			•	•	•				•	
2	30	D52XW30RQDL2-EU	51			•	•	•				•	
2	30	D52XW30WQDL5-EU	53			•	•	•				•	
2	30	D52XW30WQDL7-EU	50			•	•	•				•	
2	30	D52XW30WQL5-EU	53			•	•	•				•	
2	30	D52XW30WQL7-EU	49			•	•	•				•	
2	30	i-D61PW30K1R5	47			•	•	•			•		
2	30	i-D61PW30K1R6	46			•	•	•			•		
2	30	i-D61XWE30K1R6-EU	45			•	•	•			•		
5	9	i-C12P9 EU	87								•		
5	15	i-C71M15R-M	84	•			•	•			•		•
5	15	V41PWE15T-M	79	•						•			
5	15	V41XWE15T-M	79	•						•			
5	15	V71P15T-M	78	•			•	•		•			•
5	15	R31P15-M	75	•		•	•	•		•			•
5	15	RC31P15-M	76	•		•	•	•		•			•
5	15	R41ME15-M	73	•		•	•	•		•			•
5	15	RC41M15-M	74	•		•	•	•		•			•
5	15	i-D33P15B1R5-M	58	•		•	•	•		•	•		
5	15	i-D33P15B1R1-M	58	•		•	•	•		•	•		
5	15	i-D33XW15B1R5EU-M	57	•		•	•	•		•			
5	15	i-D33XW15B1R1EU-M	57	•		•	•	•		•			
5	15	i-D33XW15C1R5EU-M	56	•		•	•	•		•			
5	15	i-D33XW15C1R6EU-M	56	•		•	•	•		•			
5	15	D52P15RQDR5-M	54	•		•	•	•		•		•	
5	15	D52P15RQDR1-M	52	•		•	•	•		•		•	
5	15	D52P15RQDL5-M	54	•		•	•	•		•		•	
5	15	D52P15RQDL2-M	52	•		•	•	•		•		•	
5	15	D52XW15RQDR1-M	51	•		•	•	•		•		•	
5	15	D52XW15RQDL2-M	51	•		•	•	•		•		•	
5	15	D52XW15WQDR5-M	53	•		•	•	•		•		•	
5	15	D52XW15WQDR6-M	50	•		•	•	•		•		•	
5	15	D52XW15WQDL5-M	53	•		•	•	•		•		•	
5	15	D52XW15WQDL7-M	53	•		•	•	•		•		•	
5	15	D52XW15WQDR5-M	53	•		•	•	•		•		•	
5	15	D52XW15WQDR6-M	49	•		•	•	•		•		•	
5	15	D52XW15WQL5-M	53	•		•	•	•		•		•	
5	15	D52XW15WQL7-M	49	•		•	•	•		•		•	
5	15	i-D61PW15K1R5-M	47	•		•	•	•		•	•		
5	15	i-D61PW15K1R6-M	46	•		•	•	•		•	•		
5	15	i-D61XWE15K1R6-M	45	•		•	•	•		•	•		•
5	30	V12P30	82							•			
5	30	V12P30T-M	82	•						•			
5	30	V22PWE30T	81				•				•		
5	30	V22XWE30T	81								•		
5	30	D24PE30FR	60						•				
5	30	i-D33P30B1R5	58			•	•	•		•	•		•
5	30	i-D33P30B1R1	58			•	•	•		•	•		•
5	30	i-D33XW30B1R5-EU	57			•	•	•		•			•
5	30	i-D33XW30B1R1-EU	57			•	•	•		•			•
5	30	i-D33XW30C1R5-EU	56			•	•	•		•			•
5	30	i-D33XW30C1R6-EU	56			•	•	•		•			•
5	35	R71MD35-EU-M	72	•									
5	60	D52P60RQDL5	54			•	•	•				•	
5	60	D52P60RQDL2	52			•	•	•				•	
5	60	D52P60RQDX5	54			•	•	•				•	
5	60	D52P60RQDX2	52			•	•	•				•	
5	60	D52P60RQDV3	52			•	•	•				•	
5	60	D52XW60RQDL2-EU	51			•	•	•				•	
5	60	D52XW60RQDX2-EU	51			•	•	•				•	
5	60	D52XW60RQDV3-EU	51			•	•	•				•	
5	60	D52XW60WQDL5-EU	53			•	•	•				•	
5	60	D52XW60WQDL7-EU	50			•	•	•				•	
5	60	D52XW60WQDX5-EU	53			•	•	•				•	
5	60	D52XW60WQDX7-EU	50			•	•	•				•	
5	60	D52XW60WQDV8-EU	50			•	•	•				•	
5	60	D52XW60WQL5-EU	53			•	•	•				•	
5	60	D52XW60WQL7-EU	49			•	•	•				•	
5	60	D52XW60WQX5-EU	53			•	•	•				•	
5	60	D52XW60WQX7-EU	49			•	•	•				•	
5	60	D52XW60WQV8-EU	49			•	•	•				•	
5	60	i-D61PW60K1L5	47			•	•	•			•		
5	60	i-D61PW60K1L7	46			•	•	•			•		
5	60	i-D61XWE60K1L7-EU	45			•	•	•			•		

Industrial scales

Readability [d]	Capacity max (kg)	Model	Page	M* EC type approval	Internal Adjustment	Counting	RS232 Interface	USB Interface	IP 65 / 68	Rechargeable battery	Battery operation	Bluetooth / Wi-Fi	Ethernet
10	20	i-C12P20 EU	87										
10	30	i-C52M30R EU	85										
10	30	V71P30T-M	78										
10	30	R31P30-M	75										
10	30	RC31P30-M	76										
10	30	R41ME30-M	73										
10	30	RC41M30-M	74										
10	30	i-D25P30B2R1-M	108										
10	30	i-D33P30B1R5-M	58										
10	30	i-D33P30B1R1-M	58										
10	30	i-D33XW30B1R5EU-M	57										
10	30	i-D33XW30B1R1EU-M	57										
10	30	i-D33XW30C1R5EU-M	56										
10	30	i-D33XW30C1R6EU-M	56										
10	30	D52P30RQDL5-M	54										
10	30	D52P30RQDL2-M	52										
10	30	D52XW30RQDL2-M	51										
10	30	D52XW30WQDL5-M	53										
10	30	D52XW30WQDL7-M	50										
10	30	D52XW30WQL5-M	53										
10	30	D52XW30WQL7-M	49										
10	30	i-D61PW30K1R5-M	47										
10	30	i-D61PW30K1R6-M	46										
10	30	i-D61XWE30K1R6-M	45										
10	60	i-C71M60R EU	84										
10	60	i-C71M60L EU	84										
10	60	R71MD60-EU-M	72										
10	60	D24PE60FR	60										
10	60	D24PE60FL	60										
10	60	i-D33P60B1R5	58										
10	60	i-D33P60B1R1	58										
10	60	i-D33P60B1L5	58										
10	60	i-D33P60B1L2	58										
10	60	i-D33XW60B1R5-EU	57										
10	60	i-D33XW60B1R1-EU	57										
10	60	i-D33XW60C1R5-EU	56										
10	60	i-D33XW60C1R6-EU	56										
10	60	i-D33XW60C1L5-EU	56										
10	60	i-D33XW60C1L7-EU	56										
10	150	D52P150RQDX5	54										
10	150	D52P150RQDX2	52										
10	150	D52P150RQDV3	52										
10	150	D52P150RTDV3	52										
10	150	D52XW150RQDX2-EU	51										
10	150	D52XW150RQDV3-EU	51										
10	150	D52XW150RTDV3-EU	51										
10	150	D52XW150WQDX5-EU	53										
10	150	D52XW150WQDX7-EU	50										
10	150	D52XW150WQDV8-EU	50										
10	150	D52XW150WTDV8-EU	50										
10	150	D52XW150WQX5-EU	53										
10	150	D52XW150WQX7-EU	49										
10	150	D52XW150WQV8-EU	49										
10	150	i-D61PW150K1L5	47										
10	150	i-D61PW150K1L7	46										
10	150	i-D61XWE150K1L7-EU	45										
20	35	i-C31M35R EU	86										
20	50	i-C52M50R EU	85										
20	50	i-C52M50L EU	85										
20	60	i-C71M60R-M	84										
20	60	i-C71M60L-M	84										
20	60	i-D25P60B2R1-M	108										
20	60	i-D25P60B2L2-M	108										
20	60	i-D33P60B1R5-M	58										
20	60	i-D33P60B1R1-M	58										
20	60	i-D33P60B1L5-M	58										
20	60	i-D33P60B1L2-M	58										
20	60	i-D33XW60B1R5EU-M	57										
20	60	i-D33XW60B1R1EU-M	57										
20	60	i-D33XW60B1L5EU-M	57										
20	60	i-D33XW60B1L2EU-M	57										
20	60	i-D33XW60C1R5EU-M	56										
20	60	i-D33XW60C1R6EU-M	56										
20	60	i-D33XW60C1L5EU-M	56										
20	60	i-D33XW60C1L7EU-M	56										
20	60	D52P60RQDL5-M	54										
20	60	D52P60RQDL2-M	52										
20	60	D52P60RQDX5-M	54										
20	60	D52P60RQDX2-M	52										
20	60	D52P60RQDV3-M	52										
20	60	D52XW60RQDL2-M	51										
20	60	D52XW60RQDX2-M	51										
20	60	D52XW60RQDV3-M	51										
20	60	D52XW60WQDL5-M	53										
20	60	D52XW60WQDL7-M	50										
20	60	D52XW60WQDX5-M	53										
20	60	D52XW60WQDX7-M	50										
20	60	D52XW60WQDV8-M	50										
20	60	D52XW60WQL5-M	53										
20	60	D52XW60WQL7-M	49										
20	60	D52XW60WQX5-M	53										
20	60	D52XW60WQX7-M	49										
20	60	D52XW60WQV8-M	49										
20	60	i-D61PW60K1L5-M	47										
20	60	i-D61PW60K1L7-M	46										
20	60	i-D61XWE60K1L7-M	45										
20	150	i-C71M150L EU	84										
20	150	i-C71M150X EU	84										
20	150	D24PE150FX	60										
20	150	D24PE150FL	60										
20	150	i-D33P150B1L5	58										
20	150	i-D33P150B1L2	58										
20	150	i-D33P150B1X5	58										
20	150	i-D33P150B1X2	58										

Readability [d]	Capacity max (kg)	Model	Page	M* EC type approval	Internal Adjustment	Counting	RS232 Interface	USB Interface	IP 65 / 68	Rechargeable battery	Battery operation	Bluetooth / Wi-Fi	Ethernet
20	150	i-D33XW150B1L5-EU	57			•	•	•	•	•			•
20	150	i-D33XW150B1L2-EU	57			•	•	•	•	•			•
20	150	i-D33XW150B1X5-EU	57			•	•	•	•	•			•
20	150	i-D33XW150B1X2-EU	57			•	•	•	•	•			•
20	150	i-D33XW150C1L5-EU	56			•	•	•	•	•			•
20	150	i-D33XW150C1L7-EU	56			•	•	•	•	•			•
20	150	i-D33XW150C1X5-EU	56			•	•	•	•	•			•
20	150	i-D33XW150C1X7-EU	56			•	•	•	•	•			•
20	300	D52P300RQDV3	52			•	•	•	•	•		•	•
20	300	D52P300RTDV3	52			•	•	•	•	•		•	•
20	300	D52XW300RQDV3-EU	51			•	•	•	•	•		•	•
20	300	D52XW300RTDV3-EU	51			•	•	•	•	•		•	•
20	300	D52XW300WQDV8-EU	50			•	•	•	•	•		•	•
20	300	D52XW300WTDV8-EU	50			•	•	•	•	•		•	•
50	75	i-C12P75 EU	87										
50	75	i-C31M75R EU	86								•		
50	75	i-C31M75L EU	86								•		
50	100	i-C52M100L EU	85					•			•		
50	150	i-C71M150L-M	84	•							•		
50	150	i-C71M150X-M	84	•			•				•		•
50	150	i-D25P150B2L2-M	108	•									
50	150	i-D25P150B2X2-M	108	•								•	
50	150	i-D33P150B1L5-M	58		•								
50	150	i-D33P150B1L2-M	58		•								
50	150	i-D33P150B1X5-M	58		•					•			
50	150	i-D33P150B1X2-M	58		•					•	•		
50	150	i-D33XW150B1L5EU-M	57	•	•				•	•			
50	150	i-D33XW150B1L2EU-M	57	•	•				•	•			
50	150	i-D33XW150B1X5EU-M	57	•	•				•	•			
50	150	i-D33XW150B1X2EU-M	57	•	•				•	•			
50	150	i-D33XW150C1L5EU-M	56	•	•				•	•			
50	150	i-D33XW150C1L7EU-M	56	•	•				•	•			
50	150	i-D33XW150C1X5EU-M	56	•	•				•	•			
50	150	i-D33XW150C1X7EU-M	56	•	•				•	•			
50	150	D52P150RQDX5-M	54	•	•				•	•			
50	150	D52P150RQDX2-M	52	•	•				•	•			
50	150	D52P150RQDV3-M	52	•	•				•	•			
50	150	D52P150RTDV3-M	52	•	•				•	•			
50	150	D52XW150RQDX2-M	51	•	•				•	•			
50	150	D52XW150RQDV3-M	51	•	•				•	•			
50	150	D52XW150RTDV3-M	51	•	•				•	•			
50	150	D52XW150WQDX5-M	53	•	•				•	•			
50	150	D52XW150WQDX7-M	50	•	•				•	•			
50	150	D52XW150WQDV8-M	50	•	•				•	•			
50	150	D52XW150WTDV8-M	50	•	•				•	•			
50	150	D52XW150WQX5-M	53	•	•				•	•			
50	150	D52XW150WQX7-M	49	•	•				•	•			
50	150	D52XW150WQX8-M	49	•	•			•	•	•			
50	150	i-D61PW150K1L5-M	47	•	•			•	•		•		
50	150	i-D61PW150K1L7-M	46	•	•			•	•		•		
50	150	i-D61XWE150K1L7-M	45	•	•			•	•	•			
50	300	D24PE300FX	60			•	•			•			
50	300	D24PE300FV	60			•	•			•			
50	300	DF52P300F1S	38			•						•	
50	300	i-D33P300B1X5	58			•				•			
50	300	i-D33P300B1X2	58			•				•			
50	300	i-D33P300B1V3	58			•				•			
50	300	i-D33XW300B1X5-EU	57			•			•	•			
50	300	i-D33XW300B1X2-EU	57			•			•	•			
50	300	i-D33XW300B1V3-EU	57			•			•	•			
50	300	i-D33XW300C1X5-EU	56			•			•	•			
50	300	i-D33XW300C1X7-EU	56			•			•	•			
50	300	DF52P300F1R	38			•							•
50	300	DF52XW300G1S-EU	37			•			•	•			
50	300	DF52XW300G1R-EU	37			•			•	•			
50	600	D52P600RTDV3	52			•			•	•			
50	600	D52XW600RTDV3-EU	51			•			•	•			
50	600	D52XW600WTDV8-EU	50			•			•	•			
100	200	i-C31M200R EU	86								•		
100	200	i-C31M200L EU	86								•		
100	200	i-C52M200L EU	85					•			•		
100	200	i-C52M200X EU	85								•		
100	300	D24PE600FV	60			•	•			•			
100	300	i-D25P300B2X2-M	108	•									
100	300	i-D33P300B1X5-M	58			•				•		•	
100	300	i-D33P300B1X2-M	58			•				•			
100	300	i-D33P300B1V3-M	58			•				•		•	
100	300	i-D33XW300B1X5EU-M	57			•			•	•			
100	300	i-D33XW300B1X2EU-M	57			•			•	•			
100	300	i-D33XW300B1V3EU-M	57			•			•	•			
100	300	i-D33XW300C1X5EU-M	56			•			•	•			
100	300	i-D33XW300C1X7EU-M	56			•			•	•			
100	300	D52P300RQDV3-M	52			•				•			
100	300	D52P300RTDV3-M	52			•				•			
100	300	D52XW300RQDV3-M	51			•			•	•			
100	300	D52XW300RTDV3-M	51			•			•	•			
100	300	D52XW300WQDV8-M	50			•			•	•			
100	300	D52XW300WTDV8-M	50			•			•	•			
100	300	DF52P300F1S-M	38			•							
100	300	DF52P300F1R-M	38			•							
100	300	DF52XW300G1S-M	37			•			•	•			
100	300	DF52XW300G1R-M	37	•		•			•	•			
100	600	i-D33P600B1V3	58									•	
100	600	DF52P600F1S	38							•			
100	600	DF52P600F1R	38							•			
100	600	DF52P600F1LX	38							•			
100	600	DF52P600F1X	38							•			
100	600	DF52XW600G1S-EU	37						•	•			
100	600	DF52XW600G1R-EU	37						•	•			
100	600	DF52XW600G1LX-EU	37						•	•			
100	600	DF52XW600G1X-EU	37						•	•			
100	600	DFP52P600J1	43										
100	600	DFP52XW600K1-EU	43						•	•			
200	600	i-D33P600B1V3-M	58	•						•			
200	500	DF52P1500F1XV	38										
200	600	i-D33XW600R1V3-FUJ	57						•	•			

Industrial scales

Readability [d]	Capacity max (kg)	Model	Page	M* EC type approval	Internal Adjustment	Counting	RS232 Interface	USB Interface	IP 65 / 68	Rechargeable battery	Battery operation	Bluetooth / Wi-Fi	Ethernet
200	600	D52P600RTDV3-M	52	•		•	•	•		•		•	•
200	600	D52XW600RTDV3-M	51	•		•	•	•	•	•		•	•
200	600	D52XW600WTDV8-M	50	•		•	•	•		•		•	•
200	600	I-DF33P600B1R	40	•		•	•	•		•	•	•	•
200	600	I-DF33P600B1R-M	40	•		•	•	•		•	•	•	•
200	600	I-DF33XW600C1R-EU	39	•		•	•	•	•	•		•	•
200	600	I-DF33XW600C1R-EU-M	39	•		•	•	•	•	•		•	•
200	600	DF52P600F1S-M	38	•		•	•	•		•		•	•
200	600	DF52P600F1R-M	38	•		•	•	•		•		•	•
200	600	DF52P600F1LX-M	38	•		•	•	•		•		•	•
200	600	DF52P600F1X-M	38	•		•	•	•		•		•	•
200	600	DF52XW600G1S-M	37	•		•	•	•	•	•		•	•
200	600	DF52XW600G1R-M	37	•		•	•	•	•	•		•	•
200	600	DF52XW600G1LX-M	37	•		•	•	•	•	•		•	•
200	600	DF52XW600G1X-M	37	•		•	•	•	•	•		•	•
200	600	DFP52P600J1-M	43	•		•	•	•		•		•	•
200	600	DFP52XW600K1-M	43	•		•	•	•	•	•		•	•
200	1500	DF52P1500F1R	38	•		•	•	•		•		•	•
200	1500	DF52P1500F1LX	38	•		•	•	•		•		•	•
200	1500	DF52P1500F1X	38	•		•	•	•		•		•	•
200	1500	DF52XW1500G1R-EU	37	•		•	•	•	•	•		•	•
200	1500	DF52XW1500G1LX-EU	37	•		•	•	•	•	•		•	•
200	1500	DF52XW1500G1X-EU	37	•		•	•	•	•	•		•	•
200	1500	DFP52P1500J1	43	•		•	•	•		•		•	•
200	1500	DFP52XW1500K1-EU	43	•		•	•	•	•	•		•	•
500	1500	DF32M1500B5	41	•		•	•	•		•		•	•
500	1500	DF32M1500BR	41	•		•	•	•		•		•	•
500	1500	DF32M1500BL	41	•		•	•	•		•		•	•
500	1500	DF32M1500BX	41	•		•	•	•		•		•	•
500	1500	DF32M1500BLX	41	•		•	•	•		•		•	•
500	1500	I-DF33P1500B1R	40	•		•	•	•		•	•	•	•
500	1500	I-DF33P1500B1L	40	•		•	•	•		•	•	•	•
500	1500	I-DF33P1500B1X	40	•		•	•	•		•	•	•	•
500	1500	I-DF33P1500B1R-M	40	•		•	•	•		•	•	•	•
500	1500	I-DF33P1500B1L-M	40	•		•	•	•		•	•	•	•
500	1500	I-DF33P1500B1X-M	40	•		•	•	•		•	•	•	•
500	1500	I-DF33XW1500C1R-EU	39	•		•	•	•	•	•		•	•
500	1500	I-DF33XW1500C1L-EU	39	•		•	•	•	•	•		•	•
500	1500	I-DF33XW1500C1X-EU	39	•		•	•	•	•	•		•	•
500	1500	I-DF33XW1500C1R-EU-M	39	•		•	•	•	•	•		•	•
500	1500	I-DF33XW1500C1L-EU-M	39	•		•	•	•	•	•		•	•
500	1500	I-DF33XW1500C1X-EU-M	39	•		•	•	•	•	•		•	•
500	1500	I-DFP33P1500B1	44	•		•	•	•	•	•		•	•

Readability [d]	Capacity max (kg)	Model	Page	M* EC type approval	Internal Adjustment	Counting	RS232 Interface	USB Interface	IP 65 / 68	Rechargeable battery	Battery operation	Bluetooth / Wi-Fi	Ethernet
500	1500	I-DFP33P1500B1-M	44	•		•	•	•		•		•	•
500	1500	I-DFP33XW1500C1-EU	44	•		•	•	•		•		•	•
500	1500	I-DFP33XW1500C1-EU-M	44	•		•	•	•	•	•		•	•
500	1500	DF52P1500F1R-M	38	•		•	•	•		•		•	•
500	1500	DF52P1500F1LX-M	38	•		•	•	•		•		•	•
500	1500	DF52P1500F1X-M	38	•		•	•	•		•		•	•
500	1500	DF52P1500F1XV-M	38	•		•	•	•		•		•	•
500	1500	DF52XW1500G1R-M	37	•		•	•	•	•	•		•	•
500	1500	DF52XW1500G1LX-M	37	•		•	•	•	•	•		•	•
500	1500	DF52XW1500G1X-M	37	•		•	•	•	•	•		•	•
500	1500	DFP52P1500J1-M	43	•		•	•	•		•		•	•
500	1500	DFP52XW1500K1-M	43	•		•	•	•	•	•		•	•
500	3000	DF52P3000F1LX	38	•		•	•	•		•		•	•
500	3000	DF52P3000F1X	38	•		•	•	•		•		•	•
500	3000	DF52P3000F1XV	38	•		•	•	•		•		•	•
500	3000	DF52XW3000G1LX-EU	37	•		•	•	•	•	•		•	•
500	3000	DF52XW3000G1X-EU	37	•		•	•	•	•	•		•	•
1000	3000	DF32M3000BL	41	•		•	•	•		•		•	•
1000	3000	DF32M3000BX	41	•		•	•	•		•		•	•
1000	3000	DF32M3000BLX	41	•		•	•	•		•		•	•
1000	3000	I-DF33P3000B1L	40	•		•	•	•		•	•	•	•
1000	3000	I-DF33P3000B1X	40	•		•	•	•		•	•	•	•
1000	3000	I-DF33P3000B1L-M	40	•		•	•	•		•	•	•	•
1000	3000	I-DF33P3000B1X-M	40	•		•	•	•		•	•	•	•
1000	3000	I-DF33XW3000C1L-EU	39	•		•	•	•	•	•		•	•
1000	3000	I-DF33XW3000C1X-EU	39	•		•	•	•	•	•		•	•
1000	3000	I-DF33XW3000C1L-EU-M	39	•		•	•	•	•	•		•	•
1000	3000	I-DF33XW3000C1X-EU-M	39	•		•	•	•	•	•		•	•
1000	3000	DF52P3000F1LX-M	38	•		•	•	•		•		•	•
1000	3000	DF52P3000F1X-M	38	•		•	•	•		•		•	•
1000	3000	DF52P3000F1XV-M	38	•		•	•	•		•		•	•
1000	3000	DF52XW3000G1LX-M	37	•		•	•	•	•	•		•	•
1000	3000	DF52XW3000G1X-M	37	•		•	•	•	•	•		•	•

Semi-Micro, Analytical, Precision and Portable Precision balances

Readability [d]	Capacity max (g)	Model	Page	Internal Adjustment	ExCal	M* EC type approval	RS232 Interface	USB Interface	PCS	%	Battery operation**
0.01mg	52	PX125D	13	•			•		•		
0.01mg	52	PX125DM	13	•		•	•		•		
0.01mg	52	EX125D	11	•			•		•		
0.01mg	52	EX125DM	11	•			•		•		
0.01mg	82	AX85	12	•			•		•		
0.01mg	82	AX85M	12	•		•	•		•		
0.01mg/0.01mg	82/120	AX125D	12	•			•		•		
0.01mg/0.01mg	82/120	AX125DM	12	•			•		•		
0.01mg/0.01mg	102/220	AX225D	12	•			•		•		
0.01mg/0.01mg	102/220	AX25DM	12	•		•	•		•		
0.01mg	82	PX85	13	•			•		•		
0.01mg	82	PX85M	13	•		•	•		•		
0.01mg	82	PX225D	13	•			•		•		
0.01mg	82	PX225DM	13	•			•		•		
0.01mg	120	EX125	11	•			•		•		
0.01mg	120	EX125M	11	•		•	•		•		
0.01mg	120	EX225D	11	•			•		•		
0.01mg	120	EX225DM	11	•		•	•		•		
0.01mg	120	EX225D/AD	11	•			•		•		
0.01mg	120	EX225DM/AD	11	•		•	•		•		
0.01mg	220	EX225/AD	11	•			•		•		
0.01mg	220	EX225M/AD	11	•		•	•		•		
0.1mg	120	PX125D	13	•			•		•		
0.1mg	120	PX125DM	13	•			•		•		
0.1mg	120	PR124/E	18	•	•		•		•		
0.1mg	120	PR124	18	•			•		•		
0.1mg	120	PR124M	18	•			•		•		
0.1mg	120	PX124/E	17	•	•		•		•		
0.1mg	120	PX124	17	•			•		•		
0.1mg	120	PX124M	17	•		•	•		•		
0.1mg	120	AX124/E	16	•	•		•		•		
0.1mg	120	AX124	16	•			•		•		
0.1mg	120	AX124M	16	•		•	•		•		
0.1mg	120	EX124	16	•			•		•		
0.1mg	120	EX124/AD	15	•			•		•		
0.1mg	120	EX125D	11	•			•		•		
0.1mg	120	EX125DM	11	•		•	•		•		
0.1mg	220	PR224/E	18	•	•		•		•		
0.1mg	220	PR224	18	•			•		•		
0.1mg	220	PR224M	18	•			•		•		
0.1mg	220	PX224/E	17	•	•		•		•		
0.1mg	220	PX224	17	•			•		•		
0.1mg	220	PX224M	17	•		•	•		•		
0.1mg	220	PX225D	13	•			•		•		
0.1mg	220	PX225DM	13	•			•		•		
0.1mg	220	AX224/E	16	•	•		•		•		
0.1mg	220	AX224	16	•			•		•		
0.1mg	220	AX224M	16	•		•	•		•		
0.1mg	220	EX224	15	•			•		•		
0.1mg	220	EX224/AD	15	•		•	•		•		
0.1mg	220	EX224M/AD	15	•		•	•		•		
0.1mg	220	EX225D	11	•			•		•		
0.1mg	220	EX225DM	11	•		•	•		•		
0.1mg	220	EX225D/AD	11	•			•		•		
0.1mg	220	EX225DM/AD	11	•		•	•		•		
0.1mg	320	AX324	16	•			•		•		
0.1mg	320	AX324M	16	•		•	•		•		
0.1mg	320	EX324	15	•			•		•		
0.1mg	320	EX324/AD	15	•		•	•		•		
0.1mg	320	EX324M/AD	15	•		•	•		•		
0.1mg	320	EX324M/AD	15	•		•	•		•		
1mg	120	NV123	27	•			•		•		•
1mg	120	SKX123	26	•			•		•		•
1mg	120	STX123	25	•			•		•		•
1mg	220	STX223	25	•			•		•		•
1mg	220	NV223	27	•			•		•		•
1mg	220	PR223/E	23	•	•		•		•		•
1mg	220	PR223	23	•			•		•		•
1mg	220	PR223M	23	•			•		•		•
1mg	220	NV323	27	•	•		•		•		•
1mg	320	PX323/E	22	•			•		•		•
1mg	320	PX323	22	•			•		•		•
1mg	320	PX323M	22	•		•	•		•		•
1mg	220	AX223/E	21	•	•		•		•		•
1mg	220	AX223	21	•			•		•		•
1mg	220	AX223M	21	•		•	•		•		•
1mg	220	EX223/E	20	•	•		•		•		•
1mg	220	EX223	20	•			•		•		•
1mg	420	PR423/E	23	•	•		•		•		•
1mg	420	PR423	23	•			•		•		•
1mg	420	PR423M	23	•			•		•		•
1mg	420	AX423/E	21	•	•		•		•		•
1mg	420	AX423	21	•			•		•		•
1mg	420	AX423M	21	•		•	•		•		•
1mg	420	EX423/E	20	•	•		•		•		•
1mg	420	EX423	20	•			•		•		•
1mg	420	EX423M	20	•		•	•		•		•
1mg	520	PR523/E	23	•	•		•		•		•
1mg	520	PR523	23	•			•		•		•
1mg	520	PR523M	23	•			•		•		•
1mg	620	AX623/E	21	•	•		•		•		•
1mg	620	AX623	21	•			•		•		•
1mg	620	AX623M	21	•		•	•		•		•
1mg	620	PX623/E	22	•	•		•		•		•
1mg	620	PX623	22	•			•		•		•
1mg	620	PX623M	22	•		•	•		•		•
1mg	620	EX623	20	•			•		•		•
1mg	1100	EX1103	20	•			•		•		•
1mg	1100	EX1103M	20	•			•		•		•
0.01g	100	YA102	29	•			•		•		•
0.01g	220	NV222	27	•			•		•		•
0.01g	220	SKX222	26	•			•		•		•
0.01g	220	STX222	25	•			•		•		•
0.01g	420	NV422	27	•			•		•		•
0.01g	420	SKX422	26	•			•		•		•
0.01g	420	STX422	25	•			•		•		•
0.01g	420	AX422/E	21	•	•		•		•		•
0.01g	420	AX422	21	•			•		•		•
0.01g	620	NV622	27	•			•		•		•
0.01g	620	SKX622	26	•			•		•		•
0.01g	620	STX622	25	•			•		•		•
0.01g	820	AX822/E	21	•	•		•		•		•
0.01g	820	AX822	21	•			•		•		•
0.01g	1200	NV1202	27	•			•		•		•
0.01g	1200	SKX1202	26	•			•		•		•
0.01g	1200	STX1202	25	•			•		•		•
0.01g	1520	AX1502/E	21	•			•		•		•
0.01g	1520	AX1502	21	•			•		•		•

Readability [d]	Capacity max (g)	Model	Page	Internal Adjustment	ExCal	M* EC type approval	RS232 Interface	USB Interface	PCS	%	Battery operation**
0.01g	1520	AX1502M	21	•		•	•	•	•	•	
0.01g	1600	PR1602/E	23		•		•		•	•	
0.01g	1600	PR1602	23	•			•		•	•	
0.01g	1600	PR1602M	23	•			•		•	•	
0.01g	2200	NV2202	27		•		•	•	•	•	•
0.01g	2200	SKX2202	26		•		•	•	•	•	•
0.01g	2200	STX2202	25		•		•		•	•	•
0.01g	2200	PR2202/E	23				•		•	•	•
0.01g	2200	PR2202	23	•			•		•	•	•
0.01g	2200	PR2202M	23	•			•		•	•	•
0.01g	2200	AX2202/E	21		•		•		•	•	•
0.01g	2200	AX2202	21	•			•	•	•	•	•
0.01g	2200	AX2202M	21	•		•	•		•	•	•
0.01g	2200	EX2202/E	20		•		•		•	•	•
0.01g	2200	EX2202	20	•			•		•	•	•
0.01g	3200	NV3202	27		•		•	•	•	•	•
0.01g	3200	PX3202/E	22		•		•	•	•	•	•
0.01g	3200	PX3202	22	•			•		•	•	•
0.01g	3200	PX3202M	22	•		•	•		•	•	•
0.01g	4200	PR4202/E	23		•		•		•	•	•
0.01g	4200	PR4202	23	•			•		•	•	•
0.01g	4200	PR4202M	23		•		•		•	•	•
0.01g	4200	AX4202/E	21				•	•	•	•	•
0.01g	4200	AX4202	21	•		•	•		•	•	•
0.01g	4200	AX4202M	21			•	•		•	•	•
0.01g	4200	EX4202/E	20		•		•		•	•	•
0.01g	4200	EX4202	20	•			•		•	•	•
0.01g	4200	EX4202M	20	•		•	•		•	•	•
0.01g	6200	AX6202	21				•		•	•	•
0.01g	6200	AX6202M	21	•		•	•		•	•	•
0.01g	6200	AX6202/E	21	•			•		•	•	•
0.01g	6200	PX6202/E	22		•		•		•	•	•
0.01g	6200	PX6202	22	•			•		•	•	•
0.01g	6200	PX6202M	22			•	•		•	•	•
0.01g	6200	EX6202/E	20		•		•		•	•	•
0.01g	6200	EX6202	20	•			•		•	•	•
0.01g	10200	EX10202	20	•			•		•	•	•
0.01g	10200	EX10202M	20	•		•	•		•	•	•
0.05g	300	YA302	29		•						•
0.1g	220	CR221	28		•						•
0.1g	220	CX221	28		•						•
0.1g	220	NV221	27		•		•		•	•	•
0.1g	420	SKX421	26		•		•	•		•	•
0.1g	420	STX421	25		•		•	•	•	•	•
0.1g	500	YA501	29		•						•
0.1g	620	CR621	28		•						•
0.1g	620	CX621	28		•						•
0.1g	620	NV621	27		•		•		•	•	•
0.1g	620	SKX621	26		•		•	•		•	•
0.1g	620	STX621	25		•		•	•	•	•	•
0.1g	1200	CX1201	28		•		•		•	•	•
0.1g	1200	NV1201	27		•		•	•		•	•
0.1g	2200	NV2201	27		•		•		•	•	•
0.1g	2200	NVT2201	27		•		•	•	•	•	•
0.1g	2200	SKX2201	26		•		•	•	•	•	•
0.1g	2200	SKX2201	26		•		•	•	•	•	•
0.1g	2200	STX2201	25		•		•	•	•	•	•
0.1g	2200	AX2201/E	21		•		•		•	•	•
0.1g	2200	AX2201	21	•			•		•	•	•
0.1g	4200	NVT4201	27		•		•	•		•	•
0.1g	4200	PR4201/E	23		•		•		•	•	•
0.1g	4200	PR4201	23	•			•		•	•	•
0.1g	4200	PR4201M	23		•		•		•	•	•
0.1g	4200	PX4201/E	22		•		•		•	•	•
0.1g	4200	PX4201	22	•			•	•	•	•	•
0.1g	4200	PX4201M	22	•		•	•		•	•	•
0.1g	4200	AX4201/E	21		•		•		•	•	•
0.1g	4200	AX4201	21	•			•		•	•	•
0.1g	6200	NVT6201	27		•		•	•	•	•	•
0.1g	6200	SKX6201	26		•		•	•	•	•	•
0.1g	6200	STX6201	25		•		•	•	•	•	•
0.1g	6200	PR6201/E	23		•		•		•	•	•
0.1g	6200	PR6201	23	•			•		•	•	•
0.1g	6200	PR6201M	23		•		•		•	•	•
0.1g	6200	EX6201/E	20		•		•	•	•	•	•
0.1g	6200	EX6201	20	•			•		•	•	•
0.1g	8200	AX8201/E	21		•		•		•	•	•
0.1g	8200	AX8201	21	•			•		•	•	•
0.1g	8200	AX8201M	21	•		•	•		•	•	•
0.1g	12000	AX12001	21	•			•		•	•	•
0.1g	12000	AX12001M	21			•	•		•	•	•
0.1g	12000	AX12001/E	21		•		•		•	•	•
0.1g	12000	PX12001	22	•			•		•	•	•
0.1g	12000	NV12001/E	27		•		•		•	•	•
0.1g	10200	NVT10201	27				•	•	•	•	•
0.1g	10200	EX10201	20				•		•	•	•
0.1g	10200	EX10201M	20		•		•		•	•	•
0.1g	12000	EX12001	20	•			•		•	•	•
0.1g	12000	EX12001M	20	•		•	•		•	•	•
0.1g	24000	EX24001	20	•			•		•	•	•
0.1g	24000	EX24001M	20	•		•	•		•	•	•
0.1g	35000	EX35001	20				•		•	•	•
0.1g	35000	EX35001M	20	•		•	•		•	•	•
0.5g	1500	NVT1601M	27				•	•	•	•	•
1g	2200	CR2200	28		•						•
1g	2200	CX2200	28		•						•
1g	2200	NVT2200	27		•						•
1g	3200	NVT3200M	27		•		•		•	•	•
1g	5200	CR5200	28		•						•
1g	5200	CX5200	28		•						•
1g	6200	NVT6200	27		•		•		•	•	•
1g	8200	SKX8200	26		•		•	•		•	•
1g	8200	STX8200	25		•		•	•	•	•	•
1g	12000	NVT12000	27		•		•	•	•	•	•
1g	12000	NVT12000	27		•		•	•	•	•	•
2g	6400	NVT6400M	27		•		•		•	•	•
5g	16000	NVT16000M	27		•		•	•	•	•	•

Do more with OHAUS



EXPLORER™ Semi-Micro

Built for sophisticated, professional weighing



Application

Weighing, percent weighing, parts counting, check weighing, dynamic/animal weighing, filling, totalization, formulation, fill weight variation, differential weighing, density determination, peak hold, pipette adjustment, SQC

Display

Full-colour VGA graphic display (diagonal 5.7"/145mm), detachable, 4-wire resistive touch screen. QWERTY keyboard and numeric keypad to quickly input GLP and GMP data and other application data

Operation

AC Adapter (included)

Communication

Easy access communication ports including 2 units of USB port, RS232 and an optional 4th Ethernet port, GLP/GMP data output with real time clock, Direct Data Transfer, print to USB drive

Construction

Metal base, ABS top housing, stainless steel pan, anti-static glass draftshield with flip-top door, replaceable in-use cover, automatic doors on specific models

Design Features

AutoCal™ fully-automatic internal adjustment system equipped with 2 internal adjustment weights, routine test and minimum weight determination, fast stabilization time, four touchless sensors, built-in ionizer with automatic doors models, up to 14 operating languages, approved models available, save to USB function, menu lock switch, security bracket, integral weigh below hook for below balance weighing applications, removable stainless steel weighing platform, stability indicator, overload and underload indicators, auto standby



Explorer's weighing cells precision-machined from a solid metal block, together with AutoCal™ provides exceptional accuracy and durability suitable for laboratory and industrial environments.

Features such as four-level user management with password protection, and un-editable system log plus a high level of configurability make Explorer balances well suited for regulated applications.

A detachable terminal with a large colour touchscreen, programmable IR sensors, optional built-in ionizer & automatic doors, and a frameless draftshield make Explorer balances easy & convenient to use.

Model	EC Type Approved	Auto Door	Ionizer	Capacity	Readability (d)	(e)	Platform	Cal. weight, class	Item N°	Price €
EX125D				52g/120g	0,01mg/ 0,1mg		Ø 80mm	100g, E2	30218998	5865,00
EX125				120g	0,01mg		Ø 80mm	100g, E2	30218999	6070,00
EX225D				120g/220g	0,01mg/ 0,1mg		Ø 80mm	200g, E2	30219000	6365,00
EX225D/AD		•	•	120g/220g	0,01mg/ 0,1mg		Ø 80mm	200g, E2	30219001	8590,00
EX225/AD		•	•	220g	0,01mg		Ø 80mm	200g, E2	30219002	13710,00
EX125DM	•			52g/120g	0,01mg/ 0,1mg	1mg	Ø 80mm	100g, E2	30139515	6305,00
EX125M	•			120g	0,01mg	1mg	Ø 80mm	100g, E2	30139516	6520,00
EX225DM	•			120g/220g	0,01mg/ 0,1mg	1mg	Ø 80mm	200g, E2	30139517	6810,00
EX225DM/AD	•	•	•	120g/220g	0,01mg/ 0,1mg	1mg	Ø 80mm	200g, E2	30139518	9030,00
EX225M/AD	•	•	•	220g	0,01mg	1mg	Ø 80mm	200g, E2	30139519	14140,00

ADVENTURER™ Semi-Micro

Intuitive Balances Designed for Routine Weighing



**EASY OPERATION
& TIME SAVING**



Application

Weighing (17 units + custom unit), parts counting, percentage weighing, check weighing, animal/dynamic weighing, totalization/statistics, formulation, density determination, display hold, batch printing

Display

4.3" (109 mm) Full-colour VGA graphic touchscreen with user-controlled brightness, 6 mechanical buttons

Operation

AC adapter (included)

Communication

RS232, USB Host, USB Device (included), GLP/GMP compliant data output with real-time clock

Construction

Metal base, ABS top housing, stainless steel pan, glass draftshield with two piece top mounted side doors and sliding top door (1mg models only), illuminated up-front level indicator, integral weigh below hook, security bracket, adjustment lock, full housing replaceable in-use cover

Design Features

Selectable environmental filters, auto tare, user-selectable span adjustment points, software lockout and reset menu, user-selectable communication settings and data print options, user definable project and user IDs, up to 9 operating languages



Adventurer balances feature a color touchscreen, icon-based user interface, and an ergonomic design – making them easy to configure and use.

Features such as specialized weighing modes, multiple connectivity options, and AutoCal™ provide versatility and flexibility for a variety of applications.

Durable construction, large weighing surfaces, a space-saving draftshield design, and full housing in-use cover allow for use in lab, education and industrial environments.

Model	EC Type Approved	Capacity	Readability (d)	(e)	Platform	Cal. weight, class	Item N°	Price €
AX85		82g	0,01mg		Ø 80mm	50g, E2	30852308	4400,00
AX125D		82/120g	0,01mg/ 0,1mg		Ø 80mm	100g, E2	30852307	4800,00
AX225D		102/220g	0,01mg/ 0,1mg		Ø 80mm	200g, E2	30852309	5100,00
AX85M	•	85g	0,01mg	1mg	Ø 80mm	50g, E2	30852311	4660,00
AX125DM	•	82/120g	0,01mg/ 0,1mg	1mg	Ø 80mm	100g, E2	30852310	5090,00
AX225DM	•	102/220g	0,01mg/ 0,1mg	1mg	Ø 80mm	200g, E2	30852312	5400,00

PIONEER™ Semi-Micro

Combining Economy and High Performance for Essential Weighing



OPTION

**\$ECONOMY
& ACCURACY**

Application

Basic weighing, parts counting, percent weighing, animal weighing, density determination, formulation

Display

Backlit Liquid Crystal Display (LCD); Second line display for additional information or guidance (dot matrix)

Operation

AC adapter (included)

Communication

RS232 and USB device (included); GLP/GMP data output with real-time clock

Construction

Metal base, ABS top housing, removable stainless steel pan, glass draftshield, integrated weigh-below-hook, security bracket, calibration lock and replaceable in-use cover

Design Features

User-selectable environmental filters and brightness settings, auto-tare, auto-dim, user-selectable span calibration points, software lockout and reset menu, user-selectable communication settings and data print options, user-definable project and user IDs, software overload/underload indicator, stability indicator, four operating languages



The PX offers high accuracy and repeatability for essential weighing applications in laboratory, industrial and education settings at an economical price point.

Featuring a cast metal lower housing, sub-pan and stainless steel weighing pan, the PX is durably constructed for versatile, long-term use.

Pioneer features a second line display for additional information or guidance, a static removal bar for convenient grounding, and USB connectivity.

Model	EC Type Approved	AutoCal™	Capacity	Readability (d)	(e)	Platform	Cal. weight, class	Item N°	Price €
PX125D		•	52g/120g	0,01mg/0,1mg		Ø 80mm	100g, E2	30478196	3015,00
PX85		•	82g	0,01mg		Ø 80mm	50g+20g, E2	30485340	3105,00
PX225D		•	82g/220g	0,01mg/0,1mg		Ø 80mm	200g, E2	30478197	3230,00
PX125DM	•	•	52g/120g	0,01mg/0,1mg	1mg	Ø 80mm	100g, E2	30480168	3240,00
PX85M	•	•	82g	0,01mg	1mg	Ø 80mm	50g+20g, E2	30485341	3330,00
PX225DM	•	•	82g/220g	0,01mg/0,1mg	1mg	Ø 80mm	200g, E2	30480169	3455,00

Accessories for Semi-Micro Balances

Accessories	Explorer™ Semi-Micro	Adventurer™ Semi-Micro	Pioneer™ Semi-Micro	Item N°	Price €
Density determination kit	•	•	•	80253384	680,00
Sinker Glass for Density kit	•	•	•	83034024	355,00
ION-100A EU stand-alone ionizer	•	•	•	30095929	1460,00
Grid pan	•			30298582	197,00
Pipette adjustment kit	•			30252145	1230,00
Ethernet interface	•			83021082	138,00
Interface Kit, RS232-USB	•	•	•	30304101	78,25
Cable, USB A to Mini USB, Black	•	•		28120263	21,20
In-Use Cover for EX	•			83033633	31,75
In-Use Cover for AX		•		30111792	33,25
In-Use Cover for PX			•	30372546	30,00
Dust Cover for EX	•			30093334	41,00
Anti-theft device (Models over 8 kg)	•	•	•	80850000	68,50
Terminal extension cable 3m	•			83021083	22,60
Terminal extension cable 9m	•			30078078	30,50
SF40A impact printer	•	•	•	30064202	715,00
Auxiliary display, AD7-RS	•	•	•	30472064	96,00
Anti theft device	•	•	•	80850043	211,00



DENSITY DETERMINATION KIT
80253384



SINKER GLASS FOR DENSITY KIT
83034024



ION-100A EU STANDALONE IONIZER
30095929



GRID PAN
30298582



PIPETTE ADJUSTMENT KIT
30252145



ETHERNET INTERFACE
83021082



INTERFACE KIT, RS232-USB
30304101



IN-USE COVER
30372546



IN-USE COVER
83033633



DUST COVER
30093334



ANTI THEFT DEVICE
80850000



TERMINAL EXTENSION CABLE 3M
83021083



TERMINAL EXTENSION CABLE 9M
30078078




SF40A IMPACT PRINTER
30064202



AUXILIARY DISPLAY, AD7-RS
30472064

EXPLORER™ Analytical

Highest performance, most technically advanced balance from OHAUS available also with auto door feature



Application

Weighing (17 units + custom units), parts counting, percentage weighing, animal/dynamic weighing, check weighing, filling, totalization, formulation, differential weighing, density determination, peak hold, ingredient costing, pipette adjustment, fill weight variation, gross/net/tare weighing

Display

Full-colour VGA graphic display (diagonal 5.7"/145mm), detachable, 4-wire resistive touch screen. QWERTY keyboard and numeric keypad to quickly input GLP and GMP data and other application data

Operation

AC adapter (included)

Communication


Easy access communication ports including 2 units of USB port, RS232 and an optional 4th Ethernet port, GLP and GMP including date and time, Direct Data Transfer, print to USB drive

Construction

Metal base, ABS top housing, stainless steel pan, anti-static glass draftshield with flip-top door, replaceable in-use cover, automatic door on specific models

Design Features

3-doors frameless draftshield with removable panels with optional automatic draftshield door, modular display design with adjustable viewing angle (attach to base, wall or optional tower), AutoCal™ automatic internal Adjustment, level assist feature with illuminated level indicator, selectable environmental and auto-print settings, integral weigh below hook, protective in-use cover, software lock switch



Equipped with a weighing cell that's precision-machined from a solid metal block, Explorer balances are extremely accurate and durable, and are ideal for laboratory and industrial environments.

Features such as four-level user management with password protection, and un-editable system log plus a high level of configurability make Explorer balances well suited for regulated applications.

A detachable terminal with a large colour touchscreen, programmable IR sensors, automatic doors, and a frameless draftshield make Explorer analytical balances easy and convenient to use.

Model	EC Type Approved	Auto Door	Capacity	Readability (d)	(e)	Platform	Cal. weight, class	Item N°	Price €
EX124			120g	0.1mg		Ø 90mm	100g, E2	83021331	3580,00
EX224			220g	0.1mg		Ø 90mm	200g, E2	83021332	3980,00
EX324			320g	0.1mg		Ø 90mm	300g, E2	83021335	4590,00
EX124/AD		•	120g	0.1mg		Ø 90mm	100g, E2	30061976	4625,00
EX224/AD		•	220g	0.1mg		Ø 90mm	200g, E2	30061977	5015,00
EX324/AD		•	320g	0.1mg		Ø 90mm	300g, E2	30061978	5625,00
EX224M	•		220g	0.1mg	1mg	Ø 90mm	200g, E2	83021333	4240,00
EX224M/AD	•	•	220g	0.1mg	1mg	Ø 90mm	200g, E2	30061995	5280,00
EX324M	•		320g	0.1mg	1mg	Ø 90mm	300g, E2	83021336	4850,00
EX324M/AD	•	•	320g	0.1mg	1mg	Ø 90mm	300g, E2	30061996	5895,00

ADVENTURER™ Analytical

Intuitive Balances Designed for Routine Weighing



**EASY OPERATION
& TIME SAVING**

Application

Weighing (17 units + custom unit), parts counting, percentage weighing, check weighing, animal/dynamic weighing, totalization/statistics, formulation, density determination, display hold, batch printing

Display

4.3" (109 mm) Full-colour VGA graphic touchscreen with user-controlled brightness, 6 mechanical buttons

Operation

AC adapter (included)

Communication

RS232, USB Host, USB Device (included), GLP/GMP compliant data output with real-time clock

Construction

Metal base, ABS top housing, stainless steel pan, glass draftshield with two piece top mounted side doors and sliding top door (1mg models only), illuminated up-front level indicator, integral weigh below hook, security bracket, adjustment lock, full housing replaceable in-use cover

Design Features

Selectable environmental filters, auto tare, user-selectable span adjustment points, software lockout and reset menu, user-selectable communication settings and data print options, user definable project and user IDs, up to 9 operating languages



Adventurer balances feature a color touchscreen, icon-based user interface, and an ergonomic design – making them easy to configure and use.

Features such as specialized weighing modes, multiple connectivity options, and AutoCal™ provide versatility and flexibility for a variety of applications.

Durable construction, large weighing surfaces, a space-saving draftshield design, and full housing in-use cover allow for use in lab, education and industrial environments.

Model	EC Type Approved	AutoCal™	Capacity	Readability (d)	(e)	Platform	Cal. weight, class	Item N°	Price €
AX124		•	120g	0.1mg		Ø 90mm	100g, E2	30122610	2240,00
AX224		•	220g	0.1mg		Ø 90mm	200g, E2	30122611	2340,00
AX324		•	320g	0.1mg		Ø 90mm	200g+100g, E2	30122613	2410,00
AX124M	•	•	120g	0.1mg	1mg	Ø 90mm	100g, E2	30100620	2450,00
AX224M	•	•	220g	0.1mg	1mg	Ø 90mm	200g, E2	30100621	2560,00
AX324M	•	•	320g	0.1mg	1mg	Ø 90mm	200g+100g, E2	30100622	2630,00
AX124/E			120g	0.1mg		Ø 90mm	100g, E2	30122657	1970,00
AX224/E			220g	0.1mg		Ø 90mm	200g, E2	30122612	2080,00

PIONEER™ Analytical

Combining Economy and High Performance for Essential Weighing



**\$ECONOMY
& ACCURACY**

Application

Basic weighing, parts counting, percent weighing, animal weighing, density determination, formulation

Display

Backlit Liquid Crystal Display (LCD); second line display for additional information or guidance (dot matrix)

Operation

AC adapter (included)

Communication

RS232 and USB device (included); GLP/GMP data output with real-time clock

Construction

Metal base, ABS top housing, removable stainless steel pan, glass draftshield, integrated weigh-below-hook, security bracket, calibration lock and replaceable in-use cover

Design Features

User-selectable environmental filters and brightness settings, auto-tare, auto-dim, user-selectable span calibration points, software lockout and reset menu, userselectable communication settings and data print options, user-definable project and user IDs, software overload/underload indicator, stability indicator, four operating languages



The PX offers high accuracy and repeatability for essential weighing applications in laboratory, industrial and education settings at an economical price point.

Featuring a cast metal lower housing, sub-pan and stainless steel weighing pan, the PX is durably constructed for versatile, long-term use.

Pioneer features a second line display for additional information or guidance, a static removal bar for convenient grounding, and USB connectivity.

Model	EC Type Approved	Internal Adjustment	Capacity	Readability (d)	(e)	Platform	Cal. weight, class	Item N°	Price €
PX124		•	120g	0,1mg		Ø 90mm	100g, E2	30429803	1530,00
PX224		•	220g	0,1mg		Ø 90mm	200g, E2	30429804	1640,00
PX124M	•	•	120g	0,1mg	1mg	Ø 90mm	100g, E2	30430010	1745,00
PX224M	•	•	220g	0,1mg	1mg	Ø 90mm	200g, E2	30430013	1855,00
PX124/E			120g	0,1mg		Ø 90mm	100g, E2	30429810	1335,00
PX224/E			220g	0,1mg		Ø 90mm	200g, E2	30429811	1450,00

PR SERIES Analytical

Designed for Routine Weighing Applications in Your Workplace



**ROUTINE WEIGHING
FOR ANY BUDGET**

Application

Basic weighing, parts counting, percent weighing

Display

Backlit Liquid Crystal Display (LCD)

Operation

AC adapter (included)

Communication

RS232 (one port included)

Construction

ABS top housing, removable stainless steel pan, glass draftshield with sliding top door, integrated weigh-below-hook, security bracket, calibration lock, replaceable in-use cover

Design features

User-selectable environmental filters and brightness settings, auto-tare, auto-dim, user-selectable span calibration points, software lockout and reset menu, user-selectable communication settings and data print options, user-definable project and user IDs, software overload/underload indicator, stability indicator



Basic Functionality for Routine Weighing Applications.

Designed for Uncomplicated Operation with Easy-to-Use Display and Interface.

Smart Design and Durable Construction.

Model	EC Type Approved	Internal Adjustment	Capacity	Readability (d)	(e)	Platform	Cal. weight, class	Item N°	Price €
PR124		•	120g	0,1mg		Ø 90mm	100g, E2	30430067	1170,00
PR224		•	220g	0,1mg		Ø 90mm	200g, E2	30430069	1250,00
PR124M	•	•	120g	0,1mg	1mg	Ø 90mm	100g, E2	30524384	1330,00
PR224M	•	•	220g	0,1mg	1mg	Ø 90mm	200g, E2	30524385	1430,00
PR124/E			120g	0,1mg		Ø 90mm	100g, E2	30430068	998,00
PR224/E			220g	0,1mg		Ø 90mm	200g, E2	30430070	1075,00

Accessories for Analytical Balances

Accessories	Explorer™ Analytical	Adventurer™ Analytical	Pioneer™ Analytical	PR Series Analytical	Item N°	Price €
Density determination kit	•	•	•		80253384	680,00
Sinker Glass for Density kit	•	•	•		83034024	355,00
ION-100A EU stand-alone ionizer	•	•	•	•	30095929	1460,00
In-Use Cover for EX	•				83033633	31,75
Full housing in use cover for 0.1mg and 1mg model for AX		•			30111792	33,25
In-Use Cover for PX			•		30372546	30,00
In-Use Cover for PR				•	30372547	30,00
Dust Cover for EX and AX (0.1mg)	•	•			30093334	41,00
Ethernet interface	•				83021082	138,00
Cable, USB A to Mini USB, Black	•				28120263	21,20
Interface Kit, RS232-USB	•	•	•	•	30304101	78,25
Cable, RS422, Extension, 3m	•				83021083	22,60
Cable, RS422, Extension, 9m	•				30078078	30,50
Auxiliary display, AD7-RS		•	•	•	30472064	96,00
Interface Cable, USB (Type A to B)		•			83021085	22,60
Anti theft device	•	•	•	•	80850043	211,00
SF40A impact printer	•	•	•	•	30064202	715,00
Carrying Case, EX	•				30251106	630,00
Carrying Case, AX		•			30251107	555,00
Carrying Case, PX			•		30477091	424,00



DENSITY DETERMINATION
80253384



SINKER GLASS FOR DENSITY KIT
83034024



ION-100A EU STANDALONE IONIZER
30095929



IN-USE COVER EX
83033633



DUST COVER
30093334



ETHERNET INTERFACE
83021082



ANTI-THEFT DEVICE
80850000



INTERFACE KIT, RS232-USB
30304101



CABLE, RS422, EXTENSION, 3M
83021083



CABLE, RS422, EXTENSION, 9M
30078078



AUXILIARY DISPLAY, AD7-RS
30472064



FULL HOUSING IN USE COVER AX
30111792



INTERFACE CABLE, USB (TYPE A TO B)
83021085



IN-USE COVER PX
30372546



IN-USE COVER PR
30372547



SF40A IMPACT PRINTER
ER 30064202

EXPLORER™ Precision

Intelligent Performance No Matter How Difficult the Task!



INTELLIGENT PERFORMANCE

Application

Weighing (17 units + custom units), parts counting, percentage weighing, animal/dynamic weighing, check weighing, filling, totalization, formulation, differential weighing, density determination, peak hold, ingredient costing, pipette adjustment (except for Explorer High Capacity), gross/net/tare weighing

Display

High resolution 5.7" (145mm diagonal) full colour VGA touch display including QWERTY, detachable and numeric key pads plus weighing capacity indication

Operation

AC adapter (included) or rechargeable battery (High Capacity models only)

Communication

Easy access communication ports including 2 units of USB port, RS232 and an optional 4th Ethernet port, GLP and GMP including date and time, Direct Data Transfer, print to USB drive

Construction

Modular design, metal base, ABS top housing, stainless steel pan, glass draftshield with top mounted side doors and flip/sliding top door (1mg models only), illuminated up-front level indicator, integral weigh below hook, security bracket, Adjustment lock, four touchless sensors, full housing in-use cover

Design Features

Selectable environmental filters, auto tare, user-selectable span Adjustment points, software lockout and reset menu, user-selectable communication settings and data print options, user-definable project and user IDs, auto standby, up to 14 operating language



Equipped with a weighing cell that's precision-machined from a solid metal block, Explorer balances are extremely accurate and durable, and are ideal for laboratory and industrial environments.

Features such as four-level user management with password protection, and un-editable system log plus a high level of configurability make Explorer balances well suited for regulated applications.

A detachable terminal with a large colour touchscreen, programmable IR sensors, and a frameless draftshield make Explorer precision balances easy and convenient to use.

Model	EC Type Approved	AutoCal™	Capacity	Readability (d)	(e)	Platform	Cal. weight, class	Item N°	Price €
EX223		•	220g	1mg		Ø 130mm	200g, F1	83021341	2590,00
EX423		•	420g	1mg		Ø 130mm	400g, F1	83021342	3345,00
EX623		•	620g	1mg		Ø 130mm	500g, F1	83021348	3895,00
EX1103		•	1100g	1mg		Ø 130mm	1000g, F1	83021350	4865,00
EX2202		•	2200g	0.01g		190mm x 200mm	2000g, F1	83021356	2565,00
EX4202		•	4200g	0.01g		190mm x 200mm	4000g, F1	83021357	3310,00
EX6202		•	6200g	0.01g		190mm x 200mm	5000g, F1	83021360	3885,00
EX10202		•	10200g	0.01g		190mm x 200mm	10000g, F1	83021362	4545,00
EX6201		•	6200g	0.1g		190mm x 200mm	5000g, F1	83021361	2440,00
EX10201		•	10200g	0.1g		190mm x 200mm	10000g, F1	83021365	3025,00
EX423M	•	•	420g	1mg	0.01g	Ø 130mm	400g, F1	83021343	3575,00
EX1103M	•	•	1100g	1mg	0.01g	Ø 130mm	1000g, F1	83021351	5095,00
EX4202M	•	•	4200g	0.01g	0.1g	190mm x 200mm	4000g, F1	83021358	3535,00
EX10202M	•	•	10200g	0.01g	0.1g	190mm x 200mm	10000g, F1	83021363	4770,00
EX10201M	•	•	10200g	0.1g	0.1g	190mm x 200mm	10000g, F1	83021366	3250,00
EX223/E			220g	1mg		Ø 130mm	200g, F1	30061979	2290,00
EX423/E			420g	1mg		Ø 130mm	400g, F1	30061980	3045,00
EX2202/E			2200g	0.01g		190mm x 200mm	2000g, F1	30061981	2265,00
EX4202/E			4200g	0.01g		190mm x 200mm	4000g, F1	30061982	3015,00
EX6202/E			6200g	0.01g		190mm x 200mm	5000g, F1	30061983	3585,00
EX6201/E			6200g	0.1g		190mm x 200mm	5000g, F1	30061984	2150,00
EX12001	•	•	12000g	0.1g		377mm x 311mm	10kg, F1	30078072	4345,00
EX24001	•	•	24000g	0.1g		377mm x 311mm	20kg, F1	30078073	4935,00
EX35001	•	•	35000g	0.1g		377mm x 311mm	10+20kg, F1	30078074	5980,00
EX12001M	•	•	12000g	0.1g	1g	377mm x 311mm	10kg, F1	30057067	4575,00
EX24001M	•	•	24000g	0.1g	1g	377mm x 311mm	20kg, F1	30057068	5165,00
EX35001M	•	•	35000g	0.1g	1g	377mm x 311mm	10+20kg, F1	30057069	6220,00

High Capacity

ADVENTURER™ Precision

Intuitive Balances Designed for Routine Weighing



**EASY OPERATION
& TIME SAVING**

Application

Weighing (17 units + custom unit), parts counting, percentage weighing, check weighing, animal/dynamic weighing, totalization/statistics, formulation, density determination, display hold, batch printing

Display

4.3" (109 mm) Full-colour VGA graphic touchscreen with user-controlled brightness, 6 mechanical buttons

Operation

AC adapter (included)

Communication

RS232, USB Host, USB Device (included), GLP/GMP compliant data output with real-time clock

Construction

Metal base, ABS top housing, stainless steel pan, glass draftshield with two piece top mounted side doors and sliding top door (1mg models only), illuminated up-front level indicator, integral weigh below hook, security bracket, adjustment lock, full housing replaceable in-use cover

Design Features

Selectable environmental filters, auto tare, user-selectable span adjustment points, software lockout and reset menu, user-selectable communication settings and data print options, user definable project and user IDs, up to 9 operating languages



Adventurer balances feature a color touchscreen, icon-based user interface, and an ergonomic design – making them easy to configure and use.

Features such as specialized weighing modes, multiple connectivity options, and AutoCal™ provide versatility and flexibility for a variety of applications.

Durable construction, large weighing surfaces, a space-saving draftshield design, and full housing in-use cover allow for use in lab, education and industrial environments.

Model	EC Type Approved	AutoCal™	Capacity	Readability (d)	(e)	Platform	Cal. weight, class	Item N°	Price €
AX223		•	220g	1mg		Ø 130mm	200g, F1	30127960	1670,00
AX423		•	420g	1mg		Ø 130mm	200g+200g, F1	30127962	1830,00
AX623		•	620g	1mg		Ø 130mm	500g, F1	30805890	1970,00
AX422		•	420g	0,01g		175mm x 195mm	200g+200g, F1	30206079	1350,00
AX822		•	820g	0,01g		175mm x 195mm	500g+100g, F1	30206081	1510,00
AX1502		•	1520g	0,01g		175mm x 195mm	1000g+500g, F1	30127966	1550,00
AX2202		•	2200g	0,01g		175mm x 195mm	2000g, F1	30127968	1670,00
AX4202		•	4200g	0,01g		175mm x 195mm	2000g+2000g, F1	30122616	1830,00
AX6202		•	6200g	0,01g		175mm x 195mm	5000g, F1	30805892	1970,00
AX2201		•	2200g	0,1g		175mm x 195mm	2000g, F1	30206083	1300,00
AX4201		•	4200g	0,1g		175mm x 195mm	2000g+2000g, F1	30127970	1530,00
AX8201		•	8200g	0,1g		175mm x 195mm	5000g+1000g, F1	30122619	1630,00
AX223M	•	•	220g	1mg	0,01g	Ø 130mm	200g, F1	30100623	1880,00
AX423M	•	•	420g	1mg	0,01g	Ø 130mm	200g+200g, F1	30100624	2050,00
AX623M	•	•	620g	1mg	0,01g	Ø 130mm	500g, F1	30805896	2190,00
AX1502M	•	•	1520g	0,01g	0,1g	175mm x 195mm	1000g+500g, F1	30100626	1770,00
AX2202M	•	•	2200g	0,01g	0,1g	175mm x 195mm	2000g, F1	30100627	1880,00
AX4202M	•	•	4200g	0,01g	0,1g	175mm x 195mm	2000g+2000g, F1	30100628	2050,00
AX6202M	•	•	6200g	0,01g	0,1g	175mm x 195mm	5000g, F1	30805897	2200,00
AX8201M	•	•	8200g	0,1g	1g	175mm x 195mm	5000g+1000g, F1	30122651	1860,00
AX223/E			220g	1mg		Ø 130mm	200g, F1	30127961	1360,00
AX423/E			420g	1mg		Ø 130mm	200g+200g, F1	30127963	1530,00
AX623/E			620g	1mg		Ø 130mm	500g, F1	30805891	1670,00
AX422/E			420g	0,01g		175mm x 195mm	200g+200g, F1	30206080	1040,00
AX822/E			820g	0,01g		175mm x 195mm	500g+100g, F1	30206082	1210,00
AX1502/E			1520g	0,01g		175mm x 195mm	1000g+500g, F1	30127967	1230,00
AX2202/E			2200g	0,01g		175mm x 195mm	2000g, F1	30127969	1360,00
AX4202/E			4200g	0,01g		175mm x 195mm	2000g+2000g, F1	30122617	1530,00
AX6202/E			6200g	0,01g		175mm x 195mm	5000g, F1	30805893	1670,00
AX2201/E			2200g	0,1g		175mm x 195mm	2000g, F1	30206084	977,00
AX4201/E			4200g	0,1g		175mm x 195mm	2000g+2000g, F1	30127971	1210,00
AX8201/E			8200g	0,1g		175mm x 195mm	5000g+1000g, F1	30122620	1340,00
AX12001		•	12000g	0,1g		175mm x 195mm	10000g, F1	30805894	1710,00
AX12001M	•	•	12000g	0,1g	1g	175mm x 195mm	10000g, F1	30805898	1950,00
AX12001/E			12000g	0,1g		175mm x 195mm	10000g, F1	30805895	1400,00

High Capacity

PIONEER™ Precision

Combining Economy and High Performance for Essential Weighing



\$ ECONOMY & ACCURACY

Application

Basic weighing, parts counting, percent weighing, animal weighing, density determination, formulation

Display

Backlit Liquid Crystal Display (LCD); Second line display for additional information or guidance (dot matrix)

Operation

AC adapter (included)

Communication

RS232 and USB device (included); GLP/GMP data output with real-time clock

Construction

Metal base, ABS top housing, removable stainless steel pan, glass draftshield (1 mg models only), integrated weigh-below-hook, security bracket, calibration lock and replaceable in-use cover

Design Features

User-selectable environmental filters and brightness settings, auto-tare, auto-dim, user-selectable span calibration points, software lockout and reset menu, user-selectable communication settings and data print options, user-definable project and user IDs, software overload/underload indicator, stability indicator, four operating languages



The PX offers high accuracy and repeatability for essential weighing applications in laboratory, industrial and education settings at an economical price point.

Featuring a cast metal lower housing, sub-pan and stainless steel weighing pan, the PX is durably constructed for versatile, long-term use.

Pioneer features a second line display for additional information or guidance, a static removal bar for convenient grounding, and USB connectivity.

Model	EC Type Approved	Internal Adjustment	Capacity	Readability (d)	(e)	Platform	Cal. weight, class	Item N°	Price €
PX323		•	320g	1mg		Ø 120mm	200g+100g, F1	30429805	1145,00
PX623		•	620g	1mg		Ø 120mm	500g+100g, F1	30619383	1280,00
PX3202		•	3200g	0,01g		Ø 180mm	2kg+1kg, F1	30429807	1145,00
PX6202		•	6200g	0,01g		Ø 180mm	5kg+1kg, F1	30619385	1260,00
PX4201		•	4200g	0,1g		Ø 180mm	2kg+2kg, F1	30429809	1110,00
PX323M	•	•	320g	1mg	0,01g	Ø 120mm	200g+100g, F1	30430015	1315,00
PX623M	•	•	620g	1mg	0,01g	Ø 120mm	500g+100g, F1	30619386	1460,00
PX3202M	•	•	3200g	0,01g	0,1g	Ø 180mm	2kg+1kg, F1	30430021	1350,00
PX6202M	•	•	6200g	0,01g	0,1g	Ø 180mm	5kg+1kg, F1	30714491	1465,00
PX4201M	•	•	4200g	0,1g	0,1g	Ø 180mm	2kg+2kg, F1	30430024	1280,00
PX323/E			320g	1mg		Ø 120mm	200g+100g, F1	30429812	920,00
PX623/E			620g	1mg		Ø 120mm	500g+100g, F1	30619382	1050,00
PX3202/E			3200g	0,01g		Ø 180mm	2kg+1kg, F1	30429814	930,00
PX6202/E			6200g	0,01g		Ø 180mm	5kg+1kg, F1	30619384	1175,00
PX4201/E			4200g	0,1g		Ø 180mm	2kg+2kg, F1	30429816	875,00
PX120001		•	12000g	0,1g		193mm x 248mm	10kg, F1	30868072	1220,00
PX120001/E			12000g	0,1g		193mm x 248mm	10kg, F1	30868073	960,00

High Capacity

PR SERIES Precision

Designed for Routine Weighing Applications in Your Workplace



**ROUTINE WEIGHING
FOR ANY BUDGET**



Application

Basic weighing, parts counting, percent weighing

Display

Backlit Liquid Crystal Display (LCD)

Operation

AC adapter (included)

Communication

RS232 (one port included)

Construction

ABS top housing, removable stainless steel pan, glass draftshield (1mg models only) with sliding top door, integrated weigh-below-hook, security bracket, calibration lock

Design features

User-selectable environmental filters and brightness settings, auto-tare, auto-dim, user-selectable span calibration points, software lockout and reset menu, user-selectable communication settings and data print options, user-definable project and user IDs, software overload/underload indicator, stability indicator

Basic Functionality for Routine Weighing Applications.

Designed for Uncomplicated Operation with Easy-to-Use Display and Interface.

Smart Design and Durable Construction.

Model	EC Type Approved	Internal Adjustment	Capacity	Readability (d)	(e)	Platform	Cal. weight, class	Item N°	Price €
PR223		•	220g	1mg		Ø 120mm	200g, F1	30492730	869,00
PR423		•	420g	1mg		Ø 120mm	200g+200g, F1	30492731	879,00
PR523		•	520g	1mg		Ø 120mm	500g, F1	30492732	931,00
PR1602		•	1600g	0,01g		Ø 180mm	1kg+500g, F1	30492733	885,00
PR2202		•	2200g	0,01g		Ø 180mm	2kg, F1	30492734	895,00
PR4202		•	4200g	0,01g		Ø 180mm	2kg+2kg, F1	30492735	946,00
PR4201		•	4200g	0,1g		Ø 180mm	2kg+2kg, F1	30492736	833,00
PR6201		•	6200g	0,1g		Ø 180mm	5kg+1kg, F1	30498176	854,00
PR223M	•	•	220g	1mg	0,01g	Ø 120mm	200g, F1	30524386	1020,00
PR423M	•	•	420g	1mg	0,01g	Ø 120mm	200g+200g, F1	30524387	1030,00
PR523M	•	•	520g	1mg	0,01g	Ø 120mm	500g, F1	30524608	1075,00
PR1602M	•	•	1600g	0,01g	0,1g	Ø 180mm	1kg+500g, F1	30524609	1035,00
PR2202M	•	•	2200g	0,01g	0,1g	Ø 180mm	2kg, F1	30524610	1055,00
PR4202M	•	•	4200g	0,01g	0,1g	Ø 180mm	2kg+2kg, F1	30524611	1090,00
PR4201M	•	•	4200g	0,1g	0,1g	Ø 180mm	2kg+2kg, F1	30524612	972,00
PR6201M	•	•	6200g	0,1g	1g	Ø 180mm	5kg+1kg, F1	30524613	998,00
PR223/E			220g	1mg		Ø 120mm	200g, F1	30430071	710,00
PR423/E			420g	1mg		Ø 120mm	200g+200g, F1	30430072	720,00
PR523/E			520g	1mg		Ø 120mm	500g, F1	30430073	766,00
PR1602/E			1600g	0,01g		Ø 180mm	1kg+500g, F1	30430074	725,00
PR2202/E			2200g	0,01g		Ø 180mm	2kg, F1	30430075	736,00
PR4202/E			4200g	0,01g		Ø 180mm	2kg+2kg, F1	30430076	782,00
PR4201/E			4200g	0,1g		Ø 180mm	2kg+2kg, F1	30430077	664,00
PR6201/E			6200g	0,1g		Ø 180mm	5kg+1kg, F1	30430078	684,00

Accessories for Precision Balances

Accessories	Explorer™ Precision	Explorer™ Precision High Capacity	Adventurer™ Precision	Pioneer™ Precision	PR Series Precision	Item N°	Price €
Density determination kit (for 1mg only)	•		•	•		80253384	680,00
Sinker Glass for Density kit	•		•	•		83034024	355,00
Auxiliary display, AD7-RS	•	•	•	•	•	30472064	96,00
In-Use Cover for EX	•	•				83033633	31,75
Full housing in-use cover for AX (1mg models)			•			30111792	33,25
Full housing in-use cover for AX (0.01g and 0.1g models)			•			30111777	33,25
In-Use Cover for PX (0.01g and 0.1g models)				•		30759721	30,00
In-Use Cover for PX				•		30372546	30,00
In-Use Cover for PR					•	30372547	30,00
Ethernet interface	•	•				83021082	138,00
Draft shield kit for EX 0,01g and 0,1g models	•					83021084	334,00
Tower Kit	•					83021102	86,75
Tower Kit		•				30078082	86,75
Anti theft device (models over 12kg)		•				80850000	68,50
Anti theft device	•		•	•	•	80850043	211,00
Cable, USB A to Mini USB, Black	•					28120263	21,20
Interface Kit, RS232-USB	•	•	•	•	•	30304101	78,25
Cable, RS422, Extension, 3m	•	•				83021083	22,60
Cable, RS422, Extension, 9m	•	•				30078078	30,50
SF40A impact printer	•	•	•	•	•	30064202	715,00
Battery Kit		•				30041295	222,00
Adjustable rolling feet		•				30041470	320,00
Interface Cable, USB (Type A to B)			•			83021085	22,60
Carrying Case, AX			•			30251108	485,00
Carrying Case, PX				•		30477092	384,00
Pan, 12kg, 300 x 225mm, PX (only for High Capacity models)				•		30868077	35,00



DENSITY DETERMINATION KIT
80253384



SINKER GLASS FOR DENSITY KIT
83034024



AUXILIARY DISPLAY, AD7-RS
30472064



IN-USE COVER
83033633



ETHERNET INTERFACE
83021082



DRAFT SHIELD KIT
83021084



TOWER KIT
83021102



ANTI THEFT DEVICE
80850000



INTERFACE CABLE, USB (TYPE A TO B)
83021085



INTERFACE KIT, RS232-USB
30304101



CABLE, RS422, EXTENSION, 3M
83021083



FULL HOUSING IN USE COVER
30111777



SF40A IMPACT PRINTER
ER30064202



IN-USE COVER
30372546



BATTERY KIT
30041295



ADJUSTABLE ROLLING FEET
30041470

SCOUT™ STX

High-Performance Portable Balances with Intuitive Touchscreens for High Accuracy and Productivity



**HIGH PERFORMANCE
& INTUITIVE OPERATION**

Application

Weighing, percent weighing, parts counting, check weighing, dynamic/animal weighing, totalization, density determination, display hold, mole weighing

Display

Full-colour touchscreen WQVGA Graphic LCD (4.3 in. diagonal), QWERTY keyboard and numeric keypad to quickly input GLP and GMP data and other USB device, Bluetooth or Ethernet

Operation

AC Adapter (Included), batteries (not included)

Communication

Optional RS232, USB host, USB device, Bluetooth or Ethernet

Construction

ABS housing, stainless steel pan, replaceable in-use cover

Design Features

Fast stabilization time, up-front level indicator, up to 9 operating languages, save-to-USB function, menu lock switch, stackable design, security slot, integral weigh below hook for below balance weighing applications, transportation lock, removable stainless steel weighing platform, stability indicator, overload and underload indicators, low battery indicator, auto shut-off



Large, intuitive, colour touchscreen display allows for easy operation - guiding the user through every step of the weighing application.

Advanced weighing technology allows for one-second stabilization time for increased productivity.

Integrated overload protection system rated at ten times the capacity protects the balance in demanding environments.

Model	Capacity	Readability (d)	Platform	Cal. weight, class*	Item N°	Price €
STX123	120g	0,001g	Ø 93mm	100g, M1*	30253058	639,00
STX223	220g	0,001g	Ø 93mm	200g, M1*	30253059	752,00
STX222	220g	0,01g	Ø 120mm	200g, M1*	30253060	486,00
STX422	420g	0,01g	Ø 120mm	200g, M1*	30253061	613,00
STX622	620g	0,01g	Ø 120mm	300g, M1*	30253062	655,00
STX1202	1200g	0,01g	170mm x 140mm	1000g, M1	30253063	716,00
STX2202	2200g	0,01g	170mm x 140mm	2000g, M1	30253064	752,00
STX421	420g	0,1g	Ø 120mm	200g, M1	30253065	295,00
STX621	620g	0,1g	170mm x 140mm	300g, M1	30253066	336,00
STX2201	2200g	0,1g	170mm x 140mm	2000g, M1	30253067	486,00
STX6201	6200g	0,1g	170mm x 140mm	5000g, M1	30253068	629,00
STX8200	8200g	1g	170mm x 140mm	8000g, M1	30253069	506,00

* test weight included for models up to 620 g capacity

SCOUT™ SKX

The Next Generation of Portable Balance



PRECISION & SPACE SAVING

Application

Weighing, density determination, mole weighing

Display

Backlit LCD: 6-digit 7-segment with white LED backlight

Operation

AC Adapter (included), batteries (not included)

Communication

Optional RS232, USB host, USB device, Bluetooth or Ethernet

Construction

ABS housing, stainless steel pan, replaceable in-use cover

Design Features

Fast stabilization time, menu lock switch, stackable design, security slot, integral weigh below hook for below balance weighing applications, transportation lock, removable stainless steel weighing platform, stability indicator, overload and underload indicators, low battery indicator, auto shut-off



Space-saving design allows for stackable storage (excluding the 1 mg model).

Integrated overload protection system rated at ten times the capacity protects the balance in demanding classroom environments.

Large backlit display improves readability in low-light conditions and makes it easier to view results from across a desk or lab bench. Multiple connectivity options for data collection and transfer.

Model	Capacity	Readability (d)	Platform	Cal. weight, class*	Item N°	Price €
SKX123	120g	0,001g	Ø 93mm	100g, M1*	30253070	406,00
SKX222	220g	0,01g	Ø 120mm	200g, M1*	30253071	304,00
SKX422	420g	0,01g	Ø 120mm	200g, M1*	30253072	354,00
SKX622	620g	0,01g	Ø 120mm	300g, M1*	30253073	388,00
SKX1202	1200g	0,01g	170mm x 140mm	1000g, M1	30253074	458,00
SKX2202	2200g	0,01g	170mm x 140mm	2000g, M1	30253075	486,00
SKX421	420g	0,1g	Ø 120mm	200g, M1	30253076	192,00
SKX621	620g	0,1g	170mm x 140mm	300g, M1	30253077	211,00
SKX2201	2200g	0,1g	170mm x 140mm	2000g, M1	30253078	295,00
SKX6201	6200g	0,1g	170mm x 140mm	5000g, M1	30253079	394,00
SKX8200	8200g	1g	170mm x 140mm	8000g, M1	30253080	252,00

* test weight included for models up to 620 g capacity

NAVIGATOR™

Rugged, Multi-Purpose Balances for a Wide Range of General Weighing Applications

Application

Basic weighing, parts counting, percent weighing, check weighing

Display

Backlit Liquid Crystal Display (LCD)

Operation

AC adapter (included) or 4 C batteries (not included)

Communication

RS232, USB, or Ethernet (available as accessories)

Construction

ABS housing, stainless steel pan, transportation lock, slip-resistant and adjustable feet, leveling bubble

Design Features

Menu and calibration lock switch, mechanical and software overload/underload protection, stability indicator, low battery indicator, auto shut-off, user-selectable print options, user-selectable communication settings



MULTI-PURPOSE
\$ AT COMPETITIVE PRICE



Designed for Accuracy and Repeatability in Essential Weighing Applications.

Improved Productivity with Fast Stabilization Time.

Mechanical Overload Protection Safeguards the Balance and Increases Durability
OIML/EC Type Approvals Ensure Accurate Operation in Trading Applications.

Navigator™ NV

Model	Capacity	Readability (d)	Platform	Cal. weight, class	Item N°	Price €
NV123*	120g	0,001g	Ø93mm	100g, F2	30615581	283,00
NV223*	220g	0,001g	Ø93mm	200g, F2	30615582	308,00
NV323*	320g	0,001g	Ø93mm	300g, F2	30615583	323,00
NV222	220g	0,01g	Ø 145mm	200g, F2	30456382	263,00
NV422	420g	0,01g	Ø 145mm	200g, F2	30456383	277,00
NV622	620g	0,01g	Ø 145mm	300g, F2	30456384	285,00
NV1202	1200g	0,01g	190mm x 144mm	1000g, F2	30615584	292,00
NV2202	2200g	0,01g	190mm x 144mm	2000g, F2	30615585	300,00
NV3202	3200g	0,01g	190mm x 144mm	3000g, F2	30615586	323,00
NV221	220g	0,1g	190mm x 144mm	200g, F2	30456385	146,00
NV621	620g	0,1g	190mm x 144mm	600g, F2	30456386	179,00
NV1201	1200g	0,1g	190mm x 144mm	500g, F2	30456387	217,00
NV2201	2200g	0,1g	190mm x 144mm	1000g, F2	30456388	251,00

*model with draftshield

Navigator™ NVT

Model	Capacity	Readability (d)	Platform	Cal. weight, class	Item N°	Price €
NVT2201	2200g	0,1g	230mm x 174mm	1000g, F2	30456389	258,00
NVT4201	4200g	0,1g	230mm x 174mm	2000g, F2	30456390	263,00
NVT6201	6200g	0,1g	230mm x 174mm	5000g, F2	30456391	277,00
NVT10201	10200g	0,1g	230mm x 174mm	5000g, F2	30615587	289,00
NVT2200	2200g	1g	230mm x 174mm	1000g, F2	30456392	145,00
NVT6200	6200g	1g	230mm x 174mm	5000g, F2	30456393	157,00
NVT12000	12000g	1g	230mm x 174mm	5000g, F2	30456394	236,00
NVT22000	22000g	1g	230mm x 174mm	10000g, F2	30456395	258,00

Model	EC Type Approved	Capacity	Readability (d)	Verification interval (e)	Platform	Cal. weight, class	Item N°	Price €
NVT1601M	•	1600g	0,5g	0,5g	230mm x 174mm	1000g, F2	30467835	270,00
NVT3200M	•	3200g	1g	1g	230mm x 174mm	2000g, F2	30467836	270,00
NVT6400M	•	6400g	2g	2g	230mm x 174mm	5000g, F2	30467837	270,00
NVT16000M	•	16000g	5g	5g	230mm x 174mm	10000g, F2	30467838	270,00

COMPASS™ CX

Designed with Energy Efficient Weighing in Mind for Workplace & In-the-Field Weighing

Application

Basic Weighing, Display Hold

Display

Liquid Crystal Display (LCD) with backlight

Operation

3 AA batteries (included) or AC adapter (not included)

Construction

ABS housing, stainless steel pan, transportation lock, slip-resistant feet

Design Features

Stability indicator, mechanical and software overload/underload protection, low battery indicator, auto shut-off



Low Power Consumption Design Allows for Uninterrupted Operation - 1,000 hours.

Improved Functionality with Large Backlit LCD and Two-Button Operation.

Practical Space-Saving Design for Convenient Storage.

Model	Capacity	Readability (d)	Platform	Cal. weight, class	Item N°	Price €
CX221	220g	0,1g	142mm x 128mm	200g, M1	30428208	64,00
CX621	620g	0,1g	142mm x 128mm	500g, M1	30428209	64,00
CX1201	1200g	0,1g	142mm x 128mm	1000g, M1	30428210	77,00
CX2200	2200g	1g	142mm x 128mm	1000g, M1	30428211	64,00
CX5200	5200g	1g	142mm x 128mm	2000g, M1	30428212	64,00

COMPASS™ CR

Basic Weighing Capabilities at an Economical Price



**COST & SPACE
CONSCIOUS DESIGN**

Application

Basic Weighing

Display

Liquid Crystal Display (LCD)

Operation

3 AA batteries (included) or AC adapter (not included)

Construction

ABS housing and pan, transportation lock, slip-resistant feet

Design Features

Stability indicator, mechanical and software overload/underload protection, auto shut-off



Basic Weighing at a Competitive Price. Long battery life - 300h.

Space-Saving Stackable Design for Convenient Storage.

Enlarged Weighing Platform Allows for Versatile Use.

Model	Capacity	Readability (d)	Platform	Cal. weight, class	Item N°	Price €
CR221	220g	0,1g	132mm x 125mm	200g, M1	30428213	35,00
CR621	620g	0,1g	132mm x 125mm	500g, M1	30428214	39,00
CR2200	2200g	1g	132mm x 125mm	1000g, M1	30428215	35,00
CR5200	5200g	1g	132mm x 125mm	2000g, M1	30428216	39,00

YA GOLD SERIES

Reliable pocket-sized portable scale



Unique hinged cover design.

High contrast backlit LCD.

Energy saving auto shut-off.

Model	Capacity	Readability (d)	Platform	Item N°	Price €
YA102	100g	0.01g	50mm x 70mm	80251910	41,00
YA302	300g	0.05g	50mm x 70mm	80251911	37,00
YA501	500g	0.1g	50mm x 70mm	80251912	37,00

Accessories for Precision Portable Balances

Accessories	Scout™ STX	Scout™ SKX	Navigator™ NV	Navigator™ NVT	Compass™ CX	Compass™ CR	YA Gold	Item N°	Price €
RS232 Interface Kit	•	•						30268982	81,50
USB Host Interface Kit	•	•						30268983	125,00
USB Device Interface Kit	•	•						30268984	91,50
Bluetooth Interface Kit	•	•						30268985	155,00
Ethernet Interface Kit	•	•						30268986	176,00
Stacking Kit, set (6)	•	•						30268987	65,75
Stacking Kit, set (1)	•	•						30268988	17,50
In-use cover	•	•						30269022	22,75
Auxiliary display, AD7-MD	•	•						30472063	96,00
Top loading kit for Density Determination	•	•						30269020	22,75
Anti theft device	•	•	•	•				80850043	211,00
Carrying Case, Scout	•	•						30269021	125,00
RS232 interface kit			•	•				83032107	45,00
USB interface kit			•	•				83032108	45,00
Ethernet interface kit			•	•				83032109	92,00
Battery rechargeable kit			•	•				83032106	55,25
Carrying Case			•	•				30467963	96,00
Stacking Kit, x1*			•					30467917	14,75
Stacking Kit, x4*			•					30467918	59,25
Weigh below kit			•	•				30628012	10,25
Carrying case CX CR					•	•		30467763	45,00
Pan SST CX					•			30428138	21,25
Adapter & plugs (EU-US-UK-AU-KR)						•		30467987	26,50
Adapter, 12V 0.5A, fixed UK plug						•		30467994	26,70
Adapter, 5V 1A EU, fixed plug					•			30333866	26,70
Adapter, 5V 1A UK, fixed plug					•			30333867	26,70
In-use cover NV			•					30809921	22,75
In-use cover NVT				•				30809922	22,75

*To use the storage stacking cover on models with draft shield (readability 0.001g), the draft shield has to be removed.



BLUETOOTH INTERFACE KIT
30268985



RS232 INTERFACE KIT
30268982



USB DEVICE INTERFACE KIT
30268984



STACKING KIT, SET (1)
30268988



AUXILIARY DISPLAY, AD7-MD
30472063



RS232 INTERFACE KIT
83032107



USB INTERFACE KIT
83032108



ETHERNET INTERFACE KIT
83032109



PAN SST CX
30428138



BATTERY RECHARGEABLE KIT
83032106



CARRYING CASE
30467963



TOP LOADING KIT FOR DENSITY
DETERMINATION 30269020



IN-USE COVER NV
30809921



IN-USE COVER NVT
30809922

MB120

Faster, Smarter, More Efficient Moisture Analysis



Application

Percent moisture determination, percent moisture regain determination, percent solid determination, weighing

Display

4.3" TFT colour touchscreen display

Operation

Power cord (included)

Communication

RS232 with GLP/GMP data output, USB Host and USB Device (included)

Construction

Metal base, ABS top housing, halogen heat source, stainless steel pan support, stainless steel pan handler, in-use cover

Design Features

SmartGuide™, Seven shut-off criteria (manual, timed, autofree %/s, autofree mg/s, A30, A60 and A90), four drying profiles (standard, fast, ramp and step), storage up to 100 methods and 1000 test results for all methods, user guide, 40°C to 230°C heating temperature range (1°C increments), 13 operational languages

Get started quickly with SmartGuide™. This revolutionary feature automatically analyzes a sample and creates a method, thereby removing the guesswork from the initial setup.

The MB120 has the capability to produce fast and accurate results due to features such as a halogen heating system, customizable shut-off criteria and four built in drying profiles.

The MB120 allows storage of up to 100 methods and 1,000 measurement results. Stored results can be exported, and methods can be shared across multiple devices.



Model	Capacity	Readability (d)	Platform	Cal. weight, class	Item N°	Price €
MB120	120g	1mg/0.01%	Ø 90mm	50g, F1	30241169	4060,00

MB90

Professional Moisture Analyzers for Advanced Moisture Analysis



Application

Percent moisture determination, percent moisture regain determination, percent solid determination, weighing

Display

4.3" TFT colour touchscreen display

Operation

Power cord (included)

Communication

RS232 with GLP/GMP data output, USB Host and USB Device (included)

Construction

Metal base, ABS top housing, halogen heat source, stainless steel pan support, stainless steel pan handler, in-use cover

Design Features

Three shut-off criteria (manual, timed or auto), two drying profiles (standard and fast), storage up to two methods, user guide, 40°C to 200°C heating temperature range (1°C increments), 13 operational languages

Increase Efficiency in your Lab with Fast Heating Halogen Technology.

Effortless Cleaning with the Tool-Free Cleaning Structure.

Operate Intuitively Using Icon-Based Navigation on the Colour Touchscreen Display.



Model	Capacity	Readability (d)	Platform	Cal. weight, class	Item N°	Price €
MB90	90g	1mg/0.01%	Ø 90mm	50g, F1	30241168	2730,00

MB27

Outstanding performance, flexible configuration and reliable results



**FAST,
RELIABLE RESULTS**

Application

Percent moisture determination, percent solid determination, weighing

Display

Backlit liquid crystal display (LCD)

Operation

Mains power

Communication

RS232 Interface

Construction

Halogen heat source, ABS housing, metal pan support, metal pan handler, in-use cover

Design Features

3 shut off criteria (manual, timed or auto), 1 standard heating profile, 50° to 160° C heating range (5° C increments)



With 0.01% 0.001g accuracy, the MB27 Moisture Analyzer is suitable for most applications.

Halogen heating and dual weight-temperature calibration provide fast, reliable results.

User-friendly touch-button design for easy operation.

Model	Capacity	Readability (d)	Platform	Cal. weight, class	Item N°	Price €
MB27	90g	1mg/0.01%	Ø 90mm	50g, F1	30251541	1870,00

MB25

Fast, easy-to-use, compact halogen analyser



**HALOGEN HEATING
WITH RELIABLE RESULTS**

Application

Moisture determination from 50 °C to 160 °C, 5 °C increments, repeatability 0.05 % (10 g sample)

Display

Custom backlit LCD

Operation

Mains power

Communication

RS232 Interface

Construction

Halogen heat source, durable ABS

Design Features

Compact design, in-use cover, sample pan holder, halogen heating, manual, automatic or timed duration drying time



Halogen heating and dual weight-temperature calibration provide fast, reliable results.

User-friendly touch-button design for easy operation.

Easy-clean heating chamber is ideal for frequent use.

Model	Capacity	Readability (d)	Platform	Cal. weight, class	Item N°	Price €
MB25	110g	5mg/0.01% (sample <10g) 5mg/0.05% (sample >10g)	Ø 90mm	50g, M1	80252472	1760,00

MB23

Compact analyser that is easy to use and clean



**MOISTURE ANALYSIS
FOR ANY BUDGET**

Application

Moisture determination from 50 °C to 160 °C, 5 °C increments, repeatability 0.2 % (10 g sample)

Display

Custom backlit LCD

Operation

Mains operation

Communication

RS232 Interface

Construction

Durable ABS

Design Features

Compact design, in-use cover, sample pan holder, infrared heating, manual, automatic or timed duration drying time



Glass-less infrared heating increases efficiency for faster results.

Compact design suits any workspace.

Model	Capacity	Readability (d)	Platform	Cal. weight, class	Item N°	Price €
MB23	110g	10mg/0.1%	Ø 90mm	50g, M1	80252470	1250,00

Accessories for Moisture Analysers

Accessories	MB120	MB90	MB27	MB25	MB23	Item N°	Price €
Temperature Calibration Kit	•	•	•	•	•	11113857	427,00
Pan Holder	•	•				30246947	58,70
Aluminium Sample pans, 50 per pack	•	•	•	•	•	30585411	21,00
Glass fibre discs, 200 per pack	•	•	•	•	•	80850087	72,75
Sample Cage	•	•	•	•	•	80252477	50,25
Reusable sample pans, 3 per pack, 14mm	•	•	•	•	•	80252479	116,00
Reusable sample pans, 3 per pack, 7mm	•	•	•	•	•	80252478	114,00
Scoop, Accessory	•	•				30284477	38,50
In-Use-Cover, Accessory - MB90, MB120	•	•				30284478	50,25
Anti theft device	•	•	•	•	•	80850043	211,00
Interface Kit, RS232-USB	•	•	•	•	•	30304101	78,25
Pan Holder			•	•	•	80252476	58,00
In use display cover MB2x			•	•	•	80252475	54,00
Carrying Case - MB90, MB120	•	•				30298583	507,00



TEMPERATURE CALIBRATION KIT
11113857



REUSABLE SAMPLE PANS, 3 PER PACK,
14MM 80252479



ALUMINIUM SAMPLE PANS, 50 PER PACK
30585411



GLASS FIBRE DISCS, 200 PER PACK
80850087



SAMPLE CAGE
80252477



REUSABLE SAMPLE PANS, 3 PER PACK,
7MM 80252478



SCOOP, ACCESSORY
30284477



ANTI THEFT DEVICE
80850043



IN-USE-COVER, ACCESSORY - MB90,
MB120 30284478



INTERFACE KIT, RS232-USB
30304101



IN USE DISPLAY COVER, MB2X
80252475

CENT-O-GRAM™ 311-00 DIAL-O-GRAM™ 310-00

High precision and remarkable value

Model	Capacity	Readability (d)	Platform	Item N°	Price €
311-00	311g	0.01g	Ø 89mm	80000010	285,00
310-00	310g	0.01g	Ø 89mm	80000009	341,00

Metal base and beam, stainless steel pan.

Self-aligning beam design.

Vernier dial convenience on Dial-o-Gram (model 310-00).



TRIPLE BEAM PRO™

The Triple Beam balance that offers even more

Model	Capacity	Readability (d)	Platform	Item N°	Price €
TP2611	2610g	0.1g	147mm	80251370	285,00

Easy grip carrying handle.

Included weight set with integrated storage in the balance.



TRIPLE BEAM JUNIOR

OHAUS Quality at an Affordable Price

Model	Capacity	Readability (d)	Platform	Item N°	Price €
TJ2611	2610g*	0.1g	Ø 150mm	80252049	224,00
TJ611	610g	0.1g	Ø 150mm	80252052	185,00

* with included 3-piece mass set

Notched and tiered beams for easy reading.

Security Features (lockable connections).

Zero adjust compensator ensures correct zero.



TRIPLE BEAM 700

Rugged enough to deliver accurate results year after year

Model	Capacity	Readability (d)	Platform	Item N°	Price €
710-00	610g*	0.1g	Ø 152mm	80000011	271,00
750-S0	610g*	0.1g	Ø 152mm	80000012	244,00
760-00	610g*	0.1g	Ø 152mm	80000013	295,00

*2610g capacity with optional 80780108 weight set

Fixed or removable pan (710-00) to fit any application.

The outstanding quality of the beam allows for smooth positioning of the poise, resulting in superior measurement.

Easy to read deep notched and tiered beams, counterbalancing knob for quick zeroing, self-aligning beam design.



HARVARD JUNIOR

Traditional style two-pan design at a great price

Model	Capacity	Readability (d)	Platform	Item N°	Price €
HJ2001	2000g	0.5g	Ø 150mm	80252050	141,00



Easy grip carrying handle.

Lock to protect the balance during transport or storage (can be stacked).

Included 8-piece weight set and bowl that fits over the flat platform.

HARVARD TRIP

Two pan balance designed for comparative weighing

Model	Capacity	Readability (d)	Platform	Item N°	Price €
1450-SD	2000g*	0.1g	Ø 152mm	80000003	394,00
1550-SD	2000g*	0.1g	Ø 152mm	80000004	453,00
1560-SD	2000g*	0.1g	Ø 152mm	80000005	505,00

*2000g capacity with optional Weight Set, 1kg-1g (80850115)



The built-in sliding masses make weighing very simple and quick. 1560-SD with tare range of 225g.

Magnetic damping minimizes oscillation and helps speed up weighing.

A counterbalancing knob provides quick zeroing of the balance while floating agate bearings helps ensure accurate results.

Accessories for Mechanical Scales & Balances

Accessories	Cent-O-Gram™ 311-00	Triple Beam Pro™	Triple Beam Junior	Triple Beam 700	Harvard Junior	Harvard Trip	Item N°	Price €
3 pieces counter weight set increases capacity to 2610g (2x1kg and 1x500g)		•	•	•			80780108	60,75
Weight Set, 1kg-1g						•	80850115	145,00



3 PIECES COUNTER WEIGHT SET INCREASES CAPACITY TO 2610G 80780108

DF 5000 SERIES Stainless Steel Floor Scales

Prime German-Quality Floor Scales Fulfill Even the Most Demanding Industrial Applications

Application

Weighing, Parts Counting, Percent Weighing, Checkweighing, Dynamic Weighing (Display Hold), Filling

Display

Backlit Liquid Crystal Display (LCD), 256 x 80 pixels

Operation

Internal power supply with power cord (included) or rechargeable battery (sold separately as an accessory)

Communication


Easy access communication ports including RS232 and Micro SD slot (included), second RS232/RS485/USB device, Ethernet and WiFi/Bluetooth, Discrete I/O and Analog Output (sold separately as accessories); GLP/GMP data output with real-time clock

Construction

Stainless steel indicator, bead-blasted stainless steel platform, IP67 stainless steel junction box with top access, IP68 OIML R60 electro-polished stainless steel load cells, available pit frames mounting

Design Features

Structural channel bracing and pendulum design feet with removable rubber cap, 150% overload protection. Checkweighing symbols with selectable operation and audible signal settings, menu lock switch, multiple selectable operating languages, selectable environmental and auto-print settings, stability indicator, overload/underload indicators, selectable adjust contrast settings, auto-dim display, auto-off, auto tare.



EasyConnect



OHAUS Defender 5000 series Floor Scale is equipped with a new TD52XW series indicator with an easy-to-use keypad, large backlit display and plain text on-screen prompts which guide you through operation.

Four IP68 electro-polished stainless steel load cells that are laser-weld, hermetically-sealed and OIML-certified are standard on all models.

Combination floor-pit mount design provides versatility for portable or stationary placement, above or below ground installation.

Model	EC Type Approved	Capacity	Readability (d)	(e)	Platform (mm)	Item N°	Price €
DF52XW300G1S-EU		300kg	50g		800 x 800 x 90	30525406	4355,00
DF52XW600G1S-EU		600kg	100g		800 x 800 x 90	30525407	4355,00
DF52XW300G1R-EU		300kg	50g		1000 x 1000 x 90	30525408	4785,00
DF52XW600G1R-EU		600kg	100g		1000 x 1000 x 90	30525409	4745,00
DF52XW1500G1R-EU		1500kg	200g		1000 x 1000 x 90	30525410	4785,00
DF52XW600G1LX-EU		600kg	100g		1500 x 1250 x 90	30525411	5955,00
DF52XW1500G1LX-EU		1500kg	200g		1500 x 1250 x 90	30525412	5955,00
DF52XW3000G1LX-EU		3000kg	500g		1500 x 1250 x 90	30525413	5955,00
DF52XW600G1X-EU		600kg	100g		1500 x 1500 x 90	30525414	6330,00
DF52XW1500G1X-EU		1500kg	200g		1500 x 1500 x 90	30525415	6345,00
DF52XW3000G1X-EU		3000kg	500g		1500 x 1500 x 90	30525416	6345,00
DF52XW300G1S-M	•	300kg		100g	800 x 800 x 90	30525438	4465,00
DF52XW600G1S-M	•	600kg		200g	800 x 800 x 90	30525439	4465,00
DF52XW300G1R-M	•	300kg		100g	1000 x 1000 x 90	30525440	4895,00
DF52XW600G1R-M	•	600kg		200g	1000 x 1000 x 90	30525441	4855,00
DF52XW1500G1R-M	•	1500kg		500g	1000 x 1000 x 90	30525442	4895,00
DF52XW600G1LX-M	•	600kg		200g	1500 x 1250 x 90	30525443	6055,00
DF52XW1500G1LX-M	•	1500kg		500g	1500 x 1250 x 90	30525444	6055,00
DF52XW3000G1LX-M	•	3000kg		1000g	1500 x 1250 x 90	30525445	6055,00
DF52XW600G1X-M	•	600kg		200g	1500 x 1500 x 90	30525446	6420,00
DF52XW1500G1X-M	•	1500kg		500g	1500 x 1500 x 90	30525447	6440,00
DF52XW3000G1X-M	•	3000kg		1000g	1500 x 1500 x 90	30525448	6440,00

DF 5000 SERIES Floor Scales

Prime German-Quality Floor Scales Fulfill Even the Most Demanding Industrial Applications

Application

Weighing, Parts Counting, Percent Weighing, Checkweighing, Dynamic Weighing (Display Hold), Filling

Display

Backlit Liquid Crystal Display (LCD), 256 x 80 pixels

Operation

Internal power supply with power cord (included) or rechargeable battery (sold separately as an accessory)

Communication

Easy access communication ports including RS232 and Micro SD slot (included), second RS232/RS485/USB device, Ethernet and WiFi/Bluetooth, Discrete I/O and Analog Output (sold separately as accessories); GLP/GMP data output with real-time clock

Construction

ABS plastic indicator, black painted steel platform, IP67 stainless steel junction box with top access, IP68 OIML R60 electro-polished stainless steel load cells, available pit frames mounting

Design Features

Structural channel bracing and pendulum design feet with removable rubber cap, 150% overload protection. Checkweighing symbols with selectable operation and audible signal settings, menu lock switch, multiple selectable operating languages, selectable environmental and auto-print settings, stability indicator, overload/underload indicators, selectable adjust contrast settings, auto-dim display, auto-off, auto tare.



OHAUS Defender 5000 series Floor Scale is equipped with a new TD52P series indicator with an easy-to-use keypad, large backlit display and plain text on-screen prompts which guide you through operation.

Four IP68 electro-polished stainless steel load cells that are laser-weld, hermetically-sealed and OIML-certified are standard on all models.

Combination floor-pit mount design provides versatility for portable or stationary placement, above or below ground installation.

Model	EC Type Approved	Capacity	Readability (d)	(e)	Platform (mm)	Item N°	Price €
DF52P300F1S		300kg	50g		800 x 800 x 90	30525380	2920,00
DF52P600F1S		600kg	100g		800 x 800 x 90	30525381	2920,00
DF52P300F1R		300kg	50g		1000 x 1000 x 90	30525382	3075,00
DF52P600F1R		600kg	100g		1000 x 1000 x 90	30525383	3065,00
DF52P1500F1R		1500kg	200g		1000 x 1000 x 90	30525384	3100,00
DF52P600F1LX		600kg	100g		1500 x 1250 x 90	30525385	3365,00
DF52P1500F1LX		1500kg	200g		1500 x 1250 x 90	30525386	3385,00
DF52P3000F1LX		3000kg	500g		1500 x 1250 x 90	30525387	3365,00
DF52P600F1X		600kg	100g		1500 x 1500 x 90	30525388	3770,00
DF52P1500F1X		1500kg	200g		1500 x 1500 x 90	30525389	3770,00
DF52P3000F1X		3000kg	500g		1500 x 1500 x 90	30525390	3770,00
DF52P1500F1XV		1500kg	200g		2000 x 1500 x 90	30525391	4535,00
DF52P3000F1XV		3000kg	500g		2000 x 1500 x 90	30525392	4540,00
DF52P300F1S-M	•	300kg		100g	800 x 800 x 90	30525425	3020,00
DF52P600F1S-M	•	600kg		200g	800 x 800 x 90	30525426	3020,00
DF52P300F1R-M	•	300kg		100g	1000 x 1000 x 90	30525427	3180,00
DF52P600F1R-M	•	600kg		200g	1000 x 1000 x 90	30525428	3160,00
DF52P1500F1R-M	•	1500kg		500g	1000 x 1000 x 90	30525429	3200,00
DF52P600F1LX-M	•	600kg		200g	1500 x 1250 x 90	30525430	3485,00
DF52P1500F1LX-M	•	1500kg		500g	1500 x 1250 x 90	30525431	3510,00
DF52P3000F1LX-M	•	3000kg		1000g	1500 x 1250 x 90	30525432	3485,00
DF52P600F1X-M	•	600kg		200g	1500 x 1500 x 90	30525433	3885,00
DF52P1500F1X-M	•	1500kg		500g	1500 x 1500 x 90	30525434	3885,00
DF52P3000F1X-M	•	3000kg		1000g	1500 x 1500 x 90	30525435	3885,00
DF52P1500F1XV-M	•	1500kg		500g	2000 x 1500 x 90	30525436	4640,00
DF52P3000F1XV-M	•	3000kg		1000g	2000 x 1500 x 90	30525437	4645,00

DF 3000 SERIES Stainless Steel Floor Scales (i-DF33)

Floor Scales for Basic Industrial Applications

Application

Weighing, Checkweighing, Counting, Totalization

Display

7 Segment LCD with multicolor backlight

Operation

Internal power supply with power cord (included) or Rechargeable Battery (included)

Communication

Standard RS232 and optional USB Device or Ethernet

Construction

304-grade stainless steel IP66 indicator, 5mm fixed top plate, stainless steel platform, IP65 stainless steel junction box with side access, 4 x IP67 OIML R60 stainless steel load cells, available stainless steel ramp

Design Features

125% overload protection, four load cell mounted feet with self-swiveling ball-in-cup design, painted and stainless steel floor stands (available as accessories)



304 stainless steel, sand-blasted flat surface with four IP67 stainless steel load cells, fixed top plate (5mm thick) with structural bracing that ensures a stiff & solid construction.

Defender 3000 come in six models ranging from platform sizes of 1 x 1 m, 1.25 x 1.25 m and 1.5 x 1.5 m, and capacities from 600 kg to 3000 kg, allowing greater flexibility for your weighing needs.

Results are clearly shown on the large, bright backlit LCD display. In Checkweighing mode, display changes color (yellow/green/red) to alert users when within the target weight range.

Model	EC Type Approved	Capacity	Readability (d)	Readability (e)	Platform (mm)	Item N°	Price €
i-DF33XW600C1R-EU		600kg	200g		1000 x 1000	30700625	2545,00
i-DF33XW1500C1R-EU		1500kg	500g		1000 x 1000	30700626	2545,00
i-DF33XW1500C1L-EU		1500kg	500g		1250 x 1250	30700629	2960,00
i-DF33XW3000C1L-EU		3000kg	1000g		1250 x 1250	30700630	2960,00
i-DF33XW1500C1X-EU		1500kg	500g		1500 x 1500	30700633	4510,00
i-DF33XW3000C1X-EU		3000kg	1000g		1500 x 1500	30700634	4510,00
i-DF33XW600C1R-EU-M	•	600kg		200g	1000 x 1000	30764203	2630,00
i-DF33XW1500C1R-EU-M	•	1500kg		500g	1000 x 1000	30764204	2630,00
i-DF33XW1500C1L-EU-M	•	1500kg		500g	1250 x 1250	30764259	3040,00
i-DF33XW3000C1L-EU-M	•	3000kg		1000g	1250 x 1250	30764260	3040,00
i-DF33XW1500C1X-EU-M	•	1500kg		500g	1500 x 1500	30764265	4595,00
i-DF33XW3000C1X-EU-M	•	3000kg		1000g	1500 x 1500	30764266	4595,00

DF 3000 SERIES Floor Scales (i-DF33)

Floor Scales for Basic Industrial Applications

Application

Weighing, Checkweighing, Counting, Totalization

Display

7 Segment LCD with multicolor backlight

Operation

Internal power supply with power cord (included), 6C (LR14) Batteries (not included) or Rechargeable Battery (not included)

Communication


Standard RS232 and optional USB Device or Ethernet

Construction

ABS plastic indicator, 5mm fixed top plate, black treaded painted steel platform, IP65 stainless steel junction box with side access, 4 x IP67 OIML R60 alloy steel load cells, two carbon steel eye bolts, available painted steel ramp

Design Features

125% overload protection, four load cell mounted feet with self-swiveling ball-in-cup design, painted and stainless steel floor stands (available as accessories)



Painted carbon steel, raised treaded surface with four IP67 alloy steel load cells, fixed top plate (5mm thick) with structural bracing that ensures a stiff and solid construction.


Defender 3000 comes in six models ranging from platform sizes of 1 x 1 m, 1.25 x 1.25 m and 1.5 x 1.5 m, and capacities from 600 kg to 3000 kg, allowing greater flexibility for your weighing needs.

Results are clearly shown on the large, bright backlit LCD display. In Checkweighing mode, display changes color (yellow/green/red) to alert users when within the target weight range.

Model	EC Type Approved	Capacity	Readability (d)	Readability (e)	Platform (mm)	Item N°	Price €
i-DF33P600B1R		600kg	200g		1000 x 1000	30700623	1515,00
i-DF33P1500B1R		1500kg	500g		1000 x 1000	30700624	1515,00
i-DF33P1500B1L		1500kg	500g		1250 x 1250	30700627	1560,00
i-DF33P3000B1L		3000kg	1000g		1250 x 1250	30700628	1560,00
i-DF33P1500B1X		1500kg	500g		1500 x 1500	30700631	1695,00
i-DF33P3000B1X		3000kg	1000g		1500 x 1500	30700632	1695,00
i-DF33P600B1R-M	•	600kg		200g	1000 x 1000	30764201	1600,00
i-DF33P1500B1R-M	•	1500kg		500g	1000 x 1000	30764202	1600,00
i-DF33P1500B1L-M	•	1500kg		500g	1250 x 1250	30764207	1645,00
i-DF33P3000B1L-M	•	3000kg		1000g	1250 x 1250	30764258	1645,00
i-DF33P1500B1X-M	•	1500kg		500g	1500 x 1500	30764263	1775,00
i-DF33P3000B1X-M	•	3000kg		1000g	1500 x 1500	30764264	1775,00

DF 3000 SERIES Floor Scales (DF32)

Versatile, low cost floor and pit platform with T32M indicator



Application

Weighing (5 units), Parts Counting, Dynamic, Checkweighing, Weight Accumulation.

Display

Backlit LCD

Operation

AC adapter (included)

Communication

RS232 interface

Construction

Metal indicator, black painted steel platform, raised embossed treads, fixed top plate 5mm thick, 4 x H8C IP67 OIML R60 alloy steel load cells, IP65 stainless steel junction box

Design Features

Indicator with ABS front panel, painted carbon steel rear housing and mounting bracket



Indicator with intuitive usage and sturdy metal housing.

Accumulation and printing capability via RS232.

IP67 alloy steel load cells.

Model	Capacity	Readability (d)	Platform (mm)	Item N°	Price €
DF32M1500BS	1500kg	500g	800 x 800 x 110	83997598	1150,00
DF32M1500BR	1500kg	500g	1000 x 1000 x 110	83997597	1265,00
DF32M1500BL	1500kg	500g	1200 x 1200 x 110	83997596	1440,00
DF32M3000BL	3000kg	1000g	1200 x 1200 x 110	83997595	1440,00
DF32M1500BLX	1500kg	500g	1200 x 1500 x 110	83997585	1605,00
DF32M3000BLX	3000kg	1000g	1200 x 1500 x 110	83997584	1605,00
DF32M1500BX	1500kg	500g	1500 x 1500 x 110	83997594	1740,00
DF32M3000BX	3000kg	1000g	1500 x 1500 x 110	83997593	1740,00

Accessories for Floor Scales

Accessories	DF 5000 Series Stainless Steel Floor Scales	DF 5000 Series Floor Scales	DF 3000 Series Floor Scales (i-DF33)	DF 3000 Series Floor Scales (DF32)	Item N°	Price €
Ramp DF Stainless Steel G1 800mm	•				22015473	908,00
Ramp DF Stainless Steel G1 1000mm	•				22015474	1012,00
Ramp DF Stainless Steel G1 1250mm	•				22015475	1264,00
Ramp DF Stainless Steel G1 1500mm	•				22015476	1531,00
Pit Frame DF Stainless Steel G1S	•				22015493	1384,00
Pit Frame DF Stainless Steel G1R	•				22015494	1487,00
Pit Frame DF Stainless Steel G1LX	•				22015495	1711,00
Pit Frame DF Stainless Steel G1X	•				22015496	1837,00
Footplates DF Stainless Steel G1 Set (2)	•				22015581	104,00
Pit Frame DF Painted F1S		•			22015392	995,00
Pit Frame DF Painted F1R		•			22015393	1061,00
Pit Frame DF Painted F1LX		•			22015394	1181,00
Pit Frame DF Painted F1X		•			22015396	1257,00
Pit Frame DF Painted F1XV		•			22015397	1373,00
Ramp DF Painted F1 800mm		•			22015464	494,00
Ramp DF Painted F1 1000mm		•			22015465	539,00
Ramp DF Painted F1 1250mm		•			22015466	624,00
Ramp DF Painted F1 1500mm		•			22015467	722,00
Footplates DF Painted F1 Set (2)		•			22015499	84,00
Floor scale Column Kit Painted 980mm		•	•	•	30531244	127,00
Floor scale Column Kit SST 980mm	•		•	•	30539553	235,00
Ramp 800mm DF				•	80252722	235,00
Ramp 1000mm DF				•	80252723	274,00
Ramp 1200mm DF				•	80252724	320,00
Ramp 1500mm DF				•	80252725	383,00
Pit frame, 800mm for DF				•	83033708	207,00
Pit frame, 1000mm for DF				•	83033709	260,00
Pit frame, 1200mm for DF				•	83033710	297,00
Pit frame, 1200mm x 1500mm for DF				•	83033714	330,00
Pit frame, 1500mm for DF				•	83033711	352,00
Ramp, 1000mm, Painted, DF-B1			•		72198338	345,00
Ramp, 1250mm, Painted, DF-B1			•		72198339	396,00
Ramp, 1500mm, Painted, DF-B1			•		72198340	428,00
Ramp, 1000mm, SST, DF-C1			•		72198341	919,00
Ramp, 1250mm, SST, DF-C1			•		72198342	1039,00
Ramp, 1500mm, SST, DF-C1			•		72198343	1427,00
Foot plate kit, Painted, DF-B1			•		72198344	46,00
Foot plate kit, SST, DF-C1			•		72198345	99,00
Ethernet Kit i-DT33			•		30429666	111,00
USB Device Kit i-DT33			•		30699120	72,00
Dust Cover (5) i-DT33			•		30699121	37,00
Rechargeable Battery Kit i-DT33P			•		30699122	44,00
SF40A Printer	•	•	•	•	30064202	715,00
Cable RS232 i-DT33P			•		80500525	99,70
Cable RS232 i-DT33XW			•		80500552	91,80
LC Extension Cable, 3.5 m *1	•	•	•	•	30803566	78,00

1 Cannot be used for legal-for-trade applications



RRAMP DF STAINLESS STEEL G1 800MM
800MM 22015473



FOOTPLATES DF PAINTED F1 SET (2)
22015499



FLOOR SCALE COLUMN KIT PAINTED
980MM 30531244



FLOOR SCALE COLUMN KIT SST 980MM
30539553



FOOTPLATES DF STAINLESS STEEL G1 SET
(2) 22015581



RAMP 800MM DF
80252722



PIT FRAME DF STAINLESS STEEL G1R
22015494



PIT FRAME DF PAINTED F1S
22015392

DFP 5000 SERIES Pallet Scales

Prime German-Quality Pallet Scales Fulfill Even the Most Demanding Industry Applications

Application

Weighing, Parts Counting, Percent Weighing, Checkweighing, Dynamic Weighing (Display Hold), Filling

Display

Backlit Liquid Crystal Display (LCD), 256 x 80 pixels

Operation

Internal power supply with power cord (included) or rechargeable battery (sold separately as an accessory)

Communication

Easy access communication ports including RS232 and Micro SD slot (included), second RS232/RS485/USB device, Ethernet and WiFi/Bluetooth, Discrete I/O and Analog Output (sold separately as accessories); GLP/GMP data output with real-time clock

Construction DFP52XW

Stainless steel indicator IP68, bead-blasted stainless steel U-frame, IP67 stainless steel junction box, IP68 OIML R60 electro-polished stainless steel load cells

Construction DFP52P

ABS plastic indicator, black painted steel U-frame, IP67 stainless steel junction box, IP68 OIML R60 electro-polished stainless steel load cells

Design Features

Structural channel bracing and pendulum design feet with removable rubber cap, 150% overload protection. Checkweighing symbols with selectable operation and audible signal settings, menu lock switch, multiple selectable operating languages, selectable environmental and auto-print settings, stability indicator, overload/underload indicators, selectable adjust contrast settings, auto-dim display, auto-off, auto tare.



OHAUS Defender 5000 series Pallet Scale is equipped with new TD52 series indicator with an easy-to-use keypad, large backlit display and plain text on-screen prompts which guide you through operation.

The U-frame provides easy on-and-off loading of open-bottom pallets and containers, reducing risk of injury and over-exertion.

Four IP68 electro-polished stainless steel load cells that are laser-weld hermetically-sealed and OIML-certified are standard on all models.

Model	EC Type Approved	IP Dust and Water Protected	Capacity	Readability (d)	(e)	Platform	Item N°	Price €
DFP52P600J1			600kg	100g		1256mm x 840mm	30525419	2595,00
DFP52P1500J1			1500kg	200g		1256mm x 840mm	30525420	2595,00
DFP52XW600K1-EU		•	600kg	100g		1256mm x 840mm	30525423	3940,00
DFP52XW1500K1-EU		•	1500kg	200g		1256mm x 840mm	30525424	3940,00
DFP52XW600K1-M	•	•	600kg		200g	1256mm x 840mm	30525453	4040,00
DFP52XW1500K1-M	•	•	1500kg		500g	1256mm x 840mm	30525454	4040,00
DFP52P600J1-M	•		600kg		200g	1256mm x 840mm	30525451	2690,00
DFP52P1500J1-M	•		1500kg		500g	1256mm x 840mm	30525452	2690,00
Accessories							Item N°	Price €
LC Extension Cable, 3.5 m *1							30803566	78,00

*1 Cannot be used for legal-for-trade applications

DFP 3000 SERIES Pallet Scales

Pallet Scales for Basic Industrial Applications



Application

Weighing, Checkweighing, Counting, Totalization

Display

7 Segment LCD with multicolor backlight

Operation

Internal power supply with power cord (included) or Rechargeable Battery (as standard with XW model, available as accessory for P model)

Communication

Standard RS232 and optional USB Device or Ethernet

Construction

DFP33XW

304-grade stainless steel IP66 indicator, stainless steel U-frame, IP65 stainless steel junction box, 4 IP67 OIML R60 stainless steel load cells, handle and front wheels, foot plate kit (accessories sold separately)

DFP33P

ABS plastic indicator, painted carbon steel U-frame, IP65 stainless steel junction box, 4 H8C IP67 OIML R60 alloy steel load cells, handle and front wheels, foot plate kit (accessories sold separately)

Design Features

Four load cell mounted feet with self-swiveling ball-in-cup design, painted and stainless steel floor stands (available as accessories).



Available in popular 1 500 kg capacity, painted carbon steel or 304 stainless steel models with handle and rollers for convenient lifting and relocation.

Four OIML R60 approved IP67 alloy steel or stainless steel load cells with self-swiveling ball-in-cup design feet to best suit your weighing needs.

Results are clearly shown on the large, bright backlit LCD display of the i-DT33 indicator. In Checkweighing mode, the display changes color (yellow/green/red) to alert users when within the target weight range.

Stainless Steel Pallet Scale with i-DT33XW Stainless Steel Indicator

Model	EC Type Approved	Capacity	Readability (d)	Readability (e)	Platform (mm)	Item N°	Price €
i-DFP33XW1500C1-EU		1500kg	500g		1260 x 840 x 90	30712896	2545,00
i-DFP33XW1500C1-EU-M	•	1500kg		500g	1260 x 840 x 90	30764271	2630,00

Painted Steel Pallet Scale with i-DT33P ABS Indicator

Model	EC Type Approved	Capacity	Readability (d)	Readability (e)	Platform (mm)	Item N°	Price €
i-DFP33P1500B1		1500kg	500g		1260 x 840 x 90	30712894	1350,00
i-DFP33P1500B1-M	•	1500kg		500g	1260 x 840 x 90	30764269	1430,00

Accessories	Item N°	Price €
Ethernet Kit i-DT33	30429666	111,00
USB Device Kit i-DT33	30699120	72,00
Dust Cover (5) i-DT33	30699121	37,00
Rechargeable Battery Kit i-DT33P	30699122	44,00
SF40A Printer	30064202	715,00
Cable RS232 i-DT33P	80500525	99,70
Cable RS232 i-DT33XW	80500552	91,80
Floor Scale Column Kit Painted	30531244	127,00
Floor Scale Column Kit SST	30539553	235,00
LC Extension Cable, 3.5 m *1	30803566	78,00

1 Cannot be used for legal-for-trade applications



ETHERNET KIT I-DT61XW
30429666



RECHARGEABLE BATTERY KITS
30699122




SF40A IMPACT PRINTER
30064202



FLOOR SCALE COLUMN KIT SST 980MM
30539553

DEFENDER™ 6000 Stainless Steel

Advanced IP68/IP69K Washdown Bench Scale



Application

Weighing, Counting with APW optimization, Percent Weighing, Checkweighing, Checkcounting, Filling, Dynamic Weighing/Display Hold

Display

6-digit, 7-segment white LED, 20 mm high digits

Operation

100-240 VAC / 50/60 Hz universal power supply, hardwired

Communication

Easy access communication ports including RS232 and USB host (included), second RS232/RS485/USB device, Ethernet and WiFi/Bluetooth, Discrete I/O and Analog output (sold separately as accessories); GLP/GMP data output with real-time clock

Construction

316 stainless steel platform and frame, hermetically-sealed, stainless steel load cell, IP68/IP69K 316 stainless steel housing indicator and mounting bracket, knife-resistant front panel

Design Features

150% overload capacity protection, bar matrix check light, customizable print templates, menu and key lock, selectable environmental and auto-print settings, stability indicator, overload/underload indicators, auto shut-off, auto tare, adjustable non-slip rubber feet



Defender 6000 stands up to harsh cleaning in processing and packaging environments with IP68 immersion/ IP69K washdown protection, 316 stainless steel construction and hermetically-sealed load cell.

Durable Defender 6000 is NSF certified, supporting HACCP systems with a 316 stainless steel base and sturdy indicator providing corrosion protection in food, chemical processing, and packaging.

Versatile scales include fast weight display and multiple software modes to meet food, chemical processing, and packaging needs.

Model	EC Type Approved	Capacity	Readability (d)	Readability (e)	Platform	Column	Item N°	Price €
i-D61XWE3K1S6-EU		3kg	0,5g		250mm x 250mm	260mm	30575556	2420,00
i-D61XWE6K1S6-EU		6kg	0,5g		250mm x 250mm	260mm	30575557	2420,00
i-D61XWE15K1R6-EU		15kg	1g		300mm x 300mm	260mm	30575558	2620,00
i-D61XWE30K1R6-EU		30kg	2g		300mm x 300mm	260mm	30575559	2620,00
i-D61XWE60K1L7-EU		60kg	5g		500mm x 400mm	670mm	30575560	2955,00
i-D61XWE150K1L7-EU		150kg	10g		500mm x 400mm	670mm	30575561	2955,00
i-D61XWE3K1S6-M	•	3kg		1g	250mm x 250mm	260mm	30612499	2535,00
i-D61XWE6K1S6-M	•	6kg		2g	250mm x 250mm	260mm	30612500	2535,00
i-D61XWE15K1R6-M	•	15kg		5g	300mm x 300mm	260mm	30612501	2730,00
i-D61XWE30K1R6-M	•	30kg		10g	300mm x 300mm	260mm	30612502	2730,00
i-D61XWE60K1L7-M	•	60kg		20g	500mm x 400mm	670mm	30612503	3070,00
i-D61XWE150K1L7-M	•	150kg		50g	500mm x 400mm	670mm	30612504	3070,00

DEFENDER™ 6000 Hybrid

Advanced IP68/IP69K Washdown Bench Scale



OPTION
EasyConnect

Application

Weighing, Counting with APW optimization, Percent Weighing, Checkweighing, Checkcounting, Dynamic Weighing/Display Hold

Display

6-digit, 7-segment LCD display with white backlight, 45 mm high digits

Operation

6 D (LR20) Batteries (Included)

Communication

Easy access IR communication kit (available as an accessory); GLP/GMP data output with real-time clock

Construction

316 stainless steel platform and frame, hermetically-sealed stainless steel load cell, IP68/IP69K food grade polycarbonate housing indicator and 316 stainless steel mounting bracket

Design Features

150% overload capacity protection, 3-colour check light bar, menu and key lock, selectable environmental and auto-print settings, stability indicator, overload/underload indicators, low battery indicator, auto shut-off, auto tare, adjustable rubber feet



Defender 6000 stands up to harsh cleaning in processing and packaging environments with IP68 immersion/ IP69K washdown protection, 316 stainless steel construction and hermetically-sealed load cell.

Durable Defender 6000 is NSF certified, supporting HACCP systems with a 316 stainless steel base and sturdy indicator providing corrosion protection in food, chemical processing, and packaging.

Versatile scales include fast weight display and multiple software modes to meet food, chemical processing, and packaging needs.

Model	EC Type Approved	Capacity	Readability (d)	Readability (e)	Platform	Column	Item N°	Price €
i-D61PW3K1S6		3kg	0,5g		250mm x 250mm	260mm	30575550	1975,00
i-D61PW6K1S6		6kg	0,5g		250mm x 250mm	260mm	30575551	1975,00
i-D61PW15K1R6		15kg	1g		300mm x 300mm	260mm	30575552	2170,00
i-D61PW30K1R6		30kg	2g		300mm x 300mm	260mm	30575553	2170,00
i-D61PW60K1L7		60kg	5g		500mm x 400mm	670mm	30575554	2520,00
i-D61PW150K1L7		150kg	10g		500mm x 400mm	670mm	30575555	2520,00
i-D61PW3K1S6-M	•	3kg		1g	250mm x 250mm	260mm	30612473	2090,00
i-D61PW6K1S6-M	•	6kg		2g	250mm x 250mm	260mm	30612474	2090,00
i-D61PW15K1R6-M	•	15kg		5g	300mm x 300mm	260mm	30612475	2285,00
i-D61PW30K1R6-M	•	30kg		10g	300mm x 300mm	260mm	30612476	2285,00
i-D61PW60K1L7-M	•	60kg		20g	500mm x 400mm	670mm	30612477	2630,00
i-D61PW150K1L7-M	•	150kg		50g	500mm x 400mm	670mm	30612498	2630,00

DEFENDER™ 6000 Front Mount

Advanced IP68/IP69K Washdown Bench Scale



Application

Weighing, Counting with APW optimization, Percent Weighing, Checkweighing, Checkcounting, Dynamic Weighing/Display Hold

Display

6-digit, 7-segment LCD display with white backlight, 45 mm high digits

Operation

6 D (LR20) Batteries (Included)

Communication

Easy access IR communication kit (available as an accessory); GLP/GMP data output with real-time clock

Construction

316 stainless steel platform and frame, hermetically-sealed stainless steel load cell, IP68/IP69K food grade polycarbonate housing indicator and 316 stainless steel mounting bracket

Design Features

150% overload capacity protection, 3-colour check light bar, menu and key lock, selectable environmental and auto-print settings, stability indicator, overload/underload indicators, low battery indicator, auto shut-off, auto tare, adjustable rubber feet



Defender 6000 stands up to harsh cleaning in processing and packaging environments with IP68 immersion/ IP69K washdown protection, 316 stainless steel construction and hermetically-sealed load cell.

Durable Defender 6000 is NSF certified, supporting HACCP systems with a 316 stainless steel base and sturdy indicator providing corrosion protection in food, chemical processing, and packaging.

Versatile scales include fast weight display and multiple software modes to meet food, chemical processing, and packaging needs.

Model	EC Type Approved	Capacity	Readability (d)	Readability (e)	Platform	Item N°	Price €
i-D61PW3K1S5		3kg	0,5g		250mm x 250mm	30608667	1975,00
i-D61PW6K1S5		6kg	0,5g		250mm x 250mm	30608728	1975,00
i-D61PW15K1R5		15kg	1g		300mm x 300mm	30608729	2170,00
i-D61PW30K1R5		30kg	2g		300mm x 300mm	30608730	2170,00
i-D61PW60K1L5		60kg	5g		500mm x 400mm	30608731	2495,00
i-D61PW150K1L5		150kg	10g		500mm x 400mm	30608732	2495,00
i-D61PW3K1S5-M	•	3kg		1g	250mm x 250mm	30612505	2090,00
i-D61PW6K1S5-M	•	6kg		2g	250mm x 250mm	30612506	2090,00
i-D61PW15K1R5-M	•	15kg		5g	300mm x 300mm	30612507	2285,00
i-D61PW30K1R5-M	•	30kg		10g	300mm x 300mm	30612508	2285,00
i-D61PW60K1L5-M	•	60kg		20g	500mm x 400mm	30612509	2615,00
i-D61PW150K1L5-M	•	150kg		50g	500mm x 400mm	30612510	2615,00

DEFENDER™ 6000 Accessories

Accessories	Item N°	Price €
Wi-Fi/Bluetooth Dongle i-DT61XW	30412537	282,00
Analog Output kit i-DT61XW *1	30424403	183,00
2 nd RS232/RS485/USB Device kit i-DT61XW *1	30424404	156,00
LC Extension Cable, 9 meters *2	30424409	111,00
M16 Gland with Ferrite	30379716	20,00
Light Tower, 3 Colors, 243mm, IP65 (Requires 30097591)	30424021	241,00
Discrete I/O with Relay (2 In/4 Out) i-DT61XW *1	30097591	93,00
Ethernet kit i-DT61XW	30429666	111,00
External power bank holder i-DT61XW (Requires external USB power bank)	30572907	290,00
Cable RS232 i-DT61XW	80500552	91,80
Column Kit 260 mm 316 stainless steel	30572908	235,00
Column Kit 670 mm 316 stainless steel	30572909	291,00
Front Mounting kit 316 SS i-DT61PW	30626654	235,00
Pan 305 x 305 mm 304 stainless steel	30427798	156,00
Pan 254 x 254 mm 304 stainless steel	30427799	131,00
IR communication kit i-DT61PW	30572910	159,00
Wheel kit SS Defender Series	30616409	235,00
Alibi Memory i-DT61XW	80500503	52,00
LC Extension Cable, 3.5 m *2	30803566	78,00

1 Only one of them can be used at the same time

2 Cannot be used for legal-for-trade applications



WI-FI/BLUETOOTH DONGLE I-DT61XW
30412537



ANALOG OUTPUT KIT I-DT61XW
30424403



2 ND RS232/RS485/USB DEVICE KIT
I-DT61XW 30424404



LC EXTENSION CABLE, 9 METERS
30424409



M16 GLAND WITH FERRITE
30379716



LIGHT TOWER, 3 COLORS, 243MM, IP65
30424021



DISCRETE I/O WITH RELAY (2 IN/4 OUT)
I-DT61XW 30097591



ETHERNET KIT I-DT61XW
30429666



EXTERNAL POWER BANK HOLDER
I-DT61XW 30572907



CABLE RS232 I-DT61XW
80500552



COLUMN KIT STAINLESS STEEL
30572908 30572909



FRONT MOUNTING KIT 316 SS I-DT61PW
30626654



PAN STAINLESS STEEL
30427798 30427799



IR COMMUNICATION KIT I-DT61PW
30572910



WHEEL KIT SS DEFENDER SERIES
30616409



ALIBI MEMORY I-DT61XW
80500503

DEFENDER™ 5000 Stainless Steel (Single Range)

Multifunctional Stainless Steel Bench Scale for Standard Industrial Applications

Application

Weighing, parts counting, percent weighing, check weighing, dynamic weighing (display hold), Filling

Display

Backlit Liquid Crystal Display (LCD), 256 x 80 pixels

Operation

Internal power supply with power cord (included) or rechargeable battery (sold separately as an accessory)

Communication


Easy access communication ports including RS232 and Micro SD slot (included), second RS232/RS485/USB device, Ethernet and WiFi/Bluetooth, discrete I/O and analog output (sold separately as accessories); GLP/GMP data output with real-time clock

Construction

304-grade stainless steel platform and frame, stainless steel IP68 indicator, IP67 stainless steel load cell, non-slip adjustable rubber feet

Design Features

Resolution up to dual range 15.000d/3000e, 150% overload capacity, checkweighing display and beep, menu lock switch, selectable display languages, selectable filter, stability indicator, overload/underload indicators, selectable display controls, auto print and tare



User-friendly model with an easy-to-use keypad, large backlit display and plain text on-screen prompts which guide you through operation.

Durably constructed of a stainless steel platform, frame & IP68 indicator; IP67 stainless steel loadcell (OIML/EC Type single range).

Features common industrial modes and supports a 50,000-item library with 1,000 user profiles for easy data management. Includes GMP/GLP with time/date for result traceability.

Model	EC Type Approved	Capacity	Readability (d)	Readability (e)	Platform	Column	Item N°	Price €
D52XW15WQR6-EU		15kg	1g		305mm x 305mm	350mm	30471130	1425,00
D52XW15WQL7-EU		15kg	1g		400mm x 400mm	680mm	30471131	1540,00
D52XW30WQL7-EU		30kg	2g		400mm x 400mm	680mm	30471132	1540,00
D52XW60WQL7-EU		60kg	5g		400mm x 400mm	680mm	30471133	1540,00
D52XW60WQX7-EU		60kg	5g		500mm x 500mm	680mm	30471134	1775,00
D52XW150WQX7-EU		150kg	10g		500mm x 500mm	680mm	30471135	1775,00
D52XW60WQV8-EU		60kg	5g		610mm x 610mm	980mm	30471136	1955,00
D52XW150WQV8-EU		150kg	10g		610mm x 610mm	980mm	30471137	1955,00
D52XW15WQR6-M	•	15kg		5g	305mm x 305mm	350mm	30471158	1520,00
D52XW15WQL7-M	•	15kg		5g	400mm x 400mm	680mm	30471159	1630,00
D52XW30WQL7-M	•	30kg		10g	400mm x 400mm	680mm	30471160	1630,00
D52XW60WQL7-M	•	60kg		20g	400mm x 400mm	680mm	30471161	1630,00
D52XW60WQX7-M	•	60kg		20g	500mm x 500mm	680mm	30471162	1860,00
D52XW150WQX7-M	•	150kg		50g	500mm x 500mm	680mm	30471163	1860,00
D52XW60WQV8-M	•	60kg		20g	610mm x 610mm	980mm	30471164	2045,00
D52XW150WQV8-M	•	150kg		50g	610mm x 610mm	980mm	30471165	2045,00

DEFENDER™5000 Stainless Steel (Dual Range)

Multifunctional Stainless Steel Bench Scale for Standard Industrial Applications

Application

Weighing, parts counting, percent weighing, check weighing, dynamic weighing (display hold), Filling

Display

Backlit Liquid Crystal Display (LCD), 256 x 80 pixels

Operation

Internal power supply with power cord (included) or rechargeable battery (sold separately as an accessory)

Communication


Easy access communication ports including RS232 and Micro SD slot (included), second RS232/RS485/USB device, Ethernet and WiFi/Bluetooth, Discrete I/O and analog output (sold separately as accessories); GLP/GMP data output with real-time clock

Construction

304-grade stainless steel platform and frame, stainless steel IP68 indicator, IP67 aluminium load cell, non-slip adjustable rubber feet

Design Features

Resolution up to dual range 15.000d/3000e, 150% overload capacity, check weighing display and beep, menu lock switch, selectable display languages, selectable filter, stability indicator, overload/underload indicators, selectable display controls, auto print and tare



User-friendly model with an easy-to-use keypad, large backlit display and plain text on-screen prompts which guide you through operation.


Durably constructed of a stainless steel platform, frame & IP 68 indicator; IP 67 aluminum load cell (OIML/EC Type dual range).

Features common industrial modes and supports a 50,000-item library with 1,000 user profiles for easy data management. Includes GMP/GLP with time/date for result traceability.

Model	EC Type Approved	Capacity	Readability (d)	Readability (e)	Platform	Column	Item N°	Price €
D52XW6WQDR6-EU		3kg / 6kg	0.2g / 0.5g		305mm x 305mm	350mm	30444791	1245,00
D52XW15WQDR6-EU		6kg / 15kg	0.5g / 1g		305mm x 305mm	350mm	30444792	1245,00
D52XW15WQDL7-EU		6kg / 15kg	0.5g / 1g		400mm x 400mm	680mm	30444793	1345,00
D52XW30WQDL7-EU		15kg / 30kg	1g / 2g		400mm x 400mm	680mm	30444794	1345,00
D52XW60WQDL7-EU		30kg / 60kg	2g / 5g		400mm x 400mm	680mm	30444795	1345,00
D52XW60WQDX7-EU		30kg / 60kg	2g / 5g		500mm x 500mm	680mm	30444796	1535,00
D52XW150WQDX7-EU		60kg / 150kg	5g / 10g		500mm x 500mm	680mm	30444797	1535,00
D52XW60WQDV8-EU		30kg / 60kg	2g / 5g		610mm x 610mm	980mm	30444798	1720,00
D52XW150WQDV8-EU		60kg / 150kg	5g / 10g		610mm x 610mm	980mm	30444799	1720,00
D52XW300WQDV8-EU		150kg / 300kg	10g / 20g		610mm x 610mm	980mm	30444800	1720,00
D52XW150WTDV8-EU		60kg / 150kg	5g / 10g		800mm x 600mm	980mm	30444801	2125,00
D52XW300WTDV8-EU		150kg / 300kg	10g / 20g		800mm x 600mm	980mm	30444802	2125,00
D52XW600WTDV8-EU		300kg / 600kg	20g / 50g		800mm x 600mm	980mm	30444803	2125,00
D52XW6WQDR6-M	•	3kg / 6kg		1g / 2g	305mm x 305mm	350mm	30461436	1345,00
D52XW15WQDR6-M	•	6kg / 15kg		2g / 5g	305mm x 305mm	350mm	30461437	1345,00
D52XW15WQDL7-M	•	6kg / 15kg		2g / 5g	400mm x 400mm	680mm	30461568	1430,00
D52XW30WQDL7-M	•	15kg / 30kg		5g / 10g	400mm x 400mm	680mm	30461569	1430,00
D52XW60WQDL7-M	•	30kg / 60kg		10g / 20g	400mm x 400mm	680mm	30461570	1430,00
D52XW60WQDX7-M	•	30kg / 60kg		10g / 20g	500mm x 500mm	680mm	30461571	1630,00
D52XW150WQDX7-M	•	60kg / 150kg		20g / 50g	500mm x 500mm	680mm	30461572	1630,00
D52XW60WQDV8-M	•	30kg / 60kg		10g / 20g	610mm x 610mm	980mm	30461573	1805,00
D52XW150WQDV8-M	•	60kg / 150kg		20g / 50g	610mm x 610mm	980mm	30461574	1805,00
D52XW300WQDV8-M	•	150kg / 300kg		50g / 100g	610mm x 610mm	980mm	30461575	1805,00
D52XW150WTDV8-M	•	60kg / 150kg		20g / 50g	800mm x 600mm	980mm	30461576	2230,00
D52XW300WTDV8-M	•	150kg / 300kg		50g / 100g	800mm x 600mm	980mm	30461577	2230,00
D52XW600WTDV8-M	•	300kg / 600kg		100g / 200g	800mm x 600mm	980mm	30461578	2230,00

DEFENDER™ 5000 Hybrid

Multifunctional Bench Scale for Standard Industrial Applications



Application

Weighing, parts counting, percent weighing, check weighing, dynamic weighing (display hold), Filling

Display

Backlit Liquid Crystal Display (LCD), 256 x 80 pixels

Operation

Internal power supply with power cord (included) or rechargeable battery (sold separately as an accessory)

Communication

Easy access communication ports including RS232 and Micro SD slot (included), second RS232/RS485/USB device, Ethernet and WiFi/Bluetooth, Discrete I/O and Analog Output (sold separately as accessories); GLP/GMP data output with real-time clock

Construction

IP68 stainless steel indicator, 304-grade stainless steel platform with sand gray powder-coated steel frame, IP67 aluminium load cell, non-slip adjustable rubber feet

Design Features

User-configurable resolution up to dual range 15,000d/3,000e, 150% overload capacity protection, check weighing symbols with selectable operation and audible signal settings, menu lock switch, multiple selectable operating languages, selectable environmental and auto-print settings, auto tare



User-friendly model with an easy-to-use keypad, large backlit display and plain text on-screen prompts which guide you through operation.


Durable construction with a stainless steel platform, sand gray powder-coated steel frame & aluminum load cell. Defender meets OIML/EC type approvals.

Features common industrial modes and supports a 50,000-item library with 1,000 user profiles for easy data management. Includes GMP/GLP with time/date for result traceability.

Model	EC Type Approved	Capacity	Readability (d)	Readability (e)	Platform	Column	Item N°	Price €
D52XW6RQDR1-EU		3kg / 6kg	0.2g / 0.5g		305mm x 305mm	350mm	30444804	990,00
D52XW15RQDR1-EU		6kg / 15kg	0.5g / 1g		305mm x 305mm	350mm	30444805	990,00
D52XW15RQDL2-EU		6kg / 15kg	0.5g / 1g		400mm x 400mm	680mm	30444806	1095,00
D52XW30RQDL2-EU		15kg / 30kg	1g / 2g		400mm x 400mm	680mm	30444807	1095,00
D52XW60RQDL2-EU		30kg / 60kg	2g / 5g		400mm x 400mm	680mm	30444808	1095,00
D52XW60RQDX2-EU		30kg / 60kg	2g / 5g		500mm x 500mm	680mm	30444809	1375,00
D52XW150RQDX2-EU		60kg / 150kg	5g / 10g		500mm x 500mm	680mm	30444810	1375,00
D52XW60RQDV3-EU		30kg / 60kg	2g / 5g		610mm x 610mm	980mm	30444811	1480,00
D52XW150RQDV3-EU		60kg / 150kg	5g / 10g		610mm x 610mm	980mm	30444812	1480,00
D52XW300RQDV3-EU		150kg / 300kg	10g / 20g		610mm x 610mm	980mm	30444813	1480,00
D52XW150RTDV3-EU		60kg / 150kg	5g / 10g		800mm x 600mm	980mm	30444814	1775,00
D52XW300RTDV3-EU		150kg / 300kg	10g / 20g		800mm x 600mm	980mm	30444815	1775,00
D52XW600RTDV3-EU		300kg / 600kg	20g / 50g		800mm x 600mm	980mm	30444816	1775,00
D52XW6RQDR1-M	•	3kg / 6kg		1g / 2g	305mm x 305mm	350mm	30461579	1090,00
D52XW15RQDR1-M	•	6kg / 15kg		2g / 5g	305mm x 305mm	350mm	30461580	1090,00
D52XW15RQDL2-M	•	6kg / 15kg		2g / 5g	400mm x 400mm	680mm	30461581	1185,00
D52XW30RQDL2-M	•	15kg / 30kg		5g / 10g	400mm x 400mm	680mm	30461582	1185,00
D52XW60RQDL2-M	•	30kg / 60kg		10g / 20g	400mm x 400mm	680mm	30461583	1185,00
D52XW60RQDX2-M	•	30kg / 60kg		10g / 20g	500mm x 500mm	680mm	30461584	1460,00
D52XW150RQDX2-M	•	60kg / 150kg		20g / 50g	500mm x 500mm	680mm	30461585	1460,00
D52XW60RQDV3-M	•	30kg / 60kg		10g / 20g	610mm x 610mm	980mm	30461586	1580,00
D52XW150RQDV3-M	•	60kg / 150kg		20g / 50g	610mm x 610mm	980mm	30461587	1580,00
D52XW300RQDV3-M	•	150kg / 300kg		50g / 100g	610mm x 610mm	980mm	30461588	1580,00
D52XW150RTDV3-M	•	60kg / 150kg		20g / 50g	800mm x 600mm	980mm	30461589	1860,00
D52XW300RTDV3-M	•	150kg / 300kg		50g / 100g	800mm x 600mm	980mm	30461590	1860,00
D52XW600RTDV3-M	•	300kg / 600kg		100g / 200g	800mm x 600mm	980mm	30461591	1860,00

DEFENDER™5000 Standard

Multifunctional Bench Scale for Standard Industrial Applications



Application

Weighing, parts counting, percent weighing, check weighing, dynamic weighing (display hold), Filling

Display

Backlit Liquid Crystal Display (LCD), 256 x 80 pixels

Operation

Internal power supply with power cord (included) or rechargeable battery (sold separately as an accessory)

Communication

Easy access communication ports including RS232 and Micro SD slot (included), second RS232/RS485/USB device, Ethernet and WiFi/Bluetooth, discrete I/O and analog output (sold separately as accessories); GLP/GMP data output with real-time clock

Construction

ABS plastic indicator, 304-grade stainless steel platform with sand gray powder-coated steel frame, IP67 aluminium load cell, non-slip adjustable rubber feet

Design Features

User-configurable resolution up to dual range 15,000d/3,000e, 150% overload capacity protection, check weighing symbols with selectable operation and audible signal settings, menu lock switch, multiple selectable operating languages, selectable environmental and auto-print settings, auto tare



User-friendly model with an easy-to-use keypad, large backlit display and plain text on-screen prompts which guide you through operation.


Durable construction with a stainless steel platform, sand gray powder-coated steel frame & aluminum load cell. Defender meets OIML/EC type approvals.

Features common industrial modes and supports a 50,000-item library with 1,000 user profiles for easy data management. Includes GMP/GLP with time/date for result traceability.

Model	EC Type Approved	Capacity	Readability (d)	Readability (e)	Platform	Column	Item N°	Price €
D52P6RQDR1		3kg / 6kg	0.2g / 0.5g		305mm x 305mm	350mm	30424312	760,00
D52P15RQDR1		6kg / 15kg	0.5g / 1g		305mm x 305mm	350mm	30424313	760,00
D52P15RQDL2		6kg / 15kg	0.5g / 1g		400mm x 400mm	680mm	30424314	860,00
D52P30RQDL2		15kg / 30kg	1g / 2g		400mm x 400mm	680mm	30424315	860,00
D52P60RQDL2		30kg / 60kg	2g / 5g		400mm x 400mm	680mm	30424316	860,00
D52P60RQDX2		30kg / 60kg	2g / 5g		500mm x 500mm	680mm	30424317	1135,00
D52P150RQDX2		60kg / 150kg	5g / 10g		500mm x 500mm	680mm	30424318	1135,00
D52P60RQDV3		30kg / 60kg	2g / 5g		610mm x 610mm	980mm	30424319	1255,00
D52P150RQDV3		60kg / 150kg	5g / 10g		610mm x 610mm	980mm	30424320	1255,00
D52P300RQDV3		150kg / 300kg	10g / 20g		610mm x 610mm	980mm	30424321	1255,00
D52P150RTDV3		60kg / 150kg	5g / 10g		800mm x 600mm	980mm	30424322	1535,00
D52P300RTDV3		150kg / 300kg	10g / 20g		800mm x 600mm	980mm	30424323	1535,00
D52P600RTDV3		300kg / 600kg	20g / 50g		800mm x 600mm	980mm	30424324	1535,00
D52P6RQDR1-M	•	3kg / 6kg		1g / 2g	305mm x 305mm	350mm	30461545	855,00
D52P15RQDR1-M	•	6kg / 15kg		2g / 5g	305mm x 305mm	350mm	30461546	855,00
D52P15RQDL2-M	•	6kg / 15kg		2g / 5g	400mm x 400mm	680mm	30461547	945,00
D52P30RQDL2-M	•	15kg / 30kg		5g / 10g	400mm x 400mm	680mm	30461548	945,00
D52P60RQDL2-M	•	30kg / 60kg		10g / 20g	400mm x 400mm	680mm	30461549	945,00
D52P60RQDX2-M	•	30kg / 60kg		10g / 20g	500mm x 500mm	680mm	30461550	1230,00
D52P150RQDX2-M	•	60kg / 150kg		20g / 50g	500mm x 500mm	680mm	30461551	1230,00
D52P60RQDV3-M	•	30kg / 60kg		10g / 20g	610mm x 610mm	980mm	30461552	1350,00
D52P150RQDV3-M	•	60kg / 150kg		20g / 50g	610mm x 610mm	980mm	30461553	1350,00
D52P300RQDV3-M	•	150kg / 300kg		50g / 100g	610mm x 610mm	980mm	30461554	1350,00
D52P150RTDV3-M	•	60kg / 150kg		20g / 50g	800mm x 600mm	980mm	30461555	1630,00
D52P300RTDV3-M	•	150kg / 300kg		50g / 100g	800mm x 600mm	980mm	30461556	1630,00
D52P600RTDV3-M	•	300kg / 600kg		100g / 200g	800mm x 600mm	980mm	30461557	1630,00

DEFENDER™ 5000 Front Mount

Multifunctional Stainless Steel Bench Scale for Standard Industrial Applications



Application

Weighing, parts counting, percent weighing, check weighing, dynamic weighing (display hold), Filling

Display

Backlit Liquid Crystal Display (LCD), 256 x 80 pixels

Operation

Internal power supply with power cord (included) or rechargeable battery (sold separately as an accessory)

Communication

Easy access communication ports including RS232 and Micro SD slot (included), second RS232/RS485/USB device, Ethernet and WiFi/Bluetooth, discrete I/O and analog output (sold separately as accessories); GLP/GMP data output with real-time clock

Construction

304-grade stainless steel platform, frame and handle, stainless steel IP68 indicator, IP67 aluminium (for dual range models) or stainless steel load cell (for single range models), non-slip adjustable rubber feet

Design Features

Resolution up to dual range 15.000d/3000e, 150% overload capacity, check weighing display and beep, menu lock switch, selectable display languages, selectable filter, stability indicator, overload/underload indicators, selectable display controls, auto print and tare



User-friendly model with an easy-to-use keypad, large backlit display and plain text on-screen prompts which guide you through operation.

Durably constructed of a stainless steel platform, handle, frame & IP68 indicator; IP67 stainless steel load cell (OIML/EC Type single range) or aluminium load cell (OIML/EC Type dual range).

Features common industrial modes and supports a 50,000-item library with 1,000 user profiles for easy data management. Includes GMP/GLP with time/date for result traceability.

Dual range with aluminium load cell


Model	EC Type Approved	Capacity	Readability (d)	Readability (e)	Platform	Item N°	Price €
D52XW6WQDR5-EU		3kg / 6kg	0.2g / 0.5g		305mm x 305mm	30444765	1245,00
D52XW15WQDR5-EU		6kg / 15kg	0.5g / 1g		305mm x 305mm	30444766	1245,00
D52XW15WQDL5-EU		6kg / 15kg	0.5g / 1g		400mm x 400mm	30444767	1345,00
D52XW30WQDL5-EU		15kg / 30kg	1g / 2g		400mm x 400mm	30444788	1345,00
D52XW60WQDL5-EU		30kg / 60kg	2g / 5g		400mm x 400mm	30444789	1345,00
D52XW60WQDX5-EU		30kg / 60kg	2g / 5g		500mm x 500mm	30444790	1540,00
D52XW150WQDX5-EU		60kg / 150kg	5g / 10g		500mm x 500mm	30461428	1540,00
D52XW6WQDR5-M	•	3kg / 6kg		1g / 2g	305mm x 305mm	30461565	1345,00
D52XW15WQDR5-M	•	6kg / 15kg		2g / 5g	305mm x 305mm	30461566	1345,00
D52XW15WQDL5-M	•	6kg / 15kg		2g / 5g	400mm x 400mm	30461431	1430,00
D52XW30WQDL5-M	•	15kg / 30kg		5g / 10g	400mm x 400mm	30461432	1430,00
D52XW60WQDL5-M	•	30kg / 60kg		10g / 20g	400mm x 400mm	30461433	1430,00
D52XW60WQDX5-M	•	30kg / 60kg		10g / 20g	500mm x 500mm	30461434	1630,00
D52XW150WQDX5-M	•	60kg / 150kg		20g / 50g	500mm x 500mm	30461435	1630,00

Single range with stainless steel load cell

Model	EC Type Approved	Capacity	Readability (d)	Readability (e)	Platform	Item N°	Price €
D52XW15WQR5-EU		15kg	1g		305mm x 305mm	30471124	1295,00
D52XW15WQL5-EU		15kg	1g		400mm x 400mm	30471125	1425,00
D52XW30WQL5-EU		30kg	2g		400mm x 400mm	30471126	1425,00
D52XW60WQL5-EU		60kg	5g		400mm x 400mm	30471127	1425,00
D52XW60WQX5-EU		60kg	5g		500mm x 500mm	30471128	1595,00
D52XW150WQX5-EU		150kg	10g		500mm x 500mm	30471129	1595,00
D52XW15WQR5-M	•	15kg		5g	305mm x 305mm	30471152	1390,00
D52XW15WQL5-M	•	15kg		5g	400mm x 400mm	30471153	1515,00
D52XW30WQL5-M	•	30kg		10g	400mm x 400mm	30471154	1515,00
D52XW60WQL5-M	•	60kg		20g	400mm x 400mm	30471155	1515,00
D52XW60WQX5-M	•	60kg		20g	500mm x 500mm	30471156	1690,00
D52XW150WQX5-M	•	150kg		50g	500mm x 500mm	30471157	1690,00

DEFENDER™5000 Front Mount

Multifunctional Bench Scale for Standard Industrial Applications



Application

Weighing, parts counting, percent weighing, check weighing, dynamic weighing (display hold), Filling

Display

Backlit Liquid Crystal Display (LCD), 256 x 80 pixels

Operation

Internal power supply with power cord (included) or rechargeable battery (sold separately as an accessory)

Communication

Easy access communication ports including RS232 and Micro SD slot (included), second RS232/RS485/USB device, Ethernet and WiFi/Bluetooth, discrete I/O and analog output (sold separately as accessories); GLP/GMP data output with real-time clock

Construction

ABS plastic (TD52P) or stainless steel IP68(TD52XW) indicator, 304-grade stainless steel platform and handle with sand gray powder-coated steel frame, IP67 aluminium load cell, non-slip adjustable rubber feet

Design Features

User-configurable resolution up to dual range 15,000d/3,000e, 150% overload capacity protection, check weighing symbols with selectable operation and audible signal settings, menu lock switch, multiple selectable operating languages, selectable environmental and auto-print settings, auto tare



User-friendly model with an easy-to-use keypad, large backlit display and plain text on-screen prompts which guide you through operation.

Durable construction with a stainless steel platform and handle, sand gray powder-coated steel frame & aluminum load cell. Defender meets OIML/EC type approvals (depending on region).

Features common industrial modes and supports a 50,000-item library with 1,000 user profiles for easy data management. Includes GMP/GLP with time/date for result traceability.

With TD52P ABS Indicator

Model	EC Type Approved	Capacity	Readability (d)	Readability (e)	Platform	Item N°	Price €
D52P6RQDR5		3kg / 6kg	0.2g / 0.5g		305mm x 305mm	30424325	760,00
D52P15RQDR5		6kg / 15kg	0.5g / 1g		305mm x 305mm	30424326	760,00
D52P15RQDL5		6kg / 15kg	0.5g / 1g		400mm x 400mm	30424327	860,00
D52P30RQDL5		15kg / 30kg	1g / 2g		400mm x 400mm	30444761	860,00
D52P60RQDL5		30kg / 60kg	2g / 5g		400mm x 400mm	30444762	860,00
D52P60RQDX5		30kg / 60kg	2g / 5g		500mm x 500mm	30444763	1135,00
D52P150RQDX5		60kg / 150kg	5g / 10g		500mm x 500mm	30444764	1135,00
D52P6RQDR5-M	•	3kg / 6kg		1g / 2g	305mm x 305mm	30461558	855,00
D52P15RQDR5-M	•	6kg / 15kg		2g / 5g	305mm x 305mm	30461559	855,00
D52P15RQDL5-M	•	6kg / 15kg		2g / 5g	400mm x 400mm	30461560	945,00
D52P30RQDL5-M	•	15kg / 30kg		5g / 10g	400mm x 400mm	30461561	945,00
D52P60RQDL5-M	•	30kg / 60kg		10g / 20g	400mm x 400mm	30461562	945,00
D52P60RQDX5-M	•	30kg / 60kg		10g / 20g	500mm x 500mm	30461563	1230,00
D52P150RQDX5-M	•	60kg / 150kg		20g / 50g	500mm x 500mm	30461564	1230,00

With TD52XW Stainless Steel Indicator

Model	EC Type Approved	Capacity	Readability (d)	Readability (e)	Platform	Item N°	Price €
D52XW6RQDR5		3kg / 6kg	0.2g / 0.5g		305mm x 305mm	30498216	990,00
D52XW15RQDR5		6kg / 15kg	0.5g / 1g		305mm x 305mm	30498217	990,00
D52XW15RQDL5		6kg / 15kg	0.5g / 1g		400mm x 400mm	30498498	1095,00
D52XW30RQDL5		15kg / 30kg	1g / 2g		400mm x 400mm	30498499	1095,00
D52XW60RQDL5		30kg / 60kg	2g / 5g		400mm x 400mm	30498500	1095,00
D52XW60RQDX5		30kg / 60kg	2g / 5g		500mm x 500mm	30498501	1375,00
D52XW150RQDX5		60kg / 150kg	5g / 10g		500mm x 500mm	30498502	1375,00
D52XW6RQDR5-M	•	3kg / 6kg		1g / 2g	305mm x 305mm	30461610	1180,00
D52XW15RQDR5-M	•	6kg / 15kg		2g / 5g	305mm x 305mm	30461611	1180,00
D52XW15RQDL5-M	•	6kg / 15kg		2g / 5g	400mm x 400mm	30461612	1270,00
D52XW30RQDL5-M	•	15kg / 30kg		5g / 10g	400mm x 400mm	30461613	1270,00
D52XW60RQDL5-M	•	30kg / 60kg		10g / 20g	400mm x 400mm	30461614	1270,00
D52XW60RQDX5-M	•	30kg / 60kg		10g / 20g	500mm x 500mm	30461615	1565,00
D52XW150RQDX5-M	•	60kg / 150kg		20g / 50g	500mm x 500mm	30461616	1565,00

DEFENDER™5000 Accessories

Accessories	Item N°	Price €
WiFi/BT Dongle, OHAUS (Requires 30424406)	30412537	282,00
Analog Output Kit, TD52 *3	30424403	183,00
2nd RS232/RS485/USB, TD52 *3	30424404	156,00
Battery Kit, Li-ion, TD52	30424405	127,00
USB Host for Dongle, TD52 *2	30424406	111,00
Light Tower, R-Y-G, IP65, 243 mm, 24VDC (Requires 30097591)	30424021	241,00
In-use Cover Set, TD52P	30424022	20,00
In-use Cover Set, TD52XW	30424023	27,00
In-use Cover Set (10), TD52P	30469948	84,00
In-use Cover Set (10), TD52XW	30469949	149,00
Wall Mount Kit, SST, TD52XW	30424026	46,00
Wall Mount Kit, CS, TD52P	30424027	41,00
LC Extension Cable, 9 meters, D52 *1	30424409	111,00
M16 Gland with Ferrite, to TD52	30379716	20,00
Micro SD Card, 8G	30303533	35,00
I/O with Relay (2 In/4 Out) *3	30097591	93,00
Ethernet Kit, D52 *2	30429666	111,00
Carrying Handle Kit, D52xx	30424024	93,00
Printer Holder Kit, D52xx	30424025	52,00
Cable Assembly 6m to J-box DF52	30424408	91,00
Column Kit, 35 cm, CS, D52xx	30424410	74,00
Column Kit, 68 cm, CS, D52xx	30424411	93,00
Column Kit, 98 cm, CS, D52xx	30424412	104,00
Column Kit, 35 cm, SST, D52xx	30424413	99,00
Column Kit, 68 cm, SST, D52xx	30424414	165,00
Column Kit, 98 cm, SST, D52xx	30424415	189,00
Cable Kit, Extension for base, D52	30429931	46,00
Wheel Kit, Defender Series (up to 300 kg models) *4	30372146	65,00
USB Host for Flash Drive *2	30572923	128,00
Interface Kit RS232-USB *5	30304101	78,25
Interface Kit RS232-Ethernet *5	30304102	116,00
Cable, RS232, 1.5m, to connect a reference balance *6	30057595	78,40
LC Extension Cable, 3.5 m *1	30803566	78,00

1 Cannot be used for legal-for-trade applications

2 Only one of them can be used at the same time

3 Only one of them can be used at the same time

4 Not available for models with 600 kg capacity

5 Suitable for all Defender scales and Terminals equipped with RS232

6 Reference Balance automatically accepts higher-precision APW data from a scale through the included RRS232 port



ANALOG OUTPUT KIT, TD52
30424403



WIFI/BT DONGLE, OHAUS
30412537



2ND RS232/RS485/USB, TD52
30424404



BATTERY KIT, LI-ION, TD52
30424405



USB HOST FOR DONGLE, TD52
30424406



LIGHT TOWER, R-Y-G, IP65, 243 MM,
24VDC 30424021



I/O WITH RELAY (2 IN/4 OUT)
30097591



PRINTER HOLDER KIT
30424025



WHEEL KIT, DEFENDER SERIES
30372146



INTERFACE KIT RS232-USB
30304101




CABLE KIT, EXTENSION FOR BASE
30429931



COLUMN KIT, CS

DEFENDER™3000 Stainless Steel

All stainless steel modular scale for washdown environments



Application

Weighing, Checkweighing, Counting, Totalization.

Display

7 Segment LCD with multicolour backlight

Operation

Power cord (included) or Rechargeable Battery (included)

Communication

Standard RS232 and optional USB Device or Ethernet

Construction

304-grade stainless steel platform and frame, stainless steel IP66 indicator, stainless steel IP67 load cell, non-slip adjustable rubber feet

Design Features

150% overload protection



Results are clearly shown using the large, bright backlit LCD display. In Checkweighing mode display changes colour (yellow/green/red) to alert when within the target weight range.

This tough industrial scale features a 304-stainless steel platform and frame, IP67 load cell. Indicator features an IP66 washdown protected 304-stainless steel housing with mounting bracket.

Features software for Parts Counting, Checkweighing and Totalization. Includes RS232 with optional USB Device or Ethernet for printing or data connectivity. Approved models are equipped with EasyConnect™ solution.

Model	EC Type Approved	Capacity	Readability (d)	Readability (e)	Platform	Column	Item N°	Price €
i-D33XW15C1R6-EU		15kg	2g		305mm x 355mm	330mm	30684775	740,00
i-D33XW30C1R6-EU		30kg	5g		305mm x 355mm	330mm	30684776	765,00
i-D33XW60C1R6-EU		60kg	10g		305mm x 355mm	330mm	30684777	765,00
i-D33XW60C1L7-EU		60kg	10g		420mm x 550mm	680mm	30684778	1040,00
i-D33XW150C1L7-EU		150kg	20g		420mm x 550mm	680mm	30684779	1040,00
i-D33XW150C1X7-EU		150kg	20g		500mm x 650mm	680mm	30684780	1405,00
i-D33XW300C1X7-EU		300kg	50g		500mm x 650mm	680mm	30684781	1405,00
i-D33XW15C1R6EU-M	•	15kg		5g	305mm x 355mm	330mm	30685134	825,00
i-D33XW30C1R6EU-M	•	30kg		10g	305mm x 355mm	330mm	30685135	840,00
i-D33XW60C1R6EU-M	•	60kg		20g	305mm x 355mm	330mm	30685136	840,00
i-D33XW60C1L7EU-M	•	60kg		20g	420mm x 550mm	680mm	30685137	1120,00
i-D33XW150C1L7EU-M	•	150kg		50g	420mm x 550mm	680mm	30685138	1120,00
i-D33XW150C1X7EU-M	•	150kg		50g	500mm x 650mm	680mm	30685139	1480,00
i-D33XW300C1X7EU-M	•	300kg		100g	500mm x 650mm	680mm	30685140	1480,00
i-D33XW15C1R5-EU		15kg	2g		305mm x 355 mm	-	30684782	765,00
i-D33XW30C1R5-EU		30kg	5g		305mm x 355 mm	-	30684783	790,00
i-D33XW60C1R5-EU		60kg	10g		305mm x 355 mm	-	30685078	790,00
i-D33XW60C1L5-EU		60kg	10g		420mm x 550 mm	-	30685079	1070,00
i-D33XW150C1L5-EU		150kg	20g		420mm x 550 mm	-	30685080	1070,00
i-D33XW150C1X5-EU		150kg	20g		500mm x 650 mm	-	30685081	1445,00
i-D33XW300C1X5-EU		300kg	50g		500mm x 650 mm	-	30685082	1445,00
i-D33XW15C1R5EU-M	•	15kg		5g	305mm x 355 mm	-	30685141	890,00
i-D33XW30C1R5EU-M	•	30kg		10g	305mm x 355 mm	-	30685142	905,00
i-D33XW60C1R5EU-M	•	60kg		20g	305mm x 355 mm	-	30685143	905,00
i-D33XW60C1L5EU-M	•	60kg		20g	420mm x 550 mm	-	30685144	1170,00
i-D33XW150C1L5EU-M	•	150kg		50g	420mm x 550 mm	-	30685145	1170,00
i-D33XW150C1X5EU-M	•	150kg		50g	500mm x 650 mm	-	30685146	1540,00
i-D33XW300C1X5EU-M	•	300kg		100g	500mm x 650 mm	-	30685147	1540,00

Handle

DEFENDER™3000 Hybrid

The modular scale for industrial environments



OPTION OPTION OPTION

EasyConnect



Application

Weighing, Checkweighing, Counting, Totalization

Display

7 Segment LCD with multicolour backlight

Operation

Power cord (included), or Rechargeable Battery (included)

Communication

Standard RS232 and optional USB Device or Ethernet

Construction Hybrid

Stainless steel indicator, 304-stainless steel platform, powder-coated steel tubular frame, aluminum load cell, non-slip adjustable rubber feet

Design Features

150% overload protection



Results are clearly shown using the large bright backlit LCD display. In Checkweighing mode display changes colour (yellow/green/red) to alert when within the target weight range.

This industrial scale features a 304 stainless steel platform, powder-coated steel tubular frame, and IP65 aluminum load cell. Indicator features a stainless steel housing and mounting bracket.


Features software for Parts Counting, Checkweighing and Totalization. Includes RS232 with optional USB Device or Ethernet for printing or data connectivity. Approved models are equipped with EasyConnect™ solution.

Model	EC Type Approved	Capacity	Readability (d)	Readability (e)	Platform	Column	Item N°	Price €
i-D33XW15B1R1-EU		15kg	2g		305mm x 355mm	330mm	30685244	590,00
i-D33XW30B1R1-EU		30kg	5g		305mm x 355mm	330mm	30685245	605,00
i-D33XW60B1R1-EU		60kg	10g		305mm x 355mm	330mm	30685246	605,00
i-D33XW60B1L2-EU		60kg	10g		420mm x 550mm	680mm	30685247	655,00
i-D33XW150B1L2-EU		150kg	20g		420mm x 550mm	680mm	30685248	655,00
i-D33XW150B1X2-EU		150kg	20g		500mm x 650mm	680mm	30685249	785,00
i-D33XW300B1X2-EU		300kg	50g		500mm x 650mm	680mm	30685250	785,00
i-D33XW300B1V3-EU		300kg	50g		600mm x 800mm	980mm	30685251	1045,00
i-D33XW600B1V3-EU		600kg	100g		600mm x 800mm	980mm	30685252	1045,00
i-D33XW15B1R1EU-M	•	15kg		5g	305mm x 355mm	330mm	30685262	660,00
i-D33XW30B1R1EU-M	•	30kg		10g	305mm x 355mm	330mm	30685263	680,00
i-D33XW60B1R1EU-M	•	60kg		20g	305mm x 355mm	330mm	30685264	680,00
i-D33XW60B1L2EU-M	•	60kg		20g	420mm x 550mm	680mm	30685265	730,00
i-D33XW150B1L2EU-M	•	150kg		50g	420mm x 550mm	680mm	30685266	730,00
i-D33XW150B1X2EU-M	•	150kg		50g	500mm x 650mm	680mm	30685267	845,00
i-D33XW300B1X2EU-M	•	300kg		100g	500mm x 650mm	680mm	30685268	845,00
i-D33XW300B1V3EU-M	•	300kg		100g	600mm x 800mm	980mm	30685269	1095,00
i-D33XW600B1V3EU-M	•	600kg		200g	600mm x 800mm	980mm	30685270	1095,00
i-D33XW15B1R5-EU		15kg	2g		305mm x 355mm	-	30685253	585,00
i-D33XW30B1R5-EU		30kg	5g		305mm x 355mm	-	30685254	585,00
i-D33XW60B1R5-EU		60kg	10g		305mm x 355mm	-	30685255	585,00
i-D33XW60B1L5-EU		60kg	10g		420mm x 550mm	-	30685256	695,00
i-D33XW150B1L5-EU		150kg	20g		420mm x 550mm	-	30685257	695,00
i-D33XW150B1X5-EU		150kg	20g		500mm x 650mm	-	30685258	855,00
i-D33XW300B1X5-EU		300kg	50g		500mm x 650mm	-	30685259	855,00
i-D33XW15B1R5EU-M	•	15kg		5g	305mm x 355mm	-	30685271	680,00
i-D33XW30B1R5EU-M	•	30kg		10g	305mm x 355mm	-	30685272	680,00
i-D33XW60B1R5EU-M	•	60kg		20g	305mm x 355mm	-	30685273	680,00
i-D33XW60B1L5EU-M	•	60kg		20g	420mm x 550mm	-	30685274	790,00
i-D33XW150B1L5EU-M	•	150kg		50g	420mm x 550mm	-	30685275	790,00
i-D33XW150B1X5EU-M	•	150kg		50g	500mm x 650mm	-	30685276	900,00
i-D33XW300B1X5EU-M	•	300kg		100g	500mm x 650mm	-	30685277	900,00

Handle

DEFENDER™3000 Standard

The best selling modular scale for basic industrial needs



Application

Weighing, Checkweighing, Counting, Totalization.

Display

7 Segment LCD with multicolour backlight

Operation

Power cord (included), 6C (LR14) Batteries (not included) or Rechargeable Battery (not included)

Communication

Standard RS232 and optional USB Device or Ethernet

Construction

ABS plastic indicator, 304-stainless steel platform, powder-coated steel tubular frame, aluminum load cell, non-slip adjustable rubber feet

Design Features

150% overload protection



Results are clearly shown using the large, bright backlit LCD display. In Checkweighing mode display changes colour (yellow/green/red) to alert when within the target weight range.

This tough industrial scale features a 304-stainless steel platform, powder-coated steel tubular frame, and IP65 aluminum load cell. Indicator features a durable ABS housing and steel mounting bracket.

Features software for Parts Counting, Checkweighing and Totalization. Includes RS232 with optional USB Device or Ethernet for printing or data connectivity. Approved models are equipped with EasyConnect™ solution.

Model	EC Type Approved	Capacity	Readability (d)	Readability (e)	Platform	Column	Item N°	Price €
i-D33P15B1R1		15kg	2g		305mm x 355mm	330mm	30684738	435,00
i-D33P30B1R1		30kg	5g		305mm x 355mm	330mm	30684739	440,00
i-D33P60B1R1		60kg	10g		305mm x 355mm	330mm	30684740	440,00
i-D33P60B1L2		60kg	10g		420mm x 550mm	680mm	30684744	500,00
i-D33P150B1L2		150kg	20g		420mm x 550mm	680mm	30684745	500,00
i-D33P150B1X2		150kg	20g		500mm x 650mm	680mm	30684748	620,00
i-D33P300B1X2		300kg	50g		500mm x 650mm	680mm	30684749	620,00
i-D33P300B1V3		300kg	50g		600mm x 800mm	980mm	30684752	840,00
i-D33P600B1V3		600kg	100g		600mm x 800mm	980mm	30684753	840,00
i-D33P15B1R1-M	•	15kg		5g	305mm x 355mm	330mm	30685098	505,00
i-D33P30B1R1-M	•	30kg		10g	305mm x 355mm	330mm	30685099	515,00
i-D33P60B1R1-M	•	60kg		20g	305mm x 355mm	330mm	30685100	515,00
i-D33P60B1L2-M	•	60kg		20g	420mm x 550mm	680mm	30685104	580,00
i-D33P150B1L2-M	•	150kg		50g	420mm x 550mm	680mm	30685105	580,00
i-D33P150B1X2-M	•	150kg		50g	500mm x 650mm	680mm	30685108	700,00
i-D33P300B1X2-M	•	300kg		100g	500mm x 650mm	680mm	30685109	700,00
i-D33P300B1V3-M	•	300kg		100g	600mm x 800mm	980mm	30685112	925,00
i-D33P600B1V3-M	•	600kg		200g	600mm x 800mm	980mm	30685113	925,00
i-D33P15B1R5		15kg	2g		305mm x 355mm	-	30684756	435,00
i-D33P30B1R5		30kg	5g		305mm x 355mm	-	30684757	435,00
i-D33P60B1R5		60kg	10g		305mm x 355mm	-	30684758	435,00
i-D33P60B1L5		60kg	10g		420mm x 550mm	-	30684762	495,00
i-D33P150B1L5		150kg	20g		420mm x 550mm	-	30684763	515,00
i-D33P150B1X5		150kg	20g		500mm x 650mm	-	30684766	690,00
i-D33P300B1X5		300kg	50g		500mm x 650mm	-	30684767	690,00
i-D33P15B1R5-M	•	15kg		5g	305mm x 355mm	-	30685116	505,00
i-D33P30B1R5-M	•	30kg		10g	305mm x 355mm	-	30685117	505,00
i-D33P60B1R5-M	•	60kg		20g	305mm x 355mm	-	30685118	505,00
i-D33P60B1L5-M	•	60kg		20g	420mm x 550mm	-	30685122	575,00
i-D33P150B1L5-M	•	150kg		50g	420mm x 550mm	-	30685123	625,00
i-D33P150B1X5-M	•	150kg		50g	500mm x 650mm	-	30685126	730,00
i-D33P300B1X5-M	•	300kg		100g	500mm x 650mm	-	30685127	730,00

Handle

DEFENDER™3000 Accessories

Accessories	Item N°	Price €
Mounting Set Front SST DT33	30696600	96,00
Column Assembly, 330mm, SST, i-DT33	30696595	66,00
Column Assembly, 680mm, SST, i-DT33	30696596	92,00
Column Assembly, 330mm, CS, i-DT33	30696597	52,00
Column Assembly, 680mm, CS, i-DT33	30696598	75,00
Column Assembly, 980mm, CS, i-DT33	30696599	79,00
USB Device Kit i-D33	30699120	72,00
Dust Cover i-DT33 (5)	30699121	37,00
Rechargeable battery kit i-DT33P	30699122	44,00
Wheel kit, Defender Series	30372146	65,00
Wheel kit SST, Defender Series	30616409	235,00
Cable RS232 i-DT33P	80500525	99,70
Cable RS232 i-DT33XW	80500552	91,80
Ethernet Kit	30429666	111,00
LC Extension Cable, 3.5 m *1	30803566	78,00

1 Cannot be used for legal-for-trade applications



WHEEL KIT, DEFENDER SERIES
30372146



WHEEL KIT SST, DEFENDER SERIES
30616409



MOUNTING SET FRONT SST DT33
30696600



COLUMN ASSEMBLY, 330MM, SST, I-DT33
30696595



COLUMN ASSEMBLY, 680MM, SST, I-DT33
30696596



COLUMN ASSEMBLY, 330MM, CS, I-DT33
30696597



COLUMN ASSEMBLY, 680MM, CS, I-DT33
30696598



COLUMN ASSEMBLY, 980MM, CS, I-DT33
30696599



USB DEVICE KIT I-D33
30699120



IN USE COVER I-DT33 (5)
30699121



RECHARGEABLE BATTERY KITS
30699122



CABLE RS232 I-DT33XW
80500552



ETHERNET KIT
30429666



CABLE RS232 I-DT33P
80500525

DEFENDER™ 2000

The best basic modular scale for dependability and durability at affordable price

Application

Weighing, Parts Counting, Accumulation

Display

Red LED

Operation

AC adapter (included) or internal rechargeable battery

Communication


RS232 interface (accessory)

Construction

ABS indicator, 430 grade stainless steel platform with painted steel tubular frame, IP65 aluminium load cell

Design Features

Hardware configuration lockout switch, battery saving auto-off, externally visible level bubble, down-stops for overload protection



OPTION



The polished 430 stainless steel platform with painted carbon steel frame.

Plastic LED indicator (T24PE).

Internal rechargeable battery.

Model	Capacity	Readability (d)	Platform	Column	Item N°	Price €
D24PE30FR	30kg	5g	300mm x 350mm	500mm	30235922	380,00
D24PE60FR	60kg	10g	300mm x 350mm	500mm	30235927	380,00
D24PE60FL	60kg	10g	400mm x 500mm	700mm	30235932	430,00
D24PE150FL	150kg	20g	400mm x 500mm	700mm	30235937	430,00
D24PE150FX	150kg	20g	500mm x 600mm	700mm	30235942	515,00
D24PE300FX	300kg	50g	500mm x 600mm	700mm	30235947	515,00
D24PE300FV	300kg	50g	600mm x 800mm	1000mm	30235952	895,00
D24PE600FV	600kg	100g	600mm x 800mm	1000mm	30235957	895,00
Accessories					Item N°	Price €
RS232 Kit					30101019	35,00
Folding Column Kit					30101016	41,00
In-Use-Cover, i-DT33P, T24P					30101017	13,00
Load Cell Cable Adapter Kit					30101021	13,40
Wall mount kit, T24P					30387043	27,00
Column adapter Kit, D2K					30101020	41,00

T82XW Indicator

Touchscreen 8000 Series stainless steel indicator with customizable applications

Application

Accumulation and formula weighing, parts counting, weighbridges, sqc, vehicle, price computing and industrial labelling, single-product dosage systems, multi-component dosage systems

Display

Backlit graphic display with touch screen technology

Operation

AC Power

Communication

2xRS232, RS485 direct mother board connection (1xRS232 through RJ11 plug), USB, Ethernet

Construction

IP68 stainless steel housing and mounting bracket

Design Features

Check weighing lights with selectable operation and audible signal settings, multi operating languages, QWERTY, AZERTY, QWERTZ keyboard options integrated in the touch screen, real-time clock, integrated alibi memory, built-in micro SD slot for recording configuration, logos and customized screens, 4 independent analog channels, selectable environmental and auto-print settings, stability indicator, overload/underload indicators, 1 000 000 d and up to 10 000 e



Multiplied communication interfaces.

Unique applications including SQC and Vehicle Weighing.

Ability to support up to 16x350ohm load cells.

Model	Item N°	Price €
T82XW-EU	30325752	1665,00

Terminal is delivered with "accumulation and formula weighing" application software installed. Other software applications are available on <http://dmx.Ohaus.Com> for free.

* For license number, please contact SSC

Accessories	Item N°	Price €
I/O Input, Set(2), T82XW	30325822	40,00
I/O Output, Set(4), T82XW	30325823	74,00
6 In /12 Out I/O Board, T82XW	30344023	260,00
Multi L/C Input, 4 channel, T82XW	30331635	40,00
USB flash driver, T82XW	30331636	292,00
RS232-USB, T82XW	30331637	127,00
Ethernet, T82XW	30331642	253,00
Build-in optoisolated RS485, T82XW	30331643	104,00
Wall mount Kit, T82XW	30331644	111,00
In-Use-Cover Kit, T82XW	30328139	127,00
License for Display tool, T82XW	30344025	272,00
Serial to mini DIN connector, T82XW	30344026	104,00
RS232/C serial cable 1.5m, T82XW	30344027	127,00
Ethernet cable 1.5m, T82XW	30344028	127,00
RS232 Cable for balance, T82XW	30344040	117,00
Analog & 6In/12Out I/O Kit, T82XW	30422883	481,00
Analog output Board, T82XW	30344022	216,00

T72XW Indicator

Advanced performance stainless steel indicator with external device control capabilities

Application

Weighing, parts counting, check weighing, animal/dynamic weighing

Display

Transflective graphic LCD with white backlight, 240x96 pixels

Operation

AC Power (Included)

Communication

Easy access communication port including RS232 (included), Ethernet and 2nd RS232 or RS485 (accessories sold separately), GLP/GMP data output with real-time clock

Construction

IP66/NEMA 4X stainless steel housing and mounting bracket

Design Features

Check weighing lights with selectable operation and audible signal settings, sealed front panel, five operating languages (English, Spanish, French, German and Italian), selectable environmental and auto-print settings, stability indicator, overload/underload indicators, 50 000 d and up to 6 000 e



50,000 displayed resolution and the ability to support up to 10 load cells.

Two available option slots enable a wide range of connectivity interfaces.

Vivid graphic display with on-screen prompts and an alphanumeric keypad with designated keys for common functions such as tare and ID information.

Model	Item N°	Price €
T72XW-EU	30126957	790,00
Accessories	Item N°	Price €
Ethernet Kit, T72XW	30132163	93,00
USB with I/O kit, T72XW	30132164	93,00
COM2 Kit, T72XW	30132165	104,00
Floor scale Column Kit Painted	30531244	127,00
Floor scale Column Kit SST	30539553	235,00
Light Tower, R-Y-G, IP65, 243 mm, 24VDC (Requires 30097591)	30424021	241,00
Column Kit, 350mm, SST, T51/71/72XW	80500725	74,00
Column Kit, 680mm, SST, T51/71/72XW	80500726	84,00



USB WITH I/O KIT, T72XW
30132164



COM2 KIT, T72XW
30132165




LIGHT TOWER, R-Y-G, IP65, 243 MM, 24VDC
30424021



ETHERNET KIT, T72XW
30132163

i-DT61XWE Indicator

Advanced IP68/IP69K Washdown Indicator



Application

Weighing, Counting with APW optimization, Percent Weighing, Checkweighing, Checkcounting, Filling, Dynamic Weighing/Display Hold

Display

6-digit, 7-segment white LED, 20 mm high digits

Operation

100-240 VAC / 50/60 Hz universal power supply, hardwired

Communication

Easy access communication ports including RS232 and USB Host, second RS232/RS485/USB device, Ethernet and WiFi/Bluetooth, Discrete I/O and Analog output (available as accessories); GLP/GMP data output with real-time clock

Construction

IP68/IP69K 316 stainless steel housing indicator and mounting bracket, knife-resistant front panel

Design Features

3-colour check light bar, menu and key lock, selectable environmental and auto-print settings, stability indicator, overload/underload indicators, low battery indicator, auto shut-off, auto tare



Defender 6000 indicator, constructed of polycarbonate or 316 stainless steel, stands up to harsh cleaning in processing and packaging environments with IP68 immersion/IP69K washdown protection.

Large, bright display offers clear results for fast processing. Colour bar offers detailed package weight checking.

With knife-proof numeric touch keypad, i-DT61XWE indicator is suitable for harsh processing environments.

Model	Item N°	Price €
i-DT61XWE-EU	30575576	1020,00
Accessories	Item N°	Price €
Wi-Fi/Bluetooth Dongle i-DT61XW	30412537	282,00
Analog Output kit i-DT61XW *1	30424403	183,00
2nd RS232/RS485/USB Device kit i-DT61XW *1	30424404	156,00
LC Extension Cable, 9 meters.	30424409	111,00
M16 Gland with Ferrite	30379716	20,00
Light Tower, 3 Colors, 243 mm, IP65 (Requires 30097591)	30424021	241,00
Discrete I/O with Relay (2 In/4 Out) i-DT61XW *1	30097591	93,00
Ethernet kit i-DT61XW	30429666	111,00
External Power Bank Holder i-DT61XW (requires external USB power bank)	30572907	290,00
Cable RS232 i-DT61XW	80500552	91,80
Alibi Memory i-DT61XW	80500503	52,00
LC Extension Cable, 3.5 m *2	30803566	78,00

1 Only one of them can be used at the same time
2 Cannot be used for legal-for-trade applications



WI-FI/BLUETOOTH DONGLE I-DT61XW
30412537



ANALOG OUTPUT KIT I-DT61XW
30424403



2ND RS232/RS485/USB DEVICE KIT
I-DT61XW 30424404



EXTERNAL POWER BANK HOLDER I-DT61XW
30572907

i-DT61PW Indicator

Advanced IP68/IP69K Washdown Indicator



EasyConnect

Application

Weighing, Counting with APW optimization, Percent Weighing, Checkweighing, Checkcounting, Dynamic Weighing/Display Hold

Display

6-digit, 7-segment LCD display with white backlight, 45 mm high digits

Operation

6 D (LR20) Batteries (Included)

Communication

Easy access IR communication kit (available as an accessory); GLP/GMP data output with real-time clock

Construction

IP68/IP69K food grade polycarbonate housing and 316 stainless steel mounting bracket

Design Features

3-colour check light bar, menu and key lock, selectable environmental and auto-print settings, stability indicator, overload/underload indicators, low battery indicator, auto shut-off, auto tare



Defender 6000 indicator, constructed of polycarbonate or 316 stainless steel, stands up to harsh cleaning in processing and packaging environments with IP68 immersion/IP69K washdown protection.

Large, bright display offers clear results for fast processing. Colour bar offers detailed package weight checking.

Indicator operates for up to 1,600 hours on one set of batteries for months of use. Battery-powered design eliminates the need for a main connection for increased flexibility.

Model	Item N°	Price €
i-DT61PW	30575574	595,00
Accessories	Item N°	Price €
LC Extension Cable, 9 meters	30424409	111,00
IR communication kit i-DT61PW	30572910	159,00
LC Extension Cable, 3.5 m *1	30803566	78,00

*1 Cannot be used for legal-for-trade applications




LC EXTENSION CABLE, 9 METERS
30424409



IR COMMUNICATION KIT i-DT61PW
30572910

TD52XW Indicator

Multifunctional Indicator for Standard Industrial Applications



Application

Weighing, parts counting, percent weighing, check weighing, dynamic weighing (display hold), Filling

Display

Backlit Liquid Crystal Display (LCD), 256 x 80 pixels

Operation

Internal power supply with power cord (included) or rechargeable battery (sold separately as an accessory)

Communication

Easy access communication ports including RS232 and Micro SD slot (included), second RS232/RS485/USB device, Ethernet and WiFi/Bluetooth, Discrete I/O and analog output (sold separately as accessories); GLP/GMP data output with real-time clock

Construction

IP68 stainless steel housing and mounting bracket

Design Features

Check weighing symbols with selectable operation and audible signal settings, menu lock switch, multiple operating languages, selectable environmental and auto-print settings, stability indicator, overload/underload indicators, selectable adjust contrast settings, auto-dim, auto-off, auto tare, screensaver, 75 000 d and up to 10 000 e



Indicator guides the user with on-screen prompts.

Multiple connectivity options for convenient data transfer and storage.

Multiple application modes and high storage capacity for greater flexibility.

Model	Item N°	Price €
TD52XW-EU	30424066	725,00
Accessories	Item N°	Price €
WiFi/BT Dongle, OHAUS (Requires 30424406)	30412537	282,00
Analog Output Kit, TD52 *3	30424403	183,00
2nd RS232/RS485/USB, TD52 *3	30424404	156,00
Battery Kit, Li-ion, TD52	30424405	127,00
USB Host for Dongle, TD52 *2	30424406	111,00
Light Tower, R-Y-G, IP65, 243 mm, 24VDC (Requires 30097591)	30424021	241,00
In-use Cover Set, TD52XW	30424023	27,00
In-use Cover Set (10), TD52XW	30469949	149,00
Wall Mount Kit, SST, TD52XW	30424026	46,00
LC Extension Cable, 9 meters, D52 *1	30424409	111,00
M16 Gland with Ferrite, to TD52	30379716	20,00
Micro SD Card, 8 GB	30303533	35,00
I/O with Relay (2 In/4 Out) *3	30097591	93,00
Ethernet Kit, TD52 *2	30429666	111,00
Floor scale Column Kit Painted	30531244	127,00
Floor scale Column Kit SST	30539553	235,00
USB Host for Flash Drive *2	30572923	128,00
LC Extension Cable, 3.5 m *1	30803566	78,00


1 Cannot be used for legal-for-trade applications

2 Only one of them can be used at the same time

3 Only one of them can be used at the same time

TD52P Indicator

Multifunctional Indicator for Standard Industrial Applications



Application

Weighing, parts counting, percent weighing, check weighing, dynamic weighing (display hold), Filling

Display

Backlit Liquid Crystal Display (LCD), 256x80 pixels

Operation

Internal power supply with power cord (included) or rechargeable battery (sold separately as an accessory)

Communication

Easy access communication ports including RS232 and Micro SD slot (included), second RS232/RS485/USB device, Ethernet and WiFi/Bluetooth, discrete I/O and analog output (sold separately as accessories); GLP/GMP data output with real-time clock

Construction

ABS plastic housing with painted steel bracket

Design Features

Check weighing symbols with selectable operation and audible signal settings, menu lock switch, multiple operating languages, selectable environmental and auto-print settings, stability indicator, overload/underload indicators, selectable adjust contrast settings, auto-dim, auto-off, auto tare, screensaver, 75 000 d and up to 10 000 e



Indicator guides the user with on-screen prompts.

Multiple connectivity options for convenient data transfer and storage.

Multiple application modes and high storage capacity for greater flexibility.

Model	Item N°	Price €
TD52P	30424065	540,00
Accessories	Item N°	Price €
WiFi/BT Dongle, OHAUS (Requires 30424406)	30412537	282,00
Analog Output Kit, TD52 *3	30424403	183,00
2nd RS232/RS485/USB, TD52 *3	30424404	156,00
Battery Kit, Li-ion, TD52	30424405	127,00
USB Host for Dongle, TD52 *2	30424406	111,00
Light Tower, R-Y-G, IP65, 243 mm, 24VDC (Requires 30097591)	30424021	241,00
In-use Cover Set, TD52P	30424022	20,00
In-use Cover Set (10), TD52P	30469948	84,00
Wall Mount Kit, CS, TD52P	30424027	41,00
LC Extension Cable, 9 meters, D52 *1	30424409	111,00
M16 Gland with Ferrite, to TD52	30379716	20,00
Micro SD Card, 8 GB	30303533	35,00
I/O with Relay (2 In/4 Out) *3	30097591	93,00
Ethernet Kit, TD52 *2	30429666	111,00
Floor scale Column Kit Painted	30531244	127,00
Floor scale Column Kit SST	30539553	235,00
USB Host for Flash Drive *2	30572923	128,00
LC Extension Cable, 3.5 m *1	30803566	78,00

1 Cannot be used for legal-for-trade applications

2 Only one of them can be used at the same time

3 Only one of them can be used at the same time

i-DT33XW Indicator

3000 Series stainless steel indicator for washdown environments

Application

Weighing, Parts Counting, Checkweighing and Totalization

Display

7 Segment LCD with multicolour backlight

Operation

Power cord (included), Rechargeable Battery (included)

Communication


Standard RS232 and optional USB Device or Ethernet

Construction

304 stainless steel housing and mounting bracket

Design Features

Hinged front cover for easy servicing, 30 000 d and up to 10 000 e



Results are clearly shown using the large, bright backlit LCD display. In Checkweighing mode display changes colour (yellow/green/red) to alert when within the target weight range.

Features software for Parts Counting, Checkweighing and Totalization. Includes RS232 with optional USB Device or Ethernet for printing or data connectivity.

Indicator features an IP66 protected 304-stainless steel housing with bracket. Display resolution up to 1:30,000; OIML approved resolution 1:10,000.

Model	Item N°	Price €
i-DT33XW-EU	30648431	350,00
Accessories	Item N°	Price €
Cable RS232 i-DT33XW	80500552	91,80
Ethernet Kit	30429666	111,00
USB Device Kit	30699120	72,00
LC Extension Cable, 3.5 m *1	30803566	78,00

1 Cannot be used for legal-for-trade applications

i-DT33P Indicator

3000 Series basic industrial indicator

Application

Weighing, Parts Counting, Checkweighing and Totalization

Display

7 Segment LCD with multicolour backlight

Operation

Power cord (included), 6C (LR14) Batteries (not included) or Rechargeable Battery (not included)

Communication


Standard RS232 and optional USB Device or Ethernet

Construction

ABS housing and powder-coated steel mounting bracket

Design Features

Integrated key-hole mounts on back cover, 30 000 d and up to 10 000 e



Results are clearly shown using the large, bright backlit LCD display. In Checkweighing mode display changes colour (yellow/green/red) to alert when within the target weight range.

Features software for Parts Counting, Checkweighing and Totalization. Includes RS232 with optional USB Device or Ethernet for printing or data connectivity.

Indicator features a durable ABS housing and steel mounting bracket. Display resolution up to 1:30,000; OIML approved resolution 1:10,000.

Model	Item N°	Price €
i-DT33P	30648430	195,00
Accessories	Item N°	Price €
Rechargeable battery kit i-DT33P	30699122	44,00
Cable RS232 i-DT33P	80500525	99,70
Ethernet Kit	30429666	111,00
USB Device Kit	30699120	72,00
Dust Cover cover i-DT33 (5)	30699121	37,00
LC Extension Cable, 3.5 m *1	30803566	78,00

1 Cannot be used for legal-for-trade applications

T32M Indicator

Sturdy 3000 Series metal indicators

Application

Weighing (5 units), Parts Counting, Dynamic, Checkweighing, Weight Accumulation.

Display

T32MC with backlit LCD / T32ME with LED

Operation

AC adaptor (included)

Communication

RS232 Interface

Construction

ABS front panel with painted carbon steel rear housing

Design Features

Textured-surface paint minimizes scratches, 20 000d and up to 6 000e



Feature-rich, yet simple to use and easy to learn.

Intuitive usage.

Two display versions.

Model	Item N°	Price €
T32ME	83997627	215,00
T32MC	83997628	225,00

Available for all Defender Bench Scales and Terminals

Accessories	Item N°	Price €
RS232-USB (Interface Kit)	30304101	78,25
RS232-Ethernet (Interface Kit)	30304102	116,00
Floor scale Column Kit Painted	30531244	127,00
Floor scale Column Kit SST	30539553	235,00

T24PE Indicator

Dependable and durable 2000 Series indicator at affordable price

Application

Weighing (3 units), parts counting, accumulation

Display

6-digit, 7-segment red LED

Operation

AC adapter with internal rechargeable lead acid battery

Communication

Optional RS232

Construction

ABS plastic housing

Design features

Hardware configuration lockout switch, software included for parts counting mode with selectable sample sizes, selectable span adjustment points, 20 000 d



Easy-to-read red LED display.

Optional RS232C with manual/continuous/interval printing.

Versatile power capabilities including AC adapter and rechargeable battery in standard.

Model	Item N°	Price €
T24PE	30235894	175,00

Accessories	Item N°	Price €
RS232 Kit	30101019	35,00
Folding Column Kit	30101016	41,00
In-Use-Cover, T24P	30101017	13,00
Load Cell Cable Adapter Kit	30101021	13,40

Accessories	Item N°	Price €
Wall mount kit, T24P	30387043	27,00
Floor scale Column Kit Painted 980mm	30531244	127,00
Floor scale Column Kit SST 980mm	30539553	235,00

DEFENDER™ 6000 Base

Advanced IP68/IP69K Washdown Base



Construction

316 stainless steel platform and frame,
hermetically-sealed stainless steel load cell, non-slip
adjustable rubber feet, level indicator



Defender 6000 stands up to harsh cleaning in processing and packaging environments with IP68 immersion/IP69K washdown protection, 316 stainless steel base and hermetically-sealed load cell.

Durable Defender 6000 is NSF certified, supporting HACCP systems with a 316 stainless steel base providing corrosion protection in food, chemical processing, and packaging.

EasyConnect™ allows a mix of displays and bases.

Model	Capacity	Readability (d)	Platform (mm)	Item N°	Price €
i-D3K1S	3kg	0,5g	250 x 250	30554431	1190,00
i-D6K1S	6kg	0,5g	250 x 250	30554432	1190,00
i-D15K1R	15kg	1g	300 x 300	30554433	1390,00
i-D30K1R	30kg	2g	300 x 300	30554434	1390,00
i-D60K1L	60kg	5g	500 x 400	30554435	1690,00
i-D150K1L	150kg	10g	500 x 400	30554436	1690,00
Accessories				Item N°	Price €
Wheel kit SS Defender Series				30616409	235,00
LC Extension Cable, 9 meters				30424409	111,00
Column Kit 260 mm 316 stainless steel				30572908	235,00
Column Kit 670 mm 316 stainless steel				30572909	291,00
Front Mounting kit 316 SS i-DT61PW				30626654	235,00
Pan 305x305 mm in 304 stainless steel				30427798	156,00
Pan 254x254 mm 304 stainless steel				30427799	131,00
Cable Kit Extension for base				30429931	46,00
LC Extension Cable, 3.5 m				30803566	78,00



WHEEL KIT SS DEFENDER SERIES
30616409



LC EXTENSION CABLE, 9 METERS
30424409



COLUMN KIT 316 STAINLESS STEEL
30572908 30572909



CABLE KIT EXTENSION FOR BASE
30429931



FRONT MOUNTING KIT 316 SS i-DT61PW
30626654



PAN, 304 STAINLESS STEEL
30427798 30427799

DEFENDER™5000 Stainless Steel Bases

Durable Bases for the Most Demanding of Industrial Applications

Construction

Dual Range model:

304-grade stainless steel platform with sand gray powder-coated steel frame, IP67 aluminium load cell, non-slip adjustable rubber feet

Single Range model:

304-grade stainless steel platform and frame, IP67 stainless steel load cell, non-slip adjustable rubber feet

EasyConnect



Durably Constructed To Withstand Harsh Industrial Environments.

Stainless steel frame. Stainless Steel load cell for Single Range. Aluminium Load cell for Dual Range.

OIML approved dual range and single range load cell.

Model	Capacity	Platform	Item N°	Price €
D6WQDR	3kg / 6kg	305mm x 305mm	30419594	630,00
D15WQDR	6kg / 15kg	305mm x 305mm	30419595	630,00
D15WQDL	6kg / 15kg	400mm x 400mm	30419609	860,00
D30WQDL	15kg / 30kg	400mm x 400mm	30419610	860,00
D60WQDL	30kg / 60kg	400mm x 400mm	30419611	860,00
D60WQDX	30kg / 60kg	500mm x 500mm	30419614	1115,00
D150WQDX	60kg / 150kg	500mm x 500mm	30419615	1115,00
D60WQDV	30kg / 60kg	610mm x 610mm	30419619	1405,00
D150WQDV	60kg / 150kg	610mm x 610mm	30419620	1405,00
D300WQDV	150kg / 300kg	610mm x 610mm	30419621	1405,00
D150WTDV	60kg / 150kg	800mm x 600mm	30419625	2050,00
D300WTDV	150kg / 300kg	800mm x 600mm	30419626	2050,00
D600WTDV	300kg / 600kg	800mm x 600mm	30419627	2050,00
D15WQR	15kg	305mm x 305mm	30454671	650,00
D15WQL	15kg	400mm x 400mm	30454672	895,00
D30WQL	30kg	400mm x 400mm	30454673	895,00
D60WQL	60kg	400mm x 400mm	30454674	895,00
D60WQX	60kg	500mm x 500mm	30454675	1175,00
D150WQX	150kg	500mm x 500mm	30470877	1175,00
D60WQV	60kg	610mm x 610mm	30471118	1460,00
D150WQV	150kg	610mm x 610mm	30471119	1460,00

DEFENDER™3000 Stainless Steel Bases

Stainless steel platform for washdown environments

Construction

304-grade stainless steel platform, stainless steel tubular frame, IP67 stainless steel load cell



EasyConnect



This tough industrial scale features a 304-stainless steel platform, frame and IP67 load cell.

Capacities from 15 kg to 300 kg, with platforms from 305 x 355 mm to 500 x 650 mm.

Default resolution up to 7,500d; OIML approved to 3,000e for legal-for-trade performance.

Model	Capacity	Platform	Item N°	Price €
i-D15C1REZ	15kg	305mm x 355mm	30683681	405,00
i-D30C1REZ	30kg	305mm x 355mm	30683682	420,00
i-D60C1REZ	60kg	305mm x 355mm	30683683	420,00
i-D60C1LEZ	60kg	420mm x 550mm	30683686	695,00
i-D150C1LEZ	150kg	420mm x 550mm	30683687	695,00
i-D150C1XEZ	150kg	500mm x 650mm	30683690	1045,00
i-D300C1XEZ	300kg	500mm x 650mm	30683691	1045,00

DEFENDER™5000 Bases

Durable Bases for the Most Demanding of Industrial Applications

Construction

304-grade stainless steel platform, sand grey powder-coated steel frame, aluminum load cell, non-slip adjustable rubber feet, level indicator

EasyConnect



Durably Constructed To Withstand Harsh Industrial Environments.

Sand grey powder-coated steel frame.

OIML approvable dual range load cell.

Model	Capacity	Platform	Item N°	Price €
D6RQDR	3kg / 6kg	305mm x 305mm	30419592	425,00
D15RQDR	6kg / 15kg	305mm x 305mm	30419593	425,00
D15RQDL	6kg / 15kg	400mm x 400mm	30419596	565,00
D30RQDL	15kg / 30kg	400mm x 400mm	30419597	565,00
D60RQDL	30kg / 60kg	400mm x 400mm	30419608	565,00
D60RQDX	30kg / 60kg	500mm x 500mm	30419612	705,00
D150RQDX	60kg / 150kg	500mm x 500mm	30419613	705,00
D60RQDV	30kg / 60kg	610mm x 610mm	30419616	1015,00
D150RQDV	60kg / 150kg	610mm x 610mm	30419617	1015,00
D300RQDV	150kg / 300kg	610mm x 610mm	30419618	1015,00
D150RTDV	60kg / 150kg	800mm x 600mm	30419622	1405,00
D300RTDV	150kg / 300kg	800mm x 600mm	30419623	1405,00
D600RTDV	300kg / 600kg	800mm x 600mm	30419624	1405,00

DEFENDER™3000 Bases

The best platform for basic industrial needs

Construction

304-grade stainless steel platform, powder-coated steel tubular frame, aluminum load cell



EasyConnect



This tough industrial scale features a 304-stainless steel platform, powder-coated steel tubular frame and column, IP65 aluminum load cell.


Capacities from 15 kg to 600 kg, with platforms from 305 x 355 mm to 600 x 800 mm.

Default resolution up to 7,500d; OIML approved to 3,000e for legal-for-trade performance.

Model	Capacity	Platform	Item N°	Price €
i-D15B1R EZ	15kg	305mm x 355mm	30648440	220,00
i-D30B1R EZ	30kg	305mm x 355mm	30683679	225,00
i-D60B1R EZ	60kg	305mm x 355mm	30683680	225,00
i-D60B1L EZ	60kg	420mm x 550mm	30683684	300,00
i-D150B1L EZ	150kg	420mm x 550mm	30683685	300,00
i-D150B1X EZ	150kg	500mm x 650mm	30683688	410,00
i-D300B1X EZ	300kg	500mm x 650mm	30683689	410,00
i-D300B1V EZ	300kg	600mm x 800mm	30683692	655,00
i-D600B1V EZ	600kg	600mm x 800mm	30683693	655,00

RANGER™7000

Simplifying even the most complex industrial applications



Application

Weighing, parts counting, percent weighing, check weighing, dynamic weighing (display hold), filling, formulation, differential weighing, density weighing

Display

Full-colour WQVGA graphic display (4.3 in./10.9 cm diagonal); 4.3" TFT Graphic LCD

Operation

AC Power adapter (included) or rechargeable battery (available as an accessory); 100 240 VAC / 50/60 Hz universal power supply

Communication

Easy access communication ports including standard USB and RS232 and an optional third port of either RS232 or Ethernet, GLP/GMP Data Output with Real Time Clock, Data Transfer Function.

Construction

Stainless steel platform, die cast housing, column (available as an accessory)

Design Features

Removable stainless steel platform, check weighing lights with selectable operation and audible signal settings, integral weigh below hook, sealed front panel, menu lock switch, level indicator, adjustable leveling feet, several selectable operating languages, selectable environmental and auto-print settings, stability indicator, overload/underload indicators, selectable brightness settings, auto-dim, auto-sleep, auto tare, chain tare



Simplify complex applications and minimize the need for manual calculations with 10 advanced application modes, peripheral device control and an optional kit for a second scale platform.

The most precise scale in its class with one-second stabilization time, up to 75,000d or 350,000d display resolution and legal-for-trade certifications.

Extremely durable scale designed to thrive in rugged industrial environments. Features IP54-rated metal housing and sealed metal terminal.

Model	EC Type Approved	Internal Adjustment	Capacity	Readability (d)	(e)	Platform	Cal. weight, class	Item N°	Price €
R71MD3-EU			3kg	0.05g		280mm x 280mm	3kg, OIML F2 or ASTM CL 4	30088844	1090,00
R71MD6-EU			6kg	0.1g		280mm x 280mm	6kg, OIML F2 or ASTM CL 4	30088845	1090,00
R71MD15-EU			15kg	0.2g		377mm x 311mm	15kg, OIML F2 or ASTM CL 4	30088846	1215,00
R71MD35-EU			35kg	0.5g		377mm x 311mm	35kg, OIML F2 or ASTM CL 4	30088847	1215,00
R71MD60-EU			60kg	1g		377mm x 311mm	60kg, OIML F2 or ASTM CL4	30212874	1265,00
R71MHD3-EU		•	3kg	0.01g		210mm x 210mm	3kg, OIML F2 or ASTM CL 4	30088832	3080,00
R71MHD6-EU		•	6kg	0.02g		210mm x 210mm	6kg, OIML F2 or ASTM CL 4	30088833	3080,00
R71MHD15-EU		•	15kg	0.1g		377mm x 311mm	15kg, OIML F2 or ASTM CL 4	30088834	3585,00
R71MHD35-EU		•	35kg	0.1g		377mm x 311mm	35kg, OIML F2 or ASTM CL 4	30088835	3835,00
R71MHD3-EU-M	•	•	3kg		0.1g	210mm x 210mm	3kg, OIML F2 or ASTM CL 4	30104929	3145,00
R71MHD6-EU-M	•	•	6kg		0.2g	210mm x 210mm	6kg, OIML F2 or ASTM CL 4	30104940	3145,00
R71MHD15-EU-M	•	•	15kg		1g	377mm x 311mm	15kg, OIML F2 or ASTM CL 4	30104941	3655,00
R71MHD35-EU-M	•	•	35kg		1g	377mm x 311mm	35kg, OIML F2 or ASTM CL 4	30104942	3900,00
R71MD3-EU-M	•		3kg		0.5g	280mm x 280mm	3kg, OIML F2 or ASTM CL 4	30104921	1150,00
R71MD6-EU-M	•		6kg		1g	280mm x 280mm	6kg, OIML F2 or ASTM CL 4	30104922	1150,00
R71MD15-EU-M	•		15kg		2g	377mm x 311mm	15kg, OIML F2 or ASTM CL 4	30104923	1275,00
R71MD35-EU-M	•		35kg		5g	377mm x 311mm	35kg, OIML F2 or ASTM CL 4	30104924	1275,00
R71MD60-EU-M	•		60kg		10g	377mm x 311mm	60kg, OIML F2 or ASTM CL4	30238583	1335,00

RANGER™ 4000

The best value for durable industrial weighing



Application

Weighing, parts counting, check weighing, percent weighing, animal/dynamic weighing, display hold, accumulation

Display

Light-emitting diode display (LED), 3 colour check weighing LEDs

Operation

AC power (included) or rechargeable battery (included)

Communication

Easy access communication port including standard RS232 (included) and second RS232, USB or Ethernet (accessories sold separately)

Construction

Rugged cast aluminum housing, stainless steel platform

Design Features

Cast aluminum housing, stainless steel platform, integral weigh-below hook, sealed front panel, menu lock switch, up front level indicator, adjustable leveling feet, selectable environmental and auto-print settings, stability indicator, overload and underload indicators, low battery indicator, auto shut-off, auto tare



Sturdy metal housing and slip-resistant rubber feet provide the protection, stability and long product life needed for tough industrial weighing.

The most user-friendly scale on the market features Smart Text™ for easy operation and setup to ensure efficiency in the workplace.

With the largest display in its class, checkweighing LED indicators and backlit LCD display, the Ranger 4000 can function in almost any industrial environment.

Model	EC Type Approved	Capacity	Readability (d)	(e)	Platform	Item N°	Price €
R41ME3		3kg	0,1g		300mm x 225mm	30236943	600,00
R41ME6		6kg	0,2g		300mm x 225mm	30231099	600,00
R41ME15		15kg	0,5g		300mm x 225mm	30231100	600,00
R41ME30		30kg	1g		300mm x 225mm	30231101	600,00
R41ME3-M	•	3kg		1g	300mm x 225mm	30237016	665,00
R41ME6-M	•	6kg		2g	300mm x 225mm	30237017	665,00
R41ME15-M	•	15kg		5g	300mm x 225mm	30237018	665,00
R41ME30-M	•	30kg		10g	300mm x 225mm	30237019	665,00

RANGER™ 4000 COUNT

The Best Value for Durable Industrial Counting



Application

Weighing, parts counting with auto optimization, check weighing, accumulation

Display

3 window backlit liquid crystal display (LCD), 3 colour check weighing LEDs

Operation

AC power (included) or rechargeable battery (included)

Communication

Easy access communication port including standard RS232 (included) and second RS232, USB or Ethernet (accessories sold separately)

Construction

Rugged cast aluminum housing, stainless steel platform

Design Features

Cast aluminum housing, stainless steel platform, sealed front panel, menu lock switch, up front level indicator, adjustable leveling feet, selectable environmental and auto-print settings, stability indicator, overload and underload indicators, low battery indicator, auto shut-off, auto tare



Sturdy metal housing and slip-resistant rubber feet provide the protection, stability and long product life needed for tough industrial weighing.


Accurate weighing and counting performance is assured with 1:1,500,000 internal counting resolution and auto-optimization software that reduces errors.

Easy-to-operate with SmarText™ software, a centrally located count display and dedicated Sample and Tare keys.

Model	EC Type Approved	Capacity	Readability (d)	(e)	Platform	Item N°	Price €
RC41M3		3kg	0,1g		300mm x 225mm	30231102	625,00
RC41M6		6kg	0,2g		300mm x 225mm	30231103	625,00
RC41M15		15kg	0,5g		300mm x 225mm	30231104	625,00
RC41M30		30kg	1g		300mm x 225mm	30231105	625,00
RC41M3-M	•	3kg		1g	300mm x 225mm	30237020	690,00
RC41M6-M	•	6kg		2g	300mm x 225mm	30237021	690,00
RC41M15-M	•	15kg		5g	300mm x 225mm	30237022	690,00
RC41M30-M	•	30kg		10g	300mm x 225mm	30237023	690,00

RANGER™3000

The Best Compact Scale in its Class!



Application

Weighing, parts counting, check weighing, percentage weighing, dynamic weighing, display hold, accumulation

Display

Backlit liquid crystal display (LCD), 3 color checkweighing LEDs

Operation

AC power (included) or rechargeable battery (included)

Communication

Easy access communication ports including a standard RS232 and an optional second port of either RS232, USB or Ethernet

Construction

ABS housing, stainless steel platform

Design Features

210 hours of battery life, 3 LED (yellow, green, red) check weighing lights with selectable operation and audible signal settings, integral weigh-below hook, transportation handle, sealed front panel, menu lockswitch, up front level indicator, adjustable leveling feet, selectable environmental and auto-print settings, stability indicator, overload and underload indicators, low battery indicator, auto shut-off, auto tare



With the largest display in its class, checkweighing LED indicators and backlit LCD display, the Ranger 3000 can function in any industrial environment.

User-friendly scale with easy setup to ensure efficiency in the workplace.

Eliminate waiting time and increase control of your operation with accurate one-second results.

Model	EC Type Approved	Capacity	Readability (d)	(e)	Platform	Item N°	Price €
R31P1502		1.5kg	0.05g		225mm x 300mm	30031702	410,00
R31P3		3kg	0.1g		225mm x 300mm	30031703	410,00
R31P6		6kg	0.2g		225mm x 300mm	30031704	410,00
R31P15		15kg	0.5g		225mm x 300mm	30031705	410,00
R31P30		30kg	1g		225mm x 300mm	30031706	410,00
R31P1502-M	•	1.5kg		0.5g	225mm x 300mm	30060902	485,00
R31P3-M	•	3kg		1g	225mm x 300mm	30060903	485,00
R31P6-M	•	6kg		2g	225mm x 300mm	30060904	485,00
R31P15-M	•	15kg		5g	225mm x 300mm	30060905	485,00
R31P30-M	•	30kg		10g	225mm x 300mm	30060906	485,00

RANGER™ 3000 COUNT

The Best Compact Counting Scale in its Class!



Application

Weighing, parts counting with auto optimization, check weighing, accumulation

Display

3 window backlit liquid crystal display (LCD), 3 color checkweighing LEDs

Operation

AC power (included) or rechargeable battery (included)

Communication

Easy access communication ports including a standard RS232 and an optional second port of either RS232, USB or Ethernet

Construction

ABS housing, stainless steel platform

Design Features

30 item library (stores 12-digit part numbers, average piece weights, check counting or check weighing parameters, tare weights and accumulated weights or counts), 210 hours of battery life, 3 LED (yellow, green, red) check weighing lights with selectable operation and audible signal settings, auto tare, transportation handle, numeric keypad, sealed front panel, menu lockswitch, up front level indicator, adjustable leveling feet, selectable environmental and auto-print settings, stability indicator, overload and underload indicators, low battery indicator, auto shut-off, optional 2nd scale input



Accurate weighing and counting are assured with a high-resolution display, 30-item library function, and auto-optimization software that reduces errors.

Easy to use in multi-step counting operations. Features three backlit LCD displays, including a centrally located count display, and an easy-to-navigate menu.

Eliminate waiting time and increase control of your operation with accurate one-second results.

Model	EC Type Approved	Capacity	Readability (d)	(e)	Platform	Item N°	Price €
RC31P1502		1.5kg	0.05g		225mm x 300mm	30031782	430,00
RC31P3		3kg	0.1g		225mm x 300mm	30031783	430,00
RC31P6		6kg	0.2g		225mm x 300mm	30031784	430,00
RC31P15		15kg	0.5g		225mm x 300mm	30031785	430,00
RC31P30		30kg	1g		225mm x 300mm	30031786	430,00
RC31P1502-M	•	1.5kg		0.5g	225mm x 300mm	30060907	500,00
RC31P3-M	•	3kg		1g	225mm x 300mm	30060908	500,00
RC31P6-M	•	6kg		2g	225mm x 300mm	30060909	500,00
RC31P15-M	•	15kg		5g	225mm x 300mm	30060910	500,00
RC31P30-M	•	30kg		10g	225mm x 300mm	30060911	500,00

RANGER™ Accessories

Accessories	Ranger™7000	Ranger™4000	Ranger™Count 4000	Ranger™3000	Ranger™Count 3000	Item N°	Price €
PCBA Kit, Alibi Memory	•					80500503	52,00
Ethernet Kit	•	•	•	•	•	30037447	84,00
Second RS232 Kit	•	•	•	•	•	30037448	46,00
2nd platform kit, R71	•					30097590	65,00
In-use cover, R71	•					30135320	35,00
Discrete I/O Kit, 2-In 4-Out, R71	•					30097591	93,00
Rechargeable Battery, EX-HiCap R71	•					30041295	222,00
Column Kit, Painted, 400mm, R71	•					30095408	84,00
Cable, Extension, 9 m, R71	•					30101495	66,00
Light Tower, R-Y-G, IP65, 243 mm, 24VDC (Requires 30097591)	•					30424021	241,00
Cable, RS232, 1.5m, to connect a reference balance *1	•					30057595	78,40
2nd Platform, Kit		•	•		•	30037446	46,00
USB Kit		•	•	•	•	30037449	59,00
In-Use Cover 5 pcs		•	•			30656205	77,00
In-Use Cover 5 pcs				•	•	30656204	77,00
Weigh below hook					•	30037450	7,00
Carrying Case		•	•	•	•	30795772	171,00

*1 Reference Balance automatically accepts higher-precision APW data from a scale through the included RRS232 port



PCBA KIT, ALIBI MEMORY
80500503



ETHERNET KIT
30037447



SECOND RS232 KIT
30037448



2ND PLATFORM KIT, R71
30097590



IN-USE COVER, R71
30135320



DISCRETE I/O KIT, 2-IN 4-OUT, R71
30097591



RECHARGEABLE BATTERY, EX-HICAP R71
30041295



COLUMN KIT, PAINTED, 400MM, R71
30095408



WEIGH BELOW HOOK
30037450



LIGHT TOWER
30424021



2ND PLATFORM, KIT
30037446



USB KIT
30037449



IIN-USE COVER 5 PCS
30656205



IIN-USE COVER 5 PCS
30656204

VALOR™7000

The Best Compact Food Scale in its Class!



Application

Weighing, check weighing, percentage weighing, accumulation

Display

Front and rear liquid crystal display (lcd) with backlight

Operation

AC power (included) or rechargeable battery (included)

Communication

Easy access communication ports including a standard RS232 and an optional second port of either RS232, USB or Ethernet

Construction

ABS housing, stainless steel platform

Design Features

Rechargeable battery (200 hours of battery life), 3 LED (yellow, green, red) check weighing lights with selectable operation and audible signal settings, touchless sensor, transportation handle, sealed front panel, menu lockswitch, up front level indicator, adjustable levelling feet and bubble



Food-safe Valor 7000 is NSF-certified, USDA-AMS accepted, and supports HACCP-certified systems -- suitable for food applications in dry areas.


Large displays with backlight allow for easy viewing of results. Easy access communication ports for fast and useful data transfer and storage.

Valor 7000 produces repeatable and accurate results within less than one second to virtually eliminate your wait time and improve productivity.

Model	EC Type Approved	Capacity	Readability (d)	(e)	Platform	Item N°	Price €
V71P1502T		1.5kg	0.05g		225mm x 300mm	30031822	480,00
V71P3T		3kg	0.1g		225mm x 300mm	30031823	480,00
V71P6T		6kg	0.2g		225mm x 300mm	30031824	480,00
V71P15T		15kg	0.5g		225mm x 300mm	30031825	480,00
V71P30T		30kg	1g		225mm x 300mm	30031826	480,00
V71P1502T-M	•	1.5kg		0.5g	225mm x 300mm	30060912	550,00
V71P3T-M	•	3kg		1g	225mm x 300mm	30060913	550,00
V71P6T-M	•	6kg		2g	225mm x 300mm	30060914	550,00
V71P15T-M	•	15kg		5g	225mm x 300mm	30060915	550,00
V71P30T-M	•	30kg		10g	225mm x 300mm	30060916	550,00

VALOR™ 4000

Water proof IP68 table top scale with OIML approval and touchless sensor designed for wet food processing environment



Application

Weighing, percentage weighing, check weighing, accumulation

Display

Large dual red led display: 0.8 inch/20.5 mm height characters

Operation

AC adapter (included) or rechargeable battery (included)

Construction

V41PW: ABS housing with stainless steel platform; V41XW: stainless steel platform and housing

Design Features

Less than 0.5s stabilization time, easy-to-transport, power jack auto-closed design, level bubble, fixed weighing pan, overload and underload indicators, low battery indicator, AC adapter, rechargeable battery, touchless sensor for tare function.



Food-safe Valor 4000 is NSF-certified, USDA-AMS accepted, and supports HACCP-certified systems, making it suitable demanding food processing applications.

Increase productivity with its touchless Tare sensor, half-second stabilization time, checkweighing mode and LEDs, and percent weighing mode.

Unique knife-resistant keypad and IP68 protection ensures durability and safe operation in wet conditions.

Models with plastic housing

Models	EC Type Approved	Capacity	Readability (d)	(e)	Platform	Item N°	Price €
V41PWE1501T		1.5kg	0.2g		242mm x 190mm	30072307	390,00
V41PWE3T		3kg	0.5g		242mm x 190mm	30072308	390,00
V41PWE6T		6kg	1g		242mm x 190mm	30072309	390,00
V41PWE15T		15kg	2g		242mm x 190mm	30072310	390,00
V41PWE1501T-M	•	1.5kg		0.5g	242mm x 190mm	30035653	450,00
V41PWE3T-M	•	3kg		1g	242mm x 190mm	30035654	450,00
V41PWE6T-M	•	6kg		2g	242mm x 190mm	30035655	450,00
V41PWE15T-M	•	15kg		5g	242mm x 190mm	30035656	450,00

Models with stainless steel housing

Models	EC Type Approved	Capacity	Readability (d)	(e)	Platform	Item N°	Price €
V41XWE1501T		1.5kg	0.2g		242mm x 190mm	30072311	465,00
V41XWE3T		3kg	0.5g		242mm x 190mm	30072312	465,00
V41XWE6T		6kg	1g		242mm x 190mm	30072313	465,00
V41XWE15T		15kg	2g		242mm x 190mm	30072314	465,00
V41XWE1501T-M	•	1.5kg		0.5g	242mm x 190mm	30035663	535,00
V41XWE3T-M	•	3kg		1g	242mm x 190mm	30035664	535,00
V41XWE6T-M	•	6kg		2g	242mm x 190mm	30035665	535,00
V41XWE15T-M	•	15kg		5g	242mm x 190mm	30035666	535,00

Approved models need reverification once 300 x 225 mm stainless steel platforms is installed.

VALOR™3000

Full stainless steel portable scale



Application

Weighing (5 units), percentage weighing, specific gravity

Display

LCD with white LED backlight

Operation

AC adapter (included) or 4 C batteries (not included)

Construction

Stainless steel housing and stainless steel pan

Design Features

Adjustable levelling feet and bubble (only for XW models), load cell lock, user adjustable filtering and optional food accessories



High resolution models.

IP65 models.

Food safety certification (NSF) and supports an HACCP compliant quality system.

Model	IP Dust and Water Protected	Capacity	Readability (d)	Platform	Item N°	Price €
V31XH202		200g	0.01g	Ø 120mm	71210129	485,00
V31XH402		400g	0.01g	Ø 120mm	72240776	630,00
V31X501		500g	0.1g	Ø 120mm	72240398	320,00
V31XH2		2kg	0.1g	158mm x 146mm	71210130	485,00
V31X3		3kg	1g	158mm x 146mm	71210131	310,00
V31XH4		4kg	0.1g	158mm x 146mm	72240399	630,00
V31X6		6kg	1g	158mm x 146mm	71210132	345,00
V31XW301	•	300g	0.1g	158mm x 146mm	71210133	375,00
V31XW3	•	3kg	0.5g	158mm x 146mm	71210134	375,00
V31XW6	•	6kg	1g	158mm x 146mm	71210135	390,00

VALOR™2000

A rapid-response food scale that enhances productivity in harsh environments!

Application

V22PW: weighing; V22XW: weighing, check weighing

Display

Large dual red LED display: 20.5 mm height characters

Operation

AC adapter (included) or rechargeable battery (included)

Construction

V22PW: ABS housing with stainless steel platform; V22XW: stainless steel platform and housing

Design Features

Less than 0.5s stabilization time, easy-to-transport, power jack auto-closed design, level bubble, fixed weighing pan, overload and underload indicators, low battery indicator, AC adapter, rechargeable battery



The sturdy polycarbonate keypad offers protection from knives and other sharp instruments used in food preparation, processing and portioning applications.

A unique IPX8 flow-thru design ensures that the NSF-certified Valor 2000 will drain moisture in the scale when used in wet conditions.

Stabilization time is less than a half-second. Large front and rear LED displays allow for easy viewing so two operators can share one scale.

Models with plastic housing


Models	Capacity	Readability (d)	Platform	Item N°	Price €
V22PWE1501T	1.5kg	0.2g	242mm x 190mm	30072346	280,00
V22PWE3T	3kg	0.5g	242mm x 190mm	30072347	280,00
V22PWE6T	6kg	1g	242mm x 190mm	30072315	280,00
V22PWE15T	15kg	2g	242mm x 190mm	30072316	280,00
V22PWE30T	30kg	5g	242mm x 190mm	30236101	280,00

Models with stainless steel housing

Models	Capacity	Readability (d)	Platform	Item N°	Price €
V22XWE1501T	1.5kg	0.2g	242mm x 190mm	30072348	335,00
V22XWE3T	3kg	0.5g	242mm x 190mm	30072349	335,00
V22XWE6T	6kg	1g	242mm x 190mm	30072317	335,00
V22XWE15T	15kg	2g	242mm x 190mm	30072318	335,00
V22XWE30T	30kg	5g	242mm x 190mm	30236102	335,00

VALOR™ I000

Reliable Food-Safe Scale for Dry Environments



Application

Weighing

Display

LCD with white backlight (dual display rear/front only on approved models)

Operation

D-cell alkaline battery (included) or power adapter (optional), rechargeable battery (available as an accessory)

Construction

ABS housing, stainless steel pan, carry handles

Design Features

Integrated level bubble and adjustable feet, selectable weighing units



APPROVED MODELS REAR VIEW



Easy-Clean, Hygienic Design.

Low Power Consumption Design Promotes Continuous Use.

Easy Two-Button Operation.

Models	EC Type Approved	Capacity	Readability (d)	(e)	Platform	Item N°	Price €
V12P3		3kg	0.5g		245mm x 190mm	30539394	150,00
V12P6		6kg	1g		245mm x 190mm	30539395	150,00
V12P15		15kg	2g		245mm x 190mm	30539396	150,00
V12P30		30kg	5g		245mm x 190mm	30539397	150,00
V12P3T-M	•	3kg		1g	245mm x 190mm	30554438	245,00
V12P6T-M	•	6kg		2g	245mm x 190mm	30554439	245,00
V12P15T-M	•	15kg		5g	245mm x 190mm	30554440	245,00
V12P30T-M	•	30kg		10g	245mm x 190mm	30554441	245,00

VALOR Accessories

Accessories	Valor™7000	Valor™4000	Valor™3000	Valor™2000	Valor™ 1000	Item N°	Price €
In-Use Cover for V71	•					30037445	20,00
Ethernet Kit	•					30037447	84,00
2nd RS232 Kit	•					30037448	46,00
USB Kit	•					30037449	59,00
Weigh below hook	•					30037450	7,00
In-Use Cover for V71, 10pcs	•					30037469	144,00
Stainless Steel Platform, 300 x 225 mm, Valor 4000, Valor 2000		•		•		30251753	35,00
Stainless steel sandwich pan Valor 3000			•			80251147	46,00
S/Steel French fries holder Valor 3000			•			80251148	41,00
S/Steel ice cream cone holder Valor 3000			•			80251149	52,00
Carrying Case, Hard Shell Valor 3000			•			80251216	93,00
Stainless Steel platform, 300 x 225 mm, EU-UK Power Plug					•	30559107	35,00
Universal adapter for Valor 1000 (Right-Angled), Plug set (46001779) NOT included					•	46001779	13,50
Rechargeable Battery Kit V12P					•	30559104	46,00
In-use Cover (1)					•	30559105	13,00
In-use Cover (10)					•	30559106	99,00
RS232 Interface, 1.5m cable (for approved models)					•	30559103	35,00
Mildew Proof Feet for V22, V41		•		•		30572905	21,00



IN-USE COVER FOR V71
30037445



ETHERNET KIT
30037447



2ND RS232 KIT
30037448



USB KIT
30037449



WEIGH BELOW HOOK
30037450



IN-USE COVER FOR V71, 10PCS
30037469



MILDEW PROOF FEET FOR V22, V41
30572905



IN-USE COVER (10)
30559106



STAINLESS STEEL PLATFORM
30251753



STAINLESS STEEL SANDWICH PAN
80251147



S/STEEL FRENCH FRIES HOLDER
80251148



S/STEEL ICE CREAM CONE HOLDER
80251149



CARRYING CASE, HARD SHELL
80251216




STAINLESS STEEL PLATFORM
30559107



IN-USE COVER (1)
30559105

COURIER™7000

Precise Standard Scales for Dedicated Shipping Applications



Application

Weighing, Dynamic Weighing/Display Hold, Weight Alert, Totalization

Display

6-digit, 7-segment LCD display with white backlight, 28 mm high digits

Operation

Standard AC adapter with USB type C cable detachable (included) or 4 AA batteries (not included)

Communication

Easy access communication ports including standard RS232 and USB type C. Ethernet or 2nd RS232 or WiFi/Bluetooth (available as accessories sold separately)

Construction

ABS indicator, 304-stainless steel platform, powder-coated frame, aluminum load cell, non-slip adjustable rubber feet

Design Features

150% overload protection



Designed to withstand constant use from offices to distribution centers. Fast setup with connection to shipping systems like UPS & FedEx. "GeoCode" application for local calibration adjustment.

Display Hold retains weight data and presents it on a large, easily accessible backlit display. Its Weight Alert mode alerts users to packages that may be too heavy in relation to defined values.

Courier 7000 can run on USB type-C or disposable battery power. Optional Ethernet, second RS232, Bluetooth/WiFi, remote display, ball top and roller top for material handling efficiency are available.

Models	EC Type Approved	Capacity	Readability (d)	(e)	Platform	Item N°	Price €
i-C71M15R EU		15kg	2g		305mm x 355mm	30745957	375,00
i-C71M60R EU		60kg	10g		305mm x 355mm	30745958	375,00
i-C71M60L EU		60kg	10g		457mm x 457mm	30745959	630,00
i-C71M150L EU		150kg	20g		457mm x 457mm	30745960	630,00
i-C71M150X EU		150kg	20g		610mm x 610mm	30745961	935,00
i-C71M15R-M	•	15kg		5g	305mm x 355mm	30745962	450,00
i-C71M60R-M	•	60kg		20g	305mm x 355mm	30745963	450,00
i-C71M60L-M	•	60kg		20g	457mm x 457mm	30745964	710,00
i-C71M150L-M	•	150kg		50g	457mm x 457mm	30745965	710,00
i-C71M150X-M	•	150kg		50g	610mm x 610mm	30745966	1005,00

COURIER™ 5000

Low-Profile Scales Ideal for General Shipping and Animal Weighing Applications



Application

Weighing, Dynamic Weighing/Display Hold, Weight Alert, Totalization

Display

5-digit, 7-segment LCD display with white backlight, 28 mm high digits

Operation

Standard AC adapter or USB type C cable (included) or 3 AA Batteries (not Included)

Communication

Standard USB type C

Construction

ABS indicator, 304 stainless steel platform, powder coated frame

Design Features

150% overload protection



Easy-to-use scales for general shipping use with painted steel base. Stainless steel pan holds up to daily use and is easy to keep sanitary. Adjustable display for mounting to adapt to your workspace.

The Dynamic Weighing/Display Hold application is ideal for large boxes and packages as it retains weight data on the display - also ideal for animals in veterinary use.

AC adapter power included. Simple USB type C for data and power connection – includes UPS software protocol string. Can also run using disposable battery power.

Models	Capacity	Readability (d)	Platform	Item N°	Price €
i-C52M6R EU	6kg	2g	329mm x 319mm	30745890	285,00
i-C52M30R EU	30kg	10g	329mm x 319mm	30745891	285,00
i-C52M50R EU	50kg	20g	329mm x 319mm	30745892	285,00
i-C52M50L EU	50kg	20g	520mm x 400mm	30745893	435,00
i-C52M100L EU	100kg	50g	520mm x 400mm	30745894	435,00
i-C52M200L EU	200kg	100g	520mm x 400mm	30745895	435,00
i-C52M200X EU	200kg	100g	900mm x 600mm	30745896	510,00

COURIER™3000

User-Friendly Basic Low-Profile Scales for General Shipping and Veterinary Applications

Application

Weighing, Display Hold

Display

4-digit, 7-segment LCD display with white backlight, 28 mm high digits

Operation

AC Adapter (included) or 3 AA Batteries (not included)

Construction

ABS indicator, powder-coated steel platform and frame

Design Features

150% overload protection



Ideal for general shipping jobs or veterinary applications. Low-profile design with 1.5 mm thick painted steel platform to last rigorous daily use. Coiled cable ensures flexible display mounting.

The Courier 3000's Display Hold mode (for use with large boxes and packages) retains weight on the large 28 mm high backlit display, until cleared.

AC adapter power included, and disposable battery power (not included) for portable use.

Model	Capacity	Readability (d)	Platform	Item N°	Price €
i-C31M35R EU	35kg	20g	329mm x 319mm	30809911	215,00
i-C31M75R EU	75kg	50g	329mm x 319mm	30809912	215,00
i-C31M200R EU	200kg	100g	329mm x 319mm	30809913	215,00
i-C31M75L EU	75kg	50g	520mm x 400mm	30809914	245,00
i-C31M200L EU	200kg	100g	520mm x 400mm	30809915	245,00

COURIER™ I000

Compact Portable Shipping Scales for General Shipping Applications



Application

Weighing, Dynamic Weighing

Display

4-digit, 7-Segment LCD, 20 mm high digits

Operation

AC adapter (included) or 3 C (LR14) batteries (not included)

Construction

ABS housing and platform, stainless steel platform (accessory sold separately), transportation handles

Design Features

150% overload protection



Its ultra-compact ABS plastic design combines light weight and durability. Courier scales are exceptionally easy to transport thanks to their integrated handles.

Large 20 mm high display works with Dynamic Weighing mode to store data & ensure it is easily accessible when weighing packages and is also ideal for weighing small animals in veterinary use.

Battery powered (not included) for use in the office or shop. AC adapter power included.

Model	Capacity	Readability (d)	Platform	Item N°	Price €
i-C12P9 EU	9kg	5g	280mm x 316mm	30820352	135,00
i-C12P20 EU	20kg	10g	280mm x 316mm	30820353	135,00
i-C12P75 EU	75kg	50g	280mm x 316mm	30820354	135,00

COURIER™ Accessories

Accessories	Courier™7000	Courier™5000	Courier™3000	Courier™1000	Item N°	Price €
2 nd RS232 Kit i-C71	•				30745880	51,00
Ethernet Kit i-C71	•				30745882	112,00
Wi-Fi / Bluetooth Kit i-C71	•				30745881	326,00
Remote Display AD7-MD	•				30472063	96,00
Ball Top i-C71 305x355 mm *	•				30819625	195,00
Ball Top i-C71 458x458 mm *	•				30819627	256,00
Ball Top i-C71 611x611 mm *	•				30819648	368,00
Roller Top i-C71 305x355 mm *	•				30819626	288,00
Roller Top i-C71 458x458 mm *	•				30819649	332,00
Roller Top i-C71 611x611 mm *	•				30819650	491,00
Remote Column Kit 330 mm Painted Steel	•	•			30800931	82,00
Remote Column Kit 980 mm Painted Steel	•	•			30795819	145,00
Cable RS232 IBM 9P Male-to-Female	•				80500525	99,70
SF40A Printer	•				30064202	715,00
SF-110B Printer	•				30708371	805,00
Spiral cable connection base to indicator 9m max i-C51 i-C31		•	•		30809920	11,00
Rubber Mat Kit Antiskid for X base		•			30548975	74,00
Stainless Steel Weighing Pan				•	80251248	40,00

* Approved models must be recalibrated once ball top or roller top installed



2 ND RS232 KIT I-C71
30745880



ETHERNET KIT I-C71
30745882



WI-FI / BLUETOOTH KIT I-C71
30745881



REMOTE DISPLAY AD7-MD
30472063



SF40A PRINTER
30064202



SF-110B PRINTER
30708371



CABLE RS232 IBM 9P MALE-TO-FEMALE
80500525



REMOTE COLUMN KIT PAINTED STEEL
30800931



BALL TOP I-C71
30819625



ROLLER TOP I-C71
30819626



SPIRAL CABLE, I-C52 I-C31
30809920



STAINLESS STEEL WEIGHING PAN
80251248

DEFENDER™ 5000 Lift Platforms

Stainless Steel Lift Floor Platform for Hygienic Applications



Hygienic design allows for quick and easy washdown cleaning in sanitary applications. Lifting the front handles tilts the entire scale up providing complete access to the floor and underside of the scale.

Platforms are equipped with a 5 m double-shielded EasyConnect™ cable, which allows instant configuration with Defender 3000 (i-DT33), Defender 5000 (TD52) or Defender 6000 (i-DT61) indicators for quick installation.

Platform is available in 2 different models with 1000 x 1000 mm and 1250 x 1250 mm sizes and capacities of 600 to 1500 kg, allowing greater flexibility to best suit your weighing needs.

Model	Capacity	Platform (mm)	Item N°	Price €
DFL600WR EU	600kg	1000 x 1000 x 45	30540665	6710,00
DFL1500WL EU	1500kg	1250 x 1250 x 45	30540666	7810,00

DEFENDER™ 5000 Low Profile Floor Platforms

Low Profile Floor Platform for Drive-thru Applications



Ultra-low profile painted steel or stainless steel solution for weighing drums and carts by loading and driving over the platform. 35 mm height for the 600 kg models and 45 mm for the 1500 kg models.

Platforms are equipped with a 5 m double-shielded EasyConnect™ cable, which allows instant configuration with Defender 3000 (i-DT33), Defender 5000 (TD52) or Defender 6000 (i-DT61) indicators for quick installation.

Platform is available in 4 different models with 1000 x 1000 mm and 1250 x 1250 mm sizes and capacities of 600 to 1500 kg, allowing greater flexibility to best suit your weighing needs.


Model	Stainless Steel	Painted Steel	Capacity	Platform (mm)	Item N°	Price €
DFD600HR EU		•	600kg	1000 x 1000 x 35	30540661	1715,00
DFD1500HL EU		•	1500kg	1250 x 1250 x 45	30540662	1930,00
DFD600WR EU	•		600kg	1000 x 1000 x 35	30540663	4020,00
DFD1500WL EU	•		1500kg	1250 x 1250 x 45	30540664	4515,00

DF 5000 SERIES Stainless Steel Floor Scale Platforms

High quality stainless steel platform with structural bracing

Design Features

Platforms are equipped with a 5 m double-shielded EasyConnect™ cable, which allows instant configuration with Defender 3000 (i-DT33), Defender 5000 (TD52) or Defender 6000 (i-DT61) indicators for quick installation



Structural channel bracing for durability.

IP68 stainless steel load cells.

Pit mounting possibility.

Model	Capacity	Platform (mm)	Item N°	Price €
DF300G1S	300kg	800 x 800 x 90	30515665	3510,00
DF600G1S	600kg	800 x 800 x 90	30515666	3510,00
DF300G1R	300kg	1000 x 1000 x 90	30515667	3975,00
DF600G1R	600kg	1000 x 1000 x 90	30515668	3975,00
DF1500G1R	1500kg	1000 x 1000 x 90	30515669	3975,00
DF600G1LX	600kg	1500 x 1250 x 90	30515670	5095,00
DF1500G1LX	1500kg	1500 x 1250 x 90	30515671	5095,00
DF3000G1LX	3000kg	1500 x 1250 x 90	30515672	5095,00
DF600G1X	600kg	1500 x 1500 x 90	30515673	5520,00
DF1500G1X	1500kg	1500 x 1500 x 90	30515674	5520,00
DF3000G1X	3000kg	1500 x 1500 x 90	30515675	5520,00

DF 5000 SERIES Floor Scale Platforms

High quality painted steel platform

Design Features

Platforms are equipped with a 5 m double-shielded EasyConnect™ cable, which allows instant configuration with Defender 3000 (i-DT33), Defender 5000 (TD52) or Defender 6000 (i-DT61) indicators for quick installation



Structural channel bracing for durability.

Pit mounting possibility.

IP68 stainless steel load cells.

Model	Capacity	Platform (mm)	Item N°	Price €
DF300F1S	300kg	800 x 800 x 90	30515562	2355,00
DF600F1S	600kg	800 x 800 x 90	30515563	2355,00
DF300F1R	300kg	1000 x 1000 x 90	30515564	2530,00
DF600F1R	600kg	1000 x 1000 x 90	30515565	2530,00
DF1500F1R	1500kg	1000 x 1000 x 90	30515566	2530,00
DF600F1LX	600kg	1500 x 1250 x 90	30515567	2810,00
DF1500F1LX	1500kg	1500 x 1250 x 90	30515568	2810,00
DF3000F1LX	3000kg	1500 x 1250 x 90	30515569	2810,00
DF600F1X	600kg	1500 x 1500 x 90	30515660	3180,00
DF1500F1X	1500kg	1500 x 1500 x 90	30515661	3180,00
DF3000F1X	3000kg	1500 x 1500 x 90	30515662	3180,00
DF1500F1XV	1500kg	2000 x 1500 x 90	30515663	3925,00
DF3000F1XV	3000kg	2000 x 1500 x 90	30515664	3925,00

DF 3000 SERIES Floor Scale Platforms (i-DF)

Platforms for Basic Industrial Applications

Construction

i-DFxxC - 5mm fixed top plate, 304-grade stainless steel platform, IP65 stainless steel junction box with side access, 4 x IP67 OIML R60 stainless steel load cells, available stainless steel ramp

i-DFxxB - 5mm fixed top plate, black treaded painted steel platform, IP65 stainless steel junction box with side access, 4 x IP67 OIML R60 alloy steel load cells, two carbon steel eye bolts, available painted steel ramp

Design Features

125% overload protection, four load cell mounted feet with self-swiveling ball-in-cup design, painted and stainless steel floor stands (available as accessories)



Available in a painted carbon steel, raised treaded surface or a 304 stainless steel, sand-blasted flat surface – both with a fixed top plate of 5 mm thickness. Both versions feature structural bracing that ensures a stiff and solid construction.

Install easily with EasyConnect™ – Platforms are equipped with a 5 m EasyConnect™ cable, which allows instant configuration with Defender indicators for quick installation.

OIML R60 approved Load Cells – equipped with four IP67 alloy steel or stainless steel load cells with self-swiveling ball-in-cup design feet to best suit your weighing needs.

Model	Stainless	Steel	Capacity	Platform (mm)	Item N°	Price €
i-DF600C1R	•		600kg	1000 x 1000 x 90	30613550	2360,00
i-DF1500C1R	•		1500kg	1000 x 1000 x 90	30613551	2360,00
i-DF1500C1RL	•		1500kg	1000 x 1250 x 90	30786250	2585,00
i-DF1500C1L	•		1500kg	1250 x 1250 x 90	30613554	2775,00
i-DF3000C1RL	•		3000kg	1000 x 1250 x 90	30786251	2585,00
i-DF3000C1L	•		3000kg	1250 x 1250 x 90	30613555	2775,00
i-DF1500C1LX	•		1500kg	1250 x 1500 x 90	30786254	3705,00
i-DF1500C1X	•		1500kg	1500 x 1500 x 90	30613558	4325,00
i-DF3000C1LX	•		3000kg	1250 x 1500 x 90	30786255	3705,00
i-DF3000C1X	•		3000kg	1500 x 1500 x 90	30613559	4325,00
i-DF600B1R		•	600kg	1000 x 1000 x 90	30613548	1400,00
i-DF1500B1R		•	1500kg	1000 x 1000 x 90	30613549	1400,00
i-DF1500B1RL		•	1500kg	1000 x 1250 x 90	30786248	1425,00
i-DF1500B1L		•	1500kg	1250 x 1250 x 90	30613552	1445,00
i-DF3000B1RL		•	3000kg	1000 x 1250 x 90	30786249	1425,00
i-DF3000B1L		•	3000kg	1250 x 1250 x 90	30613553	1445,00
i-DF1500B1LX		•	1500kg	1250 x 1500 x 90	30786252	1470,00
i-DF1500B1X		•	1500kg	1500 x 1500 x 90	30613556	1575,00
i-DF3000B1LX		•	3000kg	1250 x 1500 x 90	30786253	1470,00
i-DF3000B1X		•	3000kg	1500 x 1500 x 90	30613557	1575,00

DF 3000 SERIES Floor Scale Platforms (DF)

Versatile, low cost painted steel floor scale platform



IP67 OIML R60 alloy steel load cells.

Platform with raised embossed treads.

Model	Capacity	Platform (mm)	Item N°	Price €
DF1500BS	1500kg	800 x 800 x 110	83032905	965,00
DF1500BR	1500kg	1000 x 1000 x 110	83032906	1085,00
DF1500BL	1500kg	1200 x 1200 x 110	83032907	1260,00
DF3000BL	3000kg	1200 x 1200 x 110	83032908	1260,00
DF1500BLX	1500kg	1200 x 1500 x 110	83033738	1420,00
DF3000BLX	3000kg	1200 x 1500 x 110	83033739	1420,00
DF1500BX	1500kg	1500 x 1500 x 110	83032909	1560,00
DF3000BX	3000kg	1500 x 1500 x 110	83032910	1560,00

Accessories for Floor Scale Platforms

Accessories	DFLxx	DFDxx	DFxxG1	DFxxF1	i-DF	DFxxB	Item N°	Price €
Ramp, DFL SS 1000 mm	•						30542186	799,00
Ramp, DFL SS 1250 mm	•						30542187	1023,00
Stop Bar, DFD DFL SS 1000 mm	•	•					30550093	172,00
Stop Bar, DFD DFL SS 1250 mm	•	•					30550094	172,00
Floor column kit, Painted 980 mm	•	•	•	•	•	•	30531244	127,00
Floor column kit, stainless steel 980 mm	•	•	•	•	•	•	30539553	235,00
Ramp, DFD CS 1000 mm		•					30478876	286,00
Ramp, DFD CS 1250 mm		•					30571736	374,00
Ramp, DFD SS 1000 mm		•					30478885	684,00
Ramp, DFD SS 1250 mm		•					30571737	908,00
Stop Bar, DFD CS 1000 mm		•					30550089	102,00
Stop Bar, DFD CS 1250 mm		•					30550090	102,00
Foot Plate, Set(2), Painted, DF-B			•				72225239	41,00
Ramp, DF Stainless Steel G1 800 mm			•				22015473	908,00
Ramp, DF Stainless Steel G1 1000 mm			•				22015474	1012,00
Ramp, DF Stainless Steel G1 1250 mm			•				22015475	1264,00
Ramp, DF Stainless Steel G1 1500 mm			•				22015476	1531,00
Pit Frame, DF Stainless Steel G1S			•				22015493	1384,00
Pit Frame, DF Stainless Steel G1R			•				22015494	1487,00
Pit Frame, DF Stainless Steel G1LX			•				22015495	1711,00
Pit Frame, DF Stainless Steel G1X			•				22015496	1837,00
Footplates, DF Stainless Steel G1 Set (2)			•				22015581	104,00
Pit Frame, DF Painted F1S				•			22015392	995,00
Pit Frame, DF Painted F1R				•			22015393	1061,00
Pit Frame, DF Painted F1LX				•			22015394	1181,00
Pit Frame, DF Painted F1X				•			22015396	1257,00
Pit Frame, DF Painted F1XV				•			22015397	1373,00
Ramp, DF Painted F1 800 mm				•			22015464	494,00
Ramp, DF Painted F1 1000 mm				•			22015465	539,00
Ramp, DF Painted F1 1250 mm				•			22015466	624,00
Ramp, DF Painted F1 1500 mm				•			22015467	722,00
Footplates, DF Painted F1 Set (2)				•			22015499	84,00
Ramp, Painted DF-B1 1000mm					•		72198338	345,00
Ramp, Painted DF-B1 1250mm					•		72198339	396,00
Ramp, Painted DF-B1 1500mm					•		72198340	428,00
Ramp, Stainless Steel DF-C1 1000mm					•		72198341	919,00
Ramp, Stainless Steel DF-C1 1250mm					•		72198342	1039,00
Ramp, Stainless Steel DF-C1 1500mm					•		72198343	1427,00
Foot plate kit, Painted DF-B1					•		72198344	46,00
Foot plate kit, Stainless Steel DF-C1					•		72198345	99,00
Ramp, 800 mm DF						•	80252722	235,00
Ramp, 1000 mm DF						•	80252723	274,00
Ramp, 1200 mm DF						•	80252724	320,00
Ramp, 1500 mm DF						•	80252725	383,00
Pit frame, 800 mm for DF						•	83033708	207,00
Pit frame, 1000 mm for DF						•	83033709	260,00
Pit frame, 1200 mm for DF						•	83033710	297,00
Pit frame, 1200 mm x 1500 mm for DF						•	83033714	330,00
Pit frame, 1500 mm for DF						•	83033711	352,00
LC Extension Cable, 3.5 m *1	•	•	•	•	•		30803566	78,00

1 Cannot be used for legal-for-trade applications



RAMP DFL SS 1000 MM
30542186



FLOOR COLUMN KIT, PAINTED 980 MM
30531244



RAMP, DFD CS 1000 MM
30478876



RAMP DF STAINLESS STEEL G1 800 MM
22015473

DFP SERIES Stainless Steel and Steel Pallet Platforms



Prime German quality stainless steel pallet platform

Design Features

Platforms are equipped with a 5 m double-shielded EasyConnect™ cable, which allows instant configuration with Defender 3000 (i-DT33), Defender 5000 (TD52) or Defender 6000 (i-DT61) indicators for quick installation



The U-frame provides easy on and off loading of open-bottom pallets and containers, reducing risk of injury and over-exertion.

Four IP68 electro-polished stainless steel load cells that are laser-weld hermetically sealed and OIML approved are standard on all models.

Structural channel bracing and pendulum design feet with removable rubber cap.

Model	IP Dust and Water Protected	Stainless Steel	Painted Steel	Capacity	Platform	Item N°	Price €
DFP600K1	•	•		600kg	1256mm x 840mm	30515680	3245,00
DFP1500K1	•	•		1500kg	1256mm x 840mm	30515681	3245,00
DFP600J1			•	600kg	1256mm x 840mm	30515678	2315,00
DFP1500J1			•	1500kg	1256mm x 840mm	30515679	2315,00

DF 3000 SERIES Stainless Steel and Steel Pallet Platforms



Platforms for Basic Industrial Applications



Defender Series Pallet Platforms are available in the popular 1 500 kg capacity as painted carbon steel or 304 stainless steel models with a handle and rollers for convenient lifting and relocation.

Install easily with EasyConnect™ – Platforms are equipped with a 5 m EasyConnect™ cable, which allows instant configuration with Defender indicators for quick installation.

OIML R60 approved Load Cells – equipped with four IP67 alloy steel or stainless steel load cells with self-swiveling ball-in-cup design feet to best suit your weighing needs.

Model	Stainless	Steel	Capacity	Platform	Item N°	Price €
i-DFP1500C1	•		1500kg	1260mm x 840mm x 90mm	30700618	2360,00
i-DFP1500B1		•	1500kg	1260mm x 840mm x 90mm	30700557	1230,00

OIML E1 WEIGHTS

Wire shape, stainless steel

Model	Item N°	Price €
Weight 1 mg, OIML E1 Wire	30517704	81,50
Weight 2 mg, OIML E1 Wire	30517705	81,50
Weight 5 mg, OIML E1 Wire	30517706	81,50
Weight 10 mg, OIML E1 Wire	30517707	81,50
Weight 20 mg, OIML E1 Wire	30517958	81,50
Weight 50 mg, OIML E1 Wire	30517959	81,50
Weight 100 mg, OIML E1 Wire	30517960	81,50
Weight 200 mg, OIML E1 Wire	30517961	81,50
Weight 500 mg, OIML E1 Wire	30517962	81,50

With NVLAP Certificate

Item N°	Price €
30517963	177,00
30517964	177,00
30517965	177,00
30517966	177,00
30517967	177,00
30517968	177,00
30517969	177,00
30517970	177,00
30517971	177,00



OIML E1 WEIGHTS

Cylindrical shape, Troemner Alloy 8

Model	Item N°	Price €
Weight 1 g, OIML E1	30517972	151,00
Weight 2 g, OIML E1	30517973	155,00
Weight 5 g, OIML E1	30517974	158,00
Weight 10 g, OIML E1	30517975	173,00
Weight 20 g, OIML E1	30517976	182,00
Weight 50 g, OIML E1	30517977	217,00
Weight 100 g, OIML E1	30517978	257,00
Weight 200 g, OIML E1	30517979	304,00
Weight 500 g, OIML E1	30517980	430,00
Weight 1 kg, OIML E1	30517981	655,00
Weight 2 kg, OIML E1	30517982	969,00
Weight 5 kg, OIML E1	30517983	1620,00
Weight 10 kg, OIML E1	30517984	2320,00
Weight 20 kg, OIML E1	30517985	4080,00

With NVLAP Certificate

Item N°	Price €
30390941	433,00
30390940	442,00
30390939	452,00
30390938	464,00
30390937	478,00
30390936	507,00
30390935	546,00
30390934	598,00
30390933	732,00
30390932	943,00
30390931	1630,00
30390930	2270,00
30390929	2980,00
30390928	5360,00



OIML E1 WEIGHTS SETS

Flat polygonal sheet, Troemner Alloy 8, stainless steel, aluminum

Model	Item N°	Price €
Weight Set (12), 500 mg-1 mg, OIML E1	30517987	1060,00
Weight Set (20), 50 g-1 mg, OIML E1	30517988	2110,00
Weight Set (21), 100 g-1 mg, OIML E1	30517989	2320,00
Weight Set (23), 2x200 g-1 mg, OIML E1	30517990	2870,00
Weight Set (24), 500 g-1 mg, OIML E1	30517991	3320,00
Weight Set (25), 1 kg-1 mg, OIML E1	30517992	3920,00
Weight Set (27), 2x2 kg-1 mg, OIML E1	30517993	5730,00
Weight Set (28), 5 kg-1 mg, OIML E1	30517994	7240,00
Weight Set (29), 10 kg-1 mg, OIML E1	30517995	9370,00
Weight Set (8), 50 g-1 g, OIML E1	30517996	1260,00
Weight Set (9), 100 g-1 g, OIML E1	30517997	1470,00
Weight Set (11), 2x200 g-1 g, OIML E1	30517998	2000,00
Weight Set (12), 500 g-1 g, OIML E1	30517999	2460,00
Weight Set (13), 1 kg-1 g, OIML E1	30518000	3030,00
Weight Set (15), 2x2 kg-1 g, OIML E1	30518001	4780,00
Weight Set (16), 5 kg-1 g, OIML E1	30518002	6350,00
Weight Set (17), 10 kg-1 g, OIML E1	30518003	8510,00

With NVLAP Certificate

Item N°	Price €
30390926	1620,00
30390924	3680,00
30390922	3770,00
30390040	4890,00
30390918	5250,00
30390916	6220,00
30390038	8740,00
30390912	11100,00
30390909	13400,00
30390923	2395,00
30390921	2640,00
30390039	3410,00
30390917	4150,00
30390915	4690,00
30390037	7200,00
30390911	9370,00
30390908	12200,00



All weights come with a plastic case.

Weights and weight sets with NVLAP certificates: In order to get a personalized certificate please contact our Sales Support Center on receipt of your order and provide the calibration certificate number

WEIGHT SET COMPOSITION - CLASS E1

Item N°	NVLAP Certificate	1 mg	2 mg*	5 mg	10 mg	20 mg*	50 mg	100 mg	200 mg*	500 mg	1 g	2 g*	5 g	10 g	20 g*	50 g	100 g	200 g*	500 g	1 kg	2 kg*	5 kg	10 kg
30517987		.	.	.	.	.	.	.	.	.													
30517988		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.							
30517989		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.						
30517990		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.					
30517991		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.				
30517992		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.		
30517993		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	
30517994		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30517995		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30517996											.	.	.	.	.	.	.	.	.	.	.	.	.
30517997											.	.	.	.	.	.	.	.	.	.	.	.	.
30517998											.	.	.	.	.	.	.	.	.	.	.	.	.
30517999											.	.	.	.	.	.	.	.	.	.	.	.	.
30518000											.	.	.	.	.	.	.	.	.	.	.	.	.
30518001											.	.	.	.	.	.	.	.	.	.	.	.	.
30518002											.	.	.	.	.	.	.	.	.	.	.	.	.
30518003											.	.	.	.	.	.	.	.	.	.	.	.	.
30390926	.	.	.	.	.	.	.	.	.	.													
30390924	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390922	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390040	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390918	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390916	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390038	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390912	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390909	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390923	.										.	.	.	.	.	.	.	.	.	.	.	.	.
30390921	.										.	.	.	.	.	.	.	.	.	.	.	.	.
30390039	.										.	.	.	.	.	.	.	.	.	.	.	.	.
30390917	.										.	.	.	.	.	.	.	.	.	.	.	.	.
30390915	.										.	.	.	.	.	.	.	.	.	.	.	.	.
30390037	.										.	.	.	.	.	.	.	.	.	.	.	.	.
30390911	.										.	.	.	.	.	.	.	.	.	.	.	.	.
30390908	.										.	.	.	.	.	.	.	.	.	.	.	.	.

*Weight set contains two weights of this value

All weights come with a master carrying case.

- Sets where the largest weight is under 1 kg come with a brush and forceps.
- Sets where the largest weight is 1 kg – 5 kg come with brush, lifters and forceps.
- Weight sets where the largest weight is 10 kg and above are supplied with a Complete Accessory Kit which includes a microfibre cloth, nylon glove, small and wide brushes, and forceps.
- Weight and weight sets with NVLAP certificates: In order to get personalized certificate, please contact Sales Support Center.

The National Voluntary Laboratory Accreditation Program (NVLAP) provides third-party accreditation to testing and calibration laboratories in response to legislative actions or requests from government agencies or private-sector organizations. NVLAP-accredited laboratories are assessed against the management and technical requirements published in the International Standard, ISO/IEC 17025:2017. NVLAP has obtained approval from the International Laboratory Accreditation Cooperation (ILAC) for the use of the ILAC MRA Mark. For further details as well as list of member organizations please visit <https://ilac.org/>.

OIML E2 WEIGHTS

Flat polygonal sheet, aluminium, stainless steel

Model	Item N°	Price €
Weight 1 mg, OIML E2	30518004	36,50
Weight 2 mg, OIML E2	30518005	36,50
Weight 5 mg, OIML E2	30518006	36,50
Weight 10 mg, OIML E2	30518007	36,50
Weight 20 mg, OIML E2	30518008	36,50
Weight 50 mg, OIML E2	30518009	36,50
Weight 100 mg, OIML E2	30518010	36,50
Weight 200 mg, OIML E2	30518011	36,50
Weight 500 mg, OIML E2	30518012	36,50

With NVLAP Certificate

Item N°	Price €
80780284	71,50
80780285	71,50
80780286	71,50
80780287	71,50
80780288	71,50
80780289	71,50
80780290	71,50
80780291	71,50
80780292	71,50



OIML E2 WEIGHTS

Cylindrical shape, Troemner Alloy 8

Model	Item N°	Price €
Weight 1 g, OIML E2	30518013	72,75
Weight 2 g, OIML E2	30518014	76,50
Weight 5 g, OIML E2	30518015	81,50
Weight 10 g, OIML E2	30518016	84,25
Weight 20 g, OIML E2	30518017	97,25
Weight 50 g, OIML E2	30518018	112,00
Weight 100 g, OIML E2	30518019	120,00
Weight 200 g, OIML E2	30518020	151,00
Weight 500 g, OIML E2	30518021	253,00
Weight 1 kg, OIML E2	30518022	367,00
Weight 2 kg, OIML E2	30518023	582,00
Weight 5 kg, OIML E2	30518024	855,00
Weight 10 kg, OIML E2	30518025	1240,00
Weight 20 kg, OIML E2	30518026	2080,00

With NVLAP Certificate

Item N°	Price €
80780293	108,00
80780294	112,00
80780295	117,00
80780296	122,00
80780297	134,00
80780298	149,00
80780299	157,00
80780300	195,00
80780301	297,00
80780302	414,00
80780303	639,00
80780304	886,00
80780305	1320,00
80780306	2170,00



OIML E2 WEIGHTS SETS

Cylindrical shape, flat polygonal sheet Troemner Alloy 8, stainless steel, aluminum

Model	Item N°	Price €
Weight Set (12), 500 mg-1 mg, OIML E2	30518028	436,00
Weight Set (20), 50 g-1 mg, OIML E2	30518029	1130,00
Weight Set (21), 100 g-1 mg, OIML E2	30518030	1250,00
Weight Set (23), 2x200 g-1 mg, OIML E2	30518031	1520,00
Weight Set (24), 500 g-1 mg, OIML E2	30518032	1750,00
Weight Set (25), 1 kg-1 mg, OIML E2	30518033	2230,00
Weight Set (27), 2x2 kg-1 mg, OIML E2	30518034	3350,00
Weight Set (28), 5 kg-1 mg, OIML E2	30518035	4150,00
Weight Set (8), 50 g-1 g, OIML E2	30518036	752,00
Weight Set (9), 100 g-1 g, OIML E2	30518037	855,00
Weight Set (11), 2x200 g-1 g, OIML E2	30518038	1140,00
Weight Set (12), 500 g-1 g, OIML E2	30518039	1390,00
Weight Set (13), 1 kg-1 g, OIML E2	30518040	1840,00
Weight Set (15), 2x2 kg-1 g, OIML E2	30518041	2980,00
Weight Set (16), 5 kg-1 g, OIML E2	30518042	3760,00

With NVLAP Certificate

Item N°	Price €
80780381	711,00
80780383	1560,00
80780385	1720,00
80780387	1910,00
80780391	2290,00
80780393	2820,00
80780395	4020,00
80780399	4870,00
80780384	917,00
80780386	1060,00
80780388	1420,00
80780392	1690,00
80780394	2160,00
80780396	3380,00
80780400	4230,00



All weights come with a plastic case.

Weights and weight sets with NVLAP certificates: In order to get a personalized certificate please contact our Sales Support Center on receipt of your order and provide the calibration certificate number

WEIGHT SET COMPOSITION - CLASS E2

Item N°	NVLAP Certificate	1 mg	2 mg*	5 mg	10 mg	20 mg*	50 mg	100 mg	200 mg*	500 mg	1 g	2 g*	5 g	10 g	20 g*	50 g	100 g	200 g*	500 g	1 kg	2 kg*	5 kg	10 kg
30518028		.	.	.	.	.	.	.	.	.													
30518029		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.							
30518030		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.						
30518031		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.					
30518032		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.				
30518033		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.		
30518034		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	
30518035		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30518036											.	.	.	.	.	.	.	.	.	.	.	.	.
30518037											.	.	.	.	.	.	.	.	.	.	.	.	.
30518038											.	.	.	.	.	.	.	.	.	.	.	.	.
30518039											.	.	.	.	.	.	.	.	.	.	.	.	.
30518040											.	.	.	.	.	.	.	.	.	.	.	.	.
30518041											.	.	.	.	.	.	.	.	.	.	.	.	.
30518042											.	.	.	.	.	.	.	.	.	.	.	.	.
80780381	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780383	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780385	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780387	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780391	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780393	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780395	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780399	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780384	.										.	.	.	.	.	.	.	.	.	.	.	.	.
80780386	.										.	.	.	.	.	.	.	.	.	.	.	.	.
80780388	.										.	.	.	.	.	.	.	.	.	.	.	.	.
80780392	.										.	.	.	.	.	.	.	.	.	.	.	.	.
80780394	.										.	.	.	.	.	.	.	.	.	.	.	.	.
80780396	.										.	.	.	.	.	.	.	.	.	.	.	.	.
80780400	.										.	.	.	.	.	.	.	.	.	.	.	.	.

*Weight set contains two weights of this value

All weights come with a master carrying case.

- Sets where the largest weight is under 1 kg come with a brush and forceps.
- Sets where the largest weight is 1 kg – 5 kg come with brush, lifters and forceps.
- Weight sets where the largest weight is 10 kg and above are supplied with a Complete Accessory Kit which includes a microfibre cloth, nylon glove, small and wide brushes, and forceps.
- Weight and weight sets with NVLAP certificates: In order to get personalized certificate, please contact Sales Support Center.

The National Voluntary Laboratory Accreditation Program (NVLAP) provides third-party accreditation to testing and calibration laboratories in response to legislative actions or requests from government agencies or private-sector organizations. NVLAP-accredited laboratories are assessed against the management and technical requirements published in the International Standard, ISO/IEC 17025:2017. NVLAP has obtained approval from the International Laboratory Accreditation Cooperation (ILAC) for the use of the ILAC MRA Mark. For further details as well as list of member organizations please visit <https://ilac.org/>.

OIML F1 WEIGHTS

Flat polygonal sheet, aluminium, stainless steel

Model	Item N°	Price €
Weight, 1 mg, OIML F1	80780307	30,00
Weight, 2 mg, OIML F1	80780308	30,00
Weight, 5 mg, OIML F1	80780309	30,00
Weight, 10 mg, OIML F1	80780310	30,00
Weight, 20 mg, OIML F1	80780311	30,00
Weight, 50 mg, OIML F1	80780312	30,00
Weight, 100 mg, OIML F1	80780313	30,00
Weight, 200 mg, OIML F1	80780314	30,00
Weight, 500 mg, OIML F1	80780315	30,00

With NVLAP Certificate

Item N°	Price €
80780331	41,75
80780332	41,75
80780333	41,75
80780334	41,75
80780335	41,75
80780336	41,75
80780337	41,75
80780338	41,75
80780339	41,75



OIML F1 WEIGHTS

Cylindrical shape, Troemner Alloy 8, stainless steel

Model	Item N°	Price €
Weight, 1 g, OIML F1	80780316	39,00
Weight, 2 g, OIML F1	80780317	41,75
Weight, 5 g, OIML F1	80780318	44,25
Weight, 10 g, OIML F1	80780319	48,25
Weight, 20 g, OIML F1	80780320	53,50
Weight, 50 g, OIML F1	80780321	59,75
Weight, 100 g, OIML F1	80780322	75,25
Weight, 200 g, OIML F1	80780323	94,50
Weight, 500 g, OIML F1	80780324	147,00
Weight, 1 kg, OIML F1	80780325	204,00
Weight, 2 kg, OIML F1	80780326	324,00
Weight, 5 kg, OIML F1	80780327	458,00
Weight, 10 kg, OIML F1	80780328	871,00
Weight, 20 kg, OIML F1	80780329	1380,00

With NVLAP Certificate

Item N°	Price €
80780340	64,75
80780341	67,50
80780342	70,00
80780343	73,75
80780344	79,25
80780345	85,50
80780346	103,00
80780347	123,00
80780348	176,00
80780349	234,00
80780350	361,00
80780351	521,00
80780352	943,00
80780353	1470,00



OIML F1 WEIGHTS SETS

Cylindrical shape, flat polygonal sheet Troemner Alloy 8, stainless steel, aluminum

Model	Item N°	Price €
Weight Set (12), 500 mg-1 mg, OIML F1	80780408	288,00
Weight Set (20), 50 g-1 mg, OIML F1	80780410	572,00
Weight Set (21), 100 g-1 mg, OIML F1	80780412	634,00
Weight Set (23), 2x200 g-1 mg, OIML F1	80780414	819,00
Weight Set (24), 500 g-1 mg, OIML F1	80780418	876,00
Weight Set (25), 1 kg-1 mg, OIML F1	80780420	1440,00
Weight Set (27), 2x2 kg-1 mg, OIML F1	80780422	1880,00
Weight Set (28), 5 kg-1 mg, OIML F1	80780426	2290,00
Weight Set (8), 50 g-1 g, OIML F1	80780411	355,00
Weight Set (9), 100 g-1 g, OIML F1	80780413	429,00
Weight Set (11), 2x200 g-1 g, OIML F1	80780415	546,00
Weight Set (12), 500 g-1 g, OIML F1	80780419	649,00
Weight Set (13), 1 kg-1 g, OIML F1	80780421	1240,00
Weight Set (15), 2x2 kg-1 g, OIML F1	80780423	1760,00
Weight Set (16), 5 kg-1 g, OIML F1	80780427	2020,00

With NVLAP Certificate

Item N°	Price €
80780435	376,00
80780437	783,00
80780439	850,00
80780441	1010,00
80780445	1110,00
80780447	1520,00
80780449	2050,00
80780453	2460,00
80780438	449,00
80780440	512,00
80780442	655,00
80780446	773,00
80780448	1360,00
80780450	1960,00
80780454	2340,00



All weights come with a plastic case.

Weights and weight sets with NVLAP certificates: In order to get a personalized certificate please contact our Sales Support Center on receipt of your order and provide the calibration certificate number

WEIGHT SET COMPOSITION - CLASS F1

Item N°	NVLAP Certificate	1 mg	2 mg*	5 mg	10 mg	20 mg*	50 mg	100 mg	200 mg*	500 mg	1 g	2 g*	5 g	10 g	20 g*	50 g	100 g	200 g*	500 g	1 kg	2 kg*	5 kg	10 kg
80780408		.	.	.	.	.	.	.	.	.													
80780410		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.							
80780412		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.						
80780414		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.					
80780418		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.				
80780420		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.		
80780422		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	
80780426		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780411											.	.	.	.	.	.	.	.	.	.	.	.	.
80780413											.	.	.	.	.	.	.	.	.	.	.	.	.
80780415											.	.	.	.	.	.	.	.	.	.	.	.	.
80780419											.	.	.	.	.	.	.	.	.	.	.	.	.
80780421											.	.	.	.	.	.	.	.	.	.	.	.	.
80780423											.	.	.	.	.	.	.	.	.	.	.	.	.
80780427											.	.	.	.	.	.	.	.	.	.	.	.	.
80780435	.	.	.	.	.	.	.	.	.	.													
80780437	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780439	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780441	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780445	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780447	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780449	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780453	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780438	.										.	.	.	.	.	.	.	.	.	.	.	.	.
80780440	.										.	.	.	.	.	.	.	.	.	.	.	.	.
80780442	.										.	.	.	.	.	.	.	.	.	.	.	.	.
80780446	.										.	.	.	.	.	.	.	.	.	.	.	.	.
80780448	.										.	.	.	.	.	.	.	.	.	.	.	.	.
80780450	.										.	.	.	.	.	.	.	.	.	.	.	.	.
80780454	.										.	.	.	.	.	.	.	.	.	.	.	.	.

*Weight set contains two weights of this value.

All weights come with a master carrying case.

- Sets where the largest weight is under 1 kg come with a brush and forceps.
- Sets where the largest weight is 1 kg – 5 kg come with brush, lifters and forceps.
- Weight sets where the largest weight is 10 kg and above are supplied with a Complete Accessory Kit which includes a microfibre cloth, nylon glove, small and wide brushes, and forceps.
- Weight and weight sets with NVLAP certificates: In order to get personalized certificate, please contact Sales Support Center.

The National Voluntary Laboratory Accreditation Program (NVLAP) provides third-party accreditation to testing and calibration laboratories in response to legislative actions or requests from government agencies or private-sector organizations. NVLAP-accredited laboratories are assessed against the management and technical requirements published in the International Standard, ISO/IEC 17025:2017. NVLAP has obtained approval from the International Laboratory Accreditation Cooperation (ILAC) for the use of the ILAC MRA Mark. For further details as well as list of member organizations please visit <https://ilac.org/>.

OIML F2 WEIGHTS

Flat polygonal sheet, aluminium, stainless steel

Model	Item N°	Price €
Weight, 1 mg, OIML F2	30390876	28,75
Weight, 2 mg, OIML F2	30390875	28,75
Weight, 5 mg, OIML F2	30390874	28,75
Weight, 10 mg, OIML F2	30390873	28,75
Weight, 20 mg, OIML F2	30390872	28,75
Weight, 50 mg, OIML F2	30390871	28,75
Weight, 100 mg, OIML F2	30390870	28,75
Weight, 200 mg, OIML F2	30390869	28,75
Weight, 500 mg, OIML F2	30390868	28,75

With NVLAP Certificate

Item N°	Price €
30390900	39,00
30390899	39,00
30390898	39,00
30390897	39,00
30390896	39,00
30390895	39,00
30390894	39,00
30390893	39,00
30390892	39,00



OIML F2 WEIGHTS

Cylindrical shape, Troemner Alloy 8, stainless steel

Model	Item N°	Price €
Weight, 1 g, OIML F2	30390867	27,50
Weight, 2 g, OIML F2	30390866	28,75
Weight, 5 g, OIML F2	30390865	30,00
Weight, 10 g, OIML F2	30390864	34,00
Weight, 20 g, OIML F2	30390863	36,50
Weight, 50 g, OIML F2	30390862	40,50
Weight, 100 g, OIML F2	30390861	47,00
Weight, 200 g, OIML F2	30390860	64,75
Weight, 500 g, OIML F2	30390859	111,00
Weight, 1 kg, OIML F2	30390858	144,00
Weight, 2 kg, OIML F2	30390857	163,00
Weight, 5 kg, OIML F2	30390856	298,00
Weight, 10 kg, OIML F2	30390855	531,00
Weight, 20 kg, OIML F2	30390854	809,00

With NVLAP Certificate

Item N°	Price €
30390891	50,75
30390890	51,75
30390889	54,50
30390888	57,00
30390887	61,00
30390886	63,50
30390885	71,50
30390884	85,50
30390883	114,00
30390882	146,00
30390881	199,00
30390880	390,00
30390879	593,00
30390878	1030,00



OIML F2 WEIGHTS SETS

Cylindrical shape, flat polygonal sheet Troemner Alloy 8, Stainless Steel, Aluminum

Model	Item N°	Price €
Weight Set (12), 500 mg-1 mg, OIML F2	30390829	293,00
Weight Set (20), 50 g-1 mg, OIML F2	30390827	531,00
Weight Set (21), 100 g-1 mg, OIML F2	30390825	582,00
Weight Set (23), 2x200 g-1 mg, OIML F2	30390032	716,00
Weight Set (24), 500 g-1 mg, OIML F2	30390821	835,00
Weight Set (25), 1 kg-1 mg, OIML F2	30390819	1280,00
Weight Set (27), 2x2 kg-1 mg, OIML F2	30390030	1630,00
Weight Set (28), 5 kg-1 mg, OIML F2	30390815	1830,00
Weight Set (8), 50 g-1 g, OIML F2	30390826	306,00
Weight Set (9), 100 g-1 g, OIML F2	30390824	365,00
Weight Set (11), 2x200 g-1 g, OIML F2	30390031	492,00
Weight Set (12), 500 g-1 g, OIML F2	30390820	582,00
Weight Set (13), 1 kg-1 g, OIML F2	30390818	1040,00
Weight Set (15), 2x2 kg-1 g, OIML F2	30390029	1440,00
Weight Set (16), 5 kg-1 g, OIML F2	30390814	1610,00

With NVLAP Certificate

Item N°	Price €
30390852	333,00
30390850	737,00
30390848	809,00
30390036	953,00
30390844	1070,00
30390842	1310,00
30390034	1750,00
30390838	2090,00
30390849	380,00
30390847	469,00
30390035	634,00
30390843	711,00
30390841	1050,00
30390033	1740,00
30390837	1760,00



All weights come with a plastic case.

Weights and weight sets with NVLAP certificates: In order to get a personalized certificate please contact our Sales Support Center on receipt of your order and provide the calibration certificate number

WEIGHT SET COMPOSITION - CLASS F2

Item N°	NVLAP Certificate	1 mg	2 mg*	5 mg	10 mg	20 mg*	50 mg	100 mg	200 mg*	500 mg	1 g	2 g*	5 g	10 g	20 g*	50 g	100 g	200 g*	500 g	1 kg	2 kg*	5 kg	10 kg
30390829		.	.	.	.	.	.	.	.	.													
30390827		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.							
30390825		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.						
30390032		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.					
30390821		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.				
30390819		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.		
30390030		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	
30390815		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390826											.	.	.	.	.	.	.	.	.	.	.	.	.
30390824											.	.	.	.	.	.	.	.	.	.	.	.	.
30390031											.	.	.	.	.	.	.	.	.	.	.	.	.
30390820											.	.	.	.	.	.	.	.	.	.	.	.	.
30390818											.	.	.	.	.	.	.	.	.	.	.	.	.
30390029											.	.	.	.	.	.	.	.	.	.	.	.	.
30390814											.	.	.	.	.	.	.	.	.	.	.	.	.
30390852	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390850	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390848	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390036	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390844	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390842	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390034	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390838	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390849	.										.	.	.	.	.	.	.	.	.	.	.	.	.
30390847	.										.	.	.	.	.	.	.	.	.	.	.	.	.
30390035	.										.	.	.	.	.	.	.	.	.	.	.	.	.
30390843	.										.	.	.	.	.	.	.	.	.	.	.	.	.
30390841	.										.	.	.	.	.	.	.	.	.	.	.	.	.
30390033	.										.	.	.	.	.	.	.	.	.	.	.	.	.
30390837	.										.	.	.	.	.	.	.	.	.	.	.	.	.

*Weight set contains two weights of this value.

All weights come with a master carrying case.

- Sets where the largest weight is under 1 kg come with a brush and forceps.
- Sets where the largest weight is 1 kg – 5 kg come with brush, lifters and forceps.
- Weight sets where the largest weight is 10 kg and above are supplied with a Complete Accessory Kit which includes a microfibre cloth, nylon glove, small and wide brushes, and forceps.
- Weight and weight sets with NVLAP certificates: In order to get personalized certificate, please contact Sales Support Center.

The National Voluntary Laboratory Accreditation Program (NVLAP) provides third-party accreditation to testing and calibration laboratories in response to legislative actions or requests from government agencies or private-sector organizations. NVLAP-accredited laboratories are assessed against the management and technical requirements published in the International Standard, ISO/IEC 17025:2017. NVLAP has obtained approval from the International Laboratory Accreditation Cooperation (ILAC) for the use of the ILAC MRA Mark. For further details as well as list of member organizations please visit <https://ilac.org/>.

OIML M1 WEIGHTS

Stainless steel

Model	Item N°	Price €
Weight 1 g OIML M1	11117141	19,75
Weight 2 g OIML M1	11117142	21,00
Weight 5 g OIML M1	11117143	22,25
Weight 10 g OIML M1	11117144	24,50
Weight 20 g OIML M1	11117145	27,00
Weight 50 g OIML M1	11117146	29,50
Weight 100 g OIML M1	11117147	35,50
Weight 200 g OIML M1	11117148	40,50
Weight 500 g OIML M1	11117149	64,75
Weight 1 kg OIML M1	11117150	103,00
Weight 2 kg OIML M1	11117151	164,00
Weight 5 kg OIML M1	11117152	361,00
Weight 10 kg OIML M1	11117153	618,00

Iron

Model	Item N°	Price €
Rectangular weight 20 kg M1	87000007	171,00
Rectangular weight 10 kg M1	87000008	136,00
Rectangular weight 5 kg M1	87000009	103,00

With certificate

Model	Item N°	Price €
Weight OH 1 g M1AC	11117220	119,00
Weight OH 2 g M1AC	11117221	120,00
Weight OH 5 g M1AC	11117222	122,00
Weight OH 10 g M1AC	11117223	124,00
Weight OH 20 g M1AC	11117224	130,00
Weight OH 50 g M1AC	11117225	137,00
Weight OH 100 g M1AC	11117226	145,00
Weight OH 200 g M1AC	11117227	149,00
Weight OH 500 g M1AC	11117228	162,00
Weight OH 1 kg M1AC	11117229	217,00
Weight OH 2 kg M1AC	11117230	284,00
Weight OH 5 kg M1AC	11117231	477,00
Weight OH 10 kg M1AC	11117232	819,00

OIML M3 WEIGHTS

Stainless steel

Model	Item N°	Price €
Weight 1 g OIML M3	11117206	19,75
Weight 2 g OIML M3	11117207	21,00
Weight 5 g OIML M3	11117208	22,25
Weight 10 g OIML M3	11117209	24,50
Weight 20 g OIML M3	11117210	27,00
Weight 50 g OIML M3	11117211	29,50
Weight 100 g OIML M3	11117212	35,50
Weight 200 g OIML M3	11117213	40,50
Weight 500 g OIML M3	11117214	57,25
Weight 1 kg OIML M3	11117215	93,75
Weight 2 kg OIML M3	11117216	147,00
Weight 5 kg OIML M3	11117217	307,00
Weight 10 kg OIML M3	11117218	618,00

With certificate

Model	Item N°	Price €
Weight 1 g OIML M3	11117234	117,00
Weight 2 g OIML M3	11117235	118,00
Weight 5 g OIML M3	11117236	119,00
Weight 10 g OIML M3	11117237	122,00
Weight 20 g OIML M3	11117238	128,00
Weight 50 g OIML M3	11117239	135,00
Weight 100 g OIML M3	11117240	143,00
Weight 200 g OIML M3	11117241	147,00
Weight 500 g OIML M3	11117242	159,00
Weight 1 kg OIML M3	11117243	196,00
Weight 2 kg OIML M3	11117244	253,00
Weight 5 kg OIML M3	11117245	430,00
Weight 10 kg OIML M3	11117246	737,00

OIML Weight sets 1 mg

Model	Item N°	Price €
M1 Weight set 1 mg - 500 mg (ALU)	87000010	138,00
M1 Weight set 1 mg - 100 g (ALU/Brass, SST)	87000011	300,00
M1 Weight set 1 mg - 200 g (ALU/Brass, SST)	87000012	363,00

With certificate

Model	Item N°	Price €
F1 Weight set 1 mg – 500 mg	87000034	470,00

OIML Weight sets 1 g

Model	Item N°	Price €
M1 Weight set 1 g - 50 g (Brass)	87000013	144,00
M1 Weight set 1 g - 500 g (SST, Brass)	87000014	316,00
M1 Weight set 1 g - 1000 g (SST, Brass)	87000015	504,00
M1 Weight set 1 g - 2000 g (SST, Brass)	87000016	881,00

With certificate

Model	Item N°	Price €
F1 Weight set 1 g – 500 g	87000035	974,00

SF-110B

A trusted printer that supports industry compliance

Application

Direct, Weighing, Multi Weighing, Tare/Gross/Net, Statistic, Sum, Mean, Formulation, Loss Weight

Display

OLED 128 x 64 Dot Matrix Display

Operation

AC adapter (included)

Communication

RS232 (Included)

Construction

ABS housing

Design Features

High-speed dot matrix printing, 40 characters per line, auto tare operation directly from the printer, multi-user access control, editable header and trailer



A login password for each user fully prevents unauthorized changes to the printer settings such as date/time.

The SF-110B provides users with full functionality from the most basic printing needs to the advanced applications, including statistics, totalization, formulation, averaging and interval printing.

Additional information like Sample name/ID can be printed before each weighing that complements the report and improves data traceability.

Model	Item N°	Price €
Printer SF-110B Impact	30708371	805,00
Accessories	Item N°	Price €
Paper Roll	12120799	12,90
Ink Ribbon, Black	30529322	12,00

OHAUS Product Test With SF-110B Printer

	Direct	Weighing	Multi weighing	Statistics	Sum	Mean	Formulation	Loss Weight	Zero/Tare
EX	YES	YES	YES	YES	YES	YES	YES	YES	YES
AX	YES	YES	YES	YES	YES	YES	YES	YES	YES
PX/PR	YES	YES	YES	YES	YES	YES	YES	YES	YES
Scout STX	YES	YES	YES	YES	YES	YES	YES	YES	YES
Scout SKX	YES	YES	YES	YES	YES	YES	YES	YES	YES
Navigator	YES	NO	NO	NO	NO	NO	NO	NO	NO
MB90/120	YES*	NO	NO	NO	NO	NO	NO	NO	NO
MB23/25/27	YES	YES	YES	YES	YES	YES	YES	YES	YES
T24P	YES	NO	NO	NO	NO	NO	NO	NO	NO
i-DT33P, i-DT33XW	YES	YES	YES	YES	YES	YES	YES	YES	YES
TD52P/TD52XW	YES	YES	YES	YES	YES	YES	YES	YES	YES
i-DT61XWE	YES	YES	YES	YES	YES	YES	YES	YES	YES
T72XW	NO	NO	NO	NO	NO	NO	NO	NO	NO
R2000/RC2000	YES	YES	YES	YES	YES	YES	YES	YES	YES
R3000/RC3000	YES	YES	YES	YES	YES	YES	YES	YES	YES
R4000/RC4000	YES	YES	YES	YES	YES	YES	YES	YES	YES
R7000	YES	YES	YES	YES	YES	YES	YES	YES	YES
V12	YES	NO	NO	NO	NO	NO	NO	NO	NO
V71	YES	YES	YES	YES	YES	YES	YES	YES	YES
Catapult 5000	YES	NO	NO	NO	NO	NO	NO	NO	NO

*automatically print when the test is finished

SF40A

OHAUS value in a full-featured, portable printer.

Application

Impact Printing, Statistics, Totalization

Operation

AC adapter (included)

Communication

RS233 (included)

Construction

ABS housing

Design Features

Status indicator, 40 characters per line, compatible with continuous labelpaper, perform zero and tare operations directly from the printer



The statistics and totalization modes and real time clock complement OHAUS weighing instruments' application modes and offer functionality that assists the reporting processes.*

SF40A helps support the efficiency of the weighing and recording processes. RS232 interface and PC configuration software make setting up and using the printer fast and simple.

With the ability to print in 13 different languages, SF40A can print in the language of the connected device to support universal operation.

* compatible with most of OHAUS products

Model	Item N°	Price €
Impact printer SF40A	30064202	715,00
Accessories	Item N°	Price €
SF40A Paper rotary axis	30063920	8,05
SF40A Paper roll (57,5mm 2pcs)	12120799	12,90
SF40A Ink ribbon cassette	30529322	12,00
Adhesive Paper Rolls (3 pcs)	11600388	18,20

OHAUS Product Test With SF40A Printer

	Print	Zero/Tare	Statistics	Totalization
Adventurer AX	YES	YES	YES	YES
Explorer EX	YES	YES	YES	YES
Navigator	YES	YES	YES	YES
PR Series	YES	YES	YES	YES
Pioneer PX/PJX/PXP	YES	YES	YES	YES
Scout STX/SKX/SJX	YES	YES	YES	YES
MB23/25	YES	YES	YES	YES
MB90/120	YES	YES	YES	YES
i-DT33P	YES	NO	NO	NO
T51P	YES	YES	YES	YES
TD52	YES	YES	YES	YES
T71P	YES	YES	YES	YES
T32M	YES	NO	NO	NO

	Print	Zero/Tare	Statistics	Totalization
i-DT33XW	YES	NO	NO	NO
T51XW	YES	YES	YES	YES
T71XW	YES	YES	YES	YES
T72XW	YES	YES	YES	YES
T82XW	YES	NO	NO	NO
CKW	YES	NO	NO	NO
Ranger 3000	YES	YES	YES*	YES
Ranger 4000	YES	YES	YES	YES
Ranger 7000	YES	YES	YES	YES
Valor 3000	YES	NO	NO	NO
Valor 7000	YES	YES	YES	YES

*not supported in Dynamic Mode

STP103

Economical Palm-Sized Thermal Printer

Design Features

Status indicator, paper feed, graphics printing capability



- Ultra compact and features low power consumption.
- 42 column printing at a speed of approximately 11.8 lines per second.
- Simple paper loading and suspends printing when printing paper becomes low.

Model	Item N°	Price €
STP103 Thermal Printer	80251993	510,00
Accessories	Item N°	Price €
Thermal Paper Roll for STP103 Printer	80251931	13,50

PRINTER CABLES

Printer Cables	Item N°	Price €	Explorer EX	Adventurer AX	Pioneer PX	PR Series	MB120, 90, 27, 25, 23	Navigator	Scout	Defender DxxP, Defender TxxP	Defender DxxXW, Defender TDxxXW	Ranger	Valor	2nd RS232 Defender	Ranger 7000	2nd RS232 Ranger 3000, 4000	2nd RS232 Valor 7000
SF40A Printer Cable RS9 M/F	30529323	106,50	•	•	•	•	•			•	•					•	
SF40A Printer Cable RS9 F/Bare Wire	80500552	91,80									•			•	•		•
SF40A Printer Cable Adapter RS9 F/M (for Scout and Navigator adapter needed)	30059316	18,80						•	•								
STP103 Printer Cable RS9 M/F + Adapter	80252581	99,70	•	•	•	•	•			•						•	
STP103 Printer Cable RS9 M/Bare Wire + Adapter	80252584	77,50									•			•	•		•
STP103 Printer Cable Adapter RS25 M/ RS 9 M(for Scout and Navigator adapter needed)	80251930	18,60						•	•								
Scout RS232 Interface Kit (Required to connect to STP103 and SF40)	30268982	81,50							•								
Navigator RS232 Interface Kit (Required to connect to STP103 and SF40)	83032107	45,00						•									
RS232 Cable for connection to Zebra printers and reference balances	30057595	78,40	•							•	•	•			•		

Items 30529323 and 30059316 are included with the SF40A printer as standard.
Please check cables compatibility with our Sales Support Center prior to placing an order
Note: Scout and Navigator require RS232 accessory for printer connection
Note: 30057595 needs one connector cutting off to connect Zebra to Defender XW

SERIAL PORT DATA COLLECTION SOFTWARE

Free of charge!



Automatic data capture.

Direct data transfer to MS Office Excel, Word, Access or text file.

Compatible with OHAUS scales and balances equipped with RS232 or Ethernet.

PC CABLES

PC Cables	Item N°	Price €	Explorer EX (Version 1 Terminal)	Explorer EX (Version 2 Terminal)	Adventurer AX	Pioneer PX	PR Series	MB90, 120	MB23, 25, 27	Navigator	Scout	Defender 2000 TxxP	Defender 3000 TxxP, TxxM	Defender 5000 TDxxP	Defender 5000 TDxxXW	Defender 5000 TDxxXW 2nd RS232	Defender 6000 TD61XW	Defender 000 TD61XW 2nd RS232	Defender 7000 T72XW	Defender 7000 T72XW 2nd RS232 (COM 2)	Valor 7000	Valor 7000 2nd RS232	Ranger 2000	Ranger 3000	Ranger 3000 2nd RS232	Ranger 4000	Ranger 4000 2nd RS232	Ranger 7000	Ranger 7000 2nd RS232	
RS232 9 Pin PC Cable RS9 F/ Bare Wire	80500552	91,80													•	•		•	•	•		•				•				•
RS232 9 Pin PC Cable RS9 M/F	410024	92,20	•	•	•	•	•	•	•			•	•	•							•		•	•		•	•	•		
USB (Type A to B) Cable	83021085	22,60	•		•	•		•						•	•		•		•							•			•	
USB A to Mini USB Cable	28120263	21,20		•																										
RS232-USB Converter	30304101	78,25	•	•	•	•	•	•	•			•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
Scout RS232 Interface Kit (Required to connect to PC)	30268982	81,50									•																			
Scout USB Interface Kit (Required to connect to PC)	30268984	91,50									•																			
Navigator RS232 Interface Kit (Required to connect to PC)	83032107	45,00									•																			
Navigator USB Interface Kit (Required to connect to PC)	83032108	45,00									•																			

Please check cables compatibility with our Sales Support Center prior to placing an order

PIPETTE ADJUSTMENT KIT

For Explorer Semi-Micro balances

Application

Pipette Adjustment for Explorer Semi-Micro balances



Model	Item N°	Price €
Pipette adjustment kit	30252145	1230,00

ION-100A Ionizer

Increasing weighing accuracy by reducing static electricity

Application

Eliminating static electricity on plastic and glass weighing vessels which are particularly prone to attracting static charges

Design Features

Fan-free technology, replaceable electrodes, adjustable angle for better alignment with the balance



Fan-free technology neutralizes static with no sample disruption.

An easily replaceable long-running motor.

A compact and versatile design.

Model	Item N°	Price €
ION-100A	30095929	1460,00
Accessories	Item N°	Price €
Replacement Electrode ION-100A	30133391	95,50

DEFENDER™ 2500

For Basic Industrial Applications



Application

Weighing, Parts Counting, Totalization, Dynamic Weighing/Display Hold

Display

6-digit, 7-segment LCD display with white backlight, 28 mm high digits

Operation

Standard AC adapter with USB type C cable, detachable (included) or 4 AA batteries (not included)

Communication

Easy access communication ports including standard RS232 and USB type C

Construction

ABS indicator, 304-stainless steel platform, powder-coated frame, aluminum load cell, non-slip adjustable rubber feet

Design Features

150% overload capacity protection, selectable environmental and auto-print settings, stability indicator, overload/underload indicators, auto-off, auto tare



Features software for most common industrial weighing modes such as Parts counting, Totalization and Dynamic Weighing/Display Hold, eliminating the need for complicated manual calculations.

Built for durability, this tough industrial scale features a 304-stainless steel platform, powder-coated steel tubular frame & IP65 aluminum load cell. Indicator features ABS housing & bracket.

Flexible power solutions to suit your operation needs – Defender 2500 can be used with an AC power adapter, USB type-C or even disposable battery power.

Model	EC Type Approved	Capacity	Readability (d)	Readability (e)	Platform	Column	Item N°	Price €
i-D25P30B2R1-M	•	30kg		10g	300mm x 350mm	500mm	30844420	399,00
i-D25P60B2R1-M	•	60kg		20g	300mm x 350mm	500mm	30844421	399,00
i-D25P60B2L2-M	•	60kg		20g	400mm x 500mm	700mm	30844422	459,00
i-D25P150B2L2-M	•	150kg		50g	400mm x 500mm	700mm	30844423	459,00
i-D25P150B2X2-M	•	150kg		50g	500mm x 600mm	700mm	30844424	555,00
i-D25P300B2X2-M	•	300kg		100g	500mm x 600mm	700mm	30848746	555,00

Accessories	Item N°	Price €
Cable RS232 IBM 9P Male-to-Female	80500525	99,70
DB9 Adapter, Male-to-Male	30692384	6,00
SF40A Printer	30064202	715,00
SF-110B Printer	30708371	805,00
Remote Display AD7-MD	30472063	96,00



DB9 ADAPTER, MALE-TO-MALE
30692384



SF40A PRINTER
30064202



SF-110B PRINTER
30708371



REMOTE DISPLAY AD7-MD
30472063

Placing an order couldn't be easier

Online Ordering via our OHAUS Partner Shop

Our user-friendly Partner Shop makes ordering and tracking as simple and convenient as possible.

- ✓ Orders are received directly into our warehouse and production site for speed and efficiency
- ✓ Once your order is placed you will receive instant order acknowledgment and a confirmed delivery date within 24 hours so that you can provide your customers with accurate information

Ordering via fax, e-mail or phone to our Sales Support Team

- ✓ Once your order is placed you will receive a confirmed delivery date within 24 hours via e-mail.

Track your order easily

Online order tracking

- ✓ If you place your order via our Partner Shop you can also track it simply and efficiently online and you can also view your previous orders

Order tracking via our Sales Support Team

- ✓ Our Sales Support Team can tell you the status of your order via phone, fax or e-mail

Terms & Conditions

Our general terms and conditions of sale can be found within our Partner Shop. Please visit www.ohaus.com and select your country. Standard delivery is free for orders over € 400.

Below € 400 a handling fee of € 22 will be charged.

Contact

Technical Support Centre

e-mail: tsc@ohaus.com

Tel: 0041 22 567 53 20

Fax: 0041 22 567 53 21

Sales Support Centre

e-mail: ssc@ohaus.com

Tel: 0041 22 567 53 19

Fax: 0041 22 567 53 21

OHAUS Europe GmbH
Heuwinkelstrasse 3,
8606 Nänikon,
Switzerland

www.ohaus.com



www.ohaus.com