



- [Home](#)
- [About us](#)
- [News](#)
- [Products](#)
- [Service](#)
- [Contact](#)
- [Products](#) > [Liquid Handling](#) > BASIC



The benefits of the OPTIFIX[®] BASIC Dispenser at a glance

[\(Art.-No. 101.080\)](#)

- Best safety**
With OPTIFIX[®] BASIC, aqueous solutions, acids and crystallizing media can be dosed safely. As there are no metal springs used in the OPTIFIX dispenser, liquid samples will not be contaminated. By using a glass cylinder one can clearly see the liquid as it is being dosed. To avoid mistakes, the finger caps are labeled in different colours. The BASIC dispenser is in blue.
- High accuracy**
OPTIFIX[®] Dispensers are tested in our calibration department in accordance with ISO standard 8655. By rotating the fine-tuning wheel the volume of liquid can be fixed precisely. The wheel will not move during multiple dosings so the user can be sure that the desired volume is maintained.
- Premium quality**
By using high quality materials such as PTFE and borosilicate glass the OPTIFIX[®] BASIC Dispenser ensures outstanding performance for many years.
- Wide volume range**
The OPTIFIX[®] BASIC Dispenser is offered in nine sizes from 0.4 ml to 500 ml with a corresponding range of bottle top adaptors. This allows a wide range of volumes to be dispensed safely and accurately directly from storage bottles

For further details about the OPTIFIX BASIC[®] dispenser please click [here](#). If you would like to order the BASIC, please visit our [Webshop](#).