

LABORHITS

Pump drive PD 5201

Precise
Dosing and Dispensing



3 year warranty



Safety for continuous operation:
the motor will be switched off if
a high thermal load situation occurs
to **prevent accidents**

High precision even for flow rates
starting at 0.005 ml/min

Optional foot-pedal remote control
via cable to operate pump drive
outside a closed fume hood
facilitate operation

All models feature a smooth start operation which
prevents spills and splashing media. The speed will
ramp up slowly until your set rpm has been reached

For the use in cell biology just use a
pump head which features convex
rollers and **does not squeeze** the
tubing as conventional rollers

► **Leading Safety Standards**

► **Superior Ease of Use**

► **Reduced Cost of Ownership**

Type	PK	Cat. No.
Pump drive PD 5201	1	9.829 244