



For dispensing
solvents

FORTUNA® Liquid Handling

FORTUNA® OPTIFIX® SOLVENT

Due to our long experience in the production of manual bottle top dispensers, our products are easy to handle, safe to use and precise. They are indispensable for daily use in the laboratory.

New manufacturing techniques guarantee an optimum price-performance ratio.

If you do work with solvents, our OPTIFIX® SOLVENT is the first choice for you.

OPTIFIX® SOLVENT Features

- **for dispensing solvents** protection against breakage
- **FI - with certificate, individually numbered**
- **manufactured to DIN EN ISO 9001**
- **using high quality materials**
- **with Dispensing Piston made of glass**
- **Dispensing Cylinder, made of glass, visible dispensing process, cylinder is coated with plastic-**
- **Finger Protection Cap with Volume Setting System**, the ergonomic finger protection cap is made of PP, graduation on the scale is lasered, for right or left-handed use, fast and easy volume setting with additional fine adjustment, easy to read scale
- **Valve Block**, high technology leads to optimal function and 100 % tightness,



no metal springs in the valves - no corrosion



- **FIX Adaptor** made of white PP up to 10 ml size, made of solid PP from 30 ml

Bottle Top Dispenser FORTUNA® OPTIFIX® SOLVENT, with GL 32 thread, comes with 2 adaptors 40 + 45, 1 PTFE intake tube, single pack

Volume	grad.	Accuracy	Precision	Art.-No.
0,4 - 2 ml	0,1 ml	≤ ± 0,6 %	≤ ± 0,2 %	101.081-27
1 - 5 ml	0,1 ml	≤ ± 0,6 %	≤ ± 0,2 %	101.081-33
2 - 10 ml	0,2 ml	≤ ± 0,6 %	≤ ± 0,2 %	101.081-38
6 - 30 ml	0,5 ml	≤ ± 0,6 %	≤ ± 0,2 %	101.081-44
10 - 50 ml	1,0 ml	≤ ± 0,6 %	≤ ± 0,2 %	101.081-48
20 - 100 ml	2,0 ml	≤ ± 0,6 %	≤ ± 0,2 %	101.081-51
40 - 200 ml	5,0 ml	≤ ± 0,6 %	≤ ± 0,2 %	101.081-53
60 - 300 ml*	5,0 ml	≤ ± 0,6 %	≤ ± 0,2 %	101.081-56
100 - 500 ml*	5,0 ml	≤ ± 0,6 %	≤ ± 0,2 %	101.081-58

additional to ISO 8655)


Poulten & Graf Ltd.




FORTUNA®

Peak Works, 1 Alfreds Way
Barking, Essex, IG11 0AS, UK
Tel.: +44 (0)20 8594 4256
Fax: +44 (0)20 8594 8419
<http://www.poulten-graf.com>
e-mail: VOLAC@poulten-graf.com

Poulten & Graf GmbH



FORTUNA®



P.O. Box 1352
97863 Wertheim / Germany
Tel.: +49 9342 9229-0
Fax: +49 9342 9229-80
<http://www.poulten-graf.com>
e-mail: FORTUNA@poulten-graf.com