

## Scraper LaboPlast®

The LaboPlast® scraper, made of polystyrene, is produced in our class 7 clean room (10,000) and individually packaged for disposable use.

Sturdy and convenient scraper for disposable use. The scraper is suitable for scraping off, smoothing or cleaning as well as removing residue.

- ▶ PS, white
- ▶ Individually packaged
- ▶ Production and packaging according to clean room class 7
- ▶ With EU foodstuffs and FDA approval



### Scraper LaboPlast®

| L x W mm | Pieces per pack | Item no.  |    |
|----------|-----------------|-----------|----|
| 200 x 80 | 10              | 5378-0030 | N4 |



## Scraper SteriPlast®

The SteriPlast® scraper, made of polystyrene, is produced in our class 7 clean room (10,000), individually packaged for disposable use and then sterilised by gamma rays.

Sturdy and convenient scraper for disposable use. The scraper is suitable for scraping off, smoothing or cleaning as well as removing residue.

- ▶ PS, white
- ▶ Individually packaged and sterilised by gamma rays
- ▶ Production and packaging according to clean room class 7
- ▶ With EU foodstuffs and FDA approval



### Scraper SteriPlast®

| L x W mm | Pieces per pack | Item no.  |    |
|----------|-----------------|-----------|----|
| 200 x 80 | 10              | 5378-1030 | N4 |

