

Sliding Plate with Combifix SM B

[No product image available](#)

Description

Sliding plates with Combifix SM B are identical to Rack systems Combifix SM however with an additional base plate. ([see Rack system Combifix SM B](#))

For mounting sliding plates on a 2-storey top frame SM a mounting set (Order No. 0052 074) is necessary. With 1 mounting set 2 sliding plates can be mounted.

Scope of Delivery:

Base plate: 560 x 400 mm

Basic Rack with 4 clamping strips h and 4 clamping strips v

Details

Order No.: 0051 485

Accessory: Mounting set

Clamping Strip h 0050 400

Clamping strip v 0050 399

Cramp strip 0050 401

Spring strip 0050 402

Test tube rack with 44 drillings à 14 mm Ø

Test tube rack with 44 drillings à 16 mm Ø

Test tube rack with 44 drillings à 18 mm Ø

Test tube rack with 14 drillings à 30 mm Ø