

### Key Features

- Very efficient heating saving time and energy
- Chemically resistant ceramic top
- “Hot” warning light for user safety



CR300

## CR300

### Hotplate, Infra Red, Ceramic

Using a very efficient infra red heater of just 900W power this hotplate will boil 1 litre of water over 30% faster than a conventional ceramic hotplate of 1200W. Ideal for heating large volumes of liquid.

### Technical Specification

Plate dimensions, mm	300 x 300
Heated area, mm	140 diameter
Heater power, W	900
Dimensions, mm (w x d x h)	320 x 370 x 105
Net weight, kg	4
Electricity supply	230V, 50-60Hz, 900W
IP Rating	31

### Ordering Information

Model	Description
CR300	Hotplate, infra red