

KERN 326-02

KERN

Individual weight, OIML Class F1, 2 g, Eco Shape, stainless steel polished



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	326

Measuring System

Test weight - Density uncertainty (kg/m ³)	140
Weight OIML class	F1
Test weight - nominal value	2 g
Permissible tolerance according to OIML (+/-)	0,12 mg
Test weight - Density (kg/m ³)	7950
OIML-compliant	✓
Test weight - Density [Min] (OIML)	4000

Construction

Shape weight	Eco Shape
Dimensions (Ø×H)	6×10 mm
Material	stainless steel polished
Adjustment method	Polishing

Services (optional)

Article number for DAkkS calibration	962-432
Article number for guaranteed plus-deviation	921-432

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	6×6×10 mm

Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

Optional cases

Case - Plastic article number	347-030-400
Case - Aluminium article number	317-020-600
Case - Wood article number	317-020-100

Product information

GTIN/EAN number	4045761223738
Trade number	84239090
REACH classification	no prohibited substances

Pictograms

STANDARD



OPTION

