

# KERN 326-07

**KERN**

Individual weight, OIML Class F1, 100 g, Eco Shape, stainless steel polished



|                              |                |
|------------------------------|----------------|
| Dimensions packaging (W×D×H) | 22×22×38 mm    |
| Shipping method              | Parcel service |
| Net weight approx.           | 0,10 kg        |
| Gross weight approx.         | 0,10 kg        |
| Shipping weight              | 0,1 kg         |

## Optional cases

|                                 |             |
|---------------------------------|-------------|
| Case - Plastic article number   | 347-070-400 |
| Case - Aluminium article number | 317-070-600 |
| Case - Wood article number      | 317-070-100 |

## Product information

|                      |                          |
|----------------------|--------------------------|
| GTIN/EAN number      | 4045761223684            |
| Trade number         | 84239090                 |
| REACH classification | no prohibited substances |

## Pictograms

### STANDARD



### OPTION



## Category

|                  |                   |
|------------------|-------------------|
| Brand            | KERN              |
| Product category | Individual weight |
| Product group    | OIML F1           |
| Product family   | 326               |

## Measuring System

|                                                        |        |
|--------------------------------------------------------|--------|
| Test weight - Density uncertainty (kg/m <sup>3</sup> ) | 140    |
| Weight OIML class                                      | F1     |
| Test weight - nominal value                            | 100 g  |
| Permissible tolerance according to OIML (+/-)          | 0,5 mg |
| Test weight - Density (kg/m <sup>3</sup> )             | 7950   |
| OIML-compliant                                         | ✓      |
| Test weight - Density [Min] (OIML)                     | 7390   |
| Test weight - Density [Max] (OIML)                     | 8730   |

## Construction

|                   |                          |
|-------------------|--------------------------|
| Shape weight      | Eco Shape                |
| Dimensions (Ø×H)  | 22×38 mm                 |
| Material          | stainless steel polished |
| Adjustment method | Adjusting cavity - Seal  |

## Services (optional)

|                                              |         |
|----------------------------------------------|---------|
| Article number for DAKKS calibration         | 962-437 |
| Article number for guaranteed plus-deviation | 921-437 |

## Packing & Shipping

|               |     |
|---------------|-----|
| Delivery time | 1 d |
|---------------|-----|