

KERN PWS-A01

Battery insert, enables mobile weighing

KERN



Pictograms

STANDARD



Category

Brand	KERN
Product category	Balance
Product group	Battery Insert

Construction

Dimension (W×D×H)	85×65×27 mm
Material	A PET

Power Supply

Supplied power supply	no battery
Battery	4×1.5 V AA
Battery connection	Pad + Spring
Battery operating time - backlight on	32 h
Battery operating time - backlight off	200 h
Battery operating time	200 h

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	200×270×30 mm
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,10 kg
Shipping weight	0,3 kg

Product information

GTIN/EAN number	4045761311282
Trade number	85073020
REACH classification	no prohibited substances