

Data Sheet

Page 1/1



25M X 0.32MM ID - BP1 0.5UM GC Column

Part Number: TJ-054065
Manufacturer: Trajan Scientific
Quantity: 1 pc

Description

25M X 0.32MM ID - BP1 0.5UM GC-Säule

Technical Information

Brand:	BP	Base Material:	Dimethylpolysiloxane
Product Type:	GC Column	Temperature Minimum [°C]:	-60
Length:	25 m	Temperature Maximum [°C]:	340
Internal Diameter:	0.32 mm	Mode of Separation:	Gas Chromatography
Film Thickness [µm]:	0.5 µm		
Chemistry:	1 (1HT, 1MS)		
Chemistry Producers Term:	BP1		
USP Class:	G1, G2, G9, G38		
Country of Origin:	Australia		